Directions for Use..... 1
Mode d'emploi..... 6
For Canada Only

90835177-01

2017-06
< Canada >

Neuroform EZ[™]
Stent System

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Stryker Neurovascular representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

Resale of this device is prohibited by US law.

DEVICE DESCRIPTION

The Neuroform EZ Stent System includes:

- A self-expanding, nitinol stent with four radiopaque markerbands on each end (distal and proximal) and 3 interconnects between the central stent segments, designed to provide support of the coil mass within the aneurysm and minimize stent deflection.
- A stent delivery wire. The stent is pre-loaded on the stent delivery wire and protected by an introducer sheath.
- An accessory pouch containing an optional guidewire introducer and an optional torque device. The physician may attach the torque device to the proximal end of the stent delivery wire, which may facilitate handling and stabilization. The stent delivery wire is not designed to be torqued. The guidewire introducer should not be used with the stent delivery wire. If the microcatheter must be repositioned with a guidewire (as described in procedure step 18), the guidewire introducer may be used to facilitate introduction into the microcatheter.

Contents

- One (1) Neuroform EZ Stent System
- One (1) Torque Device
- One (1) Guidewire Introducer


Figure 1. Neuroform EZ Stent System

Table 1. Sizing Table

Labeled Stent Diameter (mm)	Self Expanded Stent Diameter (mm)	Recommended Vessel Diameter ¹ (mm)
2.5	3.0	> 2.0 and ≤ 2.5
3.0	3.5	> 2.5 and ≤ 3.0
3.5	4.0	> 3.0 and ≤ 3.5
4.0	4.5	> 3.5 and ≤ 4.0
4.5	5.0	> 4.0 and ≤ 4.5

¹Select a stent diameter based on the sizing recommendations in Table 1 and based on the larger vessel diameter (proximal or distal reference vessel diameter).

INTENDED USE/INDICATIONS FOR USE

The Neuroform EZ Stent System is intended for use with embolic coils for the treatment of wide neck, intracranial, saccular aneurysms arising from a parent vessel with a diameter of ≥ 2 mm and ≤ 4.5 mm that are not amenable to treatment with surgical clipping. Wide neck aneurysms are defined as having a neck ≥ 4 mm or a dome-to-neck ratio of < 2.

CONTRAINDICATION

- Patients in whom antiplatelet and/or anticoagulation therapy is contraindicated.
- Patients for whom angiography and other examinations show that access to the target cerebral aneurysm by the stent may be impossible because of anatomical features, or severe vascular spasm that does not respond to pharmacological therapy.

WARNINGS

- This device should only be used by physicians who have received appropriate training in interventional neuroradiology or interventional radiology and preclinical training on the use of this device as established by Stryker Neurovascular.
- Select a stent size (length and diameter) to maintain a minimum of 4 mm on each side of the aneurysm neck along the parent vessel. An incorrectly sized stent may result in damage to the vessel or stent migration. Therefore, the stent is not designed to treat an aneurysm with a neck greater than 22 mm in length.
- If excessive resistance is encountered during the use of the Neuroform EZ Stent System or any of its components at any time during the procedure, discontinue use of the stent system. Continuing to move the stent system against resistance may result in damage to the vessel or a system component.
- The Neuroform EZ Stent System should not be used for the treatment of stenosis due to intracranial atherosclerotic disease.
- Persons allergic to nickel titanium (Nitinol) or platinum may suffer an allergic response to this stent implant.

PRECAUTIONS

- The Neuroform EZ Stent System is provided STERILE for single use only.
- Use the Neuroform EZ Stent System prior to the "Use By" date printed on the package.
- Carefully inspect the sterile package and Neuroform EZ Stent System prior to use to verify that neither has been damaged during shipment. Do not use kinked or damaged components.
- Please refer to "Magnetic Resonance Imaging (MRI) Compatibility" section for MRI Information.
- The Neuroform EZ Stent System should not be used for recapturing the stent.
- Exercise caution when crossing the deployed stent with adjunct devices.
- After deployment, the stent may foreshorten up to 1.8% in 2.5 mm stents and up to 5.4% in 4.5 mm stents.
- The safety of the Neuroform EZ Stent System in patients below the age of 18 has not been established.
- In cases where multiple aneurysms are to be treated, start at the most distal aneurysm first.
- The safety of "Y" stenting or techniques of passing a guidewire through stent interstices to access other vessels for the purpose of stenting has not been clinically established.

ADVERSE EVENTS

Potential Adverse Events

Potential complications include, but are not limited to:

- Allergic reaction including, but not limited to, contrast, Nitinol, platinum and medications
- Aneurysm perforation/rupture
- Coil herniation through stent into parent vessel
- Death

- Embolus
- Hemorrhage
- In-stent stenosis
- Infection
- Ischemia
- Neurological/intracranial sequelae
- Pseudoaneurysm
- Stent fracture
- Stent migration/embolization
- Stent misplacement
- Stent thrombosis
- Stroke
- Transient ischemic attack
- Vasospasm
- Vessel occlusion or closure
- Vessel thrombosis
- Vessel perforation/rupture, dissection, trauma or damage
- Other procedural complications including but not limited to anesthetic and contrast media risks, hypotension, hypertension, access site complications

Refer to the appropriate embolic coil DFU for other complications that may occur due to coil embolization.

Observed Adverse Events from Clinical Studies

Tables 2 and 3 identify the adverse events observed in the clinical study conducted with the Neuroform® Microdelivery Stent System.

Twenty-nine patients were implanted with the stent. The tables include all adverse events through 6 months. Of the 29 patients implanted with the stent, 17 patients had 1 or more adverse events and 5 had 1 or more serious adverse events. There were 12 serious adverse events and 21 other adverse events, all of which occurred prior to or by the time of discharge. None occurred between discharge and the 6-month timepoint. Nine patients had 1 adverse event, 4 patients had 2 adverse events, 1 patient had 3 adverse events, 2 patients had 4 adverse events, and 1 patient had 5 adverse events.

Table 2 summarizes the patient rates for observed serious adverse events. Table 3 summarizes the patient rates for all other observed adverse events.

Table 2. Serious Device or Procedure-Related Adverse Events

Serious Adverse Event ¹	n (%)
Death ²	1 (3.4%)
Aneurysm Perforation ^{2,3}	2 (6.9%)
Arterial Perforation ⁴	1 (3.4%)
Subarachnoid/Intraventricular Hemorrhage ^{2,3}	2 (6.9%)
Thromboembolic Stroke ⁴	1 (3.4%)
Intracerebral Hematoma ⁴	1 (3.4%)
Left Hemiparesis ⁴	1 (3.4%)
Intraparenchymal Bleeding ³	1 (3.4%)
Retroperitoneal Hematoma ⁵	1 (3.4%)
Confusion ⁶	1 (3.4%)

¹Five patients had these 12 serious adverse events. The "n" reflects the number of occurrences of that adverse event. The % is based on 29 patients who were assessed before or at discharge when all adverse events occurred.

²One patient had 3 serious adverse events. There was perforation of the aneurysm dome with the micro guidewire during the initial catheterization of the aneurysm resulting in subarachnoid/interventricular hemorrhage and death. Death was due to complications from aneurysm perforation leading to bleeding and pre-existing hepatitis and management of anticoagulation therapy.

³One patient had 3 serious adverse events. There was perforation of the aneurysm with the microcatheter during coil placement resulting in subarachnoid hemorrhage and subsequent intraparenchymal bleeding (from the ventricular drainage line).

⁴One patient had 4 serious adverse events. Arterial perforation occurred with the tip of the exchange length guidewire prior to stent insertion, resulting in an intracerebral hematoma. This patient also had a thromboembolic stroke that led to left hemiparesis.

⁵One patient had a retroperitoneal hematoma.

⁶One patient had confusion. Confusion was categorized by the protocol as a non-serious adverse event; however, it was determined by the clinical study investigator to be a serious adverse event because the patient required a prolonged hospital stay.

Table 3. Other Device or Procedure-Related Adverse Events

Other Adverse Event ¹	n (%)
Right Hemiparesis	1 (3.4%)
Embolic Event ²	4 (13.8%)
Vasospasm ³	5 (17.2%)
Intimal Dissection ⁴	1 (3.4%)
Seizure ⁵	1 (3.4%)
Access Site Hematoma ⁶	2 (6.9%)
Liver Failure	1 (3.4%)
Vomiting	1 (3.4%)
Headache	3 (10.3%)
Fever of Unknown Origin	1 (3.4%)
Urinary Tract Infection	1 (3.4%)

¹Fifteen patients had these 21 adverse events. The "n" reflects the number of occurrences of that adverse event. The % is based on 29 patients who are accounted for and were assessed before or at discharge when all adverse events occurred.

²Includes embolic ischemic lesion, small embolic lesion, asymptomatic microemboli to brain detected by MRI, and left prolonged reversible ischemic neurological deficit (PRIND). All embolic events resulted in mild neurological deficits. Three completely resolved, and 1 patient was discharged to a rehabilitation facility.

³Includes 4 mild and 1 moderate case. All completely resolved.

⁴Occurred during placement of the guide catheter in the cervical internal carotid prior to stent placement, not in the portion of the vessel treated with the device.

⁵One patient with a history of epilepsy experienced a seizure with no permanent sequelae while in the hospital.

⁶Includes 1 mild and 1 moderate case. Both resolved.

Tables 4 and 5 summarize the adverse events observed in the clinical study conducted with the Neuroform EZ™ Stent System.

Twenty patients were implanted with 23 stents. The tables include adverse events through 6 months. In total 71 adverse events occurred in 19 patients. Of these, 3 patients experienced a total of 3 serious adverse events (SAEs). Two SAEs were unrelated to the device or procedure including one case of drug hypersensitivity and one case of urinary tract infection. One SAE of retroperitoneal hematoma was related to the procedure. Two non-serious events were classified as having a relationship to the device including one case of retinal artery occlusion reported the day of the procedure and one case of in-stent stenosis reported at 6 month follow-up. Table 4 summarizes all observed adverse events. Table 4 summarizes adverse events classified as having a relationship to the device or procedure. Table 5 summarizes all observed adverse events.

Table 4. Adverse Events Related to Device or Procedure

Related Adverse Events	n (%)
Total	17 (85.0%)
Headache	11 (55.0%)
Nausea	4 (20.0%)
Asymptomatic cerebral infarction	7 (35.0%)
Alopecia	3 (15.0%)
Bleeding from puncture site	2 (10.0%)
Subcutaneous hematoma	2 (10.0%)
3rd cranial nerve palsy	1 (5.0%)
Cerebral artery in-stent stenosis	1 (5.0%)
Diplopia	1 (5.0%)
Hemianopsia	1 (5.0%)
Hypesthesia	1 (5.0%)
Memory impairment	1 (5.0%)
Retinal artery occlusion	1 (5.0%)
Retro-peritoneal hematoma	1 (5.0%)

Table 5. Adverse Events

Adverse Events	n (%)
Total	19 (95.0%)
Headache	16 (80.0%)
Nausea	9 (45.0%)
Asymptomatic cerebral infarction	7 (35.0%)
Alopecia	3 (15.0%)
Hepatic function abnormalities	3 (15.0%)
Bleeding from puncture site	2 (10.0%)
Gastrointestinal hemorrhage	2 (10.0%)
Subcutaneous hematoma	2 (10.0%)
Subcutaneous hemorrhage	2 (10.0%)
3rd cranial nerve palsy	1 (5.0%)
Back pain	1 (5.0%)
Blood glucose increased	1 (5.0%)
Blood pressure decreased	1 (5.0%)
Blurred vision	1 (5.0%)
Cerebral artery in-stent stenosis	1 (5.0%)
Chorioretinal atrophy	1 (5.0%)
Cystitis	1 (5.0%)
Decreased urinary volume	1 (5.0%)
Diplopia	1 (5.0%)
Dizziness	1 (5.0%)
Drug hypersensitivity	1* (5.0%)
Eye pain	1 (5.0%)
Gastrointestinal problems	1 (5.0%)
Hemianopsia	1 (5.0%)
Hypesthesia	1 (5.0%)
Hyponatremia	1 (5.0%)
Memory impairment	1 (5.0%)
Optic nerve disorders	1 (5.0%)
Retinal artery occlusion	1 (5.0%)
Retinal hemorrhage	1 (5.0%)
Retro-peritoneal hematoma	1* (5.0%)
Upper abdominal pain	1 (5.0%)
Urinary tract infection	1* (5.0%)

*Classified as serious adverse events

CLINICAL STUDIES

A Clinical study was conducted in Europe to evaluate the Neuroform Microdelivery Stent System. The patient inclusion criteria were: (1) wide neck, ruptured or unruptured, saccular, intracranial aneurysm or aneurysm on the level of the skull base, where a wide neck is defined as a dome-to-neck ratio < 2 and/or neck length of ≥ 4 mm; (2) aneurysm is in artery with diameter ≥ 1.5 mm and ≤ 5.5 mm; (3) patient is ≥ 18 years old; and (4) patient provided written informed consent.

There were 31 patients entered into the study. Five (16%) were male and 26 (84%) were female. Fifty-two percent of the patients were asymptomatic prior to treatment. Two of the 31 patients did not receive the stent because of failure to access based on anatomy. The remaining 29 patients enrolled in the study had 30 aneurysms (1 patient had 2 aneurysms that were treated with one stent). Previous attempts had been made to treat 17 of the 30 aneurysms (57%) using other devices.

Table 6 summarizes the locations of the 30 aneurysms. Table 7 summarizes the sizes of the 30 aneurysms.

Table 6. Aneurysm Location

Location	n	%
Carotid ophthalmic	7	24%
Posterior communicating artery	7	24%
Carotid cavernous	5	17%
Anterior choroidal	2	7%
Basilar tip	2	7%
Carotid bifurcation	1	3%
Middle cerebral artery	1	3%
Anterior cerebral artery	1	3%
Vertebral artery	1	3%
Posterior inferior cerebellar artery	1	3%
Basilar trunk	1	3%
Other	1	3%

Table 7. Aneurysm Size

Measurement	n	Mean	SD	Min	Max
Dome width (mm)	30	7.4	4.3	2.1	20.0
Neck length (mm)	30	4.9	1.8	2.1	11.0
Dome to neck ratio	30	1.5	0.5	0.8	2.7
Parent vessel pre-aneurysm (mm)	30	3.6	0.6	2.4	4.8
Parent vessel post-aneurysm (mm)	30	3.2	0.7	1.7	4.4
Parent vessel caliber differential (mm)	30	1.0	1.0	0.3	1.7

The 29 patients were implanted with 39 stents to treat their 30 aneurysms. Twenty (69%) patients had 1 stent, 8 (28%) patients had 2 stents, and 1 (3%) patient had 3 stents. The stents implanted ranged from 3.5 mm to 4.5 mm. One patient required a secondary endovascular procedure to place a second stent in the correct location because the original stent was inadvertently not deployed at aneurysm site; this counts for 2 of the stents. One patient had the original stent successfully deployed but was removed during the embolic coiling procedure when the clinical study investigator attempted to snare the errant coil loop and dislodged the stent. A replacement stent was implanted in its place, and this counts for 2 of the stents. For 7 patients, multiple stents were used to treat one aneurysm in cases where (1) the embolic coiling procedure left the tail of an embolic coil in the vessel or (2) the neck of aneurysm was estimated at an incorrect width and a second or third stent was necessary to cover the neck of the aneurysm.

With regard to patient accounting, 31 patients were originally entered into the study; however, 2 did not receive the stent. One patient died immediately after the procedure. There are adverse event data on 29 patients, including the one death. Therefore, there were 28 patients of 31 who were expected for evaluation through 6 months. At discharge, 28 of the expected 28 were evaluated for a follow-up rate of 100%. At 6 months, 26 of 28 patients were evaluated for a follow-up rate of 93%.

The study endpoints were (1) adverse events, (2) technical feasibility, and (3) clinical outcome. The incidence of all adverse events, device or procedure-related, were assessed. Technical feasibility was assessed by the ability to access the aneurysm and place the stent accurately across the aneurysm neck. Clinical outcome was assessed by percent angiographic aneurysm occlusion.

Adverse events were presented in Tables 2 and 3 in the **Adverse Events** section above.

Table 8 below summarizes the patient rates with regard to technical feasibility.

Table 8. Technical Feasibility

Technical Feasibility	n (%)
Ability to access aneurysm	29/31 (93.5%) patients ¹
Ability to place stent across aneurysm neck	29/29 (100%) patients ^{2,3}

¹Two patients could not be accessed based on anatomy.
²One patient required a secondary endovascular procedure to place a second stent in the correct location because the original stent was inadvertently not deployed at aneurysm site.
³There were 2 intraoperative device malfunctions involving the markerband of the 2F Stabilizer Catheter inadvertently detaching from the shaft of the 2F Stabilizer Catheter after stent deployment. In one patient, the 2F Stabilizer Catheter was inside the patient at the time of the device malfunction, and the separated markerband embolized in a small, distal intracranial artery. This patient had no adverse events from this event. In the second patient, the 2F Stabilizer Catheter was outside the patient at the time of the device malfunction. Stryker Neurovascular has since increased its markerband bond strength.

Table 9 below summarizes the patient rates with regard to clinical outcome.

Table 9. Clinical Outcome

Clinical Outcome ¹	n (%)
% occlusion at discharge ²	
100%	17 (58.6%)
95-99%	13 (44.8%)
% occlusion at 6 months	
100%	18 (69.2%)
95-99%	8 (30.8%)

¹The "n" reflects the number of occurrences. The % is based on 29 patients at discharge and 26 patients at 6 months.
²One patient had 2 aneurysms, each with different resulting % occlusion. Therefore, this patient is reported twice.

Other clinical outcomes included:

- No stent stenosis or migration.
- No emboli coil migration.
- No parent vessel thrombosis, occlusion, or dissection.
- Neurological status: Of 26 patients evaluated at 6 months, 17 (65%) had an unchanged (normal) neurological assessment as compared to baseline, 3 (16%) had an improved (from abnormal to normal) neurological assessment as compared to baseline, 5 (19%) had an unchanged (abnormal) neurological assessment as compared to baseline, and 1 (4%) had a worsened (abnormal moderate confusion to abnormal severe confusion) neurological assessment as compared to baseline.

A clinical study was conducted in Japan to evaluate the safety and effectiveness of the Neuroform EZ™ Stent System. The patient inclusion criteria were: (1) wide neck, unruptured intracranial aneurysm 10 mm or larger at the largest point that is not amenable to surgery or coil embolization alone, where a wide neck is defined as a neck diameter of ≥ 4 mm or a dome-to-neck ratio < 2 ; (2) aneurysm is in an artery with a diameter of 2.0 to 4.5 mm; (3) patient is between 20 and 80 years of age at the time of consent.

The primary endpoint was procedure success defined as meeting all of the following conditions: (1) Ability to place the Neuroform® stent at the neck of the cerebral aneurysm; (2) Ability to complete the aneurysm coil embolization and retain the coil mass within the aneurysm at end of procedure; (3) No Major Adverse Cerebrovascular Event (MACVE) related to the device or procedure within 30 days after the procedure.

Patients were assessed clinically at Baseline and at 7 days, 30 days, and 6 months after the procedure. Angiographic evaluation was conducted at Baseline, Procedure and 6 months after the procedure. An independent Data Monitor evaluated safety data. An independent review of angiographic imaging was conducted by an Image Review Board to assess the status of stent placement and coil mass retention and Independent review of clinical events was conducted to determine whether they met the definition of MACVE and to assess the relationship between the event and the device or procedure.

A total of 20 subjects with 20 aneurysms were enrolled at two sites. The average age was 59.3±10.6 years and 14 of 20 (70%) were female. All patients underwent stent placement and coil embolization. Table 10 summarizes the locations of the 20 aneurysms. Table 11 summarizes sizes of the 20 aneurysms.

Table 10. Aneurysm Location

Location	n	%
Paraclinoid segment of the carotid artery	11	55.0%
Carotid cavernous sinus	3	15.0%
Vertebral artery	3	15.0%
Posterior communicating artery	1	5.0%
Basilar trunk	1	5.0%
Petrous segment of the carotid artery	1	5.0%

Table 11. Aneurysm Size

Measurement	n	Mean	SD	Min	Med	Max
Parent artery diameter: proximal (mm)	20	3.73	0.66	2.5	3.70	4.9
Parent artery diameter: distal (mm)	20	3.31	0.54	2.1	3.35	4.2
Max dome diameter (mm)	20	12.24	3.34	8.0	11.05	20
Neck diameter (mm)	20	6.92	3.49	3.2	6.40	20
Dome-to-neck ratio	20	1.94	0.60	1.0	1.80	3.3

All 20 subjects underwent stent placement and coil embolization. Three subjects required placement of 2 stents resulting in a total of 23 stents implanted. In two cases a second stent was placed due to improper location of the first stent and in one case the tip of the coil protruded into the parent artery. In all cases satisfactory neck coverage was achieved at the end of the procedure and these procedures were determined to meet criteria for procedure success by the Image Review Board. Coil embolization was performed in all cases and the coil mass was retained within the aneurysms at the end of the procedure. No incidence of MACVE was observed in any case during the study. Thus, all subjects were rated as a "procedure success" and met the primary endpoint. Data is available for 20 patients at 30 days and 19 patients at 6 months (one patient withdrew consent). Adverse events are presented in Tables 4 and 5. A summary of procedure outcome is provided in Table 12.

Table 12. Summary of Procedure Outcome

Procedure Outcome	30 days		6 months	
	n	%	n	%
Procedure Success ¹	20	100.0%	-	-
Successful stent placement ²	20	100.0%	19 ³	95%
Successful coil embolization ²	20	100.0%	-	-
Successful retention of coil mass ²	20	100.0%	19 ³	95%

¹95% confidence interval: 83.2% - 100.0%
²Data based on Image Review Board assessment.
³One patient withdrew consent after 30 days

Cranial nerve and general neurological status evaluations were conducted and included assessment of NIH Stroke Scale (NIHSS), Barthel and modified Rankin (mRS) scores. No significant changes in neurological or cranial nerve function were reported. Table 13 summarizes neurological status evaluations.

Table 13. Summary of Neurological Status

Neurological Assessment	Status at 6 months compared to Baseline ¹	
Modified Rankin Scale	No change	89.5% (17/19)
	Worse	5.3% (1/19)
	Improved	5.3% (1/19)
Barthel	No change	100% (19/19)
	Worse	0%
	Improved	0%
NIHSS	No change	94.7% (18/19)
	Worse	5.3% (1/19) ²
	Improved	0%

¹ One patient discontinued (withdrew consent) and is not included.

² In one subject with partial hemianopsia, the visual field deteriorated from a score of 0 before the procedure to a score of 1, 7 days, 30 days and 6 months after the procedure.

PATIENT INFORMATION

You should have already provided the patient with a copy of the Patient Information Booklet so that (s)he has had adequate time to review the information and ask any questions.

Immediately after the procedure, complete the Patient Information Card, which is included in the Patient Information Booklet, and provide the card to the patient before the patient leaves the hospital. The Patient Information Card includes important information about the stent that was used and includes a statement regarding MRI information.

CONCOMITANT MEDICAL THERAPY

Typical antiplatelet and anticoagulation regimen used for interventional intracranial procedure is recommended at the discretion of the treating physician. Do not use the Neuroform EZ™ Stent System in patients in whom antiplatelet and/or anticoagulation therapy is contraindicated.

MAGNETIC RESONANCE IMAGING (MRI) COMPATIBILITY

Magnetic Resonance Conditional

Non-clinical testing and analysis have demonstrated that Neuroform® Stent is MR Conditional alone, or when overlapped with a second stent, and adjacent to a Stryker Neurovascular coil mass. A patient with Neuroform Stent can be safely scanned immediately after placement of this implant, under the following conditions:

- static magnetic field of 1.5 and 3.0 Tesla
- spatial gradient field up to 2500 Gauss/cm (25 Tesla/m)
- normal operational mode for gradients and SAR (maximum whole body averaged specific absorption rate (SAR) of lower than 2.0 W/kg and maximum head SAR of lower than 3.2 W/kg) for a total active MR scan time (with RF exposure) of 15 minutes or less per scan sequence.

In an analysis based on the temperature rises in non-clinical testing of stents and the calculated SAR in the patient during an MR scan, Neuroform Stents were determined to produce an in-vivo temperature rise of 4°C or lower for 15 minutes of MR scanning in normal operational mode in 1.5 T and 3 T MR systems. The Neuroform Stent should not migrate in this MRI environment.

MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the device. In Spin Echo and Gradient Echo sequence evaluations Neuroform stent image artifact extended approximately 2 mm from the device. Lumen of the stent was partially obscured by the artifact. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this implant.

HOW SUPPLIED

Minimum Materials Required (Not included in the Neuroform EZ Stent System package)

Quantity	Material
1	Appropriate guiding catheter(s) [0.064 in (1.63 mm) minimum ID]
3	Sterile Heparinized normal saline (HepNS) flush lines 1000 U/500 mL (500 cc)
1	0.014 in (0.36 mm) guidewire access length
1	0.014 in (0.36 mm) guidewire exchange length
1	Sterile Contrast
1	Arterial sheath and dilator set
1	20 mL (20 cc) syringe with luer-lock

- 3 Three way stopcock
 1 0.027 in (0.69 mm) ID Microcatheter, 135 or 150 cm length

Do not use if package is opened or damaged.
 Do not use if labeling is incomplete or illegible.

Handling and Storage

Store in a cool, dry, dark place.

OPERATIONAL INSTRUCTIONS

Initial Access, Angiographic Assessment and Stent Selection

1. Gain vascular access according to standard practice. Select a Boston Scientific or Stryker Neurovascular microcatheter (0.027 in / 0.69 mm ID and 135 or 150 cm length) with a neurovascular indication. Establish and maintain continuous flush with sterile heparinized saline through the microcatheter per standard vascular practice. Using angiography, determine the location of the aneurysm and the size of the aneurysm neck.
2. Navigate the microcatheter over an access length guidewire at least 1.2 cm distal to the aneurysm neck.

Note: The microcatheter tip must be placed sufficiently distal to the aneurysm neck to allow for slack to be removed from the system after the stent is advanced, while maintaining adequate stent length (approximately 4 mm) distal to the aneurysm neck. Excessive tortuosity may necessitate microcatheter tip placement more than 1.2 cm distal to the aneurysm neck.

3. Remove the guidewire.
4. Select an appropriate Neuroform EZ Stent System based on the sizing recommendations in Table 1 and based on the larger reference vessel diameter. Select a stent that is at least 8 mm longer than the aneurysm neck to maintain a minimum of 4 mm on each side of the aneurysm neck along the parent vessel.

Delivery System Preparation and Stent Transfer

5. Carefully inspect the stent system packaging for damage. Do not use if damage is found; call your Stryker Neurovascular representative.
6. Peel open the pouch using aseptic technique.
7. Carefully place the dispenser hoop into the sterile field.
8. Release the stent delivery wire from the clip on the dispenser hoop. Hold the stent delivery wire and proximal end of the introducer sheath together, then carefully remove the system from the dispenser hoop.

Note: The stent delivery wire and proximal end of the introducer sheath must be held together when removing the Neuroform EZ Stent System from the dispenser hoop to prevent stent movement and premature deployment.

Note: Confirm that the stent delivery wire does not move relative to the introducer sheath during removal of the stent system from the dispenser hoop.

9. Inspect the stent system. Confirm that the tip of the stent delivery wire is entirely within the introducer sheath. Confirm that the stent delivery wire is not kinked and that the introducer sheath tip is not damaged. Do not use if damage is found; call your Stryker Neurovascular representative.
10. Partially insert the distal end of the introducer sheath into the RHV (3.5F compatible RHV) connected to the microcatheter. Tighten the RHV. Open the y-connector valve of the RHV that is connected to the sterile heparinized saline flush and verify that fluid exits the proximal end of the introducer sheath.

Warning: Purge the system carefully to avoid the accidental introduction of air into the stent system.

Note: Partial insertion of the introducer sheath into the RHV is necessary to ensure a flow path for flush.

11. Loosen the RHV. Advance the introducer sheath until the colored tip is fully inserted into the microcatheter hub, then tighten the RHV firmly. Pull gently on the introducer sheath to ensure that it is secure. The introducer sheath tip should not move.

Warning: Confirm there are no air bubbles trapped anywhere in the stent system.

Note: After tightening the RHV firmly, the introducer sheath tip should not move when pulled gently. Failure to secure the introducer sheath may result in premature deployment of the stent within the microcatheter hub or difficulty in transferring the stent.

Note: The introducer sheath colored tip must be fully inserted into the microcatheter hub to enable the stent to move into the microcatheter. Over-tightening the RHV may crush the introducer sheath, while under-tightening the RHV may result in premature deployment of the stent.

12. Advance the stent delivery wire to transfer the stent from the introducer sheath into the microcatheter.

Note: Ensure that the introducer sheath does not move while advancing the stent delivery wire. Movement of the introducer sheath during stent advancement may indicate an inadequately tightened RHV and may result in premature deployment of the stent within the microcatheter hub.

13. Continue advancing the stent delivery wire into the microcatheter until the distal edge of the white fluoro saver mark enters the introducer sheath. The white fluoro saver mark is 135 cm from the stent delivery wire distal tip. When the white fluoro saver mark enters the introducer sheath, the stent is about 90 cm inside the microcatheter.
14. Loosen the RHV, remove the introducer sheath from the proximal end of the stent delivery wire while holding the stent delivery wire fixed in place, and set introducer sheath aside in sterile field.

Note: At this point, fluoroscopy may be used at the physician's discretion.

15. Slowly advance the stent delivery wire until the distal edge of the stent delivery wire fluoro saver mark reaches the microcatheter's RHV.

Stent Positioning and Deployment

16. Under fluoroscopy, advance the stent delivery wire until the stent's distal radiopaque markers are 1 – 2 mm from the distal tip marker of the microcatheter.

Note: Do not apply undue force if resistance is encountered at any point during stent manipulation. Withdraw the microcatheter, stent and stent delivery wire as a unit and repeat the procedure with new devices.

17. Withdraw the microcatheter slightly to remove any slack from the stent system and to position the stent for deployment by aligning the stent radiopaque markers across the target aneurysm.

Note: Maintain adequate stent length (approximately 4 mm) on each side of the aneurysm neck to ensure appropriate neck coverage.

18. Before stent deployment, if the microcatheter must be repositioned with a guidewire, the stent and the stent delivery wire must be pulled back into the introducer sheath as follows: Tighten the RHV firmly to hold the stent delivery wire and the microcatheter together, then reposition the distal flexible segment of the catheter until it is in a relatively straight segment of the artery. Load the introducer sheath over the proximal end of the stent delivery wire into the RHV attached to the microcatheter, and ensure that the colored tip is fully inserted into the microcatheter hub. Tighten the RHV firmly. Pull gently on the introducer sheath to ensure that it is secure. The introducer sheath tip should not move.

Carefully retract the stent delivery wire so that the stent is pulled back into the introducer sheath proximal to the colored tip section. Verify that the sterile heparinized saline flush exits the introducer sheath to ensure that no blood products remain in the sheath. Ensure that the distal tip of the stent delivery wire is fully inside the introducer sheath, then loosen the RHV and remove the stent system from the RHV. Once the stent system is removed from the microcatheter, place the stent system in sterile heparinized saline for reuse after microcatheter repositioning. Re-insert the access length guidewire to reposition the microcatheter, using the guidewire introducer if desired.

Note: At least 50 cm of the stent delivery wire should be proximal to the RHV in order to load the introducer sheath.

Note: Be careful that the introducer sheath is loaded in the correct orientation – the colored tip should be fully inserted into the hub of the microcatheter.

Note: After tightening the RHV firmly, the introducer sheath tip should not move when pulled gently. Failure to secure the introducer sheath may result in premature deployment of the stent within the microcatheter hub or difficulty in transferring the stent.

Note: Ensure that the introducer sheath does not move while retracting the stent delivery wire. Movement of the introducer sheath during stent retraction may indicate an inadequately tightened RHV and may result in premature deployment of the stent within the microcatheter hub.

19. If stent positioning is satisfactory, carefully retract the microcatheter in a continuous movement, while maintaining the position of the stent delivery wire to allow the stent to deploy across the neck of the aneurysm. The stent's distal markers will expand as it exits the microcatheter. Confirm deployed stent position.

Note: Do not deploy the stent if it is not properly positioned in the vessel.

Note: The physician may attach the optional torque device to the proximal end of the stent delivery wire, which may facilitate handling and stabilization. Be sure to tighten the torque device to secure the stent delivery wire. Do not use the torque device to torque the stent delivery wire as it is not designed to be torqued.

Note: Do not use the stent delivery wire to push the stent out of the microcatheter while deploying.

20. Prior to removing the stent delivery wire, position the microcatheter distal to the stent to maintain access through the stent. Remove and discard the Neuroform EZ™ Stent System.
21. Advance an exchange length guidewire through the microcatheter.
22. Remove the 0.027 in microcatheter while maintaining the position of the exchange length guidewire, and replace with an appropriate microcatheter to begin aneurysm embolization.

Aneurysm Embolization

23. The aneurysm embolization can begin immediately. Standard microcatheters accepting 0.25 mm (0.010 in), 0.36 mm (0.014 in), or 0.46 mm (0.018 in) guidewires with distal tip $\leq 2F$ may be carefully placed through the interstices of the stent to place embolic coils in the aneurysm.

Note: Carefully watch the stent markerbands when passing through the deployed stent with embolic coiling microcatheters to avoid dislodging the stent.

24. Perform a standard embolic coiling procedure using accepted embolic coiling practices.

QUESTIONS AND ANSWERS

Q: What is the optimal position of the stent with respect to the aneurysm?

A: Generally, try to position the stent so that each end of the stent is secured in relatively normal areas of the parent vessel. The stent will be more stable if each end of the stent is anchored in at least 4-6 mm of normal vessel. For example, for aneurysms located in the supraclinoid carotid, it may be better to secure the stent by deploying the distal end in the M1 (middle cerebral artery, first segment) than trying to deploy it in the few millimeters between the aneurysm and the ICA (internal carotid artery) bifurcation. When deploying the stent, care should be taken to use a view that best shows the parent vessel distal to the aneurysm, so that the distal end of the stent can be accurately deployed with respect to the aneurysm. This view may be different from the view used to advance the Neuroform EZ Stent System, or the view used as a working position for aneurysm embolization.

Deploy the distal end of the stent as precisely as possible with respect to the aneurysm neck to assure at least 4 mm of each end of the stent lies along the parent vessel. Ensure accurate measurement of the aneurysm neck so that, when the stent is properly sized per Table 1, the proximal end of the stent will deploy at the correct location, even if it is difficult to see it because of curves in the vessel.

Q: Which stent size should I choose if I intend to place the stent in a vessel that has a different diameter between the proximal and distal ends of the stent? Example: Vessel increases from 2 mm PCA (posterior communicating artery) to a 3.4 mm basilar.

A: Choose the stent sized for the larger vessel. In this example, choose the 3.5 mm stent. This stent can be deployed safely in the smaller PCA and will be well anchored in the basilar artery.

Q: Is there any problem with deploying the stent across a branch vessel? Can the stent be safely deployed across the anterior choroidal artery? What about lenticulostriate arteries or perforators arising from the basilar?

A: No adverse events resulting from branch vessel occlusion or emboli to "jailed" vessels have been observed in the limited clinical study conducted on this stent (26 patients followed through 6 months). Stents have been placed extending from the M1 (middle cerebral artery, first segment) to the ICA (internal carotid artery) without problems.

Q: A loop, or several loops, of a coil (especially a small diameter coil such as 2 mm) are protruding through the interstices of the stent, and I am unable to reposition it. What should I do?

A: If the risk of leaving part of the coil in the parent vessel is unacceptable, place a second stent inside the first stent to pin the herniated coil portion against the wall of the vessel. Three-dimensional angiography using an orthogonal view (i.e., "down the barrel") may be helpful to assess whether or not a coil loop is inside the lumen of the stent, in the parent vessel, or between the wall of the vessel and the stent. Strict attention to heparinization and antiplatelet medication is important.

Q: I have accidentally started to deploy the stent, but it is not in the location that I wanted. What should I do?

A: The safest course of action generally is not to try repositioning the stent, but to continue to deploy the stent where it is, and then deploy a second stent at the desired location. Safely deploying a stent, even in an undesired location will minimize vascular injury. Animal studies have demonstrated that the stent endothelializes in less than 30 days.

Q: I misjudged the positioning of the stent and have deployed it with one end adjacent to the aneurysm rather than in the normal part of the parent vessel. What should I do?

A: Remove the Neuroform EZ Stent System from the microcatheter while maintaining the position of the microcatheter. Insert and deploy a second stent starting from inside from the first stent to the normal portion of the parent vessel (telescoping stents). The second stent should be of the same diameter or larger than the first.

WARRANTY

Stryker Neurovascular warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond Stryker Neurovascular's control directly affect the instrument and the results obtained from its use. Stryker Neurovascular's obligation under this warranty is limited to the repair or replacement of this instrument and Stryker Neurovascular shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. Stryker Neurovascular neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **Stryker Neurovascular assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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Copyright © 2017 Stryker
Stryker Australia Pty Ltd