

**Operator's manual**

Reprocessed  
**Prevalon® MATS**

**Mobile Air Transfer System**

**Reprocessed device for single-use**

**Non-sterile (exposed to EO gas)**



**Model #R3242-RP**

- (1) 39" x 81" Reprocessed Prevalon Mobile Air Transfer Mat
- (1) M<sup>2</sup> Microclimate Body Pad

**For use with the following non-reprocessed items:**

- Model #7455 — Prevalon Air Pump
- Model #7475 — Prevalon Cart for use with model #7455
- Model #7465 — HEPA Equipped Replacement Filter for model #7455
- Model #7460 — Prevalon Hose Sleeve
- Model #3010 — Prevalon Blower
- Model #3011 — Cart for Prevalon Blower
- Model #3010-F — Replacement Filter for model #3010

If the Air Pump / Blower is not equipped with a Quick Connect Nozzle, contact Sage Products at 800 323 2220

**Reprocessed by:**  
 Stryker Sustainability Solutions  
 1810 West Drake Dr.  
 Tempe, AZ 85282  
[sustainability.stryker.com](http://sustainability.stryker.com)  
 888 888 3433

**Reprocessed by:**  
 Stryker Sustainability Solutions  
 1810 West Drake Dr.  
 Tempe, AZ 85282  
[sustainability.stryker.com](http://sustainability.stryker.com)  
 888 888 3433

# Reprocessed Prevalon® MATS

## Mobile Air Transfer System

### Uses

To assist with the lateral transfer of a patient in a horizontal position from one surface to another. To assist with repositioning of a patient.

### Intended care settings

The Reprocessed Prevalon Mobile Air Transfer System (“System”) is intended to be used in hospitals, healthcare facilities, and other environments where the horizontal lateral transfer of patients is required, and where repositioning of patients is desired. This system is only intended to be used for purposes stated in this Operator’s Manual. Any other use is prohibited.

### Indications

The Reprocessed System is for patients who are dependent and unable to participate in lateral transfer and/or repositioning. The System has a weight capacity of 1000 lbs / 454 kg. The Reprocessed Prevalon Mobile Air Transfer System Mat (“Mat”) is intended for single patient use. The M<sup>2</sup> Microclimate Body Pad is intended for single use.

### Contraindications

Patients who have thoracic, cervical, or lumbar fractures that are deemed unstable should not use the System unless a clinical decision has been made by your facility.

The System should not be used for patients that either exceed the 1000 lbs / 454 kg weight capacity or for patients whose shape does not comfortably fit within the edges of the Mat.

### Prior to use


Perform a visual inspection of the System. Examine the Mat, the Prevalon Air Pump (the “Air Pump”), and the Prevalon Blower (the “Blower”) for damage or wear that is potentially hazardous and could cause the System to malfunction.

If such damage is present, **do not use**. Examine the power cord for damage. If the power cord is damaged, **do not use** the Air Pump / Blower, as a damaged power cord may cause serious injury. Ensure all parts are present. If any parts are missing, **do not use**.


### Imaging

The Reprocessed Prevalon Mobile Air Transfer System Mat with M<sup>2</sup> Microclimate Body Pad is MR safe by rationale. The device is made from all non-metal materials; therefore, MR safety testing was not performed. Compatibility tests did not show artifacts. Based on rationale, the Reprocessed Prevalon Mobile Air Transfer System Mat with M<sup>2</sup> Microclimate Body Pad is electronically non-conductive and non-magnetic.

The Air Pump / Blower is not MR safe. Please follow your facilities protocol for using the Air Pump / Blower in a MR environment.

 **Caution** Periodically check System for signs of wear or damage. If damaged, do not use.

- **Do not** launder the Mat or the M<sup>2</sup> Microclimate Body Pad.

 **Warnings** Failure to comply with this Manual could result in injury or equipment malfunction. Alterations or modifications to the System can affect the safety and functionality of the System and are strictly prohibited.

### Read and understand the entire manual before using the System.

**To reduce the risk of fire or injury:**


1. Grounding reliability can only be achieved when the equipment is connected to an equivalent receptacle marked “hospital only” or “hospital grade.”
2. **Do not** cover or block any openings on the Air Pump / Blower.
3. **Do not** tamper with or make any adjustments to the control panel, housings, safety switches, or any other parts of the Air Pump / Blower.
4. **Do not** remove any panel or cover.
5. **Do not** store or use the Air Pump / Blower outdoors or near water.
6. The Air Pump / Blower has been designed for use exclusively with the Sage air assisted technology products.
7. **Do not** operate the Air Pump / Blower if it is not working properly, or if it has been damaged, including the cord or plug. Contact Sage Products for repair or replacement.
8. **Do not** immerse cord or plug in water.
9. Keep cord away from heated surfaces.
10. **Do not** use in the presence of flammable anesthetics and other flammable solvents and materials.
11. **Do not** hang Air Pump / Blower from the bed rails.
12. If your hospital protocol directs that the Mat is to be tucked underneath a surgical table pad, ensure the Mat does not inhibit the function of the surgical table pad.
13. **Do not** use to lift patients.
14. Always follow your facility’s safe patient handling policies and procedures.
15. The Mat is not intended to secure a patient to a Support Surface. Follow facility policies and procedures for securing patients to a Support Surface. (e.g. Trendelenburg)
16. Ensure patient is not in contact with the Quick Connect Nozzle, Quick Connect Valve or Flexible Inflation Hose.
17. Patient repositioning should always be performed using at least two caregivers.

# Reprocessed Prevalon<sup>®</sup> MATS


Mobile Air Transfer System

## Reprocessed Prevalon Mobile Air Transfer System Mat (“Mat”)

### Part identification



## Prevalon Air Pump



## Prevalon Blower and Filter (“Blower”)





# Instructions for use



Figure 1




Figure 2



Figure 3

## Set up the system

### Prior to use

1. Perform a visual inspection of the System.

Examine the Mat and the Air Pump / Blower for damage that is potentially hazardous or could cause the System to malfunction. If such damage is present, **do not use**.

2. Ensure all parts are present. If any parts are missing, **do not use**.

### Set up the Mat

**Warning:** Make sure bed rails are up and bed brakes are locked to ensure that the patient and Mat do not slide off the Support Surface, which could cause a fall risk.

3. Make sure the patient's Support Surface is horizontal and at the Caregivers' waist level. Place the folded Mat on the bed with the printed arrow pointing toward the head of the bed. (Figure 1)

4. Unfold the Mat lengthwise and then unfold widthwise. If placing the Mat beneath a patient who is already on the Support Surface, place the Mat and M<sup>2</sup> Microclimate Body Pad beneath the patient per your facility's protocol.

Ensure the patient is centered on the Mat and center the Mat on the Support Surface. (Figure 2)

**Warning:** To prevent the patient from falling, before every transfer, make sure that the patient is positioned in the center of the Mat and the Mat is centered on the Support Surface.

### Prepare to transfer the patient

**Warning:** Always use at least two caregivers to perform a lateral transfer. Exterior bed rails on both surfaces should be raised prior to transfer to prevent the patient from falling. If there are no bed rails used, the Caregivers are responsible for making sure the patient does not reach outside the boundaries of either Support Surface.

5. Place both Sending and Receiving Surfaces in a horizontal position. Set the Receiving Surface at the Receiving Caregiver's ("Receiver") waist height. Place the Sending Surface slightly above, but no more than one inch above the Receiving Surface. Raise exterior bed rails. Move Sending and Receiving Surfaces as close together as possible. If the Sending Surface needs to be moved, ensure all bed rails are raised and locked. (Figure 3)



Figure 4



Figure 5



Figure 6



Figure 7



Figure 8



Figure 9

6. Lock the brakes on both Support Surfaces. (Figure 4)

### Connect the Air Pump / Blower

7. Locate the Quick Connect Valve on the bottom corner of the Mat closest to the Air Pump / Blower. Ensure the gray sleeve and two orange loops are exposed. (Figure 5)

8. Confirm the Air Pump / Blower is equipped with a Quick Connect Nozzle. - See Page 4

**NOTE:** If the Air Pump / Blower is not equipped with a Quick Connect Nozzle, contact Sage Products at 800 323 2220.

9. Grasp the Quick Connect Nozzle with one hand. With your other hand grasp one of the orange loops on the Quick Connect Valve. (Figure 6)

10. Insert one side of the Quick Connect Nozzle into the Quick Connect Valve. (Figure 7)

11. Pull the orange loop to expand the Valve open and insert the exposed end of the Quick Connect Nozzle into the Quick Connect Valve. (Figure 8)

12. Pull on the hose to confirm the Quick Connect Nozzle is secure. (Figure 9)

**Warning:** To prevent injury or accidental inflation, ensure patient is not in contact with the Valve or Hose at all times.



Figure 10



Figure 11



Figure 12



Figure 13

## Lateral transfer

### Inflate the mat

13. Sending Caregiver ("Sender"): stand on the side of the Sending Surface.  
Receiver: stand on the side of the Receiving Surface. (Figure 10)

**Warning:** Before inflating, ensure lines and tubing are free to move with the patient and that nothing obstructs the area over which the Mat will pass. Ensure that the hose will move freely with the Mat.

14. Sender: Grasp the Mat handles with both hands so Mat does not drift out of place.
15. Receiver: Press the power switch on the Air Pump / Blower to turn it on. Return to the side of the Receiving Surface. Allow Mat to fully inflate. (Figure 11)

**Warning:** If the Mat is not fully inflated, injury to the patient or caregiver could occur, or the Mat may not perform as expected. Never leave the patient unattended while the Air Pump / Blower is powered on or the Mat is inflated.

### Transfer the patient

16. Sender: Gently push the patient toward the Receiver. (Figure 12)
17. Receiver: Grasp the handles and help glide the patient until they are centered on the Receiving Surface. (Figure 12)

### Deflate the Mat

**Warning:** Ensure the patient is centered on the Receiving Surface prior to deflating the Mat.

18. Receiver: Continue to hold the handles with both hands so Mat does not drift out of place.
  19. Sender: Press the power button on the Air Pump / Blower to turn it off.
- Note: In cases of emergency when the Mat needs to be rapidly deflated, remove the Quick Connect Nozzle from the Quick Connect Valve.
20. Allow the Mat to fully deflate and smooth out any wrinkles in the Mat and M<sup>2</sup> Microclimate Body Pad. (Figure 13)
  21. Slightly separate the Receiving and Sending Surfaces until the rail on the Receiving Surface can be accessed.

22. Raise all rails.
23. Disconnect the Quick Connect Nozzle from the Quick Connect Valve.
24. Unplug the Air Pump / Blower and roll up the power cord prior to moving or when not in use. If using the Cart, unlock the wheels before moving. Make sure the wheels move freely when moving the Cart.

## Repositioning the patient

**Warning:** To avoid caregiver or patient injury, make sure the bed rails are raised and the bed brakes are locked.

1. Caregivers stand on opposite sides of the bed. Lock the bed brakes and raise the bed rails. Place the bed in a horizontal position and at waist level.
2. Ensure the patient is centered on the Mat and center the Mat on the Support Surface. With power off, connect the Quick Connect Nozzle on the Air Pump / Blower to the Mat by inserting the Quick Connect Nozzle into the Quick Connect Valve. Ensure Quick Connect Nozzle is securely connected.

**Warning:** Before inflating, ensure lines and tubing are free and that nothing obstructs the area over which the linens will pass.

**Warning:** Never leave the patient unattended while the Air Pump / Blower is powered on or the Mat is inflated.

3. While closely observing the patient, press the power switch on the Air Pump / Blower to turn it on.
4. Allow the Mat to fully inflate.
5. Use the handles to gently glide the patient to where the hips are aligned with the hinge point of the bed.
6. While one Caregiver holds the mat in place, the second Caregiver turns the Air Pump / Blower off. Disconnect the Quick Connect Nozzle from Quick Connect Valve.
7. Allow Mat to fully deflate and smooth out any wrinkles in the Mat and M<sup>2</sup> Microclimate Body Pad.

### Changing the Prevalon M<sup>2</sup> Microclimate Body Pad

Replace only with Prevalon M<sup>2</sup> Microclimate Body Pads. Align the edges of the M<sup>2</sup> Microclimate Body Pad with the M<sup>2</sup> Microclimate Body Pad Indicator Lines. Dispose of the soiled M<sup>2</sup> Microclimate Body Pad per your facility's waste management protocol.

## Changing linens

### Removing linens

**Warning:** To avoid caregiver or patient injury, make sure the bed rails are raised and the bed brakes are locked.

1. Caregivers stand on opposite sides of the bed. Lock the bed brakes and raise the bed rails. Place the bed in a horizontal position and at waist level.





Figure 14



Figure 15



Figure 16



Figure 17



Figure 18

2. Ensure the patient is centered on the Mat and center the Mat on the Support Surface. With power off, connect the Quick Connect Nozzle on the Air Pump / Blower to the Mat by inserting the Quick Connect Nozzle into the Quick Connect Valve. Ensure Quick Connect Nozzle is securely connected.

**Warning:** Before inflating, ensure lines and tubing are free and that nothing obstructs the area over which the linens will pass.

**Warning:** Never leave the patient unattended while the Air Pump / Blower is powered on or the Mat is inflated.

3. While closely observing the patient, press the power switch on the Air Pump / Blower to turn it on.
4. Allow the Mat to fully inflate.
5. Starting at the head of the bed, lift the corners of the linen off the mattress. (Figure 14)
6. Both Caregivers hold the Mat with one hand and gently pull the linen toward the foot of the bed, allowing it to pass below the inflated Mat. (Figure 15)
7. While one Caregiver ("Caregiver 1") holds the mat in place, the second Caregiver ("Caregiver 2") removes the linen from the bed. Follow your facility's protocol for handling soiled linens. (Figure 16)
8. Caregiver 2 turn the Air Pump / Blower off. Disconnect Quick Connect Nozzle from Quick Connect Valve.
9. Allow the Mat to fully deflate and smooth out any wrinkles in the Mat and M<sup>2</sup> Microclimate Body Pad. (Figure 17)

## Placing linens

**Warning:** To avoid caregiver or patient injury, make sure the bed rails are raised and the bed brakes are locked.

10. Caregivers stand on opposite sides of the bed. Lock the bed brakes and raise the bed rails. Place the bed in a horizontal position at waist level.
11. Ensure the patient is centered on the Mat and center the Mat on the Support Surface. With power off, connect the Quick Connect Nozzle on the Air Pump / Blower to the Mat by inserting Quick Connect Nozzle into the Quick Connect Valve. Ensure Quick Connect Nozzle is securely connected.
12. Tuck the bottom corners of the linen beneath the patient's pillow, Mat, and shoulders.
13. Anchor the top corners of the linen to the mattress. (Figure 18)

**Warning:** Before inflating, ensure lines and tubing are free and that nothing obstructs the area over which the linens will pass.

**Warning:** Never leave the patient unattended while the Air Pump / Blower is powered on or the Mat is inflated.



Figure 19



Figure 20

14. Caregiver 1 hold the Mat in place, Caregiver 2 turn the Air Pump / Blower on.
15. Allow the Mat to fully inflate. (Figure 19)
16. Both Caregivers hold the Mat with one hand and gently pull the linen toward the foot of the bed, allowing it to pass below the inflated Mat. (Figure 20)
17. Once the linen is at the bottom of the bed, Caregiver 1 hold the Mat in place and Caregiver 2 turn the Air Pump / Blower off. Disconnect the Quick Connect Nozzle from the Quick Connect Valve.
18. Allow the Mat to fully deflate. Anchor the bottom corners of the linen around the mattress and smooth out any wrinkles in the Mat and M<sup>2</sup> Microclimate Body Pad.

## Cleaning instructions

### Cleaning and disinfecting

Always follow your facility's infection control policies and procedures regarding the cleaning of equipment and the use of disinfectants on facility surfaces in contact with patients' skin.

### Mat cleaning instructions

The Reprocessed Prevalon Mobile Air System Mat is for single patient use. The M<sup>2</sup> Microclimate Body Pad is for single use only.

To clean the Reprocessed Prevalon Mobile Air Transfer System Mat, wipe all surfaces of the Mat with a damp cloth using soap and water.

### Air Pump / Blower cleaning instructions

#### Air Pump / Blower general cleaning:

1. Prior to cleaning, unplug the Power Cord from the wall outlet.
2. Prepare a mild detergent solution according to the detergent manufacturer's instructions.
3. Wipe the outside of the Air Pump / Blower with a large cloth, dampened with the solution. Ensure the cloth is not so wet as to drip liquid or cause liquid to pool on the Air Pump / Blower. To prevent damage to the operating parts inside the Air Pump / Blower, do not allow liquid to seep into openings of the Air Pump / Blower.
4. Stretch the Hose in sections to clean between coils. Continue until the entire Hose is clean. Do not allow solution to seep inside the Hose.
5. Wipe all excess liquid from the unit and allow to air dry. Allow the external surfaces of the Air Pump / Blower and Hose to dry thoroughly before use.

#### Air Pump / Blower disinfecting:

Use a hospital grade disinfectant. Always follow the disinfectant manufacturer's instructions for use. Prior to disinfecting, follow the general cleaning instructions described above to remove any visible debris or soiling.

1. Use rubber gloves and eye protection as recommended by the disinfectant manufacturer.



Figure 21




Figure 22




Figure 23




Figure 24

2. Prior to disinfecting, unplug the power cord from the wall outlet.
3. Prepare disinfectant solution (if applicable) according to the disinfectant manufacturer's instructions for use.
4. Thoroughly wipe down the outside of the Air Pump / Blower with the disinfectant solution. Ensure cloth is not so wet as to drip liquid or cause liquid to pool on the Air Pump / Blower. To prevent damage to the operating parts inside the Air Pump / Blower, do not allow liquid to seep into the openings.
5. Stretch the Hose in sections to clean between coils. Continue until the entire Hose is clean. Do not allow solution to seep inside the Hose.
6. Allow Air Pump / Blower and Hose to air dry thoroughly before use.

## Filter replacement

### Air Pump Filter removal and replacement

**Warning:** To reduce the risk of fire, electrical shock, or injury, unplug Air Pump from electrical outlet before servicing.

1. Unplug the Air Pump / Blower from the wall. (Figure 21)
2. Rotate Filter counterclockwise until released from housing. (Figure 22)
3. Pull Filter rearward to remove. Dispose of filter according to your facility's policy. (Figure 23)
4. To install a new Filter, insert into opening. Rotate clockwise and push Filter until you feel slight resistance. (Figure 24)
5. Rotate Filter clockwise 1/3 turn more.



Figure 25



Figure 26



Figure 27



Figure 28




Figure 29



Figure 30

## Blower Filter replacement

**Warning:** To reduce the risk of fire, electrical shock, or injury, unplug Blower from electrical outlet before servicing.

1. Unplug Blower from wall. Remove the hose by twisting the gray fitting on top of the pump 1/4 turn counterclockwise. Then pull to remove. Ensure that the teeth of the gray fitting are in alignment with the grooves in the opening at the top of the motorhead so it can be removed properly. (Figure 25)
2. Remove the power cord by lifting the retention tab upward and pulling cord from receptacle. (Figure 26)
3. Release the power cord bundle by pulling the tab to release the retention strap. (Figure 27)
4. Turn the pump on its side and remove the thumbscrew on the bottom of the filter in a counterclockwise direction. (Figure 28)
5. Remove the filter canister and remove the filter by pulling it off of the motor housing. (Dispose of the filter per your facility's policy.) (Figure 29)
6. Place the top of the motorhead into the filter canister to support it during the filter replacement procedure. (Figure 30)






Figure 31

7. Place the new filter over the end of the motor housing and ensure that the threaded rod passes through the end of the filter. (Figure 31)
8. Pull the elastic end of the filter all the way down so that it covers the yellow band. (Figure 32)
9. Remove the motor head from the filter canister and place it on its side. Then slide the filter canister over the motor head aligning the rod to pass through the bottom hole.
10. Replace the thumbscrew back on the rod and tighten by turning in a clockwise direction until snug.

**Caution:** If overtightened, the filter will not properly function.




Figure 32

11. Stand the pump up.
12. Lift the gray retention tab on the motorhead and fully insert the plug into the receptacle opening.
13. Push the retention tab downward to lock the tab into place. (Figure 33)
14. Gather cord and place the gathered cord in the retention strap on the pump and snap the strap secure. (Figure 34)
15. Grasp the gray fitting end of the hose and place the end in the opening at the top of the motorhead.
16. Ensure that the teeth on the fitting end are aligned in the grooves in the opening.
17. Push in the fitting and twist 1/4 turn in a clockwise direction until it snaps into place.
18. Secure the Quick Connect Nozzle of the hose by placing the fitting in the retention strap on the pump and snap the strap to secure the Quick Connect Nozzle in place. To keep the hose off of the floor, wrap the hose around the top of the Blower. (Figure 35)



Figure 33




Figure 34



Figure 35

## Repair

Contact Sage Products at 800 323 2220 for repair or replacement.

# Hose Protection Sleeve (HPS) for use with Prevalon Air Pump

## Part identification



## Installing the Hose Protection Sleeve (HPS)

1. Remove the HPS from the package. (Figure 36)
2. Place the HPS over the Quick Connect Nozzle. The HPS core with round pegs should be facing towards the base of the Air Pump. (Figure 37)
3. Grasp the rubber band and ensure it is seated within the Nozzle Groove. (Figure 38)



Figure 36



Figure 37



Figure 38



# Air Pump/ Blower technical specifications




Figure 39

4. Holding the sleeve at the Quick Connect Nozzle end, pull the HPS core towards the base of the Air Pump. (Figure 39)



Figure 40

5. Once the HPS core is at the base of the Air Pump, locate the round pegs and align them with the openings on the Air Pump. (Figure 40)




Figure 41

6. Insert the pegs into the opening and rotate clockwise to lock. (Figure 41)

## Removal of the Hose Protection Sleeve (HPS)




Figure 42

7. Rotate the HPS core counterclockwise to unlock. (Figure 42)



Figure 43

8. While holding the Hose near the base of the pump, pull the HPS core towards the Quick Connect Nozzle. (Figure 43)



Figure 44

9. While removing the HPS, the sleeve will turn inside out. (Figure 44)

10. Continue pulling the HPS until it is completely off the Hose and Quick Connect Nozzle. (Figure 45)

11. Discard the HPS per your facility's protocol. (Figure 46)



Figure 45



Figure 46







| Product Code   | 7455 (US Configuration)  | 3010 (US Configuration)  |
|--|--|--|
| Input Rating   | 120 V AC / 60 Hz / 10 AMP  | 120 V AC / 60 Hz / 10 AMP  |
| Ambient Operating Range  | 50°F (10°C) - 86°F (30°C)  | 50°F (10°C) - 86°F (30°C)  |
| Relative Humidity Operating Range                                    | 10 - 70% Non-Condensing  | 10 - 70% Non-Condensing  |
| Operating Pressure   | 70 kPa - 106 kPa   | 70 kPa - 106 kPa   |
| Outer Dimensions   | 10.24"H x 26.58"L x 11.80"W (260mm x 675mm x 300mm)  | 10"H x 10"W x 18"D (254mm x 254mm x 451mm)   |
| Gross Weight (G.W.)  | G.W.: 14 lbs. (6.35kg)   | G.W.: 14 lbs. (6.35kg)   |
| Max Cart Loading   | NA   | NA   |
| Classification   | Class II With Functional Earth   | Class II With Functional Earth   |
| Degree of Protection Against Ingress of Water                        | IPX0   | IPX0   |
| Mode of Operation-Duty Cycle   | 30 sec on 1 min off 5 Cycles w/30min rest  | 30 sec on 1 min off 5 Cycles w/30min rest  |
| Storage / Transport Temperature                                      | 24.8 - 158°F (-4 - 70°C)   | 24.8 - 158°F (-4 - 70°C)   |
| Storage / Transport Relative Humidity                                | 10% to 90%   | 10% to 90%   |
| Storage / Transport Atmospheric Pressure Range                       | 500 to 1060 hPa  | 500 to 1060 hPa  |
| Standards - UL / EN / CSA / IEC 60601-1                              | UL 60601-1 Medical Electrical Equipment, 2nd Edition -3rd Edition  | UL 60601-1 Medical Electrical Equipment, 2nd Edition -3rd Edition  |
| UL Classification  | IEC 60601-1: 2012, Edition 3.1   | IEC 60601-1: 2012, Edition 3.1   |
| Third Party Evaluation and Testing                                   | IEC 60601-1: 2012, Edition 3.1   | IEC 60601-1: 2012, Edition 3.1   |
| Characteristics-Construction   | Plastic Housing  | Plastic Housing  |
| Characteristics-Motor  | Two state, single speed, 120 volts, double ball bearings, low-noise bypass discharge, 40 peak horsepower | Two state, single speed, 120 volts, double ball bearings, low-noise bypass discharge, 40 peak horsepower |
| Characteristics-Airflow  | 5,300 ft/min / 1,615 m/min   | 5,900 ft/min / 1,798 m/min   |
| Characteristics-Hose   | 8 ft x 1.5 in / 3.048 m x 3.81 cm commercial strength flexible hose                                      | 8 ft x 1.5 in / 3.048 m x 3.81 cm commercial strength flexible hose                                      |
| Characteristics-Nozzle   | Custom heavy duty right-angle fast connector   | Custom heavy duty right-angle fast connector   |
| Characteristics-Cord / Connector                                     | Heavy duty medical grade - grounded IEC female   | Heavy duty medical grade - grounded IEC female   |
| Characteristics-IEC Power Entry Module with twin UL rated 10 A Fuses | NA   | NA   |
| Fuse Type  | NA   | NA   |
| Disposal   | Follow national requirements   | Follow national requirements   |
| Consult Instructions for Use   | Refer to accompanying documents  | Refer to accompanying documents  |
| Protective Earth (Ground)  | NA   | NA   |
| Applied Parts  | Air Assist Mat is Type B Applied Part  | Air Assist Mat is Type B Applied Part  |
| Caution  | DO NOT use in the presence of flammable anesthetics and other flammable gases or vapors.                 | DO NOT use in the presence of flammable anesthetics and other flammable gases or vapors.                 |
| Filter   | Glass Fiber Tested per IEST-RP-CC001.5   | ePTFE Filter Media Lot Tested per ASTM D 2986  |
|  |  | 5 micro poly filter  |











Medical equipment with respect to electric shock, fire and mechanical hazards only, in accordance with ANSI/AAMI ES60601-1 and 1 2012, CAN/CSA C22.2 NO. 60601-1-114, ANSI/AAMI/IEC 60601-2.

E472680

Note: The **emissions** characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 Class A). If used in a residential environment (for which CISPR 11 Class B is normally required), this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

|   |   |  |   |
|---|---|--|---|
|  | 1. Use two caregivers.                  |  | 4. Using head and foot indicators, center patient.  |
|  | 2. Lock bed brakes.                     |   | 5. Raise exterior bed rails.<br>Inflate Mat. Do not leave patient unattended on inflated Mat.<br>6. Transfer patient. |
|  | 3. Ensure patient is lying flat on Mat. |  | 7. Deflate Mat. Raise bed rails when finished.  |

|   |  |
|---|--|
|   | Caution: Plug directly into wall outlets only. Do not plug into bed outlets. |
|  | Consult instructions for use   |
|  | Type B Applied Part  |
|  | Double wire insulated  |

|   |                                    |
|---|------------------------------------|
|  | Date of reprocessing               |
| <b>REF</b>  | Product code                       |
|  | Do not reuse                       |
|  | Non-sterile                        |
|  | Not made with natural rubber latex |

## Warranty

### Reprocessed products

Stryker warrants all reprocessed products, subject to the exceptions provided herein, to be free from defects in reprocessing and to substantially conform to the product specifications contained in the documentation provided by Stryker with the products for one use in accordance with the instructions for use of such product.

**STRYKER SHALL NOT BE LIABLE FOR ANY DAMAGES TO THE EXTENT CAUSED BY ANY DEFECT IN MATERIAL, WORKMANSHIP OR DESIGN BY THE ORIGINAL MANUFACTURER OF THE PRODUCT OR ANY ACT OR OMISSION OF THE ORIGINAL MANUFACTURER OF THE PRODUCT.**

### Products for which Stryker is the original manufacturer

Stryker warrants all products for which it is the original manufacturer, subject to the exceptions provided herein, to be free from defects in design, materials and workmanship and to substantially conform to the product specifications contained in the documentation provided by Stryker with the products for a period of one year from the date of purchase.

### General warranty terms applicable to all products

**TO THE FULLEST EXTENT PERMITTED BY LAW, THE EXPRESS WARRANTY SET FORTH HEREIN IS THE ONLY WARRANTY APPLICABLE TO THE PRODUCTS AND IS EXPRESSLY IN LIEU OF ANY OTHER WARRANTY BY STRYKER, EXPRESSED OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OR MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL STRYKER'S LIABILITY ARISING IN CONNECTION WITH THE SALE OF THE PRODUCT (WHETHER UNDER THE THEORIES OF BREACH OF CONTRACT, TORT, MISREPRESENTATION, FRAUD, WARRANTY, NEGLIGENCE, STRICT LIABILITY OR ANY OTHER THEORY OF LAW) EXCEED THE PURCHASE PRICE, CURRENT MARKET VALUE OR RESIDUAL VALUE OF THE PRODUCTS, WHICHEVER IS LESS. STRYKER SHALL NOT BE LIABLE FOR INDIRECT, SPECIAL, INCIDENTAL, PUNITIVE OR CONSEQUENTIAL DAMAGES RESULTING FROM ANY BREACH OF WARRANTY OR UNDER ANY OTHER LEGAL THEORY.**

This warranty shall apply only to the original end-user purchaser of products directly from Stryker or a Stryker authorized distributor. This warranty may not be transferred or assigned without the express written consent of Stryker.

This warranty does not apply to: (1) products that have been misused, neglected, modified, altered, adjusted, tampered with, improperly installed or refurbished; (2) products that have been repaired by any person other than Stryker personnel without the prior written consent of Stryker; (3) products that have been subjected to unusual stress or have not been maintained in accordance with the instructions in the user manual or as demonstrated by a Stryker representative; (4) products on which any original serial numbers or other identification marks have been removed or destroyed; or (5) products that have been repaired with any unauthorized or non-Stryker components.

If a valid warranty claim is received within thirty (30) days of the expiration of the applicable warranty period, Stryker will, in its sole discretion: (1) replace the product at no charge with a product that is at least functionally equivalent to the original product or (2) refund the purchase price of the product. If a refund is provided by Stryker, the product for which the refund is provided must be returned to Stryker and will become Stryker's property. In any event, Stryker's liability for breach of warranty shall be limited to the replacement value of the defective or non-conforming part or component.

If Stryker determines in its reasonable discretion that the claimed defect or non-conformance in the product is excluded from warranty coverage as described hereunder, it will notify the customer of such determination and will provide an estimate of the cost of repair of the product. In such an event, any repair would be performed at Stryker's standard rates.

Products and product components repaired or replaced under this warranty continue to be warranted as described herein during the initial applicable warranty period or, if the initial warranty period has expired by the time the product is repaired or replaced, for thirty (30) days after delivery of the repaired or replaced product. When a product or component is replaced, the item provided in replacement will be the customer's property and the replaced item will be Stryker's property. If a refund is provided by Stryker, the product for which the refund is provided must be returned to Stryker and will become Stryker's property.

The OEM information listed on the label is provided as device ID prior to reprocessing and may contain the trademarks of unrelated third parties that do not sponsor this device.

This product and its packaging have been exposed to ethylene oxide gas (EtO). Even though the product then is processed in compliance with all applicable laws and regulations relating to EtO exposure, Proposition 65, a State of California voter initiative, requires the following notice:

Warning: This product and its packaging have been exposed to ethylene oxide. The packaging may expose you to ethylene oxide, a chemical know to the State of California to cause cancer or birth defects or other reproductive harm.

Prevalon is a registered trademark of Sage Products LLC.

MAT EL10075 Rev. A – 10/2018 RM702507