# IntelliVue MX40 Single-Patient-Use Cables

## *Instructions for Use*


### ECG Only

AAMI 5-lead grabber

**REF** 989803172031

IEC 5-lead grabber

**REF** 989803172131


### ECG + SpO<sub>2</sub>

AAMI 5-lead grabber plus SpO<sub>2</sub>

**REF** 989803172051

IEC 5-lead grabber plus SpO<sub>2</sub>

**REF** 989803172151

# PHILIPS

# About This Edition

## Notice

### Copyright © 2011

Koninklijke Philips Electronics N.V.  
All rights are reserved.


## Printing History

### Current Edition/Revision:

Printed in Germany October 2011  
Part Number 453564230191

## Contents

Instructions for Use.....	1	Návod na obsluhu.....	43
Manuel d'utilisation.....	4	Navodila za uporabo.....	46
Gebrauchsanweisung.....	7	Instrucțiuni de utilizare.....	49
Gebruiksaanwijzing.....	10	Указания за употреба.....	52
Instrucciones de uso.....	13	Kasutusjuhend.....	55
Instruções de Uso.....	16	Naudojimo instrukcijos.....	58
Istruzioni d'uso.....	19	Lietošanas pamācība.....	61
Bruksanvisning.....	22	Инструкции по эксплуатации..	64
Bruksanvisning.....	25	Οδηγίες Χρήσης.....	67
Käyttöopas.....	28	Kullanım Talimatları.....	70
Brugervejledning.....	31	Пайдалану нұсқаулары.....	73
Használati útmutató.....	34	ユーザーズガイド.....	76
Instrukcja obsługi.....	37	使用説明.....	79
Návod k obsluze.....	40	使用说明书.....	82

United States law restricts this device to sale by or on the order of a physician.

---

# Product Overview

## Intended Use

The Single-Patient-Use MX40 Patient Cables are intended to provide ambulatory and bedside monitoring of ECG/Impedance Respiration, or ECG/Impedance Respiration plus SpO<sub>2</sub> parameters, for adults who weigh greater than 40kg and pediatrics who weigh greater than 20kg. The cables are for use with IntelliVue MX40 patient wearable monitors, and are intended to be used by trained healthcare personnel in a professional healthcare facility.

## Cable Configurations

The MX40 Single-Patient-Use Cables are available in ECG only or ECG/SpO<sub>2</sub> configurations (see below).

5-Lead* ECG Only		5-Lead* ECG with attached SpO <sub>2</sub> Sensor**		
REF	ECG Color Coding	REF	ECG Color Coding	SpO <sub>2</sub> Sensor Type
989803172031	AAMI color coded grabber connectors	989803172051	AAMI color coded grabber connectors	Finger wrap sensor
989803172131	IEC color coded grabber connectors	989803172151	IEC color coded grabber connectors	Finger wrap sensor

\*All ECG leads are 0.85m (34 in) long.

\*\*The cable for the attached SpO<sub>2</sub> sensor is 2m (6.5 ft) long. The preferred SpO<sub>2</sub> patient application site is any finger except a thumb.

---

## Warnings

- These patient cables are for single patient use only. They should never be reused on a different patient due to risk of cross-infection.
- To avoid patient injury, assure that leadset is not positioned where leads could become entangled around patient, or cause choking, strangulation, or inhibit circulation in extremities.
- Protect the patient cable trunk connector, ECG lead connectors, and SpO<sub>2</sub> sensor from becoming soaked with any liquid.

---

## Patient Application

Visually inspect the ECG leadset (and SpO<sub>2</sub> sensor/cable when present) before applying the patient cable to a patient.

- If any signs of damage or deterioration are seen, immediately discard the patient cable following approved medical waste disposal methods specified by local regulations or your patient care facility.
- If no visual damage is seen, apply the ECG leads to the patient following standard AAMI or IEC practices (refer to the *IntelliVue MX40 Instructions for Use*).
- If the patient cable includes a SpO<sub>2</sub> sensor, follow the sensor application instructions supplied in the *IntelliVue MX40 Instructions for Use*. If the SpO<sub>2</sub> sensor is not being used, be sure to securely store the sensor inside the telemetry pouch.
- Recommended SpO<sub>2</sub> sensor application site is any finger except the thumb.



---

## Additional Information

For more detailed information describing patient cable use with the MX40, including warnings, patient safety, and ECG/SpO<sub>2</sub> patient application instructions, refer to the *IntelliVue MX40 Instructions for Use*.

# Product Symbols

Product packaging symbols are defined below:

 <p>Refer to instruction manual/booklet</p>	 <p>Temperature limitation</p>	 <p>Contains no latex</p>	 <p>Humidity limitation</p>	 <p>Non sterile</p>
 <p>Do not reuse</p>	 <p>Quantity per box</p>	 <p>Prescription use only</p>	 <p>Manufacturer</p>	 <p>Catalog number</p>
 <p>Storage pressure range</p>	 <p>Adult patient</p>	 <p>Adult finger</p>	 <p>Pediatric patient</p>	 <p>Pediatric finger</p>
 <p>European community authorized representative</p>				

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For more information, please call the local Philips sales office listed in your telephone directory; or contact one of the regional sales offices listed below for the location of your nearest Philips sales office.

**United States:**

Philips Medical Systems  
North America Corporation  
3000 Minuteman Road  
Andover, MA 01810-1099  
+1 800 225-0230 toll free

**Canada:**

Philips Medical Systems Canada  
281 Hillmount Road  
Markham, ON  
L6C 2S3  
(800) 291-6743

**Latin America Headquarters:**

Philips Healthcare Brazil  
Rua Verbo Divino 1400 - 7° Andar  
Chacara Santo Antonio  
04719-002 Sao Paulo, SP  
Brazil  
Tel: +55 11 2125-0744  
Fax: +55 11 2125-0761

**Europe, Middle East and Africa:**

Philips Medizin Systeme  
Böblingen GmbH  
Hewlett-Packard Str. 2  
71034 Böblingen, Germany  
Tel: (+49) 7031 463 2254  
Fax: (+49) 7031 463 1552

**Asia Pacific Headquarters:**

Philips Electronics Singapore Private  
Limited  
Philips Healthcare  
620A Lorong 1 Toa Payoh  
Singapore 319762  
Tel: 1800-PHILIPS  
Fax: +65 6255 4853

This product complies with the requirements of the Council Directive 93/42/EEC June 1993 (Medical Device Directive).

Contains no Latex

US Status: Distributed by  
Philips Medical Systems  
Made in Germany

Manufacturing Address:  
Philips Medical Systems  
3000 Minuteman Road  
Andover, MA 01810-1099 USA  
Tel: (800) 225-0230

<http://philips.com/healthcarestore>

Subject to modification

Copyright © 2011  
Koninklijke Philips Electronics N.V.  
Printed in Germany October 2011

Part No. 453564230191