

Reprocessing **overview**

Save Simply

Sage Prevalon® Mobile Air Transfer System (MATS)



1 Product availability:
Our blended MATS are packaged in case quantities of 10 units per cardboard box. Each box contains a blend of both new and reprocessed product.



4 Decontamination and cleaning: After the inspection is complete, all MATS are unfolded and individually spot cleaned by hand, utilizing a scrub brush and validated cleaning agents.



2 Initial inspection:
Once MATS arrive to our reprocessing facility, the teal bags are removed from their shipping boxes and scanned through a metal detector to confirm no needles or sutures are present.




5 Visual inspection:
All MATS are then inspected to confirm stains, debris and contamination has been removed, and to assess the overall device integrity. Any device that does not meet our strict quality standards is rejected and recycled appropriately.



3 Receiving and login:
Each box is scanned in to identify which facilities they came from. They are individually tallied to calculate total quantity per hospital.



6 Performance testing:
All MATS are connected to air pumps to ensure proper inflation and functionality. Once inflated, each handle is pulled individually to test attachment and confirm durability.



7 Device tracking: Once all visual and function tests are passed, a cycle count mark is placed underneath the handle to track reprocessing cycles. Any device that has reached its maximum number of reprocessing cycles is rejected.



9 Packaging: An M² Microclimate Body Pad is added, and the MATS are folded and placed in their packaging.



8 Ethylene oxide (EO) exposure: All MATS are exposed to EO which is validated utilizing ANSI/AAMI/ISO 11135 as a guideline for this bioburden reduction process. EO decontamination validation demonstrates the process and equipment are capable of consistently and reliably achieving a minimum six spore log reduction (SLR). EO residuals do not exceed maximum allowable limits of ANSI/AAMI/ISO 10993-7.



10 Shipping: Finished reprocessed MATS are shipped to Stryker Sage where they are placed into a blended box of both reprocessed and new MATS.

Stryker's Sustainability Solutions

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. We do not dispense medical advice and recommend that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate Stryker's products. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

MKT10204 REV B | Copyright © 2023 Stryker

Stryker
Stryker's Sustainability
Solutions
1810 West Drake Drive
Tempe, AZ 85283
stryker.com