# INSTRUCTIONS FOR USE EXETER FEMORAL AND PROXIMAL FEMORAL STEM SEALS

Howmedica Osteonics Corp.  
325 Corporate Drive  
Mahwah, NJ 07430, USA  
A subsidiary of Stryker Corporation

Telephone #: +1 201-831-5000

CE 0086

©2016 Howmedica Osteonics Corp.

## 0090-9-705 Rev AA

Refer to product label for CE Mark Status and Legal Manufacturer. The CE mark is only valid if also found on the product label.

## Labeling Symbols

Attention, see Instructions for Use

Do Not Reuse

Sterilized using Irradiation

Sterilized using Hydrogen Peroxide

Sterilized using Ethylene Oxide

Use by Date

Date of Manufacture

Legal Manufacturer

Authorized Representative in the European Community

Catalog Number





Batch Code

Serial Number

## CONTENTS

English .....	1
Français .....	3
Deutsch .....	5
Italiano .....	7
Español.....	9
Svenska .....	11
日本語.....	13
Nederlands .....	15
Suomi .....	17
Dansk .....	19
Norsk .....	21
Português.....	23
Ελληνικά.....	25
Polski .....	27
한국어 .....	29
中文.....	31

## **Instructions for Use**

### **Exeter Femoral and Proximal Femoral Stem Seals**

#### **Description**

Exeter Femoral and Proximal Femoral Stem Seals are designed for use in Femoral Cementing with the Exeter Total Hip femoral stems and SIMPLEX Surgical Bone Cement.

These instructions for use should be used in conjunction with the Exeter Operative Technique Brochure.

#### **Procedure**

##### **Preparation**

The proximal femoral seal and backing plate are positioned onto the cement gun nozzle. These are placed over the proximal end of the femur to ensure a good fit (**Fig. 1**). *If a tight fit cannot be achieved*, the other size of seal should be tried, or the sealing area of the bone at the femoral slot altered to ensure a tight seal. If a tight fit cannot be obtained after any of these measures, the use of the gun should be abandoned and the 'suck-down' technique should be used with doughy cement and vigorous finger packing.

##### **Cement Preparation**

Surgical SIMPLEX P Radiopaque or Antibiotic SIMPLEX bone cement should be mixed in a mixing bowl for 1 minute and then poured into the cement barrel which should be left to stand for 30 seconds. Three mixes may be necessary in large femora. The nozzle is fixed to the barrel, the femoral seal and backing plate positioned, and the cement gun is primed.

##### **Cementing**

Cement is introduced in a retrograde (i.e. from distal to proximal) fashion (**Fig. 2**). Typically, using Surgical SIMPLEX P Radiopaque or Antibiotic SIMPLEX bone cement and with a theatre temperature of 21°C this is usually 2 ½ to 3 minutes after commencement of mixing. As the cement reaches the upper part of the canal, the nozzle of the gun is cut off distal to the femoral seal. The cement remaining in the nozzle is ejected into the femur using the obturator.

##### **Use of Femoral Seal and Cement Pressurisation**

The proximal femoral seal is then impacted forcibly onto the upper end of the femur, and the remaining cement is slowly pumped into the femur (**Fig. 2B**). This action will maintain a steady pressurisation of the cement to overcome the bleeding pressure inside the femur.

When pressurisation is correctly performed, a steady extrusion of fat through the walls of the upper femur will be seen.

Cement injection and pressurisation is continued until the viscosity of the cement starts rising (judged by a small sample held in the hand of the surgeon). Typically, using Surgical SIMPLEX P Radiopaque or Antibiotic SIMPLEX bone cement, and with a theatre temperature of 21°C this is seldom less than 6 minutes from commencement of mixing. The femoral stem is then inserted.

### **Insertion of Stem**

The aim should be to delay stem insertion as long as possible remembering that during stem insertion interface pressures in the canal are directly related to the viscosity of the cement.

Following stem insertion, and when the final stem position has been reached, and the stem introducer removed; the stem seal and backing plate are positioned and firm pressure maintained on the top of the cement until it polymerisation is completed. (Fig. 3A, 3B, 3C). Ensure that the stem does not back out during cement polymerisation.

- Sterility cannot be assured for resterilized devices, as cleaning and re-sterilization have not been verified.

**CAUTION:** Federal law (U.S.A.) restricts this device to sale by or on the order of a licensed physician.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademark(s) or service mark(s): Exeter, Howmedica, Osteonics, Simplex, Stryker. All other trademarks are trademarks of their respective owners or holders.

Refer to product label for CE Mark Status and Legal Manufacturer. The CE mark is only valid if also found on the product label.

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