Howmedica Osteonics Corp.
Sterile Instrumentation
English

Howmedica Osteonics Corp.'s Sterile Instrumentation consists of manual surgical instruments intended for use during hip, knee, shoulder, elbow or wrist arthroplasty, or trauma surgeries.

Utilization

- Before clinical use, the surgeon should thoroughly understand all aspects of the surgical procedure and the limitations of the instrumentation.
- Howmedica Osteonics Corp.'s trial components should be used for size determination, preparation evaluation, trial reduction and range of motion evaluation, thus preserving the integrity of the actual implants and their sterile packaging.
- Radiographic templates are available to assist in the preoperative prediction of component size and style.
- Howmedica Osteonics Corp.'s Surgical Protocols provide additional procedural information.

Warnings

- Care should be taken not to cut through surgical gloves when handling any sharp-edged orthopaedic device.
- Due to different manufacturers employing different design parameters, varying tolerances, different materials and manufacturing specifications, Howmedica Osteonics Corp.'s Sterile Instrumentation should not be used to implant any other manufacturer's components with the exception of the Pelvic Alignment Level device which can be used to implant other manufacturers' hip stems. Any such use will negate the responsibility of Howmedica Osteonics Corp. for the performance of the resulting implant.
- Instruments made of non-metallic material(s) and fragments thereof may not be visible using certain forms of external imaging (e.g. x-ray) unless otherwise specified, such as radiopaque femoral head trials that are visible.
- Stryker Orthopaedics instruments are not intended to be used in or present within the MR environment. Since they are not intended to enter the MR environment, testing for safety in the MR environment has not been undertaken, and thus any potential hazards from exposure to such an environment is unknown.

Adverse Effects

- Malfunctions of instruments, including, but not limited to: weld fracture, thread gauging and material fatigue, have occurred in a small percentage of cases.
- Serious complications may be associated with any joint replacement or trauma surgery. These complications include, but are not limited to: genitourinary disorders; gastrointestinal disorders; vascular disorders, including thrombus; bronchopulmonary disorders, including emboli; myocardial infarction or death.
- Peripheral neuropathies, nerve damage, circulatory compromise and heterotopic bone formation may occur.
- Intraoperative fissure, fracture, or perforation of the bone can occur due to numerous factors including the presence of defects, poor bone stock, the use of the surgical instruments and the impact of the component into the preparation.

Sterilization

- This instrument has been sterilized by gamma radiation.
- The packaging of all sterile products should be inspected for flaws in the sterile barrier before opening. In the presence of such a flaw, the product must be assumed non-sterile.
- Care should be taken to prevent contamination of the instrument. In the event of contamination, this product must be discarded.
- If the package is opened, but the product is not used, the instrument **must not** be resterilized and must be discarded or returned to the supplier.
- Sterility cannot be assured for reesterilized devices, as cleaning and re-sterilization have not been verified.

DISPOSAL

- If a device is being returned for evaluation, please contact your local Stryker representative for shipping/handling information. If the device is not being returned to Stryker, instruments are to be disposed of in accordance with applicable laws, rules, and regulations for the disposal of biohazardous waste. Follow all guidelines for biohazardous waste in accordance with the Centers for Disease Control and Prevention guidelines as well as applicable federal/national, state, and local regulations.

CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a licensed physician.

Stryker Corporation or its divisions or its corporate affiliated entities, own, use, or have applied for the following trademark(s) or service mark(s): **Howmedica, Osteonics, Stryker.** All other trademarks are trademarks of their respective owners or holders.

Refer to product label for CE Mark Status and Legal Manufacturer. The CE mark is only valid if also found on the product label.

The following table contains a list of abbreviations that are used on Howmedica Osteonics Corp. product labeling:

Term	Abbreviation	Term	Abbreviation
Alpha Code	ALPH CDE	Neck	NK
Angle	ANG	Offset	OFFST
Degree	DEG or °	Outer Diameter	OD
Diameter	DIA	Right	RT ▶
Extra Deep	XDP	Screw Holes	SCR HLS
Extra Large	XLGE	Side	SDE
Extra Small	XSM	Size	SZE
Head	HD	Small	SM
Height	HT	Standard	STD
Inner Diameter	ID	Taper	TPR
Insert	INSR	Thickness	THKNS
Large	LG	Type	TYP
Left	◄ LFT	With	W/
Length	LNTH	Without	W/O
Medium	MED		

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Poniższe tabele zawiera listę składowych używanych w etykietkach produktów firmy Howmedica Osteonics Corp.