TRAUSON

English

Trauson (China) Medical Instrument Co., Ltd.

9 Longmen Road, Changzhou, Jiangsu, China Zip Code: 213164 Tel: +86-519-86387075 Fax: +86-519-86387379

	Medical device		Quantity
	Manufacturer		Date of manufacture
	Batch code		Do not use if package is damaged
	Do not re-use		Humidity limitation
	Non-sterile		Authorized representative in the European Community
	Consult instructions for use		Material
	Catalogue number		Caution

DESCRIPTION
The Bone Screw System is intended for use in conjunction with Trauson's products (e.g. bone plate system) for adjunct fixation. The Bone Screw System is available in various length and designs. The Bone Screw system is made from titanium alloy conforming to ISO 5832-3.

INTENDED USE
The bone screw system is indicated for fr. stabilize the fracture site, in order to maintain proper occlusion during the fracture healing. Some types of the bone screw can be used together with the bone plates.

INDICATIONS FOR USE
The Bone Screw System including Corical Screws, Cancellous Screws, Cannulated Screws, Washer, Locking Screws, and Screw Hole Insertor, is designed for fixation for fractures, fusions or osteotomies. In general, these screws are indicated for fractures requiring additional stability.

Indications include:

- Corical screw bone.
- Fracture of cortical bone.
- Fracture of cancellous bone.
- Cannulated screw and washer.
- Fracture of metaphysis of limbs and fracture of hand, foot, pelvis and sacrum.

Locking screw and filling screw

- Fracture of limbs and pelvis.

CONTRAINdicATIONS
The physician's education, training and professional judgment must be relied upon to choose the most appropriate device and treatment. Conditions presenting an increased risk of failure include:
1.Any acute or suspected latent infection or marked local inflammation in or about the affected area.

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2. Kompromize vaskülarizasyon ya da medifin adekvate blood supply to the fracture or operative site.
3. Bone stock compromised by disease, infection or prior instrumentation that can not provide adequate support and/or fixation of the fracture.
4. Malaligned sensitivity, documented or suspected.
5. Obesity. An overweight or obese patient can produce loads on the implant that can lead to failure of the fixation of the device or failure of the device itself.
6. Inadequate soft tissue coverage over the operative site.
7. Implant utilization that would interfere with anatomical structures or physiological performance.
8. Any mental or neuromuscular disorder which would create an unacceptable risk of fixation failure or complications in postoperative care.
9. Other medical or surgical conditions which would preclude the potential benefit of surgery.

INTENDED USERS
The Bone Screws System is intended to be used by licensed healthcare professionals. The healthcare professionals should be fully aware of the intended use of the products and the applicable surgical techniques and should be qualified by appropriate training methods (for example, relevant surgical residency programs). Additional user groups include nurses and reprocessing staff performing cleaning and sterilization.
PATIENT TARGET GROUP
The target group for the application of the Bone Screws System is skeletally mature patients (according to the judgment of the licensed healthcare professional) with utilization of the appropriate anatomical structures as defined in the indications.

POTENTIAL ADVERSE EVENTS AND COMPLICATIONS
In many instances, adverse results may be clinically related rather than device related. This list may not include all complications caused by the surgical procedure itself.
1.Delayed union or non-union of the fracture site.
2.These devices can break when subjected to the increased loading associated with delayed unions and/or non-unions. Internal fixation devices are load sharing devices which are intended to hold fractured bone surfaces in apposition to facilitate healing. In the event of a device failure, the fracture may be exposed and the device produced by load bearing and the patient's activity level will dictate the longevity of the device.
3.Conditions attributable to non-union, osteoporosis, osteomyelitis, diabetes, inhibited revascularization and poor bone formation can cause loosening, bending, cracking or fractures into the screw shafts or stems as these may produce premature failure of the component.
4.Improper alignment can cause a mal-union of the bone and/or bending, cracking or even breakage of the device.
5.Inadequate tissue response around the fracture site due to unstable consolidation or fixation.
6.Inadequate fixation, both deep and/or superficial.
7.Deep venous thrombosis.
8.Avascular necrosis.
9.Screw loosening due to an affected bonefracture site.
10.Substantial nerve damage may possibly occur as a result of the surgical trauma.
11.Material sensitivity reactions in patients following surgical implantation have rarely been reported, however their significance awaits further clinical evaluation.

WARNINGS AND CAUTIONS
1.**NON-STERILE PRODUCTS:** The implants and the instruments to be used must be cleaned carefully in a clean environment and sterilized thoroughly with appropriate temperature and pressure before they are used.
2.**IMPLANT SELECTION AND SIZING:** The correct application of the fracture fixation equipment is extremely important. Failure to use the appropriate appliance for the fracture condition may accelerate clinical failure. Failure to use the proper component to maintain adequate blood supply and provide rigid fixation may result in loosening, bending, cracking or fracture of the device and/or bone. The correct implant size should be determined by evaluating the patient's height, weight, height, weight, functional demands and anatomy. Every implant must be used in the correct anatomic location, consistent with accepted standards of internal fixation.
3.**PATIENT SELECTION:** Proper patient selection is critical to the success of the procedure. Only patients who satisfy the criteria set forth under the INDICATIONS section of this document AND who do not have any of the conditions set forth under the CONTRAINDICATIONS section of this document should be considered for surgery using the Bone Screw System.
4.**HANDLING:** Implant components should be handled and stored appropriately to protect them from unintended damage. The surgeon should avoid introducing air or moisture into the screw shafts or stems as these may produce premature failure of the component.
5.**PATIENT EDUCATION:** Preoperative instructions to the patient are essential. The patient should be made aware of the limitations of the implant and potential risks of the surgery. The patient should be instructed to limit postoperative activity, as this will reduce the risk of bent, broken or loose implant components. The patient must be made aware that implant components may bend, break or loosen even though restrictions in activity are followed.
6.**SINGLE USE ONLY:** Reuse of a single use device that comes in contact with blood, tissue, or other body fluids may lead to patient or user injury. Possible risks associated with reuse of a single use device include, but are not limited to, mechanical failure, material degradation, potential leachables, and transmission of infectious agents. Reesterilization may result in the loss of efficacy of the device and/or disposal.

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