

# Prime TC<sup>®</sup>: an ergonomic solution

## Situation

A recent study has shown that work-related musculoskeletal injuries are particularly prevalent in the healthcare profession.<sup>1</sup>

**The highest incidence of injuries are in caregivers who manually handle patients, including those who perform patient transport tasks.**<sup>2</sup> Traditionally, a heavy emphasis has been placed on designing transport devices that maximize patient comfort and minimize costs. However, little consideration has been given to design concepts that minimize caregiver burden.<sup>3</sup>

## Study

The purpose of this study was to compare differences in operator trunk and upper extremity muscle activity and joint angles when operating the ergonomically designed Prime TC transport chair and two other seated transport devices.

**This study was completed at the VA Pittsburgh Healthcare System in conjunction with the Human Engineering Research Laboratories and Department of Rehabilitation Science and Technology**

## Why is wrist flexor activity important?

- These muscles are present in the forearm and are often prone to injury and strain as their function is to facilitate movement of the wrist.
- May reduce musculoskeletal strain on forearm muscles and protect the wrists from developing cumulative trauma injuries.<sup>4</sup>

## Stryker solutions: Ergonomic features specific to Prime TC Transport Chair include:

Adjustable armrests and footrests for ease in patient egress/ingress

Highly maneuverable frame with anti-tip wheels



Vertically oriented push handles that position caregivers elbows at an ergonomic 90-degree angle

One touch central brake pedal that eliminates bending when initiating wheel locks

**Findings**

The traditional depot-style wheelchair required more effort to push, increasing strain on the caregivers during mobility tasks. Additionally, subjects used significantly less wrist flexor activity when using the Prime TC transport chair compared to the standard transport chair.<sup>4</sup>

These results demonstrate the importance and efficacy of ergonomic design features in transport devices in promoting better caregiver postures and helping to reduce operator musculoskeletal strain<sup>4</sup>. Additionally, by customizing Prime TC transport chair to caregivers, a more enjoyable transport experience can be provided to both the caregiver and the patient.

**In conclusion, the implementation of ergonomic seated transport devices in hospitals and clinics has the potential to mitigate incidence of work-related musculoskeletal injury.<sup>4</sup>**



**References**

1. Oranye O, Wallis B, Roer K. Do Personal Factors or Types of Physical Tasks Predict Workplace Injury? Workplace Health Safety 64: 141-151. April 2016.
2. Kothiyal K, Yuen TW. Muscle strain and perceived exertion in a patient handling with and without a transferring aid. Occupational Ergonomics (4): 185-197. 2004.
3. Lee SY, Kim SC. Comparison of Shoulder and Back Muscle Activation in Caregivers According to Various Handle Heights. J Phys Ther Sci: 1231-1233. Oct 2013.
4. Chemini A, Kulich H, Bass S, Vijayvargiya A, Koontz A. Muscle activation for three different patient transport chairs on ramps and flat ground. VA Pittsburgh Healthcare System. 30th Annual International Occupational Ergonomics and Safety Conference; July 2018.