


# Vision

## Clinical Dashboard



# Vision Clinical Dashboard

## Personalized fall prevention

**\$32,215**

A fall in the hospital can cost up to \$32,215 per injury.<sup>2</sup>

**Nearly 1,000,000**

patients in the United States experience a fall while being treated in the hospital.<sup>1</sup>


### Advanced fall prevention to help drive better outcomes

We know keeping patients safe is your #1 responsibility. As a trusted partner in fall prevention, we're focused on developing innovations that can help reduce falls and injuries from falls.

### Introducing Vision

Vision is a patient-centric clinical dashboard that gives caregivers increased visibility to bed compliance, falls protocols and bed exit alarm activity. It takes a patient's fall risk score from the EHR\* and associates it with appropriate protocols to help ensure a safer patient experience and an easier workflow for caregivers.

\*No PHI information is transmitted.





### **It's all about the patient**

Each patient's fall risk score is different and can change throughout his or her hospital stay. That's why Vision's dynamic and personalized advanced fall prevention monitoring is so critical. Featuring configurable rules for safe bed protocol compliance, Vision displays near real-time, active changes based on a patient's fall-risk – allowing you to provide enhanced care as a patient's fall risk increases or decreases.

### **Wireless bed notifications**

Violations are generated by Vision when bed protocols aren't being met (e.g., low height, siderails up or down, brake set, bed exit on, and bed exit zone) independently of the nurse call cable.<sup>3</sup>

### **Safe bed monitoring**

Hospital-specific safe bed protocol compliance automatically adjusts as the patient's EHR fall risk status changes. Vision will also display the beds head of bed compliance, if the mattress is unplugged, a not in bed timer and when nurse call has been activated.

### **Customizable display**

Customize colors, icons, violation alerts and text based on hospital preference.

### **Color coded tiles**

Enabling easy recognition of fall risk status.





## Bed exit notifications for personalized care

Hospitals can set pre-determined safe bed configuration requirements that align with fall protocols. If the safe bed configuration isn't met (e.g., the bed exit alarm is not armed for a high risk patient), the notification will be displayed on the Vision dashboard.

# Vision customization template

## Vision\* customization template - Beds

**Low fall risk (level 1)**

Included in dashboard | Yes ☐ No ☐

Bed parameters to be monitored

Score range:

☐ Head siderails in the up position

☐ Right rail

☐ Left rail

☐ Both

☐ Either

☐ Foot siderails in the up position

☐ Right rail

☐ Left rail

☐ Both

☐ Either

☐ Bed

☐ Bed Exit armed

☐ Zone 1 - Loose restrictive

☐ Zone 2 - Medium restrictive

☐ Zone 3 - Most restrictive

☐ Low bed height

☐ Brake on

☐ Nurse call push\*\*

☐ HOB Not 30 (degrees symbol)

☐ Not in bed timer

☐ Mattress Unplugged (InTouch or ZMX only)

**Medium fall risk (level 2)**

Included in dashboard | Yes ☐ No ☐

Bed parameters to be monitored

Score range:

☐ Head siderails in the up position

☐ Right rail

☐ Left rail

☐ Both

☐ Either

☐ Foot siderails in the up position

☐ Right rail

☐ Left rail

☐ Both

☐ Either

☐ Bed

☐ Bed Exit armed

☐ Zone 1 - Loose restrictive

☐ Zone 2 - Medium restrictive

☐ Zone 3 - Most restrictive

☐ Low bed height

☐ Brake on

☐ Nurse call push\*\*

☐ HOB Not 30 (degrees symbol)

☐ Not in bed timer

☐ Mattress Unplugged (InTouch or ZMX only)

**High fall risk (level 3)**

Included in dashboard | Yes ☐ No ☐

Bed parameters to be monitored

Score range:

☐ Head siderails in the up position

☐ Right rail

☐ Left rail

☐ Both

☐ Either

☐ Foot siderails in the up position

☐ Right rail

☐ Left rail

☐ Both

☐ Either

☐ Bed

☐ Bed Exit armed

☐ Zone 1 - Loose restrictive

☐ Zone 2 - Medium restrictive

☐ Zone 3 - Most restrictive

☐ Low bed height

☐ Brake on

☐ Nurse call push\*\*

☐ HOB Not 30 (degrees symbol)

☐ Not in bed timer

☐ Mattress Unplugged (InTouch or ZMX only)

**Choose color indicators from the table.**

#####	#DCE6B0	#A0C8CC	#757678	#2A2A2A	#000000
#A02087	#FF0000	#FF0000	#000000	#000000	#000000
#000000	#000000	#000000	#000000	#000000	#000000
#000000	#000000	#000000	#000000	#000000	#000000
#000000	#000000	#000000	#000000	#000000	#000000
#000000	#000000	#000000	#000000	#000000	#000000

**Dashboard visualizations**

Options

Tile configuration

☐ Words

☐ Icon and words

Violation configuration

☐ General violation

☐ Specific violation

Default configuration

If patient has not had a fall risk assessment

When do you want to default to:

☐ Fall risk incomplete

☐ High fall risk rules

Tile enhancements

Display maintenance icon

☐ Yes

☐ No

Display nurse call push\*\*\*

☐ Yes

☐ No

**Tile configurations**

Options

Tile configuration

☐ Words

☐ Icon and words

Violation configuration

☐ General violation

☐ Specific violation

Default configuration

If patient has not had a fall risk assessment

When do you want to default to:

☐ Fall risk incomplete

☐ High fall risk rules

Tile enhancements

Display maintenance icon

☐ Yes

☐ No

Display nurse call push\*\*\*

☐ Yes

☐ No

**Violation configurations**

Options

Tile configuration

☐ Words

☐ Icon and words

Violation configuration

☐ General violation

☐ Specific violation

Default configuration

If patient has not had a fall risk assessment

When do you want to default to:

☐ Fall risk incomplete

☐ High fall risk rules

Tile enhancements

Display maintenance icon

☐ Yes

☐ No

Display nurse call push\*\*\*

☐ Yes

☐ No

**Specific violations**

Prioritize violation importance

1=highest priority - 4=lowest priority.

If the tile is not desired

☐ Brake not set

☐ Bed not lowered

☐ Siderails not raised

☐ Zone incorrect

☐ Mattress unplugged

☐ Not in bed timer

**Optional alternative violation language**

Bed not lowered \_\_\_\_\_ Fall risk incomplete \_\_\_\_\_ EHR \_\_\_\_\_

Brake not set \_\_\_\_\_ Not in bed \_\_\_\_\_ ADT \_\_\_\_\_

Siderails not raised \_\_\_\_\_ Low \_\_\_\_\_ Facility name \_\_\_\_\_


Zone incorrect \_\_\_\_\_ Medium \_\_\_\_\_ Signature \_\_\_\_\_

Unplugged/no connection \_\_\_\_\_ High \_\_\_\_\_ Date \_\_\_\_\_

# Stryker's Financial Guarantee

## 40% fall reduction. Guaranteed.

Through Stryker's Financial Guarantee program, we guarantee a 40% reduction in bed-related falls occurring on Stryker beds equipped with iBed Wireless Technology and Vision\*.



\*Facility must meet eligibility requirements as outlined in the iBed Wireless Fall Reduction Contract. If a reduction of 40% of bed-related falls is not recognized after one full year (based on delivery date) of product use by the facility, Stryker will refund a portion of the license fee for Vision.

## Vision

Compatible with wireless ProCuity models, S3 MedSurg and InTouch Critical Care beds.



## Customer Support Services

### Technical Support

Our Technical Support comprises a team of professionals available to help with your product needs. Contact via phone at 1 800 STRYKER or email at [medicaltechnicalsupport@stryker.com](mailto:medicaltechnicalsupport@stryker.com).

### Stryker's ProCare Services

Every day, you count on your medical equipment to perform at its best. With ProCare Services, our people help to ensure your equipment is ready to perform when it's needed and make it easier to get the most from your investment. When an issue arises, we promise that we'll work to solve it — performing repairs quickly and correctly in accordance with the terms of your service agreement.

ProCare isn't just a service program. ProCare is a partnership you can count on to give you one less thing to worry about, and one more reason to feel confident you're doing all you can for your clinicians, staff and patients.

All ProCare agreements provide:

- Stryker-authorized service representative
- Stryker-direct factory parts
- Two-hour callback response time
- Fixed service costs up front
- Increased uptime
- 24-72 hour equipment turnaround time\*\*

\*\*Based on the provisions of the Service Agreement and the location of the product.

### Flex Financial Program

Our financial programs provide a range of smart alternatives designed to fit your organization's needs. We offer flexibility beyond a cash purchase with payment structures that can be customized to meet budgetary needs and help to build long-term financial stability. Contact your account manager for more information.

### References

1. Health Research & Educational Trust. Preventing patient falls: A systematic approach from the Joint Commission Center for Transforming Healthcare project. 2016 Oct. Retrieved from <http://www.hpoe.org/Reports-HPOE/2016/preventing-patient-falls.pdf>
2. Dykes, Patricia C., et al. "Evaluation of a Patient-Centered Fall-Prevention Tool Kit to Reduce Falls and Injuries." JAMA Network Open, vol. 3, no. 11, 2020, <https://doi.org/10.1001/jamanetworkopen.2020.25889>.
3. Per UL1069 hospitals must still provide a wired connection for patients to send nurse call request.

Vision 3.5. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: ProCuity, InTouch, S3, iBed, Flex Financial, ProCare, Stryker, Vision. All other trademarks are trademarks of their respective owners or holder.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

AC-VISION-BROC-1330607\_REV-0\_en\_us © 2024 Stryker

3800 E. Centre Avenue  
Portage, MI 49002 USA  
t: 269 329 2100  
toll free: 800 327 0770

[stryker.com](https://www.stryker.com)