

LIFEPAK® CR2 AED

Establish and confirm a Wi-Fi connection

How to perform two-button push (device must be within Wi-Fi network range):

1. Open the lid
2. After the voice prompts start, immediately press and hold the **LANGUAGE** and **CHILD MODE** buttons simultaneously until you hear DEVICE READY
3. The device will audibly give you the status and initiate a Wi-Fi connection. Immediately close the lid. Listen for the following voice prompts:
 - Wi-Fi connection initiated, followed by "CONNECTION IN PROGRESS."
There will be a pause while the AED connects.
 - "CONNECTION ESTABLISHED"
 - Once a Wi-Fi connection is established, the AED will update your online account with its current status. The AED beeps every few seconds during the update process.
 - When all updates are complete, you will hear "COMMUNICATION COMPLETE"
 - You will hear "POWERING OFF"



Actions after deployment

1. Perform two-button push (Wi-Fi location)
Flip card over for instructions
*If CPR data file is not uploaded and AED is used again, the case data file will not be accessible.
2. Order replacement pads (automated when two-button push is completed, or contact your Stryker AED program administrator)
3. Optional: user can remove battery while waiting for pad replenishment
4. Install pad tray upon receipt of replacement pads
5. If battery was removed in step 3, please remember to reinsert battery after replenishment pad installation is complete
6. Complete two-button push BEFORE placing AED back into service

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: LIFEPAK, Stryker. All other trademarks are trademarks of their respective owners or holders. The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

M0000010096 REV AA

Copyright © 2022 Stryker

Support

If your device was provided as part of the multi-state AED initiative funded by the Helmsley Charitable Trust, contact your Stryker AED program administrator for connectivity, troubleshooting support, or to order your post-event replenishment supplies.

Phone: 800 834 9743

(follow prompts to connect to your Helmsley State Stryker AED program administrator)

Email: helmsleyaedsupport@stryker.com

To order all other supplies, contact your local Stryker account manager.