

# LUCAS<sup>®</sup> 3

## chest compression system

### Suggestions for performance evaluation



Name: \_\_\_\_\_ Unit: \_\_\_\_\_

Reviewer: \_\_\_\_\_ Date: \_\_\_\_\_

**This performance evaluation is a suggested basic assessment of a trainee’s ability to operate the LUCAS 3 chest compression system. This is a limited evaluation and does not cover all the information and skills required to operate this device safely.** It is designed to be completed after viewing the appropriate in-service video or observing an equipment demonstration given by a qualified instructor. Refer to the LUCAS 3 device Instructions for Use for complete information.

Performance criteria	Complete	Incomplete	Comments
1. Opens case; presses ON/OFF for 1 second to start self-test and power up the LUCAS 3 device.			
2. Removes back plate from the case.			
3. Verbalizes “stop manual CPR” to other rescuer (instructor may “act” as other rescuer).			
4. Places back plate under patient, immediately below the patient’s armpits.			
5. Verbalizes “resume manual CPR” to other rescuer.			
6. Lifts upper part of the LUCAS 3 device from case and pulls once on the release rings to check that the claw locks are open, then lets go of rings.			
<b>Instructor is prepared and starts timer when manual compressions stop. Maintain manual CPR as long as possible.</b>			
7. Connects the upper part to the back plate, starting on the side closest to user. Listens for click.			
8. Pulls up once to ensure attachment.			
9. Positions the suction cup immediately above the end of the sternum in the center of the chest. Adjusts if necessary.			
10. Pushes the suction cup down using two fingers with the LUCAS 3 device in the ADJUST mode until the pressure pad inside suction cup touches the patient’s chest without compressing the chest.			
11. Pushes PAUSE to lock start position.			
12. Presses ACTIVE (30:2) or (continuous). <b>Instructor stops timer: _____ seconds (goal: ≤ 20 seconds).</b>			
13. Fully extends the stabilization strap at the buckles.			
14. Places support cushion under patient’s neck.			
15. Connects buckles on support cushion straps to device straps.			
16. Tightens support cushion straps firmly.			
17. Visually checks for proper position of suction cup, adjusting if needed.			

## General performance guidelines

### Consider the following issues when determining competency:

- Was application and operation performed in a safe manner?
- Were compressions started within 20 seconds from the time manual compressions were stopped to connect the upper part of the LUCAS device?
- Able to verbalize contraindications and side effects?

### Suggested advanced performance scenarios

- Change battery during patient use

### Optional activities

- Apply defibrillation pads, perform defibrillation
- Perform ventilations with bag, valve and mask
- Perform ventilations with advanced airway

## Evaluation and action plan

[Empty box for evaluation and action plan]

Refer to operating instructions for complete directions for use indications, contraindications, warnings, cautions, and potential adverse events.

For further information, please contact Stryker at 800 442 1142 (U.S.), 800 668 8323 (Canada) or visit our website at [stryker.com](http://stryker.com)

## Emergency Care

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: LUCAS, Stryker. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

M0000009036 REV AA  
Copyright © 2022 Stryker



### Manufactured by:

Jolife AB  
Scheelevägen 17  
Ideon Science Park  
SE-223 70 Lund  
Sweden

### Distributed in Canada by:

Stryker Canada  
2 Medicorum Place  
Waterdown, Ontario  
L8B 1W2  
Canada  
Toll free 800 668 8323