stryker

LIFEPAK[®] CR2 AED Trainer

Data sheet

Training improves readiness.

- Dynamic and versatile training tool
- Same intuitive, simple operation as the LIFEPAK CR2 automated external defibrillator (AED)
- Fully customizable scenarios
- No remote to manage
- ClearVoice[™] technology

The LIFEPAK CR2 AED Trainer provides realistic training in the use of the LIFEPAK CR2 AED without the actual charge and discharge of electrical energy. Designed to simulate the appearance of the LIFEPAK CR2, the Trainer simulates its actual operation using the same simple user interface. Pre-set demos can be programmed for fully automatic or semi-automatic mode to enable a quick demonstration of the key LIFEPAK CR2 features, as well as the benefits of cprINSIGHT^{**}analysis technology.

Device realism

The Trainer guides users through simulated analysis, energy delivery and prompted CPR intervals with CPR coaching, even automatically delivering email notifications via Wi-Fi® through LIFELINKcentral[™] AED program manager when the device is in use and the pads are placed. Like the LIFEPAK CR2, the Trainer uses ClearVoice[™] technology to enable voice prompts to be heard clearly even in noisy environments for clear instruction.

Programmable scenarios

All of the functions of the Trainer are programmable to simulate the full range of scenarios users might



encounter during "real life" use of the LIFEPAK CR2 AED. The Trainer provides two fully configurable scenarios that are completely customized using a setup screen and can be updated on the go for maximum flexibility.

Email notifications

The LIFEPAK CR2 Trainer demo app can be used to set up automatic email notifications or send manual email notifications via Bluetooth® from any Windows®, Android® or IOS® mobile device. This app simulates the automatic email notifications provided when a LIFEPAK CR2 is managed by LIFELINKcentral AED program manager.

Simple and affordable

The Trainer detects when pressure is placed on the OUIK-STEP[™] training electrodes and automatically advances through the scenario, eliminating the need for a remote control. A pause button on the unit enables stopping an in-progress scenario as needed. Plus, the trainer is cost-effective. Not only is it powered by four standard D-cell batteries, but the electrodes are affordable and reusable, simulating the rapid-access OUIK-STEP defibrillation ECG electrodes unique to the LIFEPAK CR2 AED.

Specifications

General

For use in temperature-controlled areas protected from precipitation. Capable of withstanding heavy usage and prolonged periods of operation.

Physical characteristics

Height: 10.7 in (27.2 cm), with handle Width: 8.6 in (21.8 cm) Depth: 3.4 in (8.6 cm) Weight: 3.9 lb (1.8 kg), with batteries

Batteries

6 volts/4 cells, range 3.3–6.6 Vdc. Four nonchargeable D-cell alkaline batteries. New batteries will provide no less than 100 hours of continuous operation at volume 4.

Buttons

Language, power, child mode, shock

Display

Power indicator: Off when device is off; On when device is in use or setup mode; flashes when device is in standby mode (lid is closed)

Scenario options

Device type: Semi, fully-automatic Number of events: 1, 2, 3, 4, repeat Shock advised: Yes, no CPR metronome: Hands-only, 30:2 CPR time: Short, 60s, 120s, 180s cprINSIGHT technology: On, Off No CPR detected: On, Off Motion detected: On, Off Bilingual: On, Off Second language: All available languages Check breathing: On, Off

Physio-Control is now part of Stryker.

For further information, please contact Stryker at 800.442.1142 or visit our website at www.strykeremergencycare.com

Physio-Control Manufacturing 11811 Willows Road NE Redmond, WA 98052 www.physio-control.com

Customer Support

P. O. Box 97006 Redmond, WA 98073 Toll free 800 442 1142 Fax 800 426 8049

Physio-Control, Inc., 11811 Willows Road NE, Redmond, WA 98052 USA