

LIFENET[®]

Alert

Incoming alert flashing?

Double-click flashing line to view

To view:



Select report and either double click or click "View"

- Use mouse wheel to zoom
- Close viewer when done

To forward:



Select report and click "Forward"

- Add comments if needed
- Select receiving target(s)

To E-mail:



Select report and click "Send E-mail"

- Confirm attachment
- Add comments if needed
- Select recipient(s)

To Print:



Select report and click "Print"

Please refer to online Help for complete instructions for use.

IUSA Rx Only

Physio Control is now part of Stryker
For further information, please contact
Physio-Control at 800.732.3081 or visit
our website at www.lifenetsystems.com

©2018 Physio-Control, Inc. Redmond, WA
GDR 3305591_C