

LIFEPAK® 20e

defibrillator/monitor



USA Rx Only

Copyright © 2020 Stryker
GDR 3302911_C

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo. Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: LIFEPAK, QUIK-COMBO, Stryker. All other trademarks are trademarks of their respective owners or holders.

For a brief summary of important disclosure and safety information, please visit our website at strykeremergencycare.com/disclosure-and-safety-information/

stryker


LIFEPAK® 20e

defibrillator/monitor

Quick Reference for clinical use


Manual operation

Defibrillation

1. Press **ON**.*
2. Apply QUIK-COMBO® therapy electrodes or hard paddles with conductive gel.
3. Select **ENERGY SELECT**. Adjust up or down as needed.
4. Press **CHARGE**.
5. Stand clear.
6. Press **SHOCK**  to deliver energy.

*To change from AED operation to manual operation press gray manual tab and open door or press **ENERGY SELECT**.

Synchronized cardioversion


1. Press **ON**.
2. Attach ECG cable and apply ECG electrodes.
3. Apply QUIK-COMBO therapy electrodes or hard paddles with conductive gel.
4. Select lead II or the lead with the greatest QRS complex amplitude.
5. Press **SYNC**. Confirm placement of triangle sense marker ▼ appears near the middle of each QRS complex.
6. Press **ENERGY SELECT**. Adjust up or down as needed.
7. Press **CHARGE**.
8. Stand clear.
9. Press and hold **SHOCK**  until energy is delivered.

Pacing operation

1. Press **ON**.
2. Attach ECG cable and apply ECG electrodes (ensure adequate spacing between ECG and QUIK-COMBO electrodes).
3. Apply QUIK-COMBO therapy electrodes.
4. Press **PACER**.
5. Press **RATE** to adjust up or down as needed.
6. Press **CURRENT** and adjust up until electrical capture is achieved.
7. Confirm mechanical capture.
8. As needed, press **PAUSE** and hold to confirm underlying rhythm.

AED operation

Confirm therapy cable is not connected to test plug before pressing **ON** button.

1. Press **ON**. Follow audio/visual instructions.
2. Press **ANALYZE** and stand clear.
3. Press **SHOCK**  when directed to deliver energy.

CPR Metronome

The CPR Metronome automatically provides audible compression “tocks” and ventilation prompts or tones during CPR.

AED mode: The Metronome turns on automatically if configured on. The default C:V ratio is Adult-no airway 30:2 because most patients in cardiac arrest are adults who have an initially unsecured airway.

Manual mode: To activate, select the CPR metronome icon using the SPEED DIAL.

Select the desired Age-Airway setting using the SPEED DIAL.

Adult—No Airway 30:2 (default)

Adult—Airway 10:1

Youth—No Airway 15:2

Youth—Airway 10:1

To stop the Metronome, select Stop Metronome in the CPR Metronome menu.

See device Operating Instructions for complete directions for use, indications, contraindications, warnings, precautions and potential adverse events.