®

stryker

HemoPore®

Hemostatic Fragmentable  
Nasal Dressing

REF


5400-020-208

5400-222-208S

Instructions For Use

R<sub>x</sub> ONLY

STERILE EO

## Introduction

This instructions for use manual is the most comprehensive source of information for the safe and effective use of your product. This manual may be used by in-service trainers, physicians, nurses, surgical technologists, and biomedical equipment technicians. Keep and consult this reference manual during the life of the product.

The following conventions are used in this manual:

- A **WARNING** highlights a safety-related issue. ALWAYS comply with this information to prevent patient and/or healthcare staff injury.
- A **CAUTION** highlights a product reliability issue. ALWAYS comply with this information to prevent product damage.
- A **NOTE** supplements and/or clarifies procedural information.

For additional information, especially safety information, or in-service training, contact your Stryker sales representative or call Stryker customer service. Outside the US, contact your nearest Stryker subsidiary.

## Indications For Use

HemoPore is intended for use in patients undergoing nasal/sinus surgery as a temporary wound dressing. HemoPore functions as a topical hemostatic aid to control mild bleeding by tamponade effect, blood absorption, platelet activation and aggregation. It acts as an adjunct to aid in the natural healing process as a space occupying stent to separate and support tissues. It prevents adhesions and minimizes edema. HemoPore is indicated for use as a nasal packing to treat epistaxis.

## Description

The dressing is a biodegradable, polyurethane foam that fragments within several days. The nasal dressing drains from the nasal cavity via natural mucus flow and/or irrigation.

## Contraindications

HemoPore should not be used on patients who have known allergies to shellfish.

## Adverse Reactions

Possible adverse reactions for the nasal dressing include, but are not limited to, the following:

- Wound infection
- Allergy

Device dislocation

- Recurrent or persistent bleeding
- Toxic shock syndrome (rare)

## User/Patient Safety

### WARNINGS:

- Only trained and experienced healthcare professionals should use this product.
- Before using this product, read and understand the instructions. Pay particular attention to WARNING information. Become familiar with the product prior to use.
- Sterile only if the package is unopened and undamaged. DO NOT use if the package is open or damaged.
- DO NOT use this product after its use-by date.
- DO NOT reuse, reprocess, or re-package this product. This product is intended for a single use only. This product may not withstand chemical, chemical vapor, or high-temperature sterilization reprocessing. Design features may make cleaning difficult. Reuse may create a serious risk of contamination and may compromise the structural integrity of the product resulting in operational failure. Critical product information may be lost if the product is re-packaged. Failure to comply may lead to infection or cross-infection and result in patient and/or healthcare staff injury.
- DO NOT use this product with ointments or creams.
- HemoPore contains material derived from shellfish.
- Avoid soaking or impregnation with solutions prior to insertion.

### CAUTIONS:

- Extended exposure to temperatures above 40 °C [104 °F] may damage the product.
- DO NOT expose this product to organic solvents (for example, acetone).

## Definitions

SYMBOL	DEFINITION
R <sub>x</sub> ONLY	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.
	Consult instructions for use
	General warning sign
	Catalog number
	Batch code
	Sterilized using ethylene oxide
	Do not use if package is damaged
	Do not re-use
	Use-by date
	Keep dry
	Upper limit of temperature
	Manufacturer
	Do not re-sterilize
	Keep away from sunlight

## Instructions

1. Applying sterile protocol, open the blister pack and remove the dressing.
2. Cut the dressing to the desired size, if required.
3. Using forceps, gently insert the nasal dressing into the nasal cavity until it is completely inside the cavity.
4. Position the dressing within the cavity as required.

### NOTES:

- Maximum use: two devices per procedure.
- The dressing will absorb surrounding fluids and/or blood.
- Stryker recommends that the patient starts saline irrigation (nose and sinus) 2-3 times a day (or as prescribed) 48 hours after insertion/surgery.

## Storage and Handling

Store in a dark, dry place at or below 8 °C [46 °F].

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