®

stryker


NasoPore® Fragmentable Nasal Dressing

REF	5400-010-004
	5400-010-008
	5400-020-004
	5400-020-008
	5400-030-004
	5400-030-008
	5400-212-004S
	5400-222-004S

Instructions For Use

R_x ONLY

STERILE EO

Introduction

This instructions for use manual is the most comprehensive source of information for the safe and effective use of your product. This manual may be used by in-service trainers, physicians, nurses, surgical technologists, and biomedical equipment technicians. Keep and consult this reference manual during the life of the product.

The following conventions are used in this manual:

- A **WARNING** highlights a safety-related issue. ALWAYS comply with this information to prevent patient and/or healthcare staff injury.
- A **CAUTION** highlights a product reliability issue. ALWAYS comply with this information to prevent product damage.
- A **NOTE** supplements and/or clarifies procedural information.

For additional information, especially safety information, or in-service training, contact your Stryker sales representative or call Stryker customer service. Outside the US, contact your nearest Stryker subsidiary.

Indications For Use

Nasopore is a fragmentable nasal dressing and is indicated for use in patients undergoing nasal/sinus surgery as a space occupying stent to separate and prevent adhesions between mucosal surfaces; to help control minimal bleeding following surgery or nasal trauma by tamponade effect and blood absorption.

Description

The dressing is a biodegradable, polyurethane foam that fragments within several days. The nasal dressing drains from the nasal cavity via natural mucus flow.

Contraindications

This product is contraindicated for surgical implantation.

Adverse Reactions

Possible adverse reactions for the nasal dressing include, but are not limited to, the following:

- Wound Infection
- Toxic Shock Syndrome
- Allergy

User/Patient Safety

WARNINGS:

- Only trained and experienced healthcare professionals should use this product.
- Before using this product, read and understand the instructions. Pay particular attention to WARNING information. Become familiar with the product prior to use.
- Sterile only if the package is unopened and undamaged. DO NOT use if the package is open or damaged.
- DO NOT use this product after its use-by date.
- DO NOT reuse, reprocess, or repackage this product. This product is intended for a single use only. This product may not withstand chemical, chemical vapor, or high-temperature sterilization reprocessing. Design features may make cleaning difficult. Reuse may create a serious risk of contamination and may compromise the structural integrity of the product resulting in operational failure. Critical product information may be lost if the product is repackaged. Failure to comply may lead to infection or cross-infection and result in patient and/or healthcare staff injury.
- DO NOT use this product with ointments or creams.
- This product is manufactured in a facility that also processes modified materials originating from crustaceans.

CAUTIONS:

- Extended exposure to temperatures above 40 °C [104 °F] may damage the product.
- DO NOT expose this product to organic solvents (for example, acetone).

Definitions

SYMBOL	DEFINITION
R_x ONLY	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.
	Consult instructions for use
	General warning sign
REF	Catalog number
LOT	Batch code
STERILE EO	Sterilized using ethylene oxide
	Do not use if package is damaged
	Do not re-use
	Use-by date
	Keep dry
	Upper limit of temperature
	Manufacturer
	Do not resterilize

Instructions

1. Applying sterile protocol, open the blister pack and remove the dressing.
 2. Cut the dressing to the desired size, if required.
- NOTE:** For patients with dry sinuses, you may moisten the nasal dressing with saline before insertion.
3. Using forceps, gently insert the nasal dressing into the nasal cavity until it is completely inside the cavity.
 4. Position the dressing within the cavity as required.

NOTES:

- The dressing will absorb surrounding fluids and/or blood.
- Stryker recommends that the patient starts saline irrigation (nose and sinus) 2-3 times a day (or as prescribed) 48 hours after insertion/surgery.

Storage and Handling

Store in a dark, dry place at or below 4 °C [39 °F].

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