®

stryker


OtoPore® Fragmentable Ear Dressing

REF 5400-010-000
5400-020-000
5400-020-100

Instructions For Use

R_x ONLY

STERILE EO

ENGLISH (EN)
FRANÇAIS (FR)

Introduction

This instructions for use manual is the most comprehensive source of information for the safe and effective use of your product. This manual may be used by in-service trainers, physicians, nurses, surgical technologists, and biomedical equipment technicians. Keep and consult this reference manual during the life of the product.

The following conventions are used in this manual:

- A **WARNING** highlights a safety-related issue. ALWAYS comply with this information to prevent patient and/or healthcare staff injury.
- A **CAUTION** highlights a product reliability issue. ALWAYS comply with this information to prevent product damage.
- A **NOTE** supplements and/or clarifies procedural information.

For additional information, especially safety information, or in-service training, contact your Stryker sales representative or call Stryker customer service. Outside the US, contact your nearest Stryker subsidiary.

Indications For Use

Otopore is a fragmentable ear packing and is indicated for use in patients undergoing ear surgery as a space occupying stent to separate and prevent adhesions between mucosal surfaces; to help control minimal bleeding following ear surgery by tamponade effect and blood absorption.

Description

The dressing is a biodegradable, polyurethane foam that fragments within several days. The ear dressing falls out of the outer ear or leaves the body by the natural mucus flow.

Contraindications

This product is contraindicated for surgical implantation.

Adverse Reactions

Possible adverse reactions for the ear dressing include, but are not limited to, the following:

- Wound Infection
- Allergy

User/Patient Safety

WARNINGS:

- Only trained and experienced healthcare professionals should use this product.
- Before using this product, read and understand the instructions. Pay particular attention to WARNING information. Become familiar with the product prior to use.
- Sterile only if the package is unopened and undamaged. DO NOT use if the package is open or damaged.
- DO NOT use this product after its use-by date.
- DO NOT reuse, reprocess, or repackage this product. This product is intended for a single use only. This product may not withstand chemical, chemical vapor, or high-temperature sterilization reprocessing. Design features may make cleaning difficult. Reuse may create a serious risk of contamination and may compromise the structural integrity of the product resulting in operational failure. Critical product information may be lost if the product is repackaged. Failure to comply may lead to infection or cross-infection and result in patient and/or healthcare staff injury.
- This product is manufactured in a facility that also processes modified materials originating from crustaceans.

CAUTIONS:

- Extended exposure to temperatures above 40 °C [104 °F] may damage the product.
- DO NOT expose this product to organic solvents (for example, acetone).

Definitions

| SYMBOL | DEFINITION |
|---------------------|---|
| R _x ONLY | Caution: Federal law (USA) restricts this device to sale by or on the order of a physician. |
| | Consult instructions for use |
| | General warning sign |
| | Catalog number |
| | Batch code |
| | Sterilized using ethylene oxide |
| | Do not use if package is damaged |
| | Do not re-use |
| | Use-by date |
| | Keep dry |
| | Upper limit of temperature |
| | Manufacturer |
| | Do not re-sterilize |

Instructions

1. Applying sterile protocol, open the blister pack and remove the dressing.
2. Cut the dressing to the desired size, if required.
3. Using forceps, gently insert the dressing into the outer or middle ear.
4. Position the dressing as required.

NOTE: The dressing will absorb surrounding fluids and/or blood.

Storage and Handling

Store in a dark, dry place at or below 4 °C [39 °F].

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