
Stryker ENT Navigation System
 Stryker ØNH-navigationsssystem
 Stryker KNO-navigatiesysteem
 Stryker KNK-navigointijärjestelmä
 Système de navigation ORL Stryker
 Stryker ENT-Navigationsssystem
 Σύστημα πλοήγησης ΩΡΛ της Stryker
 Sistema di navigazione otorinolaringoiatrico Stryker
 Stryker ENT 내비게이션 시스템
 Stryker navigasjonssystem
 Stryker ENT Navigation System
 Sistema de navegação ORL Stryker
 Sistema de navegación para ORL de Stryker
 Stryker ÖNH-navigationsssystem
 Stryker KBB Navigasyon Sistemi

TGS® Guidewire
 TGS-Guidewire
 TGS-voerdraad
 TGS Guidewire
 Fils-guides TGS
 TGS-Guidewire
 Οδηγό σύρμα TGS
 Filo guida TGS
 TGS 가이드와이어
 TGS Guidewire
 TGS Guidewire
 Fio-guia TGS
 Aguja guía de TGS
 TGS-ledare
 TGS Guidewire

REF 8000-060-009

INSTRUCTIONS FOR USE	2
BRUGSANVISNING	3
GEBRUIKSAANWIJZING	4
KÄYTTÖOHJEET	5
MODE D'EMPLOI	6
GEBRAUCHSANWEISUNG	7
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	8
ISTRUZIONI PER L'USO	9
사용 지침	10
BRUKSANVISNING	11
INSTRUKCJA UŻYCIA	12
INSTRUÇÕES DE UTILIZAÇÃO	13
INSTRUCCIONES DE USO	14
BRUKSANVISNING	15
KULLANIM TALİMATLARI	16

- ENGLISH
- DANSK
- NEDELANDS
- SUOMI
- FRANÇAIS
- DEUTSCH
- ΕΛΛΗΝΙΚΑ
- ITALIANO
- 한국어
- NORSK
- POLSKI
- PORTUGUÊS
- ESPAÑOL
- SVENSKA
- TÜRKÇE

Important Information

Instructions for Use

These instructions are the most comprehensive source of information for the safe, effective, and compliant use of the product. Read and understand these instructions as well as the respective system user manual before using the product or any component compatible with the product. Contact Stryker for training as needed.

The following signal words may be used throughout this document:

WARNING - Highlights a safety-related issue. Always comply with this information to prevent patient and/or healthcare staff injury.

CAUTION – Highlights a product reliability issue. Always comply with this information to prevent product damage.

Note – Supplements and/or clarifies information.

Intended Use

The TGS Guidewire is intended as an aid for precisely locating anatomical structures in either open or percutaneous procedures. It is intended for use with the Stryker ENT Navigation System and the XprESS LoProfile ENT Dilation System during balloon sinus dilation procedures.

WARNING – The product may only be used for its intended purpose and in accordance with the following instructions, as well as the applicable, current version of the Stryker ENT Navigation System manual. The manual is part of the product and must therefore be accessible to personnel at all times. It must be handed over to subsequent owners or users.

Indications For Use

The TGS Guidewire is indicated for any medical condition in which the use of stereotactic surgery may be appropriate, and where reference to a rigid anatomical structure in the field of ENT surgery, such as the paranasal sinuses, mastoid anatomy, can be identified relative to a CT- or MR-based model of the anatomy.

Example procedures include, but are not limited to the following ENT procedures:

- Transsphenoidal access procedures
- Intranasal procedures
- Sinus procedures, such as maxillary antrostomies, ethmoidectomies, sphenoidotomies/ sphenoid explorations, turbinate resections, and frontal sinusotomies
- ENT-related anterior skull base procedures

Contraindications

The instrument must not be exposed to MRI or used in a Magnetic Resonance Environment. The MRI exposure might magnetize the sensor.

Safety Directives

WARNING

- Do not use breached or damaged packages, since the sterility and functionality of the device may be compromised.
- The TGS Guidewire is provided sterile and intended for single use only. Do not re-sterilize and/or reuse, as it may result in compromised device performance and risk improper sterilization and cross-contamination.
- Unauthorized modifications of the product are forbidden for safety reasons.
- Due to its composition, the product is not to be used together with other magnetically sensitive medical products, devices or instruments.

PRECAUTIONS

- Store the TGS Guidewire in a cool dry place.
- Never use a device that is beyond its expiration date.
- Before use confirm that the tip of the TGS Guidewire is located at or near the ball tip of the XprESS LoProfile device.
- Be sure to fully engage the bayonet connector with the XprESS LoProfile device.
- Unbend XprESS LoProfile device from the maxillary bend configuration prior to inserting or removing the TGS Guidewire.
- Consider using a new device if cross-contamination between sinuses is a concern.

Using Your Product

Compatible Devices

The TGS Guidewire can only be used as a navigation instrument when used with the Stryker ENT Navigation System.

The TGS Guidewire is compatible with the XprESS LoProfile ENT Dilation System. It can be inserted into the working lumen and used to navigate the tip of the XprESS device.

Device Description

The TGS Guidewire is a sterile, single use, disposable instrument intended to be used with the Stryker ENT Navigation System and the XprESS LoProfile ENT Dilation System. The instrument is an electromagnetically navigated device which consists of a sensor in a protective sheath that can be inserted into the working lumen of the XprESS device. A bayonet connector secures the device to the XprESS luer and a plug connects the device to the navigation system. The navigation system displays the position of the instrument in the preoperative scans.


Figure 1 TGS Guidewire

How Supplied

The TGS Guidewire is provided sterile and is intended for single-use only. Do not resterilize and/or reuse, as it may result in compromised device performance and risk improper sterilization and cross-contamination. Do not use breached or damaged packages, since the sterility and functionality of the device may be compromised.

Preparing for Navigation

1. Remove the TGS Guidewire from the sterile packaging. Remove the silicone cable tie from the guidewire.
2. Load the guidewire into the working lumen of the XprESS device (Figure 2).
3. Attach the guidewire bayonet connector to the luer fitting of the XprESS device (Figure 3).
4. Connect the plug to the navigation system (Figure 4).
5. Confirm accuracy by touching the tip of the device to a known rigid anatomical reference point on the patient.
6. After completing the entire procedure, dispose of the devices and all waste products according to appropriate environmental health safety guidelines.


Figure 2

Figure 3


Figure 4

Instructions for Application

For instructions and procedures of navigated surgery, refer to the Stryker ENT Navigation System instructions for use.

For information and precautions regarding electromagnetic compatibility (EMC), refer to Section 11 of the Stryker ENT Navigation System instructions for use.

For instructions and procedures of transnasal balloon dilation, refer to the XprESS ENT Dilation System instructions for use.

Not made with natural rubber latex. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: the Stryker logo and XprESS. All other trademarks are trademarks of their respective owners or holders.

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Entellus Medical Inc.
ICON (LR) Limited
Entellus Medical, Inc., a Stryker Company