Calibration Body Electromagnetic

REF 8000-060-002

Instructions for Use

Rx Only

EN Instructions for use	1
DA Brugsanvisning	19
DE Gebrauchsanweisung	37
EL Οδηγίες χρήσης	55
ES Instrucciones de uso	74
FI Käyttöohjeet	92
FR Notice d'utilisation	110
IT Istruzioni per l'uso	128
JA 使用説明書	146
KO 사용 설명서	164
NL Gebruiksaanwijzing	182
NO Bruksanvisning	200
PL Instrukcja użycia	218
PT Instruções de utilização	236
RO Instrucțiuni de utilizare	254
SV Bruksanvisning	272
TR Kullanma Talimatı	290
ZH 使用说明	308

1. How to Use this Document

This manual is the most comprehensive source of information for the safe, effective, and compliant use and/or maintenance of the product. Read and understand this manual as well as the respective system user manual before using the product or any component compatible with the product. When combined with other medical devices, the user manual of these devices is to be considered as well. This manual is a permanent part of the product. Keep this manual for future reference.

The following conventions are used in this document:

The signal word **WARNING** highlights a safety-related issue. Comply with this information to prevent patient and/or medical staff injury.

The signal word **CAUTION** highlights a product reliability issue. Comply with this information to prevent product damage.

Supplements or clarifies information.

1.1. Symbol Definition

EN ISO 7010 Graphical symbols – Safety colors and safety signs – Registered safety signs

Symbol/Number Name: Definition

W001

General warning sign: To signify a general warning.

M002

Refer to instruction manual/booklet: To signify that the user instruction manual/booklet must be read.

EN ISO 15223-1 Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied – Part 1 General requirements

Symbol/Number Name: Definition

5.1.1

Manufacturer: Indicates the medical device manufacturer as defined in the European Union harmonization legislation.

5.1.3

Date of manufacture: Indicates the date the medical device was manufactured.

5.1.5

Batch code: Indicates the manufacturer's batch code so that the batch or lot can be identified.

5.1.6

Catalog number: Indicates the manufacturer's catalog number so that the medical device can be identified.

5.1.7

Serial number: Indicates the manufacturer's serial number so that the medical device can be identified.

5.2.7

Non-Sterile: Indicates a medical device that has not been subjected to a sterilization process.

Symbol/Number Name: Definition

Keep away from sunlight: Indicates a medical device that needs protection from light sources.

Keep dry: Indicates a medical device that needs to be protected from moisture.

Temperature limit: Indicates the temperature limits to which the medical device can be safely exposed.

Humidity limitation: Indicates the range of humidity to which the medical device can be safely exposed.

Atmospheric pressure limit: Indicates the range of atmospheric pressure to which the medical device can be safely exposed.

Symbol/Number

Name: Definition

5.4.3

Consult instructions for use: Indicates the need to consult the instructions for use.

Indicates a medical device according to European Union harmonization legislation.

Product-Specific Symbols

Symbol Name: Definition

Quantity: Indicates the number of medical devices in the packaging.

Note: Used to supplement or clarify information.

GTIN

Global Trade Item Number.

Directive 2012/19/EU on waste electrical and electronic equipment (WEEE)

Symbol Name: Definition

Indicates that the product must be collected separately and must not be disposed of as unsorted municipal waste.

81 FR 38911 FDA Final rule for the use of symbols in labeling

Symbol	Name: Definition
--------	------------------

Rx Only	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.
---------	--

Regulatory marks and logos

Symbol	Definition
--------	------------

	Indicates that a device is in conformity with the applicable requirements set out in applicable European Union harmonization legislation providing for its affixing.
--	--

2. Safety Information

WARNING

- The product is delivered in a non-sterile condition. Before first use, as well as before each following use, the product must be processed according to a validated procedure.
- Unauthorized modifications of the product are forbidden for safety reasons.
- Before using a medical product, assure yourself of its functional safety and proper condition via a visual inspection. Do not use the product if any defects, such as breaks, cracks, deformations or wear, are detected. Particularly important parts such as tips, notches and all moveable components should be checked with extra care.
- Check the dot pattern before using the product. If the dot pattern is damaged, the product cannot be used.
- Due to its composition, the product is not to be used together with other magnetically sensitive medical products, devices or instruments (e.g. MRI).
- The product is a high-precision instrument. Avoid subjecting it to serious strains, such as heavy impacts. The product cannot be used if there are visible defects. After a heavy impact, the product must be checked for defects.

- The healthcare provider performing any procedure is responsible for determining the appropriateness of using the product and for the specific technique for each patient. Stryker, as a manufacturer, does not recommend a specific surgical procedure.
-

The user and/or patient should report any serious product-related incident(s) to both the manufacturer and the national competent authority where the user and/or patient is established.

3. Product Information

WARNING

- The product may only be used for its intended purpose and in accordance with this document and all current versions of relevant system and software application documentation.
 - This document is part of the product and must be accessible to personnel at all times. It must be provided to subsequent owners or users.
-

3.1. Intended Use

The Calibration Body Electromagnetic is an accessory to the Electromagnetic Navigation Unit and intended for intraoperatively calibrating navigated instruments.

3.2. Indications for Use

The Stryker ENT Navigation System is indicated for any medical condition in which the use of stereotactic surgery may be appropriate, and where reference to a rigid anatomical structure in the field of ENT surgery, such as the paranasal sinuses, mastoid anatomy, can be identified relative to a CT- or MR-based model of the anatomy.

Example procedures include, but are not limited to the following ENT procedures:

- Transsphenoidal access procedures
- Intranasal procedures
- Sinus procedures, such as maxillary antrostomies, ethmoidectomies, sphenoidotomies/ sphenoid explorations, turbinate resections, and frontal sinusotomies
- ENT-related anterior skull base procedures

3.3. Contraindications

None known.

3.4. User Group

Healthcare professionals (surgeon/resident, nurse/professional caregiver) educated in computer-assisted surgery and thoroughly familiar with the instructions for use and with the operation of this product.

To request an additional in-service instruction, contact Stryker.

4. Product Overview


Figure 1: Calibration Body Electromagnetic

A Dot Pattern

For instructions on how to use the product for a surgical procedure, refer to the instructions for use supplied with the respective software application.

5. For Use with

WARNING

Use only Stryker-approved products, unless otherwise specified.

For information related to compatible software applications, refer to the user manual supplied with the respective software application.

6. Product Lifetime

The instrument is validated for ten reprocessing cycles. Thus, the instrument contains a usage counter that disables the instrument after running the tenth usage cycle in surgery. The usage counter is updated as soon as the instrument is connected during surgery and shown in the software.

To maintain accurate use count, connect the instrument to the electromagnetic port at least once during surgery. This also applies to replacement instruments if they have been taken from sterile packaging.

7. Reprocessing

WARNING

Do not wrap the cable around the calibration pattern during the reprocessing phase. Otherwise, the calibration pattern could be damaged.

For reprocessing instructions, safety directives, and reprocessing equipment, refer to the Stryker ENT Navigation System Reprocessing Instructions (TD8000010706).

8. Maintenance

Maintenance and repair may only be conducted by the manufacturer or authorized partners. The product components may only be sent back to the manufacturer in a cleaned, disinfected, and sterilized condition. Sharp or pointed components need to be sent back in a protected state.

9. Disposal

Risk of injury or infection can be avoided by safe disposal of the product components. Sharp and pointed components need to be collected and locked in a tight and break-proof container. They must be stored in such a way that they are protected from unauthorized use. Contaminated products are to be supplied to a hazardous waste site and handled in a way that contamination of third parties is excluded.

- After using the electromagnetic system, make sure that instruments that cannot be reprocessed again are disposed of.
 - In accordance with European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) product should be collected separately for recycling.
 - Do not dispose of as unsorted municipal waste. Contact local distributor for disposal information. Ensure infected equipment is decontaminated prior to recycling.
 - The plug used in the electronic products contain the following substance: Lead, CAS No. 7439-92-1 (1907/2006 REACH). Handling Instructions: No specific precautions are required for handling items manufactured from alloys containing lead in the supplied condition.
 - The silicone glue used in this product contains the following substances: Decamethylcyclopentasiloxane, CAS No. 541-02-6 and Dodecamethylcyclohexasiloxane, CAS No. 540-97-6 (1907/2006 REACH).
-

10. Technical Specifications

Environmental Limitations	Operation	Storage and Transportation
Temperature		
Relative humidity		
Atmospheric air pressure		-

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-

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

U.S. Patents: www.stryker.com/patents

Stryker Corporation or its divisions or other affiliated entities own, use or have applied for the following trademarks or service marks: Leibinger, Stryker. All other trademarks are trademarks of their respective owners or holders.

Copyright © 2019-2021 Stryker

Stryker Leibinger GmbH & Co. KG
Bötzingen Straße 41
79111 Freiburg (Germany)
t: +49 761 4512 0 (Germany)
t: +1 269 323 7700 (USA)