Field Generator Mounting Arm

REF 8000-010-005

Instructions for Use

Rx Only

EN Instructions for use	1
DA Brugsanvisning	21
DE Gebrauchsanweisung	41
EL Οδηγίες χρήσης	61
ES Instrucciones de uso	81
FI Käyttöohjeet	101
FR Notice d'utilisation	121
IT Istruzioni per l'uso	141
JA 使用説明書	161
KO 사용 설명서	181
NL Gebruiksaanwijzing	201
NO Bruksanvisning	221
PL Instrukcja użycia	241
PT Instruções de utilização	261
RO Instrucțiuni de utilizare	281
RU Инструкции по применению	301
SV Bruksanvisning	321
TR Kullanım talimatları	341
ZH 使用说明	361

1. How to Use this Document

1.1. About this Document

This manual is the most comprehensive source of information for the safe, effective, and compliant use and/or maintenance of the product. Read and understand this manual as well as the instructions for use supplied with the respective software application before using the product or any component compatible with the product. When combined with other medical devices, the instructions for use of these devices are to be considered as well. This manual is a permanent part of the product. Keep this manual for future reference.

The following conventions are used in this document:

The signal word **WARNING** highlights a safety-related issue. Comply with this information to prevent patient or medical staff injury.

The signal word **CAUTION** highlights a product reliability issue. Comply with this information to prevent product damage.

Supplements or clarifies information.

1.2. Symbol Definitions

EN ISO 7010 Graphical symbols – Safety colors and safety signs – Registered safety signs:

Symbol	Definition
	General warning sign (W001): Signifies a general warning.

ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer - Part 1 General requirements:

Symbol	Definition
	Manufacturer (5.1.1): Indicates the medical device manufacturer as defined in the European Union harmonization legislation.
	Date of manufacture (5.1.3): Indicates the date when the medical device was manufactured.
	Batch code (5.1.5): Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog number (5.1.6): Indicates the manufacturer's catalog number so that the medical device can be identified.

Symbol	Definition
	Keep away from sunlight (5.3.2): Indicates a medical device that needs protection from light sources.
	Keep dry (5.3.4): Indicates a medical device that needs to be protected from moisture.
	Temperature limit (5.3.7): Indicates the temperature limits to which the medical device can be safely exposed.
	Consult instructions for use (5.4.3): Indicates the need to consult the instructions for use.
	Caution (5.4.4): Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.

81 FR 38911 FDA Final rule for the use of symbols in labeling

Symbol	Definition
Rx Only	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Product-Specific Symbols

Symbol	Definition
	Quantity: Indicates the number of products in the packaging.
	Note: Used to supplement or clarify information
GTIN	Global Trade Item Number.
	Indicates a medical device according to European Union harmonization legislation.

Regulatory marks and logos

Symbol	Definition
	Indicates that a device is in conformity with the applicable requirements set out in applicable European Union harmonization legislation providing for its affixing.

2. Safety Information

WARNING

- Unauthorized modifications of the product are forbidden for safety reasons.
 - Due to its composition, the product is not to be used together with other magnetically sensitive medical products, devices or instruments (e.g. MRI).
 - Before using a medical product, assure yourself of its functional safety and proper condition via a visual inspection. Do not use the product if any defects, such as breaks, cracks, deformations or wear, are detected. Particularly important parts such as tips, notches and all moveable components should be checked with extra care
 - Avoid direct contact of the product with the patient.
 - The healthcare provider performing any procedure is responsible for determining the appropriateness of using the product and for the specific technique for each patient. Stryker, as a manufacturer, does not recommend a specific surgical procedure.
-

CAUTION

- Avoid subjecting the product to serious strains, such as heavy impacts. The product cannot be used if there are visible defects. After a heavy impact, the product must be checked for defects.
-

The user and/or patient should report any serious product-related incident to both the manufacturer and the national competent authority where the user and/or patient is established.

3. Product Information

WARNING

- The product may only be used for its intended purpose and in accordance with this document and all current versions of relevant system and software application documentation.
 - This document is part of the product and must be accessible to personnel at all times. It must be provided to subsequent owners or users.
-

3.1. Stryker ENT Navigation System

3.1.1. United States / Rest of World*

Intended Use

The Field Generator Mounting Arm is an accessory to the Electromagnetic Navigation Unit and is intended for positioning the Field Generator.

Indications for Use

The Stryker ENT Navigation System is indicated for any medical condition in which the use of stereotactic surgery may be appropriate, and where reference to a rigid anatomical structure in the field of ENT surgery, such as the paranasal sinuses, mastoid anatomy, can be identified relative to a CT- or MR-based model of the anatomy.

Example procedures include, but are not limited to the following ENT procedures:

- Transsphenoidal access procedures
- Intranasal procedures
- Sinus procedures, such as maxillary antrostomies, ethmoidectomies, sphenoidotomies/ sphenoid explorations, turbinate resections, and frontal sinusotomies
- ENT-related anterior skull base procedures

Contraindications

None known.

*Note: Product may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Contact your Stryker representative for product availability.

3.2. Cranial Guidance System

3.2.1. United States

Refer to the user manual supplied with the Cranial Guidance Software (REF 6000-670-000) for system indications and contraindications.

3.3. User Group

Healthcare professionals (surgeon/resident, nurse/professional caregiver) educated in computer-assisted surgery and thoroughly familiar with the instructions for use and with the operation of this product.
To request an additional in-service instruction, contact Stryker.

4. Product Overview


Figure 1: Field Generator Mounting Arm

- A Upper part of the articulated arm
- B Central clamp
- C Screw
- D Rail of the operating table
- E Clamp
- F Star handle screw

5. Assembly

Refer to Figure 1.

How to attach the Field Generator Mounting Arm to a standard rail

1. Loosen the screw (C) counterclockwise.
2. Place the clamp (E) on the standard rail. Ensure that the notch of the clamp is hooked to the rail (D).
3. Fasten the screw (C) of the clamp clockwise.

How to position the Field Generator Mounting Arm

When operating the central clamp (B), always secure the upper part of the articulated arm (A) with the other hand.

1. To loosen the fixation of the Field Generator Mounting Arm, open the central clamp (B) counterclockwise.
2. Position the Field Generator Mounting Arm as desired.
3. Rigidly attach the Field Generator Mounting Arm by tightening the central clamp (B) clockwise.

For instructions on how to use the product for a surgical procedure, refer to the instructions for use supplied with the respective software application.

6. For Use With

WARNING

Use only Stryker-approved products, unless otherwise specified.

For information related to product-specific compatibility, refer to the table below.

Description	REF
Field Generator	8000-010-004

7. Reprocessing

WARNING

- The product requires at minimum an intermediate level of disinfection for safe use. Before first use, as well as before each following use, the product must be cleaned and disinfected according to a validated procedure.
- Only use products that are suited for wipe disinfection of medical devices.
- The product cannot be sterilized. If sterility is required for the specific surgical procedure, the product must be draped.

CAUTION

- Do not sterilize the product. Failure to comply may lead to product damage.
- Do not subject the product to temperatures of more than 60°C.
- Avoid contact between the mounting arm and salt solutions (e.g. physiological NaCl), as these could corrode the metal surface.

Necessary steps in the reprocessing of medical products include preparation, cleaning, disinfection, testing, maintenance, inspection, and storage. Details on how to carry out these steps for the listed products will be provided below. Also consult the CDC Guidance Document Guideline for Disinfection and Sterilization in Healthcare Facilities, 2008. Also refer to the following for the purpose of additional education but not in lieu of validated reprocessing instructions: professional organizations' clinical practice guidelines or clinical guidelines of the CDC.

- Reprocess devices within two hours after use.
- Reprocess devices before you return devices for service.

Personal safety

Use the following personal protective equipment:

- (single-use) shirt
- (single-use) gloves
- safety glasses
- anti-splash mask

Refer to the user instructions of the cleaning agent and disinfectant for additional protective equipment.

7.1. Manual Cleaning and Disinfection

The following manual cleaning and disinfection process has been validated by Stryker. If other reprocessing methods are used, make sure the reprocessing cycle is validated.

Required equipment:

- Low-lint cloth, swab, cotton swabs
- CaviWipes from Metrex Research
- Freshly prepared purified water, highly purified water or sterile water with less than 10 cfu/ml and 0.25 EU/ml. In the U.S., refer to AAMI TIR 34 for guidance on water quality.
- Medical grade compressed air

CAUTION

- Pay attention to the instructions provided by the cleaning and disinfection aid manufacturer.
 - Do not use an ultrasonic bath to clean the product.
 - For cleaning, only use soft, clean cloths and/or swabs that are suitable for the purpose. Never use any sharp objects such as metal brushes or wire wool.
 - The mounting arm may not be immersed in liquid. The mounting arm's central clamp must be open during cleaning/disinfection.
-

7.1.1. Prepare the equipment and the devices

1. Dismount the Field Generator from the mounting arm.
2. Dismount the mounting arm from the rail of the operating table.
3. Open the mounting arm's central clamp (Figure 1, B).

7.1.2. Cleaning

1. Wipe all external surfaces of the equipment with CaviWipes wipes.
2. Keep the central clamp open while cleaning around articulated parts.
3. Use a new clean CaviWipes wipe when the wipe is contaminated.
4. Clean until all visible soil is removed.
5. Follow also the instructions for use of the manufacturer cleaning aid.

7.1.3. Disinfection

1. Using a new clean CaviWipes wipe, wipe all surfaces.
2. Adjust the mounting arm to different positions to reach all external surfaces.
3. Make sure the surfaces remain visibly wet at room temperature for at least 3 minutes. Follow also the instructions for use supplied by the disinfectant manufacturer.
4. Use a new clean CaviWipes wipe when the wipe is contaminated.
5. Remove any excess disinfectant solution using a soft, low-lint cloth slightly moistened with purified water.

7.1.4. Drying and inspection

1. Dry the external surface by letting it dry on a low-lint cloth, wipe it with a dry low-lint cloth until the external surfaces are dry or use medical grade compressed air.
2. Visually inspect the product for any damages, wear or remaining soil or disinfectant solution. If soil or disinfectant solution remains, repeat the cleaning and disinfection procedure using fresh CaviWipes wipes.

8. Maintenance

Maintenance and repair may only be conducted by the manufacturer or authorized partners. The product components may only be sent back to the manufacturer in a cleaned and disinfected condition. Sharp or pointed components need to be sent back in a protected state.

9. Disposal

Risk of injury or infection can be avoided by safe disposal of the product components. They must be stored in such a way that they are protected from unauthorized use. Contaminated products are to be supplied to a hazardous waste site and handled in a way that contamination of third parties is excluded.


10. Technical Specifications

Environmental conditions for transport and storage

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Stryker Corporation or its divisions or other affiliated entities own, use or have applied for the following trademarks or service marks: Leibinger, Stryker. All other trademarks are trademarks of their respective owners or holders.
Copyright © 2022 Stryker
Stryker Leibinger GmbH & Co. KG