# Navigation Tool Extension Cable

REF 8000-050-011

Instructions for Use

Rx Only

EN   Instructions for use	1
DA   Brugsanvisning	15
DE   Gebrauchsanweisung	29
EL   Οδηγίες χρήσης	44
ES   Instrucciones de uso	59
FI   Käyttöohjeet	74
FR   Notice d'utilisation	88
IT   Istruzioni per l'uso	102
JA   使用説明書	116
KO   사용 설명서	130
NL   Gebruiksaanwijzing	144
NO   Bruksanvisning	158
PL   Instrukcja użycia	172
PT   Instruções de utilização	187
RO   Instrucțiuni de utilizare	202
SV   Bruksanvisning	216
TR   Kullanma Talimatı	230
ZH   使用说明	244

# 1. How to Use this Document

## 1.1. About this Document

This manual is the most comprehensive source of information for the safe, effective, and compliant use and/or maintenance of the product. Read and understand this manual as well as the respective system user manual before using the product or any component compatible with the product. When combined with other medical devices, the user manual of these devices is to be considered as well. This manual is a permanent part of the product. Keep this manual for future reference.

The following conventions are used in this document:

The signal word **WARNING** highlights a safety-related issue. Comply with this information to prevent patient and/or healthcare staff injury.

The signal word **CAUTION** highlights a product reliability issue. Comply with this information to prevent product damage.

Supplements or clarifies information.

---

## 1.2. Symbol Definitions

### EN ISO 7010 Graphical symbols – Safety colors and safety signs – Registered safety signs

Symbol	Definition
	<b>General warning sign (W001):</b> to signify a general warning.
	<b>Refer to instruction manual/booklet (M002):</b> to signify that the user instruction manual/booklet must be read.

### EN ISO 15223-1 Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied - Part 1 General requirements

Symbol	Definition
	<b>Manufacturer (5.1.1):</b> indicates the medical device manufacturer as defined in the European Union harmonization legislation.
	<b>Date of manufacture (5.1.3):</b> indicates the date when the medical device was manufactured.

Symbol	Definition
	<b>Batch code (5.1.5):</b> indicates the manufacturer's batch code so that the batch or lot can be identified.
	<b>Catalog number (5.1.6):</b> indicates the manufacturer's catalog number so that the medical device can be identified.
	<b>Serial number (5.1.7):</b> indicates the manufacturer's serial number so that a specific medical device can be identified.
	<b>Non-Sterile (5.2.7):</b> indicates a medical device that has not been subjected to a sterilization process.
	<b>Keep away from sunlight (5.3.2):</b> indicates a medical device that needs protection from light sources.
	<b>Keep dry (5.3.4):</b> indicates a medical device that needs to be protected from moisture.
	<b>Temperature limit (5.3.7):</b> indicates the temperature limits to which the medical device can be safely exposed.

Symbol	Definition
	<b>Humidity limitation (5.3.8):</b> indicates the range of humidity to which the medical device can be safely exposed.
	<b>Atmospheric pressure limitation (5.3.9):</b> indicates the range of atmospheric pressure to which the medical device can be safely exposed.
	<b>Consult instructions for use (5.4.3):</b> indicates the need for the use to consult the instructions for use.
	Indicates a medical device according to European Union harmonization legislation.

#### IEC 60417 Graphical symbols for use on equipment

Symbol	Definition
	<b>Type BF applied part (5333):</b> to identify a type BF applied part complying with IEC 60601-1.

## Product-Specific Symbols

Symbol	Definition
	<b>Quantity:</b> indicates the number of products in the packaging
	<b>Note:</b> used to supplement or clarify information
<b>GTIN</b>	Global Trade Item Number

## 81 FR 38911 FDA Final rule for the use of symbols in labeling

Symbol	Definition
<b>Rx Only</b>	<b>Caution:</b> Federal law (USA) restricts this device to sale by or on the order of a physician.

## Directive 2012/19/EU on waste electrical and electronic equipment (WEEE)

Symbol

Definition

Indicates that the product must be collected separately and must not be disposed of as unsorted municipal waste.

## Regulatory marks and logos

Symbol

Definition

Indicates that a device is in conformity with the applicable requirements set out in applicable European Union harmonization legislation providing for its affixing.

## 2. Safety Information

### WARNING

- The product is delivered in a non-sterile condition. Before first use, as well as before each following use, the product must be processed according to a validated procedure.
- Unauthorized modifications of the product are forbidden for safety reasons.
- Due to its composition, the product is not to be used together with other magnetically sensitive medical products, devices or instruments (e.g. MRI).
- Before using a medical product, assure yourself of its functional safety and proper condition via a visual inspection. Do not use the product if any defects, such as breaks, cracks, deformations or wear, are detected. Particularly important parts such as tips, notches and all moveable components should be checked with extra care.
- The product is a high precision instrument. Avoid subjecting the product to serious strains, such as heavy impacts. The product cannot be used if there are visible defects. After a heavy impact, the product must be checked for defects.
- The healthcare provider performing any procedure is responsible for determining the appropriateness of using the product and for the specific technique for each patient. Stryker, as a manufacturer, does not recommend a specific surgical procedure.

The user or patient should report any serious product-related incident to both the manufacturer and the national competent authority where the user or patient is established.

### 3. Product Information

#### WARNING

- The product may only be used for its intended purpose and in accordance with this document and all current versions of relevant system and software application documentation.
  - This document is part of the product and must be accessible to personnel at all times. It must be provided to subsequent owners or users.
- 

#### 3.1. Intended Use

The Navigation Tool Extension Cable is an accessory to the Electromagnetic Navigation Unit and is intended for sterile handling of compatible navigated instruments.

#### 3.2. Indications for Use

The Stryker ENT Navigation System is indicated for any medical condition in which the use of stereotactic surgery may be appropriate, and where reference to a rigid anatomical structure in the field of ENT surgery, such as the paranasal sinuses, mastoid anatomy, can be identified relative to a CT- or MR-based model of the anatomy.

Example procedures include, but are not limited to the following ENT procedures:

- Transsphenoidal access procedures
- Intranasal procedures
- Sinus procedures, such as maxillary antrostomies, ethmoidectomies, sphenoidotomies/sphenoid explorations, turbinate resections, and frontal sinusotomies
- ENT-related anterior skull base procedures

### **3.3. Contraindications**

None known.

### **3.4. User Group**

Healthcare professionals (surgeon/resident, nurse/professional caregiver) educated in computer-assisted surgery and thoroughly familiar with the instructions for use and with the operation of this product.

To request an additional in-service instruction, contact Stryker.

## 4. Product Overview

Figure 1: Navigation Tool Extension Cable

## 4.1. Instructions for Application

### WARNING

- Ensure the instrument connection socket is kept at least 30 cm away from the Field Generator. Otherwise the transmitted signal may be lost or the accuracy of the instrument may deteriorate.
  - Do not connect multiple Navigation Tool Extension Cables to one another. Failure to comply may compromise navigation accuracy.
- 

For instructions on how to use the product for a surgical procedure, refer to the instructions for use supplied with the respective software application.

## 5. For Use With

### WARNING

Use only Stryker-approved products, unless otherwise specified.

---

The Navigation Tool Extension Cable can be used in conjunction with navigated electromagnetic instruments. For information related to compatible software applications, refer to the instructions for use supplied with the respective software application.

## **6. Disassembly**

Detach the Navigation Tool Extension Cable from the Electromagnetic Navigation Unit and from the navigated electromagnetic instrument.

## **7. Reprocessing**

For reprocessing instructions, safety directives, and reprocessing equipment, refer to the Stryker ENT Navigation System Reprocessing Instructions (TD8000010706).

## **8. Maintenance**

Maintenance and repair may only be conducted by the manufacturer or authorized partners. The product components may only be sent back to the manufacturer in a cleaned, disinfected, and sterilized condition. Sharp or pointed components need to be sent back in a protected state.

## 9. Disposal

Risk of injury or infection can be avoided by safe disposal of the product components. Sharp and pointed components need to be collected and locked in a tight and break-proof container. They must be stored in such a way that they are protected from unauthorized use. Contaminated products are to be supplied to a hazardous waste site and handled in a way that contamination of third parties is excluded.

- In accordance with European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) product should be collected separately for recycling. Do not dispose of as unsorted municipal waste. Contact local distributor for disposal information. Ensure infected equipment is decontaminated prior to recycling.
  - The plugs used in the electronic products contain the following substance: Lead, CAS No. 7439-92-1 (1907/2006 REACH).  
Handling Instructions: No specific precautions are required for handling items manufactured from alloys containing lead in the supplied condition.
  - The silicone glue used in the products contains the following substances: Decamethylcyclopentasiloxane, CAS No. 541-02-6 and Dodecamethylcyclohexasiloxane, CAS No. 540-97-6 (1907/2006 REACH).
-

## 10. Technical Specifications

Environmental Limitations	Operation	Storage and Transportation
Temperature	 <p>10°C — 30°C</p>	 <p>-10°C — 50°C</p>
Relative humidity	 <p>30% — 75%</p>	
Atmospheric air pressure	 <p>80 kPa — 106 kPa</p>	—

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This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

---

This page was intentionally left blank for your notes.

U.S. Patents: [www.stryker.com/patents](http://www.stryker.com/patents)

Stryker Corporation or its divisions or other affiliated entities own, use or have applied for the following trademarks or service marks: Leibinger, Stryker. All other trademarks are trademarks of their respective owners or holders.

Copyright © 2019-2021 Stryker

Stryker Leibinger GmbH & Co. KG  
Bötzingen Straße 41  
79111 Freiburg (Germany)  
t: +49 761 4512 0 (Germany)  
t: +1 269 323 7700 (USA)