# Field Generator

REF 8000-010-004

Instructions For Use

EN   Instructions for use.....	1
DA   Brugsanvisning.....	18
DE   Gebrauchsanweisung .....	37
EL   Οδηγίες χρήσης .....	56
ES   Instrucciones de uso .....	75
FI   Käyttöohjeet .....	94
FR   Mode d'emploi.....	112
IT   Istruzioni per l'uso.....	131
JA   使用説明書 .....	150
KO   사용 설명서 .....	168
NL   Gebruiksaanwijzing .....	186
NO   Bruksanvisning .....	205
PL   Instrukcja użycia .....	224
PT   Instruções de utilização.....	243
RO   Instrucțiuni de utilizare .....	262
RU   Инструкции по применению .....	281
SV   Bruksanvisning .....	300
TR   Kullanım talimatları.....	319
ZH   使用说明 .....	337

Table of Contents

- 1. How to Use this Document ..... 1
  - 1.1. About this Document ..... 1**
  - 1.2. Symbol Definition..... 1**
  
- 2. Safety Information ..... 3
  
- 3. Product Information..... 4
  - 3.1. Product Description ..... 4**
  - 3.2. User Group ..... 4**
  
- 4. For Use with..... 5
  
- 5. Getting Started ..... 6
  - 5.1. Unpacking ..... 6**
  - 5.2. Positioning..... 6**
  
- 6. Using Your Product ..... 8
  - 6.1. Instructions for Application..... 8**
  - 6.2. Connecting the Field Generator ..... 8**
  
- 7. Disassembly..... 10
  
- 8. Reprocessing ..... 11
  - 8.1. Manual Cleaning and Disinfection ..... 11**
    - 8.1.1. Prepare the Equipment and the Devices..... 12**
    - 8.1.2. Cleaning ..... 12**
    - 8.1.3. Disinfection ..... 12**
    - 8.1.4. Drying and Inspection ..... 13**
  
- 9. Maintenance..... 14

10. Transport, Storage and Disposal .....	15
<b>10.1. Disposal.....</b>	<b>15</b>
11. Technical Specifications.....	16
12. Electromagnetic Compatibility.....	17

# 1. How to Use this Document

## 1.1. About this Document

This manual is the most comprehensive source of information for the safe, effective, and compliant use and/or maintenance of the product. Read and understand this manual as well as the instructions for use supplied with the respective software application before using the product or any component compatible with the product. When combined with other medical devices, the instructions for use of these devices are to be considered as well.

This manual is a permanent part of the product. Keep this manual for future reference.

The following conventions are used throughout this document:

The signal word **WARNING** highlights a safety-related issue. Comply with this information to prevent patient or medical staff injury.

The signal word **CAUTION** highlights a product reliability issue. Comply with this information to prevent product damage.

 Supplements or clarifies information.

## 1.2. Symbol Definition

EN ISO 7010 Graphical symbols – Safety colors and safety signs – Registered safety signs

Symbol	Name: Definition
 W001	<b>General warning sign:</b> To signify a general warning.
 M002	<b>Refer to instruction manual/booklet:</b> To signify that the user instruction manual/booklet must be read.
 P007	<b>No access for people with active implanted cardiac devices:</b> To prohibit people with active implanted cardiac devices from entering a designated area.

EN ISO 15223-1 Medical devices — Symbols to be used with medical device labels, labeling and information to be supplied - Part 1 General requirements

Symbol/ number	Name: Definition
 5.1.3	<b>Date of manufacture:</b> Indicates the date when the medical device was manufactured.
 5.1.5	<b>Batch code:</b> Indicates the manufacturer's batch code so that the batch or lot can be identified.
 5.1.6	<b>Catalog number:</b> Indicates the manufacturer's catalog number so that the medical device can be identified.
 5.1.7	<b>Serial number:</b> Indicates the manufacturer's serial number so that the medical device can be identified.
 5.3.1	<b>Fragile, handle with care:</b> Indicates a medical device that can be broken or damaged if not handled carefully.
 5.3.2	<b>Keep away from sunlight:</b> Indicates a medical device that needs protection from light sources.

Symbol/ number	Name: Definition
 5.3.4	<b>Keep dry:</b> Indicates a medical device that needs to be protected from moisture.
 5.3.7	<b>Temperature limit:</b> Indicates the temperature limits to which the medical device can be safely exposed.
 5.3.8	<b>Humidity limitation:</b> Indicates the range of humidity to which the medical device can be safely exposed.
 5.3.9	<b>Atmospheric pressure limit:</b> Indicates the range of atmospheric pressure to which the medical device can be safely exposed.
 5.4.3	<b>Consult instructions for use:</b> Indicates the need for the use to consult the instructions for use.
 5.4.4	<b>Caution:</b> Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.

#### Product-Specific Symbols

Symbol	Name: Definition
	<b>Quantity:</b> Indicates the number of medical devices in the packaging.
	<b>Note symbol:</b> It is used to supplement or clarify information.
<b>GTIN</b>	Global Trade Item Number.

#### Directive 2012/19/EU on waste electrical and electronic equipment (WEEE)

Symbol	Name: Definition
	Indicates that the product must be collected separately and must not be disposed of as unsorted municipal waste.

#### Regulatory marks and logos

Symbol	Definition
	Indicates that a device is in conformity with the applicable requirements set out in applicable European Union harmonization legislation providing for its affixing.
	UK product marking that is used for goods being placed on the market in Great Britain.

## 2. Safety Information

### WARNING

- Before using the product, assure yourself of its functional safety and proper condition via a visual inspection. Do not use the product if any defects, such as breaks, cracks, deformations or wear, are detected.
  - Unauthorized modifications of the product are forbidden for safety reasons.
  - Do not bend or kink cables, or use cables that are damaged. Position measurements from a system with damaged cables may result in possible personal injury.
  - Due to its composition, the product is not to be used together with other magnetically sensitive medical products, devices or instruments (e.g. MRI).
  - Do not use mobile phones in the vicinity of medical devices. Mobile phones and other devices emitting large amounts of electromagnetic radiation may interfere with the functioning of the product. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the navigation system, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
  - Do not operate the Field Generator within 200 mm of an installed pacemaker. The magnetic field produced by the Field Generator may interfere with the operation of the pacemaker. This interference may result in personal injury.
  - Do not disconnect the Field Generator from the system during navigation. Failure to comply may lead to sparks being generated, and to personal injury.
  - The healthcare provider performing any procedure is responsible for determining the appropriateness of using the product and for the specific technique for each patient. Stryker, as a manufacturer, does not recommend a specific surgical procedure.
- 

**i** The user and/or patient should report any serious product-related incident to both the manufacturer and the national competent authority where the user and/or patient is established.

---

### 3. Product Information

#### **WARNING**

- The product may only be used for its intended purpose and in accordance with this document and all current versions of relevant system and software application documentation.
  - This document is part of the product and must be accessible to personnel at all times. It must be provided to subsequent owners or users.
- 

#### **3.1. Product Description**

The Field Generator is used for the electromagnetic position measurement of patient and navigated instruments during surgical procedures in which compatible clinical navigation systems are used.

#### **3.2. User Group**

Healthcare professionals (surgeon/resident, nurse/professional caregiver) educated in computer-assisted surgery and thoroughly familiar with the instructions for use and with the operation of this product.

To request an additional in-service instruction, contact Stryker.

#### 4. For Use with

#### **WARNING**

Use only Stryker-approved accessories, unless otherwise specified.

---

For information related to compatible software applications, refer to the instructions for use supplied with the respective software application. For information related to product-specific compatibility, refer to the table below.

<b>Description</b>	<b>REF</b>
Electromagnetic Navigation Unit	8000-010-003
Field Generator Mounting Arm	8000-010-005
EM Instrument Adapter Box	8900-200-100

## 5. Getting Started

### 5.1. Unpacking

After opening the packaging, check that the delivery is complete and undamaged. If the delivery provides you with any reason to make a claim, get in contact with the delivery service immediately.

The original packaging should be retained, as it can be useful in case you need to transport the product at a later point.

### 5.2. Positioning

#### WARNING

- The product cannot be sterilized. If sterility is required for the specific surgical procedure, the product must be draped.
- Avoid direct contact of the product with the patient.

- For the positioning of the Field Generator, the Field Generator Mounting Arm REF 8000-010-005 is available. The Field Generator Mounting Arm is attached to the standard rail and allows the flexible positioning of the Field Generator.

1. Ensure that the Field Generator Mounting Arm is securely fixed to the standard rail of the operating table. Refer to the Field Generator Mounting Arm Instructions for Use.
2. Position the Field Generator Mounting Arm with the teeth (Figure 1, C) facing up and the upper part of the articulated arm of the Field Generator Mounting Arm (Figure 1, A) parallel to the operating table.
3. Place the Field Generator flush on the teeth of the Field Generator Mounting Arm (Figure 1, C).


Figure 1 – Positioning the Field Generator on the Field Generator Mounting Arm (Lateral View)

- A Upper part of the articulated arm
- B Star handle screw
- C Teeth

4. Fixate the Field Generator from above with one hand, insert the star handle screw (Figure 1, B) straight into the thread of the Field Generator and fasten the star handle screw with the other hand from below.
5. Position the Field Generator Mounting Arm so that the distance between the Field Generator and the operating area of the patient is 50 mm. Refer to Figure 2.


Figure 2 – Positioning of the Field Generator (Top View)

- A Tracking Volume of the Field Generator
- B Operating Area
- C Field Generator
- D Field Generator Cable

**i** For tracking volume limitations refer to the instructions for use supplied with the respective software application.

## 6. Using Your Product

### 6.1. Instructions for Application

For instructions on how to use the product for a surgical procedure, refer to the instructions for use supplied with the respective software application.

### 6.2. Connecting the Field Generator

#### WARNING

Do not place the Field Generator cable inside the tracking volume or wrap it around the Field Generator, as it may create magnetic interference. This interference can contribute to inaccurate position measurements and possible personal injury.

#### CAUTION

Never use force to put the plug into the plug socket. Do not bend or strain the cable.

Connect the plug of the Field Generator cable to the Field Generator socket on the front side of the Electromagnetic Navigation Unit (Figure 3) or on the underside of the EM Instrument Adapter Box (Figure 4). The red marks on the plugs of the Field Generator Cable and on the connection socket must be pointing at each other. The plug must click into place.


Figure 3 – Field Generator Connection to Electromagnetic Navigation Unit

---

Figure 4 – Field Generator Connection to EM Instrument Adapter Box

---

## 7. Disassembly

Remove the Field Generator as follows:

1. Position the Field Generator Mounting Arm in such way that the the upper part of the articulated arm is parallel to the operating table and the star handle screw of the Field Generator Mounting Arm faces down. Refer to Figure 1.
2. Secure the Field Generator with one hand and loosen the star handle screw of the Field Generator Mounting Arm with the other hand.
3. Once the star handle screw is loose, remove the Field Generator.

## 8. Reprocessing

Reprocess devices within two hours after use.

Reprocess devices before you return devices for service.

### WARNING

- Before cleaning and maintenance, unplug the device from the power supply and other devices connected to the mains.
- Do not expose or immerse the electrical components of the navigation system to liquids. Do not place liquids on or above the devices. Liquids or incoming parts may result in equipment damage, produce a fire or shock hazard.
- The product requires at minimum an intermediate level of disinfection for safe use. Before first use, as well as before each following use, the product must be cleaned and disinfected according to a validated procedure.
- Only use products that are suited for wipe disinfection of medical devices.

### CAUTION

Do not sterilize the product. Failure to comply may lead to product damage.

Necessary steps in the reprocessing of medical products include preparation, cleaning, disinfection, testing, maintenance, inspection, and storage. Details on how to carry out these steps for the listed products will be provided below. Also consult the CDC Guidance Document Guideline for Disinfection and Sterilization in Healthcare Facilities, 2008. Also refer to the following for the purpose of additional education but not in lieu of validated reprocessing instructions: professional organizations' clinical practice guidelines or clinical guidelines of the CDC.

### Personal safety

Use the following personal protective equipment:

- (single-use) shirt
- (single-use) gloves
- safety glasses
- anti-splash mask

Refer to the user instructions of the cleaning agent and disinfectant for additional protective equipment.

### 8.1. Manual Cleaning and Disinfection

The following manual cleaning and disinfection process has been validated by Stryker. If other reprocessing methods are used, make sure the reprocessing cycle is validated.

Required equipment:

- Low-lint cloth, swab, cotton swabs
- CaviWipes wipes from Metrex Research
- Freshly prepared purified water, highly purified water or sterile water with less than 10 cfu/ml and 0.25 EU/ml. In the U. S., refer to AAMI TIR 34 for guidance on water quality.
- Medical grade compressed air

## CAUTION

- Pay attention to the instructions provided by the cleaning and disinfection aid manufacturer.
  - Do not use an ultrasonic bath to clean the product.
  - For cleaning, only use soft, clean cloths and/or swabs that are suitable for the purpose. Never use any sharp objects such as metal brushes or wire wool.
- 

### **8.1.1. Prepare the Equipment and the Devices**

1. Disconnect the power supply cable of the Field Generator.
2. Disassemble the Field Generator from the connected devices.

### **8.1.2. Cleaning**

1. Wipe all external surfaces of the equipment with CaviWipes wipes.
2. Pay special attention to the power supply cable and the hard to reach areas of the device.
3. Use a new clean CaviWipes wipe when the wipe is contaminated.
4. Clean until all visible soil is removed.
5. Follow also the instructions for use of the manufacturer cleaning aid.

### **8.1.3. Disinfection**

1. Using a new clean CaviWipes wipe, wipe all surfaces.
2. Pay special attention to the power supply cable and the hard to reach areas of the device.
3. Make sure the surfaces remain visibly wet at room temperature for at least 3 minutes. Follow also the instructions for use supplied by the disinfectant manufacturer.
4. Use a new clean CaviWipes wipe, when the wipe is contaminated.
5. Remove any excess disinfectant solution using a soft, low-lint cloth slightly moistened with purified water.

#### **8.1.4. Drying and Inspection**

1. Dry the external surface by letting it dry on a low-lint cloth, wipe it with a dry low-lint cloth until the external surfaces are dry or use medical grade compressed air.
2. Pay special attention regarding the dryness to the plug of the power supply cable and the hard to reach areas of the device.
3. Visually inspect the product for any remaining soil or disinfectant solution. If soil or disinfectant solution remains, repeat the cleaning and disinfection procedure using fresh CaviWipes wipes.
4. Check device functionality and examine the product for:
  - Loose parts
  - Surface damages
  - Cable damages
  - Other damage
5. Dispose of any damaged instruments.

## 9. Maintenance

### **WARNING**

Do not service the product. It does not contain any parts the user can service. If service is required, contact your Stryker sales representative.

---

Maintenance and repair may only be conducted by the manufacturer or authorized partners. The product components may only be sent back to the manufacturer in a cleaned and disinfected condition. Sharp or pointed components need to be sent back in a protected state.

## 10. Transport, Storage and Disposal

The product components may only be sent back to the manufacturer in a cleaned and disinfected condition. Sharp or pointed components need to be sent back in a protected state.

Transport and shipping can lead to product damage. If possible, use the original packaging for packing the product. In case this is not available, pad all empty spaces carefully (to the sides, above, and below). Please also pad other boxes, in order to prevent product damage during transport.

For storage and transportation environmental limitations, see Technical Specifications.

### **10.1. Disposal**

Risk of injury or infection can be avoided by safe disposal of product components. They must be stored so that they are protected from unauthorized use. Contaminated products are to be supplied to hazardous waste disposal.

In accordance with European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) the product should be collected separately for recycling. Do not dispose of as unsorted municipal waste. Contact your local distributor for disposal information. Ensure infected equipment is decontaminated prior to recycling.

## 11. Technical Specifications

### WARNING

The device needs a warm-up period of 10 minutes after storage, plus 1 minute for each 1 °C difference between storage and operating temperature. Do not turn the device on until after the warm-up time has passed. If the temperature difference is less than 3 °C, no warm-up period is necessary.

### Specifications

Model	Field Generator (REF 8000-010-004)
Dimensions	<ul style="list-style-type: none"> <li>- Size: 200 mm x 200 mm x 70 mm (W x H x D)</li> <li>- Weight: 3.2 kg</li> <li>- Measuring Range: between 0.05 m and 0.55 m</li> <li>- Cable Length: 4.5 m</li> </ul>
Materials/Substances	<p>The PCB (Printed Circuit Board) within this product contains diodes that may contain the following substances:</p> <ul style="list-style-type: none"> <li>- Diboron trioxide, CAS No. 1303-86-2</li> <li>- Lead monoxide, CAS No. 1317-36-8</li> <li>- Lead titanium trioxide, CAS No. 12060-00-3</li> <li>- Methylhexahydromethylphthalic anhydride, CAS No. 25550-51-0</li> <li>- 141 4,4 - Isopropylidenediphenol, CAS No. 80-05-7 (1907/2006 REACH)</li> </ul>

Environmental Limitations	Operation	Storage and Transportation
Temperature		
Relative Humidity		
Atmospheric Air Pressure		

## 12. Electromagnetic Compatibility

### WARNING

- Do not operate the Field Generator within 10 m of another Field Generator. To do so may contribute to inaccurate position measurement and possible personal injury.
- Do not drop the Field Generator or subject it to impact. Physical damage to the Field Generator may alter its calibration and contribute to inaccurate position measurement and possible personal injury.
- Do not place cables of electromagnetic instruments within 30 mm of the Field Generator cable. Placing cables that close to the Field Generator cable — particularly if the cables are parallel to each other — can lead to electromagnetic interferences and compromise navigation accuracy.
- Do not coil the Field Generator cable, as it produces enough electric current that a magnetic field will be created when the cable is placed in a circular formation. This magnetic field may disturb the Field Generators' magnetic field, contributing to inaccurate position measurements and possible personal injury.

Electric devices intended to be used with medical electrical equipment are subject to particular precautionary requirements regarding their electromagnetic compatibility (EMC). Instructions regarding EMC are provided in the installation instructions of the navigation system.

Electromagnetic measurement systems are susceptible to other electromagnetic fields and to metals, especially ferromagnetic substances. Electromagnetic measurement systems provide accurate navigational information only if the following conditions are met:

- There are no disturbing substances either within 80 cm of the Field Generator or between the Field Generator and measuring sensors built into the electromagnetic instruments and the patient tracker. The following materials produce interferences:
  - Steel
  - Ferromagnetic stainless steel
  - Iron
  - Aluminum
  - Copper
- The following materials have only a small impact on the measurement accuracy and may be placed near the Field Generator if necessary:
  - Stainless steel (material number: 1.4301, 1.4305, 1.4401, 1.4404)
  - Titanium Ti6Al4V
  - Cobalt chromium steel
- There are no disturbing items within 80 cm of the Field Generator. Typical objects that cause interferences are:
  - Medical devices that require creation of electromagnetic fields to carry out their function
  - Electrically driven motors, from drills and shavers
  - Navigation system
  - Surgical cart with navigation and endoscopy equipment
  - Monitors
  - Keys, watches, jewelry
  - Personal electronics, such as phones

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Stryker Corporation or its divisions or other affiliated entities own, use or have applied for the following trademarks or service marks: Leibinger, Stryker. All other trademarks are trademarks of their respective owners or holders.

UK Importer  
Stryker UK Ltd. | Stryker House  
Hambridge Road | Newbury  
Berkshire | RG14 5AW

U.S. Patents: [www.stryker.com/patents](http://www.stryker.com/patents)

Copyright © 2022 Stryker

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Stryker Leibinger GmbH & Co. KG