## Patient Tracker Electromagnetic– 10 uses

**REF** 8000-040-002

Instructions For Use

Rx Only

TD8000040701 Rev AD

2022-06-01

EN   Instructions for use	1
DA   Brugsanvisning	18
DE   Gebrauchsanweisung	35
EL   Οδηγίες χρήσης	52
ES   Instrucciones de uso	69
FI   Käyttöohjeet	86
FR   Notice d'utilisation	103
IT   Istruzioni per l'uso	120
JA   使用説明書	137
KO   사용 설명서	154
NL   Gebruiksaanwijzing	171
NO   Bruksanvisning	188
PL   Instrukcja użycia	205
PT   Instruções de utilização	222
RO   Instrucțiuni de utilizare	239
RU   Инструкции по применению	256
SV   Bruksanvisning	273
TR   Kullanım talimatları	290
ZH   使用说明	307

# 1. How to Use this Document

## 1.1. About this Document

This manual is the most comprehensive source of information for the safe, effective, and compliant use and/or maintenance of the product. Read and understand this manual as well as the respective system user manual before using the product or any component compatible with the product. When combined with other medical devices, the instructions for use of these devices are to be considered as well. This manual is a permanent part of the product. Keep this manual for future reference.

The following conventions are used in this document:

The signal word **WARNING** highlights a safety-related issue. Comply with this information to prevent patient and medical staff injury.

The signal word **CAUTION** highlights a product reliability issue. Comply with this information to prevent product damage.

 Supplements or clarifies information.

## 1.2. Symbol Definitions

### EN ISO 7010 Graphical symbols – Safety colors and safety signs – Registered safety signs

Symbol	Definition
	General warning sign (W001): To signify a general warning
	Refer to instruction manual/booklet (M002): To signify that the user instruction manual/booklet must be read.

### ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements

Symbol	Definition
	Manufacturer (5.1.1): Indicates the medical device manufacturer as defined in the European Union harmonization legislation
	Date of manufacture (5.1.3): Indicates the date when the medical device was manufactured.

Symbol	Definition
	<b>Batch code (5.1.5):</b> Indicates the manufacturer's batch code so that the batch or lot can be identified.
	<b>Catalog number (5.1.6):</b> Indicates the manufacturer's catalog number so that the medical device can be identified.
	<b>Serial number (5.1.7):</b> Indicates the manufacturer's serial number so that a specific medical device can be identified.
	<b>Non-sterile (5.2.7):</b> Indicates a medical device that has not been subjected to a sterilization process.
	<b>Keep away from sunlight (5.3.2):</b> Indicates a medical device that needs protection from light sources.
	<b>Keep dry (5.3.4):</b> Indicates a medical device that needs to be protected from moisture.
	<b>Temperature limit (5.3.7):</b> Indicates the temperature limits to which the medical device can be safely exposed.
	<b>Humidity limitation (5.3.8):</b> Indicates the range of humidity to which the medical device can be safely exposed.

Symbol	Definition
	<b>Atmospheric pressure limitation (5.3.9):</b> Indicates the range of atmospheric pressure to which the medical device can be safely exposed.
	<b>Consult instructions for use (5.4.3):</b> Indicates the need to consult the instructions for use.
	<b>Medical device (5.7.7):</b> Indicates the item is a medical device.

#### Product-specific symbols

Symbol	Definition
	<b>Quantity:</b> Indicates the number of products in the packaging.
	<b>Note:</b> Used to supplement or clarify information.
<b>GTIN</b>	Global Trade Item Number.

## IEC 60417 Graphical symbols for use on equipment

Symbol	Definition
	<b>Type BF applied part (5333):</b> To identify a type BF applied part complying with IEC 60601-1.

## 81 FR 38911 FDA Final rule for the use of symbols in labeling

Symbol	Definition
<b>Rx Only</b>	<b>Caution:</b> Federal law (USA) restricts this device to sale by or on the order of a physician.

## Directive 2012/19/EU on waste electrical and electronic equipment (WEEE)

Symbol	Definition
	Indicates that the product must be collected separately and must not be disposed of as unsorted municipal waste.

## Regulatory marks and logos

Symbol	Definition
	Indicates that a device is in conformity with the applicable requirements set out in applicable European Union harmonization legislation providing for its affixing.

## 2. Safety Information

---

### WARNING

- The product is delivered non-sterile. Before first use, as well as before each following use, the product must be processed according to a validated procedure.
- Before using a medical product, assure yourself of its functional safety and proper condition via a visual inspection. Do not use the product if any defects, such as breaks, cracks, deformations or wear, are detected. Particularly important parts such as tips, notches and all moveable components should be checked with extra care.
- Unauthorized modifications of the product are forbidden for safety reasons.
- Due to its composition, the product is not to be used together with other magnetically sensitive medical products, devices or instruments (e.g. MRI).

- The product is a high-precision instrument. Avoid subjecting it to serious strains, such as heavy impacts. The product cannot be used if there are visible defects. After a heavy impact, the product must be checked for defects.
  - The healthcare provider performing any procedure is responsible for determining the appropriateness of using the product and for the specific technique for each patient. Stryker, as a manufacturer, does not recommend a specific surgical procedure.
- 

 The user and/or patient should report any serious product-related incident to both the manufacturer and the national competent authority where the user and/or patient is established.

## 3. Product Information

---

### WARNING

- The product may only be used for its intended purpose and in accordance with this document and all current versions of relevant system and software application documentation.
  - This document is part of the product and must be accessible to personnel at all times. It must be provided to subsequent owners or users.
- 

### 3.1. Stryker ENT Navigation System

#### 3.1.1. United States / Rest of World\*

##### **Intended Use**

The Patient Tracker Electromagnetic – 10 Uses is an accessory to the Electromagnetic Navigation Unit and is intended for locating the patient.

##### **Indications for Use**

The Stryker ENT Navigation System is indicated for any medical condition in which the use of stereotactic surgery may be appropriate, and where reference to a rigid anatomical structure in the field of ENT surgery, such as the paranasal sinuses, mastoid anatomy, can be identified relative to a CT- or MR-based model of the anatomy.

Example procedures include, but are not limited to the following ENT procedures:

- Transsphenoidal access procedures
- Intranasal procedures
- Sinus procedures, such as maxillary antrostomies, ethmoidectomies, sphenoidotomies/ sphenoid explorations, turbinate resections, and frontal sinusotomies
- ENT-related anterior skull base procedures

### **Contraindications**

None known.

\*Note: Product may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Contact your Stryker representative for product availability.

## **3.2. Cranial Guidance System**

### **3.2.1. United States**

Refer to the user manual supplied with the Cranial Guidance Software (REF 6000-670-000) for system indications and contraindications.

### **3.3. User Group**

Healthcare professionals (surgeon/resident, nurse/professional caregiver) educated in computer-assisted surgery and thoroughly familiar with the instructions for use and with the operation of this product.

To request an additional in-service instruction, contact Stryker.

## 4. Product Overview

Figure 1: Patient Tracker Electromagnetic-10 uses

## 4.1. Instructions for Application

The patient tracker provides a reference for the patient during surgery. For this purpose, the tracker is securely fixed onto the patient with appropriate accessories as near as possible to the operation area.

---

### **WARNING**

- Prior to and during the surgical procedure, regularly verify that the patient tracker is firmly attached to the patient. The patient tracker must not move with respect to the patient. If the patient tracker moves, the system displays inaccurate positional information of the instruments. If you suspect that the patient tracker has moved, check landmarks to verify navigation accuracy.
- 

For instructions on how to use the product for a surgical procedure, refer to the instructions for use supplied with the respective software application.

## 5. For Use With

---

 **WARNING**

- Use only Stryker-approved products, unless otherwise specified.
- 

For information related to compatible software applications, refer to the user manual supplied with the respective software application. For information related to product-specific compatibility, refer to the table below.

Description	REF
Patient Tracker Tabs	8000-100-001

## 6. Disassembly

1. Carefully remove the patient tracker from the patient skin.
2. Remove the Patient Tracker Tab from the patient tracker, and dispose of the Patient Tracker Tab.

## 7. Product Lifetime

The instrument is validated for ten reprocessing cycles. Thus, the instrument contains a usage counter that disables the instrument after running the tenth usage cycle in surgery. The usage counter is updated as soon as the instrument is connected during surgery and shown in the software.

To maintain accurate use count, connect the instrument to the electromagnetic port at least once during surgery. This also applies to replacement instruments if they have been taken from sterile packaging.

## 8. Reprocessing

If you use the dedicated EM Instruments Insert Tray REF 8000-820-000\* for cleaning, disinfection and sterilization, refer to the Guide for Cleaning, Disinfection and Steam-Based Sterilization (TD6000005750) for reprocessing instructions, safety directives, and reprocessing equipment. The Patient Tracker Electromagnetic– 10 uses is designed for machine cleaning and disinfection like instruments without electronics. The cleaning group of the Patient Tracker Electromagnetic– 10 uses is: V.

If you use sterilization pouches for sterilization, refer to the Stryker ENT Navigation System Reprocessing Instructions (TD8000010706) for reprocessing instructions, safety directives, and reprocessing equipment.

## 9. Maintenance

Maintenance and repair may only be conducted by the manufacturer or authorized partners. The product components may only be sent back to the manufacturer in a cleaned, disinfected, and sterilized condition.

---

\* Note: Product may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Contact your Stryker representative for product availability.


## 10. Disposal

Risk of injury or infection can be avoided by safe disposal of the product components. They must be stored in such a way that they are protected from unauthorized use. Contaminated products are to be supplied to a hazardous waste site and handled in a way that contamination of third parties is excluded.

- After using the electromagnetic system, make sure that instruments that cannot be reprocessed again are disposed of.
- In accordance with European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) product should be collected separately for recycling. Do not dispose of as unsorted municipal waste. Contact local distributor for disposal information. Ensure infected equipment is decontaminated prior to recycling.
- The plugs used in the electronic products contain the following substance: Lead, CAS No. 7439-92-1 (1907/2006 REACH). Handling Instructions: No specific precautions are required for handling items manufactured from alloys containing lead in the supplied condition.
- The silicone glue used in the products contains the following substances: Decamethylcyclopentasiloxane, CAS No. 541-02-6 and Dodecamethylcyclohexasiloxane, CAS No. 540-97-6 (1907/2006 REACH).

# 11. Technical Specifications

Environmental Limitations	Operation	Storage and Transportation
Temperature	 A diagram showing a thermometer with a scale. The minimum value is labeled 10°C and the maximum value is labeled 30°C.	 A diagram showing a thermometer with a scale. The minimum value is labeled -10°C and the maximum value is labeled 50°C.
Relative humidity	 A circular icon containing a percentage sign and wavy lines representing water. A line points to the top of the circle labeled 75%, and another line points to the bottom of the wavy lines labeled 30%.	 An icon of an open umbrella with rain falling from it.
Atmospheric air pressure	 A circular icon with two arrows pointing towards a central dot. A line points to the top of the circle labeled 106 kPa, and another line points to the bottom of the circle labeled 80 kPa.	-

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Ограничения по условиям окружающей среды	Эксплуатация	Хранение и транспортировка
Температура воздуха	 A line graph showing a temperature range from 10°C to 30°C. A vertical line at 10°C has a small black circle at its base. A vertical line at 30°C has a small black circle at its top. A horizontal line connects the two vertical lines at the top, and a diagonal line connects the top of the 10°C line to the top of the 30°C line.	 A line graph showing a temperature range from -10°C to 50°C. A vertical line at -10°C has a small black circle at its base. A vertical line at 50°C has a small black circle at its top. A horizontal line connects the two vertical lines at the top, and a diagonal line connects the top of the -10°C line to the top of the 50°C line.
Относительная влажность	 A circular icon containing a percentage sign (%). A wavy line at the bottom represents humidity. A line points from the text '30%' to the wavy line. Another line points from the text '75%' to the top of the circle.	 An icon of an open umbrella with several raindrops falling from it.
Атмосферное давление	 A circular icon containing a central dot and two horizontal arrows pointing towards each other. A line points from the text '80 kPa' to the left side of the circle. Another line points from the text '106 kPa' to the right side of the circle.	-

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This page was intentionally left blank for your notes.

# stryker

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Leibinger, Stryker.

All other trademarks are trademarks of their respective owners or holders.

U.S. Patents: [www.stryker.com/patents](http://www.stryker.com/patents)

Copyright © 2022 Stryker

Stryker Leibinger GmbH & Co. KG

Bötzingen Straße 41

79111 Freiburg (Germany)

t: +49 761 4512 0 (Germany)

t: +1 269 323 7700 (USA)