# **Universal Tracker Electromagnetic**

REF 8000-060-006

Instructions for Use

Rx Only

EN   Instructions for use	1
DA   Brugsanvisning	18
DE   Gebrauchsanweisung	35
EL   Οδηγίες χρήσης	53
ES   Instrucciones de uso	71
FI   Käyttöohjeet	89
FR   Notice d'utilisation	106
IT   Istruzioni per l'uso	123
JA   使用説明書	141
KO   사용 설명서	158
NL   Gebruiksaanwijzing	175
NO   Bruksanvisning	193
PL   Instrukcja użycia	210
PT   Instruções de utilização	228
RO   Instrucțiuni de utilizare	246
SV   Bruksanvisning	264
TR   Kullanma Talimatı	282
ZH   使用说明	299

## 1. How to Use this Document

This manual is the most comprehensive source of information for the safe, effective, and compliant use and/or maintenance of the product. Read and understand this manual as well as the respective system user manual before using the product or any component compatible with the product. When combined with other medical devices, the user manual of these devices is to be considered as well. This manual is a permanent part of the product. Keep this manual for future reference.

The following conventions are used in this document:

The signal word **WARNING** highlights a safety-related issue. Comply with this information to prevent patient and medical staff injury.

The signal word **CAUTION** highlights a product reliability issue. Comply with this information to prevent product damage.

Supplements or clarifies information.

---

## 1.1. Symbol Definition

### EN ISO 7010 Graphical symbols – Safety colors and safety signs – Registered safety signs

Symbol

Name: Definition

W001

**General warning sign:** To signify a general warning.

M002

**Refer to instruction manual/booklet:** To signify that the user instruction manual/booklet must be read.

### EN ISO 15223-1 Medical devices — Symbols to be used with medical device labels, labeling and information to be supplied - Part 1 General requirements

Symbol

Name: Definition

5.1.1

**Manufacturer:** Indicates the medical device manufacturer as defined in the European Union harmonization legislation.

Symbol

Name: Definition

5.1.3

**Date of manufacture:** Indicates the date when the medical device was manufactured.

**LOT**

5.1.5

**Batch code:** Indicates the manufacturer's batch code so that the batch or lot can be identified.

**REF**

5.1.6

**Catalog number:** Indicates the manufacturer's catalog number so that the medical device can be identified.

**SN**

5.1.7

**Serial number:** Indicates the manufacturer's serial number so that the medical device can be identified.

5.2.7

**Non-Sterile:** Indicates a medical device that has not been subjected to a sterilization process.

5.3.2

**Keep away from sunlight:** Indicates a medical device that needs protection from light sources.

Symbol

Name: Definition

5.3.4

**Keep dry:** Indicates a medical device that needs to be protected from moisture.

5.3.7

**Temperature limit:** Indicates the temperature limits to which the medical device can be safely exposed.

5.3.8

**Humidity limitation:** Indicates the range of humidity to which the medical device can be safely exposed.

5.3.9

**Atmospheric pressure limit:** Indicates the range of atmospheric pressure to which the medical device can be safely exposed.

5.4.3

**Consult instructions for use:** Indicates the need to consult the instructions for use.

Indicates a medical device according to European Union harmonization legislation.

## Product-Specific Symbols

Symbol	Name: Definition
	<b>Quantity:</b> Indicates the number of medical devices in the packaging.
	<b>Note:</b> Used to supplement or clarify information.
<b>GTIN</b>	Global Trade Item Number.

## 81 FR 38911 FDA Final rule for the use of symbols in labeling

Symbol	Name: Definition
--------	------------------

<b>Rx Only</b>	<b>Caution:</b> Federal law (USA) restricts this device to sale by or on the order of a physician.
----------------	--

## Directive 2012/19/EU on waste electrical and electronic equipment (WEEE)

Symbol	Name: Definition
--------	------------------

	Indicates that the product must be collected separately and must not be disposed of as unsorted municipal waste.
--	--

## Regulatory marks and logos

Symbol	Definition
--------	------------

	Indicates that a device is in conformity with the applicable requirements set out in applicable European Union harmonization legislation providing for its affixing.
--	--

## 2. Safety Information

### WARNING

- The product is delivered in a non-sterile condition. Before first use, as well as before each following use, the product must be processed according to a validated procedure.
- Unauthorized modifications of the product are forbidden for safety reasons.
- Due to its composition, the product is not to be used together with other magnetically sensitive medical products, devices or instruments (e.g. MRI).
- Before using a medical product, assure yourself of its functional safety and proper condition via a visual inspection. Do not use the product if any defects, such as breaks, cracks, deformations or wear, are detected. Particularly important parts such as tips, notches and all moveable components should be checked with extra care.
- The product is a high-precision instrument. Avoid subjecting it to serious strains, such as heavy impacts. The product cannot be used if there are visible defects. After a heavy impact, the product must be checked for defects.
- The healthcare provider performing any procedure is responsible for determining the appropriateness of using the product and for the specific technique for each patient. Stryker, as a manufacturer, does not recommend a specific surgical procedure.

The user and/or patient should report any serious product-related incident to both the manufacturer and the national competent authority where the user and/or patient is established.

### 3. Product Information

#### WARNING

- The product may only be used for its intended purpose and in accordance with this document and all current versions of relevant system and software application documentation.
  - This document is part of the product and must be accessible to personnel at all times. It must be provided to subsequent owners or users.
- 

#### 3.1. Intended Use

The Universal Tracker Electromagnetic is an accessory to the Electromagnetic Navigation Unit and is intended for navigating conventional surgical instruments.

#### 3.2. Indications for Use

The Stryker ENT Navigation System is indicated for any medical condition in which the use of stereotactic surgery may be appropriate, and where reference to a rigid anatomical structure in the field of ENT surgery, such as the paranasal sinuses, mastoid anatomy, can be identified relative to a CT- or MR-based model of the anatomy.

Example procedures include, but are not limited to the following ENT procedures:

- Transsphenoidal access procedures
- Intranasal procedures
- Sinus procedures, such as maxillary antrostomies, ethmoidectomies, sphenoidotomies/sphenoid explorations, turbinate resections, and frontal sinusotomies
- ENT-related anterior skull base procedures

### **3.3. Contraindications**

None known.

### **3.4. User Group**

Healthcare professionals (surgeon/resident, nurse/professional caregiver) educated in computer-assisted surgery and thoroughly familiar with the instructions for use and with the operation of this product.

To request an additional in-service instruction, contact Stryker.

## 4. Product Overview


Figure 1: Universal Tracker Electromagnetic  
A Knurled Screw

## 4.1. Instructions for Application

With the help of an instrument clamp, the Universal Tracker Electromagnetic can be attached to surgical instruments that are not originally intended for navigation.

To mount the Universal Tracker Electromagnetic on an instrument clamp, follow the steps below.

1. Ensure that the instrument clamp is securely attached to the surgical instrument.
2. Insert the knurled screw (Figure 1, A) into the tracker and turn the knurled screw clockwise until the knurled screw is loosely connected with the tracker.
3. Attach the tracker to the instrument clamp.
4. Turn the knurled screw clockwise until the tracker is rigidly attached to the instrument clamp (Figure 2).


Figure 2: Mounting Universal Tracker Electromagnetic on Instrument Clamp

5. Calibrate the surgical instrument. Refer to the instructions for use supplied with the respective software application.

For instructions on how to use the product for a surgical procedure, refer to the instructions for use supplied with the respective software application.

## 5. For Use With

### WARNING

Use only Stryker-approved products, unless otherwise specified.

---

For information related to compatible software applications, refer to the instructions for use supplied with the respective software application. For information related to product-specific compatibility, refer to the table below.

Description	REF
Instrument Clamp, Forceps	8000-060-010
Instrument Clamp, 2-6 mm	8000-060-011
Instrument Clamp, 6-10 mm	8000-060-012
Instrument Clamp, 10-16 mm	8000-060-013

## 6. Disassembly

1. Unscrew and remove the knurled screw (Figure 1, A).
2. Detach the tracker from the instrument clamp.

## 7. Product Lifetime

The instrument is validated for ten reprocessing cycles. Thus, the instrument contains a usage counter that disables the instrument after running the tenth usage cycle in surgery. The usage counter is updated as soon as the instrument is connected during surgery and shown in the software.

To maintain accurate use count, connect the instrument to the electromagnetic port at least once during surgery. This also applies to replacement instruments if they have been taken from sterile packaging.

## **8. Reprocessing**

For reprocessing instructions, safety directives, and reprocessing equipment, refer to the Stryker ENT Navigation System Reprocessing Instructions (TD8000010706).

## **9. Maintenance**

Maintenance and repair may only be conducted by the manufacturer or authorized partners. The product components may only be sent back to the manufacturer in a cleaned, disinfected, and sterilized condition. Sharp or pointed components need to be sent back in a protected state.


## 10. Disposal

Risk of injury or infection can be avoided by safe disposal of the product components. Sharp and pointed components need to be collected and locked in a tight and break-proof container. They must be stored in such a way that they are protected from unauthorized use. Contaminated products are to be supplied to a hazardous waste site and handled in a way that contamination of third parties is excluded.

- After using the electromagnetic system, make sure that instruments that cannot be reprocessed again are disposed of.
  - In accordance with European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) product should be collected separately for recycling.
  - Do not dispose of as unsorted municipal waste. Contact local distributor for disposal information. Ensure infected equipment is decontaminated prior to recycling.
  - The plug used in the electronic products contains the following substance: Lead, CAS No. 7439-92-1 (1907/2006 REACH). Handling Instructions: No specific precautions are required for handling items manufactured from alloys containing lead in the supplied condition.
  - The silicone glue used in this product contains the following substances: Decamethylcyclopentasiloxane, CAS No. 541-02-6 and Dodecamethylcyclohexasiloxane, CAS No. 540-97-6 (1907/2006 REACH).
-

## 11. Technical Specifications

Environmental Limitations	Operation	Storage and Transportation
Temperature		
Relative humidity		
Atmospheric air pressure		-

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Stryker Corporation or its divisions or other affiliated entities own, use or have applied for the following trademarks or service marks: Leibinger, Stryker. All other trademarks are trademarks of their respective owners or holders.
Copyright © 2019-2021 Stryker
Stryker Leibinger GmbH & Co. KG