## Antibiotic Simplex® with Tobramycin

Howmedica Osteonics Corp.  
325 Corporate Drive  
Mahwah, NJ 07430, USA  
A subsidiary of Stryker Corporation

Telephone #: +1 201-831-5000

CE 0086

©2016 Howmedica Osteonics Corp.

0700-7-012 Rev AA

0700-7-012, Ver. AA ECN Release Date: Jan 10, 2017  
Print Date: Feb 22, 2017 16:50:18 GMT

## Labeling Symbols

Attention, see Instructions for Use

Do not Reuse

Sterilized using Irradiation

Sterilized using Hydrogen Peroxide

Sterilized using Ethylene Oxide

Sterilized using Aseptic Fill

Use by Date

Date of Manufacture

Legal Manufacturer

Authorized Representative in the European Community

Catalog Number

Batch Code

Serial Number

## CONTENTS

English .....	1
Français .....	6
Deutsch .....	11
Italiano .....	16
Español.....	20
Svenska .....	25
Nederlands.....	29
Suomi .....	34
Dansk .....	39
Português.....	43
Ελληνικά.....	48
Русский .....	53
Norsk .....	58
Íslenska .....	62
Polski .....	66
한국어 .....	71
中文.....	76
Türkçe.....	80

0700-7-012, Ver. AA ECN Release Date: Jan 10, 2017  
Print Date: Feb 22, 2017 16:50:18 GMT

## Instructions for Use

### ANTIBIOTIC SIMPLEX

#### Bone Cement, Radiopaque with Tobramycin

##### Active ingredient

Tobramycin (as sulphate)

##### DESCRIPTION

Antibiotic SIMPLEX Bone Cement with Tobramycin is packaged in two sterile components.

One component is an ampoule containing 10 or 20 ml of a colourless, flammable liquid monomer that has a sweet slightly acrid odour and has the following composition:

	Half dose	Full dose
Methyl methacrylate (monomer)	9.75 ml	19.5 ml
N, N-dimethyl para toluidine	0.25 ml	0.5 ml
Hydroquinone	0.75 mg	1.5 mg

Hydroquinone is added to help prevent premature polymerisation which may occur under certain conditions; e.g. exposure to light, elevated temperatures. N, N-dimethyl para toluidine, is added to help promote cold curing of the finished therapeutic compound.

The individual liquid component is sterilized by membrane filtration. Final packaged liquid is sterilized by ethylene oxide.

The other component is a packet of 20.5 or 41g of finely divided powder of the following composition:

	Half dose	Full dose
Methyl methacrylate-styrene copolymer (Contains 1.7 % Benzoyl Peroxide)	15.0g	30.0g
Polymethyl methacrylate	3.0g	6.0g
Barium Sulphate USP and EP	2.0g	4.0g
Tobramycin (as sulphate) USP	0.5g	1.0g

The individual powder component is sterilized by gamma irradiation. The antibiotic, prepared sterile, is blended with the cement component under sterile mixing and packing conditions.

The barium sulphate incorporated in Antibiotic SIMPLEX Bone Cement with Tobramycin acts as a contrast medium for X-ray examination.

At the time of use the powder and liquid are mixed. This mixture is designed to result in the exothermic polymeric formation of a soft, pliable, dough-like mass. It is designed so that as the reaction progresses, within a few minutes a hard, cement-like complex is formed.

##### INDICATIONS

Antibiotic SIMPLEX Bone Cement with Tobramycin is indicated as a bone cement in prosthetic surgery.


##### DOSAGE

As required for adequate fixation of the prosthesis.

##### PREPARATION AND ADMINISTRATION

The liquid monomer, the ampoule itself and the ampoule package have been pre-sterilized. If any packaging is damaged, do not use this product.

Before mixing Antibiotic SIMPLEX Bone Cement with Tobramycin, care should be taken to allow the cement components to reach operating room temperature. By referring to the following setting time chart, a reasonable indication of the cement's mixing and handling characteristics for the particular operating room temperature can be determined. The surgeon must take this information into account when using the product.

To mix, empty the entire contents of the packet containing the powder component into an appropriate sterile, inert mixing device. Add the entire contents of the ampoule containing the liquid component. Add the liquid component to the powder, not the powder to the liquid. Mix with a suitable inert device with a stirring action until the powder is completely saturated with the liquid.

For cement gun administration - stir thoroughly for a total of 1 - 2½ minutes. Introduce the liquid cement into a suitable sterile injection device and apply product to bone.

For manual administration, continue stirring until a dough-like mass is formed that does not stick to the surgical gloves of the operator. The dough-like mass is then ready for manipulation. The mixing and manipulation process should be at least 4 minutes in duration. The mixing and kneading time required, to obtain a product of proper consistency for application to the bone, is altered by the room temperature and humidity. The mixed and kneaded product is then ready to be applied to the bone manually.

The correct working consistency of the cement for application to bone is best determined by the experience of the surgeon.

After application of the product, position of the prosthesis must be maintained securely without movement until the cement is hard and the prosthesis is firmly fixed in position. Excess cement must be removed before it is totally hardened.

#### **CONTRAINdicATIONS**

The use of Antibiotic SIMPLEX Bone Cement with Tobramycin is contraindicated in patients who are allergic to any of its components including Tobramycin Sulphate. It is contraindicated in the younger patient with a normal life expectancy because of lack of information on the long-term effect. It is also contraindicated where loss of musculature or neuromuscular compromise in the affected limb would render the procedure unjustifiable.

The concomitant use of tobramycin with certain ototoxic loop diuretics including furosemide and ethacrynic acid is contraindicated.

Tobramycin is contraindicated in patients with myasthenia gravis.

#### **WARNINGS, PRECAUTIONS AND INTERACTIONS**

For safe and efficacious use of Antibiotic SIMPLEX Bone Cement with Tobramycin, the surgeon must have specific training and experience to be thoroughly familiar with the properties, handling characteristics and application of the product.

Care must be taken in the mixing of the liquid and powder components that the entire contents of the ampoule and pouch are used. The mixing of the liquid monomer and the powder component must be thorough and at least four minutes in duration. Data from clinical studies have shown by gas chromatography determinations, that the mixing and kneading of the liquid and powder components for four or more minutes depressed the level of monomer in the circulating blood and also in the urine.

Caution, however, must be taken to avoid kneading of the product too long to avoid progression of the polymerisation process to the point that the product is not adequately soft and pliable to obtain good filling of the bone cavities and fitting to the prosthesis.

The completion of polymerisation occurs in the patient and is an exothermic reaction with considerable liberation of heat. Temperatures occurring during the polymerisation have been reported as high as 110°C. Registry data for cemented hip and knee replacement have not indicated an adverse long-term effect of the heat produced along with the resulting tissue damage.

Special precautions must be taken to detect and correct the transitory fall in blood pressure that may occur when the product is implanted into the bone.

As the liquid monomer is highly volatile and inflammable, the operating room must be provided with adequate ventilation so as to eliminate the maximum amount of monomer vapour. Caution must be exercised during the mixing of the two components to prevent excessive exposure to the concentrated vapours of the monomer, which may produce irritation of the respiratory tract, eyes and possibly the liver.

The liquid component is a powerful lipid solvent. It has caused contact dermatitis in susceptible individuals. Wearing a second pair of surgical gloves and strict adherence to the mixing instructions may diminish the possibility of hypersensitivity reactions. The compound must not be allowed to come into contact with sensitive tissues or to be absorbed by the body.

It has been recommended by manufacturers of soft contact lenses that such lenses should be removed "in the presence of noxious and irritating vapors." Since soft contact lenses are quite permeable, they should not be worn in an operating room where Methyl methacrylate is being mixed.

If tobramycin, gentamicin or other aminoglycosides have been administered to the patient prior to surgery, monitoring of trough serum concentrations should be performed on the day before operation. If serum concentrations exceed 1 µg/ml, tobramycin, gentamicin or other aminoglycosides, Antibiotic SIMPLEX Bone Cement with Tobramycin must not be used.

The ototoxicity of tobramycin is enhanced by the use of certain loop diuretics such as frusemide and ethacrynic acid.

The concomitant use of tobramycin and neuromuscular blocking agents may lead to neuromuscular blockage and respiratory paralysis. This effect may be reversed by calcium salts.

Since the concentration of tobramycin reaching the VIII nerve and the kidney will be quite low, ototoxic and nephrotoxic reactions are unlikely to arise from the use of Antibiotic SIMPLEX Bone Cement with Tobramycin.

Hypersensitivity to tobramycin is rare.

Use in pregnancy: Reproduction studies have not been performed in animals. There is no adequate information on whether this product may affect fertility in human males or females or have teratogenic potential or other adverse effect on the fetus.

#### **INFORMATION FOR THE PATIENT**

- Surgeons must advise patients of the limitations of the reconstruction and the need to protect the implant from full weight bearing until adequate healing has occurred.
- Surgeons must advise patients that the product has a finite service life and may need to be replaced in the future.
- Surgeons must warn the patient of surgical risks and possible adverse effects.
- Dental procedures, endoscopic examinations and other minor surgical procedures have been associated with transient bacteraemia. Instruct the patient to inform their doctors that they have an artificial hip replacement, so that their doctors can decide whether to use antibiotic prophylaxis for such procedures.

## **ADVERSE REACTIONS**

The most serious adverse reactions, some with fatal outcome, reported with use of acrylic bone cements are:

- Cardiac arrest
- Myocardial infarction
- Pulmonary embolism
- Cerebrovascular accident

The most frequent adverse reactions reported are:

- Transitory fall in blood pressure
- Thrombophlebitis
- Haemorrhage and haematoma
- Loosening or displacement of the prosthesis
- Surgical wound infection
- Deep wound infection
- Trochanteric bursitis
- Trochanteric separation

Other adverse reactions reported are:

- Heterotopic new bone
- Short-term irregularities in cardiac conduction.

## **IMPORTANT PHYSICIAN INFORMATION**

Adverse reactions affecting the cardiovascular system appear to be related to the liquid monomer rather than the powder polymer component of the final product. It appears the more monomer available on the surface of the final product for systemic absorption, the greater the risk of hypotension and other cardiovascular reactions. The hypotension seen appears to be the result of peripheral vasodilation. The reported fat, bone marrow and air emboli appear to be the result of the direct pressure from the forcing of the cement into the "raw" medullary cavity. The aetiology of cardiac arrest is unclear, but may well be either a direct embolic effect or secondary to hypoxia produced by pulmonary embolic phenomenon.

The degree of hypotension seen appears to be more marked in patients with elevated or high normal blood pressures, in hypovolaemic states, and in those patients with pre-existing cardiovascular abnormalities. Age, sex, and anaesthetic agent or technique do not appear, according to the available literature, to be significant in the aetiology of the hypotensive reaction. The degree of hypotension reported in the literature varies from an average of 7-20% or more of pre-operative level; duration reported varies from 30 seconds to 5-6 minutes; onset following application of the bone cement varies from 10-165 seconds.

## **PHARMACEUTICAL PRECAUTIONS**

Store in dark at temperature below 25°C.

Do not use this product after its expiry date.

## **HOW SUPPLIED**

Individual unit or ten-unit dispenser carton, each unit containing:

### **Half Dose**

- 1 STERILE PACKET CONTAINING 20.5g of sterile powder.
  - 1 STERILE AMPOULE CONTAINING 10 ml of sterile liquid monomer.
- Full Dose**
- 1 STERILE PACKET CONTAINING 41g of sterile powder
  - 1 STERILE AMPOULE CONTAINING 20 ml of sterile liquid monomer.

## **SAFE DISPOSAL**

Mixed cement should be allowed to set before disposal with other clinical waste. For disposal of liquid or powder separately, contact your local disposal authority.

Product should not be used after the expiry date displayed on the label as packaging had not been validated beyond this date.

**CAUTION**

Federal Law in the USA restricts this device to sale by or on the order of a physician.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademark(s) or service mark(s): Howmedica, Osteonics, Simplex, Stryker. All other trademarks are trademarks of their respective owners or holders.

Refer to product label for CE Mark Status and Legal Manufacturer. The CE mark is only valid if also found on the product label.

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