stryker®
Howmedica
OSTEONICS

**TOTAL HIP JOINT REPLACEMENT PROSTHESES
FOR CEMENTLESS & CEMENTED APPLICATIONS**

Howmedica Osteonics Corp.
325 Corporate Drive
Mahwah, NJ 07430, USA
A subsidiary of Stryker Corporation

Telephone #: +1 201-831-5000

CE 0086

©2016 Howmedica Osteonics Corp.

96E112 Rev. AA

Refer to product label for CE Mark Status and Legal Manufacturer.
The CE mark is only valid if also found on the product label.

Labeling Symbols

Attention, See Instructions for Use

Do Not Reuse

Sterilized using Irradiation

Sterilized using Hydrogen Peroxide

Sterilized using Ethylene Oxide

Use by Date

Date of Manufacture

Legal Manufacturer

Authorized Representative in the European Community

Catalog Number

Batch Code

Serial Number

TOTAL HIP JOINT REPLACEMENT	
PROSTHESES FOR CEMENTLESS & CEMENTED APPLICATIONS.....	1
ARTHROPLASTIE TOTALE DE LA HANCHE	
PROTHÈSES POUR DES APPLICATIONS CIMENTÉES ET NON CIMENTÉES	8
HÜFTGELENKTOTALERSATZ	
FÜR ZEMENTFREIE UND ZEMENTIERTE ANWENDUNGEN	16
ARTROPLASTICA TOTALE DELL'ARTICOLAZIONE DELL'ANCA	
PROTESI PER APPLICAZIONI CON E SENZA CEMENTO.....	24
REEMPLAZO TOTAL DE LA ARTICULACIÓN DE LA CADERA	
PRÓTESIS PARA APPLICACIONES SIN Y CON CEMENTO	32
PRÓTESE PARA SUBSTITUIÇÃO TOTAL DA ARTICULAÇÃO DA ANCA	
PRÓTESES PARA APPLICAÇÕES COM E SEM CIMENTO.....	40
TOTAL HÖFTLEDSPROTES	
FÖR CEMENTFRÄA OCH CEMENTERADE APPLIKATIONER.....	48
SEMENTÖITÄVÄT JA ILMAN SEMENTTÄÄ KIINNITETÄVÄT	
LONKAN KOKOTEKONIVELET.....	55
TOTAL HOFTLEDSALLOPLASTIKPROTESER	
TIL ANVENDELSE MED OG UDEN CEMENT.....	63
TOTALE HEUPVERVANGING PROTHESEN	
VOOR ONGECEMENTEERDE & GECEMENTEERDE TOEPASSING.....	70
INNSATTING AV TOTAL HOFTLEDSPROTESE	
FOR SEMENTLØS OG SEMENTERT ANVENDELSE	78
CAŁKOWITE ENDOPROTEZY STAWU BIODROWEGO	
DO OSADZANIA METODĄ BEZCEMENTOWĄ I CEMENTOWĄ.....	85
ΠΡΟΘΕΣΕΙΣ ΟΛΙΚΗΣ ΑΡΘΡΟΠΛΑΣΤΙΚΗΣ ΙΞΙΟΥ	
ΓΙΑ ΕΦΑΡΜΟΓΕΣ ΜΕ Ή ΧΩΡΙΣ ΤΣΙΜΕΝΤΟ	94
无骨水泥和骨水泥固定型全髋关节置换假体应用	103
무시멘트 및 시멘트 접합을 위한 완전 고관절 대치 보철	111
セメント無使用およびセメント使用のための	
全股関節置換プロテーゼ	120

English

TOTAL HIP JOINT REPLACEMENT PROSTHESES FOR CEMENTLESS & CEMENTED APPLICATIONS

DESCRIPTION

Howmedica Osteonics Corp. manufactures a wide variety of reconstructive total hip systems to satisfy anatomical requirements and surgeon/patient needs. Porous coated implants are generally intended for cementless use in primary applications where there is adequate bone stock, and for cemented use in revision procedures. Hydroxyapatite-coated implants are intended for cementless use. See package label for specific product features and the appropriate hip system surgical protocol for additional procedural information and product information.

Femoral Stems. Femoral stems are available in straight or anatomic designs, collared or collarless; with and without coatings; and with or without cement mantles

Femoral Heads. Femoral heads are available in a wide range of outside diameters and several tapers to mate with varying acetabular components.

Acetabular Cups. Acetabular cups are available in a wide range of outer diameter sizes and several inner diameters. Cups include one-piece all polyethylene (may include an x-ray location wire); one-piece metal backed/polyethylene; and two-piece modular cups consisting of a metal shell requiring a polyethylene or ceramic insert. Cups are available with and without screw holes and coatings. Some acetabular cups are available with polymethylmethacrylate (PMMA) spacers.

Acetabular Inserts. The inserts are available with varying inner diameters to help accommodate the mating femoral head. The inserts have a snap-lock mechanism to securely lock into the metal shell. The UHMWPE insert may contain an x-ray location wire. Alumina inserts are also available.

Surgical Mesh. Surgical Mesh is intended to be used to reinforce bone or tissue in any situation where additional strengthening and support is required due to poor bone/tissue quality.

Hip system implants are manufactured from materials that meet the following ASTM and/or ISO Standards:

Material	ASTM Standard / ISO Standard	Implant or Implant Component
Vitallium Alloys		
Cast	F75 / 5832-4	Stems, Heads, Cups, Beaded coating
Forged	F799	Stems, Modular Necks, Heads, Cups
Wrought	F1537 / 5832-12	Stems, Cups. X-Ray Wire
Gas Atomized Dispersion Strengthened (GADS)	F1537	Stems, Modular Necks, Cups
Forged Titanium Alloy (Ti 6Al-4V)	F136/ 5832-3	Cups, Stems, X-ray location wire
TMZF (Ti 12Mo-6Zr-2Fe)	F1813	Stems, Modular Stems
Alumina (Al ₂ O ₃)**	F603 / 6474	Heads
Hydroxyapatite (powder)	F1185	Coating
PMMA (polymethylmethacrylate)	D788/ 8257-1	Stems (mantles), Cups (spacers)
Stainless Steel Alloys		
Stainless Steel	F138 / 5832-1	Cups, Heads, Stems, Surgical Mesh
ORTHINOX Stainless Steel	5832-9	Cups, Heads, Stems, X-ray location wire
Tantalum (beads)	F560	RSA*** spikes
Titanium Commercially Pure (CP)	F1580 / 5832-2	Coating
UHMWPE (polyethylene)	F648/ 5834-2	Cups, Acetabular Inserts

COMPATABILITY

Howmedica Osteonics Corp. manufactures femoral heads with several tapers and diameters in both metal and ceramic**. Care must be taken to use stems and heads with the same taper angle and diameter. See product label for the specific taper information, such as angle and diameter.

Ceramic* inserts are to be used exclusively with HOWMEDICA OSTEONICS Alumina Ceramic** Heads.

Based upon laboratory testing, the +8mm and larger modular heads are not recommended for use on all HOWMEDICA OSTEONICS Hip Stems. See additional labeling for the specific stems that are contraindicated for use with these heads.

See maximum femoral head offset warning on package label and HIPSTAR trunnion.

HOWMEDICA OSTEONICS Ceramic** Femoral Heads are intended to mate only with HOWMEDICA OSTEONICS Femoral Stem Components labeled for use with ceramic** heads and with the correct taper angle and diameter. HOWMEDICA OSTEONICS Ceramic** Heads must not be used with other manufacturers' femoral stem components. HOWMEDICA OSTEONICS Ceramic** Femoral Heads should not be used in procedures where only the femoral head component of a total hip replacement is being revised.

Screws/spikes are manufactured from materials to match that of the acetabular cups, and screw/cup materials should not be mixed.

HOWMEDICA OSTEONICS inserts and acetabular cups made of UHMWPE are compatible (except when specifically mentioned) with all HOWMEDICA OSTEONICS femoral heads of the same diameter.

INDICATIONS

The indications for use for total hip and hemi hip arthroplasty include:

1. noninflammatory degenerative joint disease including osteoarthritis and avascular necrosis;
2. rheumatoid arthritis;
3. correction of functional deformity;
4. revision procedures where other treatments or devices have failed; and,
5. treatment of nonunion, femoral neck and trochanteric fractures of the proximal femur with head involvement that are unmanageable using other techniques.

INDICATIONS FOR RESTORATION DLS ARE LIMITED TO:

The indications for use for total hip and hemi hip arthroplasty include:

1. revision procedures where other treatments or devices have failed

See constrained acetabular package insert (QIN 4357) for indications for constrained acetabular inserts.

HIPSTAR Femoral Stems Cups are intended for cementless use only.

In the USA, when mated with a constrained acetabular liner the HIPSTAR Hip Stem is indicated for use in primary and revision total hip arthroplasty for patients at high risk of hip dislocation due to a history of prior dislocations, bone loss, joint or soft tissue laxity, neuromuscular disease or intraoperative instability.

CONTRAINDICATIONS

1. active infection or suspected latent infection in or about the hip joint;
2. bone stock that is inadequate for support or fixation of the prosthesis
3. skeletal immature
4. any mental or neuromuscular disorder that would create an unacceptable risk of prosthesis instability, prosthesis fixation failure, or complications in postoperative care

Additional conditions presenting increased risk of failure pertaining to the X-Change Mesh include:

1. Compromised vascularity that would inhibit adequate blood supply to the operative site.
2. Malignancy in the surgical area

WARNINGS

In using this system, the surgeon should be aware of the following:

- In selecting patients for total joint replacements, the following factor is of extreme importance to the eventual success of the procedure: The patient's weight. The heavier the patient, the greater the load on the prosthesis. As the loads on the prosthesis increase, the chance a patient will suffer adverse reactions increases, including but not limited to failure of fixation, loosening, fracture and dislocation of the device and can lead to a decreased service life. The effect of these loads will be accentuated when a small sized prosthesis is used in larger patients. Overweight or obese patients impose greater loads on the prosthesis. As obesity is a clinical diagnosis, we leave it to the surgeon to make the diagnosis based on his/her own clinical judgment. However, the World Health Organization (WHO) defines "overweight" as a BMI equal to or more than 25, and "obesity" as a BMI equal to or more than 30.
- Do not substitute another manufacturer's device, except if identified in Compatibility section above, for any component of the HOWMEDICA OSTEONICS Total Hip System. Any such use will negate the responsibility of Howmedica Osteonics Corp. for the performance of the resulting mixed component implant
- **DO NOT IMPACT** directly on the integral cement spacer.
- Acetabular bone screws—If bone screws or spikes are used:
 - Select screw/spike length and placement carefully to avoid damaging underlying soft tissue structures

Perforation of the pelvic wall can cause internal bleeding and possible damage to vital organs

- Ensure that bone screws/spikes are fully seated to avoid interference with the liner component
- Modular Acetabular Shell/Liner. Clear all surgical debris from inside the shell component before assembling the liner. Debris may prevent proper liner assembly, and cause dissociation of the liner from the shell.
- Modular Junctions: Mate modular components firmly to prevent dissociation. Machined taper surfaces must be clean, dry and firmly mated to ensure proper seating and assembly. Repeated assembly/disassembly or failure to clean, dry and firmly mate the components could compromise the taper lock and contribute to wear/corrosion and significant clinical consequences to the patient. (See Adverse Effects below).
- Skirted and 22.2 mm femoral heads may result in reduced range of motion which may increase risk of dislocation and disassociation.
- Do not impact (hammer) directly on the ceramic** head. HOWMEDICA OSTEONICS Ceramic** Heads should not be used if dropped or damaged during preparation.
- Patient post-operative pain. Inherent to all joint replacement is the risk that a patient will develop post-operative pain; pain is a commonly reported symptom regardless of the device implanted. The clinical literature reveals numerous potential causes of pain not directly related to the implant performance including, but not limited to, prior history of trauma and natural disease progression.
- For patients who present with pain following implantation of any orthopedic implant system, the physician should consider all potential causes of the symptoms identified in the clinical literature, including infection, soft tissue impingement, and possible adverse local tissue reactions associated with wear debris, metal ions or corrosion. Accurate diagnosis of the source of pain and directed, timely intervention is essential to ensuring effective treatment of pain.

PRECAUTIONS

- Before clinical use, the surgeon must thoroughly understand all aspects of the surgical procedure, limitations of the device, instruments and implant characteristics, prior to performing surgery.
- Appropriate selection, placement, positioning, and fixation of the total hip components are critical factors which affect implant service life. Proper implant selection must consider design, fixation, and environmental variables including: patient weight, age, bone quality and size, activity level and pre-operative level of health,

as well as the surgeon's experience and familiarity with the device. Strict adherence to the indications, contraindications, precautions and warnings for this product is essential to potentially maximize service life.

- Periodic, long-term follow-up is recommended to monitor the position and condition of the prosthetic components, as well as the condition of the adjoining bone.
- The smaller, sized femoral stem implants are intended for use in patients with smaller intramedullary femoral canals. Their geometry has been reduced to accommodate the anatomy of the smaller intramedullary femoral canal, which thereby decreases their fatigue-strength and load-bearing characteristics. Therefore, patients with high physical activity levels, poor bone quality, or who are overweight are not candidates for the smaller femoral implant stem.
- Patients with high-activity level and/or higher weight patients are at greater risk for implant complications or failures. For patients with poor proximal bone quality, the use of supplemental adjunctive proximal fixation/support is advised for implant stability.
- Single use devices cannot be explanted and subsequently reimplanted as the physical forces exerted by these actions may compromise the physical integrity, dimensions and/or surface finishes of the devices. Also, sterility cannot be assured for reused devices as cleaning and re-sterilization procedures have not been verified.
- Do not allow coated surfaces to contact cloth or other fiber-releasing materials.
- Protect polished bearing areas and machined taper surfaces from contact with hard or abrasive surfaces.
- Avoid handling hydroxyapatite coated regions, as it may compromise the effectiveness of the device.
- Radiographic templates are available to assist in the preoperative prediction of component size and style.
- To preserve the integrity of the actual implants and their sterile packaging, use the recommended trial components for size determination, trial reduction and range-of-motion evaluation.
- Surgeons should warn patients with metallic implants of the potential risks of undergoing a Magnetic Resonance Imaging (MRI) scan. The electromagnetic field created by an MRI scanner can interact with the metallic implant, resulting in displacement of the implant, heating of the tissue near the implant, implant damage or malfunction, or other undesirable effects. In addition, the presence of a metallic implant can produce an image artifact that may appear as a void region or geometric distortion of the true image. If the image artifact is near the area of interest, it may make the MRI scan uninformative or may lead to inaccurate clinical diagnosis or treatment.

INTERACTION with MAGNETIC RESONANCE IMAGING

The Stryker Hip Systems have not been evaluated for safety and compatibility in the MR environment. The Stryker Hip Systems have not been tested for heating or migration in the MR environment.

ADVERSE EFFECTS

1. Although rare, sensitivity/allergic reactions to the materials in the implant have occurred in patients following joint replacement. Implantation of foreign material in tissues can result in immune responses and in histological reactions involving macrophages and fibroblasts.
2. Peripheral neuropathies, nerve damage, circulatory compromise, and heterotopic bone formation may occur.
3. Dislocation of the hip prosthesis can occur due to inappropriate patient activity, trauma or other biomechanical consideration. Muscle and fibrous tissue laxity can also contribute to these conditions.
4. Loosening of total hip components can occur. Early mechanical loosening may result from inadequate initial fixation, latent infection, premature loading of the prosthesis or trauma. Late loosening may result from trauma, infection, biological complications, including osteolysis, or mechanical problems, with the subsequent possibility of bone erosion and/or pain.
5. Undesirable shortening or lengthening of the limb.
6. Infection can lead to failure of the joint replacement.
7. Fatigue fracture of femoral stems has occurred in a small percentage of cases. Stem fracture is more likely to occur in heavy, physically active patients, or when contralateral joint disability results in a disproportionate distribution of weight on the reconstructed joint.

8. Intraoperative fissure, fracture, or perforation of the femur, acetabulum, or trochanter can occur due to impaction of the component into the prepared femoral canal or acetabulum. Postoperative femoral or acetabular fracture can occur due to trauma, the presence of defects, or poor bone stock.
9. Polyethylene particles and metal particles from mechanisms other than wear. Very small particles from metal and polyethylene components can be shed from non-articulating surfaces during normal use and over time. Although most of these particles stay in the relevant joint (i.e. contained in the synovium) or are trapped by surrounding scar tissue, microscopic particles can migrate throughout the body and on occasions have been described as accumulating in lymph nodes and other parts of the body. Although no significant medical complications have been reported as a result of these particles, their migration and/or accumulation in the body have been described in the literature. The long-term effects, if any, from these particles, are unknown. The long-term effects have been theorized to include:
 - Cancer: There is presently no scientific evidence that links metallic or polyethylene particles with cancer. However, the possibility cannot be ruled out.
 - Lymphadenopathy and Accumulation in Other Tissues/Organs: There have been a few reports of the accumulation of particles in lymph nodes (proximal and distal). Although no medical complications or disease process has been reported as stemming from these accumulations, their existence should be recognized to facilitate diagnosis and avoid confusion with suspicious lesions, cancerous or otherwise.
 - Systemic Disease: It is possible that some long-term effects may be demonstrated at some point in the future, but because there is very little scientific data suggesting association between migration of particles and systemic disease, it is believed that the benefits of these devices clearly outweigh the potential risks for any such theoretical long-term effect.
10. Non-metallic wear debris. Wear debris is generated by interaction between components, as well as between components and bone, primarily through wear mechanisms of adhesion, abrasion and fatigue. Secondarily, particulate can also be generated by third-body wear. With all implant devices, asymptomatic, localized progressive bone resorption (osteolysis) may occur around the prosthetic components as a consequence of foreign-body reaction to the particulate matter of cement, ultra-high molecular weight polyethylene (UHMWPE) and/or ceramic. Osteolysis can lead to future complications, including loosening, necessitating the removal and replacement of prosthetic components.
11. Metal wear debris. Metal wear debris, metal ions and corrosion of metal implants. Generation of metal wear debris, metal ions and/or corrosion occurs whenever two surfaces are in contact, and at least one surface is metal. There have been reports in the literature of cases of adverse local tissue reactions associated with wear and/or corrosion at modular junctions formed by modular heads (stem/head interface). Local joint chemistry and/or other patient-specific conditions including, but not limited to infection may affect the potential for *in vivo* corrosion and its possible clinical consequences.
12. Metal debris and corrosion by-products. Modular junctions may release metal debris and/or metal ions due to fretting, galvanic corrosion, crevice corrosion, or other processes. There are several factors involved in these processes, including forces across a junction, which are not fully understood. These corrosion products or metal debris may affect the tissues surrounding the implant and adversely affect the duration of service life. There have been reports of patients developing adverse local tissue reactions (including, but not limited to, tissue necrosis, pseudotumors, cysts and fluid accumulations, metallosis and aseptic lymphocyte dominated vasculitis associated lesions), elevated metal ion levels in the blood and/or wear and hypersensitivity/allergic reactions associated with corrosion and /or wear-related debris in the implant vicinity. Affected patients may present with symptoms similar to those associated with infection, including pain (most likely during weight-bearing) and swelling at the local joint site. These reactions should be carefully monitored and may result in early revision surgery. Medical literature describes systemic reactions to by-products arising from contemporary metal on metal articulating bearing surfaces. While there are no metal on metal bearing interfaces in this system, it could be theorized that similar systemic reactions may arise from fretting and corrosion arising from any metal interface.
13. Serious complications may be associated with any total joint replacement surgery. These complications include, but are not limited to: genitourinary disorders, gastrointestinal disorders, vascular disorders including thrombus, bronchopulmonary disorders including emboli, myocardial infarction, or death.
14. Adverse effects may require medical intervention including reoperation, revision, arthrodesis of the involved joint, Girdlestone, or amputation of the limb

INFORMATION FOR PATIENTS

- The surgeon must advise patients of both the limitations of the reconstruction and the need for protection of the implant from full weight bearing until adequate fixation and healing have occurred. Excessive activity and trauma affecting the joint replacement have been implicated in failure of the reconstruction by loosening, fracture and/or wear of the prosthetic implants. Loosening of the components can result in increased production of wear particles, as well as damage to the bone, making successful revision surgery more difficult.
- The surgeon must caution patients to limit activities and protect the replaced joint from unreasonable stresses, and to follow the instructions of the physician with respect to follow-up care and treatment.
- The surgeon must warn the patient of surgical risks and all of the possible adverse effects detailed in the "Adverse Effects" section above. The patient should be informed of the limitations of biomaterials and how they will react in the body. It is understood that metal components can wear and/or corrode when in contact with hard surfaces. While the body's reaction to such biomaterials is different for every individual, there exists potential for local and systemic effects.
- The surgeon must warn patients that the device does not replicate the flexibility, strength, reliability, or durability of a normal healthy joint, that the implant can break or become damaged as a result of strenuous activity or trauma, and that the device has a finite service life and may need to be replaced in the future.
- Transient bacteremia can occur in daily life. Dental manipulation, endoscopic examination and other minor surgical procedures have also been associated with transient bacteremia. To prevent infection at the implant site, it may be advisable to use antibiotic prophylaxis before and after such procedures.
- If the patient is involved in an occupation or activity which includes substantial walking, running, lifting, or muscle strain, the resultant forces can cause failure of the fixation, the device, or both. The prosthesis will not restore function to the level expected with normal healthy bone. The surgeon must counsel the patient against having unrealistic functional expectations.

HOW SUPPLIED

These components have been sterilized by gamma radiation or hydrogen peroxide. Refer to the package label for the specific sterilization method.

- Do NOT resterilize.
- Take care to prevent contamination of any components.
- Inspect the packaging of sterile products for flaws before opening. In the presence of any flaws, assume that the product is nonsterile.
- Discard ALL nonsterile or contaminated products.
- Device should not be used after the expiry date displayed on the label as packaging had not been validated beyond this date.

Caution: Federal Law in the USA restricts this device to sale by or on the order of a physician.

* Consult the TRIDENT Acetabular Component System IFU (QIN 4365) for indications, contraindications, warnings, and precautions on the TRIDENT Alumina ceramic insert.

** Consult the TRIDENT Acetabular Component System IFU (QIN 4350) for indications, contraindications, warnings, and precautions on the ceramic femoral heads.

*** Roentgen Stereophotogrammetric Analysis (RSA) spikes are not available for sale in the United States.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademark(s) or service mark(s): ABG, Hipstar, Howmedica, Orthinox, Osteonics, Stryker, Stryker Orthopaedics, TMZF, Trident, Vitallium. All other trademarks are trademarks of their respective owners or holders.

BIOLOX® delta is a registered trademark of CeramTec AG.

Refer to product label for CE Mark Status and Legal Manufacturer. The CE mark is only valid if also found on the product label.

The following table contains a list of abbreviations that are used on Howmedica Osteonics Corp. product labeling:

Term	Abbreviation	Term	Abbreviation
Alpha Code	ALPH CDE	Neck	NK
Angle	ANG	Offset	OFFST
Degree	DEG or °	Outer Diameter	OD
Diameter	DIA	Right	RT ►
Extra Deep	XDP	Screw Holes	SCR HLS
Extra Large	XLGE	Side	SDE
Extra Small	XSM	Size	SZE
Head	HD	Small	SM
Height	HT	Standard	STD
Inner Diameter	ID	Taper	TPR
Insert	INSR	Thickness	THKNS
Large	LGE	Type	TYP
Left	◀ LFT	With	W/
Length	LNTH	Without	W/O
Medium	MED		

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