Total Hip Replacement Acetabular Component

Howmedica Osteonics Corp.
325 Corporate Drive
Mahwah, NJ 07430, USA
A subsidiary of Stryker Corporation

Telephone #: +1 201-831-5000

CE 0086

©2016 Howmedica Osteonics Corp.

QIN4301 Rev. AA

Refer to product label for CE mark status and Legal Manufacturer.
The CE mark is only valid if also found on the product label.

Labeling Symbols

Attention, see Instructions for Use

Do not Reuse

Sterilized using Irradiation

Sterilized using Hydrogen Peroxide

Sterilized using Ethylene Oxide

Use By Date

Date of Manufacture

Legal Manufacturer

Authorized Representative in the European Community

Catalog Number

Batch Code

Serial Number

CONTENTS

English	1
Français	9
Deutsch	17
Italiano	26
Español	34
Nederlands	42
Svenska	51
日本語	59
Suomi	67
Dansk	75
Português	83
Ελληνικά	91
Norsk	100
Polski	108
한국어	117
中文	125

English

TOTAL HIP REPLACEMENT ACETABULAR COMPONENT

The advancement of partial and total hip replacement has provided the surgeon with a means of restoring mobility and reducing pain with the use of implanted prosthetic devices. While these devices have proven to be largely successful in obtaining these goals, they are manufactured from metal, plastic, or other biomaterials. Any partial or total hip replacement system, therefore, cannot be expected to withstand the same activity and loads of normal healthy bone. The system will not be as strong, reliable or durable as a natural human hip joint and does not have an infinite lifetime. The surgeon must warn patients about device limitations.

In using total joint implants, the surgeon should be aware of the following:

- A. The correct selection of the implant is extremely important. The potential for success in total joint replacement is increased by selection of the proper size, shape, and design of the implant. Total joint prostheses require careful seating and adequate bone support. As in the case of all prosthetic implants, the durability of these components is affected by numerous biologic, biomechanic, and other extrinsic factors, which limit their service life. Accordingly, strict adherence to indications, contraindications, precautions and warnings for this product is essential to potentially maximize service level.
- B. In selecting patients for total joint replacements, the following factors can be of extreme importance to the eventual success of the procedure:
 1. A condition of senility, mental illness, chemical dependence or alcoholism. These conditions, among others, may cause the patient to ignore certain necessary limitations and precautions in the use of the implant, leading to failure or other complications.
 2. Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made prior to material selection or implantation.

Description

STRYKER Total Hip Replacement Acetabular Components are intended for cemented fixation within the prepared acetabulum. The acetabular components are used in conjunction with any appropriately sized STRYKER Stems and femoral heads of compatible head size, to help achieve total reconstructive replacement of the hip joint.

THR Acetabular Cups: The total hip replacement acetabular cups are single piece components available in all-polyethylene designs. The all-polyethylene components include the OSTEONICS UHMWPE Acetabular Cups, EXETER X3 RimFit Cups and the TRIDENT All Poly Cups. TRIDENT All Poly Cups have factory-assembled cement spacers which provide a minimum 2mm cement mantle upon cemented implantation. EXETER X3 RimFit Cups have an outer flange and 4 factory-assembled cement spacers which provide a minimum of 2mm or 3mm cement mantle upon cemented implantation. Refer to product label for cement mantle sizing. STRYKER THR Acetabular Cups are compatible (except specific mention) with all STRYKER Heads of the same diameter.

These cups are manufactured from materials certified to ASTM and/or ISO Standards.

Materials:

• ASTM F 648 (UHMWPE)	TRIDENT Crossfire All Poly Cup
• ASTM F 90 (Cobalt Chromium Alloy X-Ray Wire) • ASTM D 788 & ISO 8257-1 (PMMA Cement Spacers)	EXETER X3 RimFit All Poly Cup
• ASTM F 648 (UHMWPE) • ISO 5832-9 & ASTM F 1586 (Stainless Steel X-Ray wire) • ASTM D 788 & ISO 8257-1 (PMMA Cement Spacers)	Osteonics Polyethylene Acetabular Cup Osteonics Crossfire Acetabular Cup

Label Information

The product label provides information regarding the characteristics specific to each device (including the specific material(s) from which the product is manufactured). See product label for information regarding the specific product referenced in this package insert.

Indications

The indications for use for total hip arthroplasty include:

- Painful, disabling joint disease of the hip resulting from: degenerative arthritis, rheumatoid arthritis, post-traumatic arthritis or late stage avascular necrosis.
- Revision of previous unsuccessful femoral head replacement, cup arthroplasty or other procedure.
- Clinical management problems where arthrodesis or alternative reconstructive techniques are less likely to achieve satisfactory results.
- Where bone stock is of poor quality or inadequate for other reconstructive techniques, such as cementless fixation, as indicated by deficiencies of the acetabulum.

Stryker's Exeter X3 RimFit Acetabular Cup is intended for cemented use only. Stryker THR All-Polyethylene Acetabular Cups are intended for cemented use only.

Contraindications

- Any active or suspected latent infection in or about the hip joint.
- Any mental or neuromuscular disorder which would create an unacceptable risk of prosthesis instability, prosthesis fixation failure, or complications in postoperative care.
- Bone stock compromised by disease, infection or prior implantation which cannot provide adequate support and/or fixation to the cement mantle around the prosthesis.
- Skeletal immaturity.
- Obesity. An overweight or obese patient can produce loads on the prosthesis which can lead to failure of the fixation of the device or to failure of the device itself.

Precautions

- Before clinical use, the surgeon must thoroughly understand all aspects of the surgical procedure and limitations of the device.
- Appropriate selection, placement and fixation of the femoral stem and/or acetabular components are critical factors which affect implant service life. As in the case of all prosthetic implants, the durability of these components is affected by numerous biologic, biomechanic and other extrinsic factors, which limit their service life. Accordingly, strict adherence to the indications, contraindications, precautions and warnings for this product is essential to potentially maximize service life.
- Care must be taken to protect the components from being marred, nicked or notched as a result of contact with metal or abrasive objects.
- Surgeons should warn patients with metallic implants of the potential risks of undergoing a Magnetic Resonance Imaging (MRI) scan. The electromagnetic field created by an MRI scanner can interact with the metallic implant, resulting in displacement of the implant, heating of the tissue near the implant, implant damage or malfunction, or other undesirable effects. In addition, the presence of a metallic implant can produce an image artifact that may appear as a void region or geometric distortion of the true image. If the image artifact is near the area of interest, it may make the MRI scan uninformative or may lead to inaccurate clinical diagnosis or treatment.

MRI Safety Information

The Stryker THR Acetabular System has not been evaluated for safety and compatibility in the MR environment. The acetabular component system has not been tested for heating, migration or image artifact in the MR environment. The safety of Stryker THR Acetabular System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Utilization and Implantation

- The recommended gauge and trial components should be used for size determination, trial reduction and range of motion evaluation, thus preserving the integrity of the actual implants and their sterile packaging.
- Radiographic templates are available to assist in the preoperative prediction of component size and style.
- Howmedica Osteonics Corp.'s Surgical Protocols provide additional procedural information.
- All HOWMEDICA OSTEONICS Femoral Heads are compatible with all HOWMEDICA OSTEONICS Acetabular Liners. The flange of the EXETER X3 RimFit cup may be trimmed by the surgeon in accordance with instructions mentioned in the surgical protocol.

Warnings

- In selecting patients for total joint replacements, the following factors can be of extreme importance to the eventual success of the procedure: the patient's weight. The heavier the patient, the greater the load on the prosthesis. As the loads on the prosthesis increase, the chance a patient will suffer adverse reactions increases, including but not limited to a failure of fixation, loosening, fracture and dislocation of the device and can lead to a decreased service life. The effect of these loads will be accentuated when a small sized prosthesis is used in larger patients. Overweight or obese patients impose greater loads on the prosthesis. As obesity is a clinical diagnosis, we leave it to the surgeon to make the diagnosis based on his/her own clinical judgment. However, the World Health Organization (WHO) defines "overweight" as a BMI equal to or more than 25, and "obesity" as a BMI equal to or more than 30.
- Discard all damaged or mishandled implants.
- Never reuse an implant, even though it may appear undamaged.
- Bearing areas must always be clean and free of debris prior to assembly.
- Contouring or bending of an implant may reduce its fatigue strength and cause failure under load.
- Care should be taken not to cut through surgical gloves when handling any sharp-edged orthopaedic device.
- Howmedica Osteonics Corp. strongly advises against the use of another manufacturer's femoral head or stem component with any HOWMEDICA OSTEONICS Acetabular System Component. Any such use will negate the responsibility of Howmedica Osteonics Corp. for the performance of the resulting mixed component implant.
- Intentional removal of an acetabular component can be accomplished by careful use of cutting burrs, thin and narrow osteotomes and cautious extraction forces.
- Return all packages with flaws in the sterile barrier to the supplier. **Do not resterilize.**
- Improper selection, placement, positioning, and fixation of the implant components may result in unusual stress conditions and subsequent reduction in service life of the prosthetic implant. The surgeon must be thoroughly familiar with the surgical procedure, instruments, and implant characteristics, prior to performing surgery. Periodic, long-term follow-up is recommended to monitor the position and condition of the prosthetic components, as well as the condition of the adjoining bone.
- Patient post-operative pain. Inherent to all joint replacement is the risk that a patient will develop post-operative pain; pain is a commonly reported symptom regardless of the device implanted. The clinical literature reveals numerous potential causes of pain not directly related to the implant performance, including, but not limited to, prior history of trauma and natural disease progression. For patients who present with pain following implantation of any orthopedic implant system, the physician should consider all potential causes of the symptoms identified in the clinical literature, including infection, soft tissue impingement, and possible adverse local tissue reactions associated with wear debris, metal ions or corrosion. Accurate diagnosis of the source of pain and directed, timely intervention is essential to ensuring effective treatment of pain.

- **Ultra-High-Molecular-Weight Polyethylene (UHMWPE) Polyethylene Wear.** As would be expected, wear of the polyethylene articulating surfaces of acetabular components has been reported following total hip replacement. Higher rates of wear may be initiated by particles of cement, metal, or other debris which can cause abrasion of the articulating surfaces. Higher rates of wear may shorten the useful life of the prosthesis, and lead to relatively early revision surgery to replace the worn prosthetic components. **Care must be taken if the flange of the EXETER X3 RimFit cup is trimmed. Refer to surgical protocol for appropriate method. Never implant a cup that has been incorrectly trimmed.**
- Avoid excessive verticalization, which may accelerate bearing wear.

Adverse Effects

- While the expected life of total hip replacement components is difficult to estimate, it is finite. These components are made of foreign materials which are placed within the body for the potential restoration of mobility or reduction of pain. However, due to the many biological, mechanical and physicochemical factors which affect these devices but cannot be evaluated *in vivo*, the components cannot be expected to indefinitely withstand the activity level and loads of normal healthy bone.
- Dislocation of the hip prosthesis can occur due to inappropriate patient activity, trauma or other biomechanical considerations. Muscle and fibrous tissue laxity can also contribute to these conditions.
- Loosening of total hip components can occur. Early mechanical loosening may result from inadequate initial fixation, latent infection, premature loading of the prosthesis or trauma. Late loosening may result from trauma, infection, biological complications, including osteolysis, or mechanical problems, with the subsequent possibility of bone erosion and/or pain.
- Peripheral neuropathies, nerve damage, circulatory compromise and heterotopic bone formation may occur.
- Serious complications may be associated with any total joint replacement surgery. These complications include, but are not limited to: genitourinary disorders; gastrointestinal disorders; vascular disorders, including thrombus; bronchopulmonary disorders, including emboli; myocardial infarction or death.
- Acetabular pain may occur due to loosening of the implant.
- Wear of polyethylene components has occurred and literature reports have associated its occurrence with bone resorption, loosening and infection.
- Intraoperative fissure, fracture, or perforation of the femur, acetabulum or trochanter can occur due to impaction of the component into the prepared femoral canal or acetabulum. Postoperative femoral or acetabular fracture can occur due to trauma, the presence of defects, or poor bone stock.
- Although rare, sensitivity/allergic reactions to the materials in the implant have occurred in patients following joint replacement. Implantation of foreign material in tissues can result in immune responses and in histological reactions involving macrophages and fibroblasts.
- Adverse effects may necessitate reoperation, revision, arthrodesis of the involved joint, Girdlestone and/or amputation of the limb.

- Polyethylene particles and metal particles from mechanisms other than wear. Very small particles from metal and polyethylene components can be shed from non-articulating surfaces during normal use and over time. Although most of these particles stay in the relevant joint (i.e. contained in the synovium) or are trapped by surrounding scar tissue, microscopic particles can migrate throughout the body and on occasions have been described as accumulating in lymph nodes and other parts of the body. Although no significant medical complications have been reported as a result of these particles, their migration and/or accumulation in the body have been described in the literature. The long-term effects, if any, from these particles, are unknown. The long-term effects have been theorized to include:
 - Cancer: There is presently no scientific evidence that links metallic or polyethylene particles with cancer. However, the possibility cannot be ruled out.
 - Lymphadenopathy and Accumulation in Other Tissues/Organs: There have been a few reports of the accumulation of particles in lymph nodes (proximal and distal). Although no medical complications or disease process has been reported as stemming from these accumulations, their existence should be recognized to facilitate diagnosis and avoid confusion with suspicious lesions, cancerous or otherwise.
 - Systemic Disease: It is possible that some long-term effects may be demonstrated at some point in the future, but because there is very little scientific data suggesting association between migration of particles and systemic disease, it is believed that the benefits of these devices clearly outweigh the potential risks for any such theoretical long-term effect.

Information for patients

- The surgeon must warn patients of surgical risks, and inform them of possible adverse effects. The surgeon must warn patients that the implant does not replicate the flexibility, strength, reliability, or durability of a normal healthy joint, that the implant can break or become damaged for numerous reasons, including as a result of strenuous activity or trauma, and that the implant has a finite service life and may need to be replaced in the future.
- The surgeon must warn patients of the limitations of the reconstruction and the need to protect the implant from full weight bearing until adequate fixation and healing have occurred. The surgeon must advise the patient to limit activities and protect the implant from strenuous activity, trauma or impact loading, and to follow the surgeon's instructions regarding activity level, follow-up care, and treatment.
- The surgeon must advise patients that the implant cannot be expected to withstand the same activity levels and loads as a normal healthy joint, and that the implant will not restore function to the level expected with normal healthy bone. If the patient is involved in an occupation or activity which includes substantial walking, running, lifting, or muscle strain, the resultant forces can cause failure of the fixation, the implant, or both. The surgeon must advise the patient against having unrealistic functional expectations.
- The surgeon must warn patients that strenuous activity, trauma or impact loading affecting the implant have been implicated in failure of the implant by loosening, fracture and/or wear of the implants. Many factors, including loosening of the implant components can result in increased production of wear particles, as well as damage to the bone, making successful revision surgery more difficult.

- Transient bacteremia can occur in daily life. Dental manipulation, endoscopic examination and other surgical procedures have also been associated with transient bacteremia. To help minimize the risk of infection at the implant site, it may be advisable to use antibiotic prophylaxis before and after such procedures. Surgeons should advise the patient to inform their doctors/dentists if they have an artificial joint replacement so that a decision can be made regarding antibiotic prophylaxis for such procedures.

Sterilization

- This acetabular component has been sterilized by gamma radiation or hydrogen peroxide. Refer to product label for sterilization method.
- The packaging of all sterile products should be inspected for flaws in the sterile barrier before opening. In the presence of such a flaw, the product must be assumed nonsterile. Special trial prostheses are available to avoid having to open any aspect of the sterile package prior to component use.
- If the package is opened, but the product is not used, the component must not be resterilized and must be discarded or returned to the supplier.
- Device should not be used after the expiry date displayed on the label as packaging has not been validated beyond this date.
- Single use devices cannot be explanted and subsequently reimplanted as the physical forces exerted by these actions may compromise the physical integrity, dimensions and/or surface finishes of the devices. Also, sterility cannot be assured for reused devices as cleaning and re-sterilization procedures have not been verified.

CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a licensed physician.

CAUTION: THIS DEVICE CANNOT BE USED WITHOUT BONE CEMENT.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademark(s) or service mark(s): Howmedica, Exeter, Osteonics, Stryker, Trident, UHR, V40, X3. All other trademarks are trademarks of their respective owners or holders.

Refer to product label for CE Mark Status and Legal Manufacturer. The CE mark is only valid if also found on the product label.

The following table contains a list of abbreviations that are used on Howmedica Osteonics Corp. product labeling:

Term	Abbreviation	Term	Abbreviation
Alpha Code	ALPH CDE	Neck	NK
Angle	ANG	Offset	OFFST
Degree	DEG or °	Outer Diameter	OD
Diameter	DIA	Right	RT ►
Extra Deep	XDP	Screw Holes	SCR HLS
Extra Large	XLGE	Side	SDE
Extra Small	XSM	Size	SZE
Head	HD	Small	SM
Height	HT	Standard	STD
Inner Diameter	ID	Taper	TPR
Insert	INSR	Thickness	THKNS
Large	LGE	Type	TYP
Left	◀ LFT	With	W/
Length	LNTH	Without	W/O
Medium	MED		

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