Restoration ADM and MDM Liner

Howmedica Osteonics Corp.
325 Corporate Drive
Mahwah, NJ 07430, USA
A subsidiary of Stryker Corporation
Telephone #: +1 201-831-5000

CE 2797

©2019 Howmedica Osteonics Corp.

QIN4400 Rev. AB

Refer to product label for CE Mark Status, Notified Body Number and Legal Manufacturer. The CE mark is only valid if also found on the product label.

Labeling Symbols

Attention, see Instructions for Use

Do not Reuse

Sterilized using Irradiation

Sterilized using Hydrogen Peroxide

Sterilized using Ethylene Oxide

Use by Date

Date of Manufacture

Legal Manufacturer

Authorized Representative in the European Community

Catalog Number

Batch Code

Serial Number

MR Safe

MR Conditional

MR Unsafe

CONTENTS

| | |
|-----------------|-----|
| English | 1 |
| Français | 8 |
| Deutsch | 15 |
| Italiano | 22 |
| Español..... | 29 |
| Nederlands..... | 36 |
| Svenska | 43 |
| 日本語..... | 50 |
| Suomi | 57 |
| Dansk..... | 64 |
| Português..... | 71 |
| Ελληνικά..... | 78 |
| Norsk | 86 |
| Polski | 93 |
| 한국어 | 100 |
| 中文..... | 107 |

Restoration ADM and MDM Liner

ATTENTION OPERATING SURGEON

The advancement of partial and total hip replacement has provided the surgeon with a means of restoring mobility and reducing pain with the use of implanted prosthetic devices. While these devices have proven to be largely successful in obtaining these goals, they are manufactured from metal, plastic, or other biomaterials. Any partial or total hip replacement system, therefore, cannot be expected to withstand the same activity and loads of normal healthy bone.

The system will not be as strong, reliable or durable as a natural human hip joint and does not have an infinite lifetime. The surgeon must warn patients about device limitations.

In using total joint implants, the surgeon should be aware of the following:

- A. The correct selection of the implant is extremely important. The potential for success in total joint replacement is increased by selection of the proper size, shape, and design of the implant. Total joint prostheses require careful seating and adequate bone support.
- B. In selecting patients for total joint replacements, the following factors can be of extreme importance to the eventual success of the procedure:
 1. The patient's occupation or activity. If the patient is involved in an occupation or activity which includes significant impact loads (walking, running, lifting, or twisting), the resultant forces can cause failure of the fixation, the device, or both. High levels of physical activity over the years can also accentuate the normal wear process that occurs with prosthetic joints. The prosthesis will not restore function to the level expected with normal healthy bone, and the patient should not have, and should be disabused of, unrealistic functional expectations. (See PRECAUTIONS section for more information).
 2. A condition of senility, mental illness, chemical dependence or alcoholism. These conditions, among others, may cause the patient to ignore certain necessary limitations and precautions in the use of the implant, leading to failure or other complications.
 3. Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made prior to material selection or implantation.

DESCRIPTION

Stryker Orthopaedics manufactures a wide variety of reconstructive total hip systems to help satisfy anatomical requirements and surgeon/patient needs. Please refer to the applicable, Surgical Protocol for further product specific detail.

ADM Acetabular Cups

Acetabular cups are available in a range of outer diameter sizes. The Restoration ADM system is a two-piece system consisting of a metal HA coated shell and a polyethylene insert.

The outer acetabular cup is manufactured from cobalt chrome alloy and coated with CP Ti, and hydroxylapatite. The cobalt chrome cups are manufactured from materials certified to ASTM Standard ASTM F1537.

MDM Liner

The MDM (Modular Dual Mobility) Liner is a highly polished modular liner with a Trident locking mechanism. This will allow for compatibility with acetabular cups containing the same locking mechanism. Refer to the surgical protocol for outer acetabular shell/ dual mobility insert compatibility. The Restoration ADM Inserts articulate within the MDM Liner. The MDM Liner is manufactured from cobalt chrome alloy certified to ASTM Standard ASTM F1537.

Polyethylene Insert

The inner acetabular polyethylene inserts are available with an inner diameter that accommodates the mating femoral head.

The inserts used with the cups are designed to move freely in the cup. The Ultra High Molecular Weight Polyethylene (UHMWPE) inserts are manufactured from materials certified to ASTM Standard F648 and/or ISO 5834-2.

HOWMEDICA OSTEONICS Inserts and Acetabular Cups made of polyethylene (UHMWPE) are compatible (except specific mention) with all HOWMEDICA OSTEONICS Heads of the same diameter.

Coatings

The coatings incorporated on the Restoration ADM HA metal shells consist of titanium plasma spray, and the titanium bond layer overlaid with a layer of hydroxylapatite. These materials are certified to ASTM and/or ISO Standards:

| Material: | ASTM Standard / ISO Standard |
|---------------------------------|-------------------------------------|
| Titanium Commercially Pure (CP) | F1580 / 5832-2 |
| Hydroxylapatite (powder) | F1185 |

Label Information

The product label provides information regarding specific material(s) from which the product is manufactured.

INDICATIONS

The indications for use for total hip arthroplasty include:

1. Noninflammatory degenerative joint disease including osteoarthritis and avascular necrosis;
2. Rheumatoid arthritis;
3. Correction of functional deformity;
4. Revision procedures where other treatments or devices have failed; and,
5. Treatment of nonunion, femoral neck and trochanteric fractures of the proximal femur with head involvement that are unmanageable using other techniques.
6. Dislocation risks

RESTORATION ADM HA Cups and MDM Liners are intended for cementless use only.

CONTRAINdications

- 1) Overt infection;
- 2) Distant foci of infections (which may cause hematogenous spread to the implant site);
- 3) Rapid disease progression as manifested by joint destruction or bone resorption apparent on roentgenogram;
- 4) Skeletally immature patients; and
- 5) Cases where there is a loss of abductor musculature, poor bone stock, or poor skin coverage around the hip joint which would make the procedure unjustifiable.

Conditions presenting increased risk of failure include:

- 1) Uncooperative patient or patient with neurologic disorders, incapable of following instructions;
- 2) Osteoporosis;
- 3) Metabolic disorders which may impair bone formation; and
- 4) Osteomalacia.

WARNINGS

In using this system, the surgeon should be aware of the following:

In selecting patients for total joint replacements, the following factor is of extreme importance to the eventual success of the procedure: The patient's weight. The heavier the patient, the greater the load on the prosthesis. As the loads on the prosthesis increase, the chance a patient will suffer adverse reactions increases, including but not limited to failure of fixation, loosening, fracture and dislocation of the device and can lead to a decreased service life. The effect of these loads will be accentuated when a small sized prosthesis is used in larger patients. Overweight or obese patients impose greater loads on the prosthesis. As obesity is a clinical diagnosis, we leave it to the surgeon to make the diagnosis based on his/her own clinical judgment. However, the World Health Organization (WHO) defines "overweight" as a BMI equal to or more than 25, and "obesity" as a BMI equal to or more than 30.

Improper selection, placement, positioning, and fixation of the implant components may result in unusual stress conditions and subsequent reduction in service life of the prosthetic implant. The surgeon must be thoroughly familiar with the surgical procedure, instruments, and implant characteristics, prior to performing surgery. Periodic, long-term follow-up is recommended to monitor the position and condition of the prosthetic components, as well as the condition of the adjoining bone.

Do not substitute another manufacturer's device for any component of the Howmedica Osteonics Total Hip System. Any such use will negate the responsibility of Howmedica Osteonics Corp. for the performance of the resulting mixed component implant.

Restoration ADM HA Cup is intended for Press-Fit Application only. Secure fixation at the time of surgery is critical to the success of the procedure. Bone stock must be adequate to support the device. Press-fit acetabular components require a precise operative technique and use of specified instruments. Adequate acetabular preparation is important for the success of the procedure.

Hydroxylapatite coated implants are not intended for use with bone cement.

Modular Acetabular Shell/Liner. Fixation devices, when used, should be fully seated to help assure stable fixation of the shell, and to help avoid interference with the liner component. Prior to seating the liner component into the shell component, surgical debris must be cleaned from the interior of the shell. Debris may inhibit the liner from functioning correctly and may increase polyethylene wear.

Patient post-operative pain. Inherent to all joint replacement is the risk that a patient will develop post-operative pain; pain is a commonly reported symptom regardless of the device implanted. The clinical literature reveals numerous potential causes of pain not directly related to the implant performance including, but not limited to, prior history of trauma and natural disease progression.

For patients who present with pain following implantation of any orthopedic implant system, the physician should consider all potential causes of the symptoms identified in the clinical literature, including infection, soft tissue impingement, and possible adverse local tissue reactions associated with wear debris, metal ions or corrosion. Accurate diagnosis of the source of pain and directed, timely intervention is essential to ensuring effective treatment of pain.

- **Modular Junctions:** Mate modular components firmly to prevent dissociation. Machined taper surfaces must be clean, dry and firmly mated to ensure proper seating and assembly. Repeated assembly/disassembly or failure to clean, dry and firmly mate the components could compromise the taper lock and contribute to wear/corrosion and significant clinical consequences to the patient. (See Adverse Effects below).

Information for patients. The surgeon must warn patients of surgical risks, and inform them of possible adverse effects. The surgeon must warn patients that the implant does not replicate the flexibility, strength, reliability, or durability of a normal healthy joint, that the implant can break or become damaged for numerous reasons, including as a result of strenuous activity or trauma, and that the implant has a finite service life and may need to be replaced in the future.

The surgeon must warn patients of the limitations of the reconstruction and the need to protect the implant from full weight bearing until adequate fixation and healing have occurred. The surgeon must advise the patient to limit activities and protect the implant from strenuous activity, trauma or impact loading, and to follow the surgeon's instructions regarding activity level, follow-up care, and treatment.

The surgeon must advise patients that the implant cannot be expected to withstand the same activity levels and loads as a normal healthy joint, and that the implant will not restore function to the level expected with normal healthy bone. If the patient is involved in an occupation or activity which includes substantial walking, running, lifting, or muscle strain, the resultant forces can cause failure of the fixation, the implant, or both. The surgeon must advise the patient against having unrealistic functional expectations.

The surgeon must warn patients that strenuous activity, trauma or impact loading affecting the implant have been implicated in failure of the implant by loosening, fracture and/or wear of the implants. Many factors, including loosening of the implant components can result in increased production of wear particles, as well as damage to the bone, making successful revision surgery more difficult.

Transient bacteremia can occur in daily life. Dental manipulation, endoscopic examination and other surgical procedures have also been associated with transient bacteremia. To help minimize the risk of infection at the implant site, it may be advisable to use antibiotic prophylaxis before and after such procedures. Surgeons should advise the patient to inform their doctors/dentists if they have an artificial joint replacement so that a decision can be made regarding antibiotic prophylaxis for such procedures.

PRECAUTIONS

Instruments. Specialized instruments are available and must be used to help assure the accurate implantation of prosthetic components.

While rare, intraoperative fracture or breaking of instruments can occur. Instruments which have experienced extensive use or excessive force are more susceptible to fracture. Instruments must be examined for wear, or damage, prior to surgery.

Re-use. An implant must never be reused. While it may appear undamaged, a used implant may have acquired blemishes or latent compromise of its integrity which would reduce its service life.

Handling. Proper handling of implants is important. The highly polished portion of the implant should not come in contact with hard surfaces.

Surgeons should warn patients with metallic implants of the potential risks of undergoing a Magnetic Resonance Imaging (MRI) scan. The electromagnetic field created by an MRI scanner can interact with the metallic implant, resulting in displacement of the implant, heating of the tissue near the implant, implant damage or malfunction, or other undesirable effects. In addition, the presence of a metallic implant can produce an image artifact that may appear as a void region or geometric distortion of the true image. If the image artifact is near the area of interest, it may make the MRI scan uninformative or may lead to inaccurate clinical diagnosis or treatment.

MRI Safety Information

Restoration ADM Acetabular Shell, MDM Liners & MDM/ADM Inserts have been evaluated for safety and compatibility in the MR environment.

Non-clinical testing has demonstrated that the devices listed above are MR Conditional. A patient with these listed devices can be safely scanned in an MR scanner meeting the following conditions:

- Static magnetic field of 1.5-Tesla and 3.0-Tesla, only
- Maximum spatial gradient magnetic field of 2,310 gauss/cm (23 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system.
- Evaluation was performed using a quadrature body coil only

Under the scan conditions defined above, these devices are expected to produce a temperature rise of less than 4.1 °C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 84mm from the device when imaged with a gradient echo pulse sequence using a 3.0 T/128 MHz MRI system.

This MRI information is also available at <https://www.stryker.com/us/en/joint-replacement/MRI.html>.

ADVERSE EFFECTS

- 1) Although rare, sensitivity/allergic reactions to the materials in the implant have occurred in patients following joint replacement. Implantation of foreign material in tissues can result in immune responses and in histological reactions involving macrophages and fibroblasts.
- 2) Peripheral neuropathies have been reported following total joint surgery. Subclinical nerve damage has been reported, and may occur as the result of surgical trauma.
- 3) Dislocation of the hip prosthesis can occur due to inappropriate patient activity, trauma or other biomechanical considerations. Muscle and fibrous tissue laxity can also contribute to these conditions.
- 4) Implants can loosen or migrate due to trauma or loss of fixation.
- 5) Infection can lead to failure of the joint replacement.
- 6) Fatigue fracture of the implant can occur as a result of trauma, strenuous activity, improper alignment, and/or duration of service, singularly or in combination.
- 7) Polyethylene particles and metal particles from mechanisms other than wear. Very small particles from metal and polyethylene components can be shed from non-articulating surfaces during normal use and over time. Although most of these particles stay in the relevant joint (i.e. contained in the synovium) or are trapped by surrounding scar tissue, microscopic particles can migrate throughout the body and on occasions have been described as accumulating in lymph nodes and other parts of the body. Although no significant medical complications have been reported as a result of these particles, their migration and/or accumulation in the body have been described in the literature. The long-term effects, if any, from these particles, are unknown. The long-term effects have been theorized to include:
 - Cancer: There is presently no scientific evidence that links metallic or polyethylene particles with cancer. However, the possibility cannot be ruled out.
 - Lymphadenopathy and Accumulation in Other Tissues/Organs: There have been a few reports of the accumulation of particles in lymph nodes (proximal and distal). Although no medical complications or disease process has been reported as stemming from these accumulations, their existence should be recognized to facilitate diagnosis and avoid confusion with suspicious lesions, cancerous or otherwise.
 - Systemic Disease: It is possible that some long-term effects may be demonstrated at some point in the future, but because there is very little scientific data suggesting association between migration of particles and systemic disease, it is believed that the benefits of these devices clearly outweigh the potential risks for any such theoretical long-term effect.
- 8) Non-metallic wear debris. Wear debris is generated by interaction between components, as well as between components and bone, primarily through wear mechanisms of adhesion, abrasion and fatigue. Secondarily, particulate can also be generated by third-body wear. With all implant devices, asymptomatic, localized progressive bone resorption (osteolysis) may occur around the prosthetic components as a consequence of foreign-body reaction to the particulate matter of cement, ultra-high molecular weight polyethylene (UHMWPE) and/or ceramic. Osteolysis can lead to future complications, including loosening, necessitating the removal and replacement of prosthetic components.
- 9) Metal wear debris. Metal wear debris, metal ions and corrosion of metal implants. Generation of metal wear debris, metal ions and/or corrosion occurs whenever two surfaces are in contact, and at least one surface is metal. There have been reports in the literature of cases of adverse local tissue reactions associated with wear and/or corrosion at modular junctions formed by modular heads (stem/ head interface). Local joint chemistry and/or other patient-specific conditions including, but not limited to infection may affect the potential for in vivo corrosion and its possible clinical consequences.
- 10) Metal debris and corrosion by-products. Modular junctions may release metal debris and/or metal ions due to fretting, galvanic corrosion, crevice corrosion, or other processes. There are several factors involved in these processes, including forces across a junction, which are not fully understood. These corrosion products or metal debris may affect the tissues surrounding the implant and adversely affect the duration of service life. There have been reports of patients developing adverse local tissue

reactions (including, but not limited to, tissue necrosis, pseudotumors, cysts and fluid accumulations, metallosis and aseptic lymphocyte dominated vasculitis associated lesions), elevated metal ion levels in the blood and/or urine and hypersensitivity/allergic reactions associated with corrosion and/or wear-related debris in the implant vicinity. Affected patients may present with symptoms similar to those associated with infection, including pain (most likely during weight-bearing) and swelling at the local joint site. These reactions should be carefully monitored and may result in early revision surgery. Medical literature describes systemic reactions to by-products arising from contemporary metal on metal articulating bearing surfaces. While there are no metal on metal bearing interfaces in this system, it could be theorized that similar systemic reactions may arise from fretting and corrosion arising from any metal interface.

Intraoperative and early postoperative complications can include:

- 1) Femoral or acetabular perforation, or fracture;
- 2) Femoral fracture while seating the device;
- 3) Damage to blood vessels;
- 4) Temporary or permanent nerve damage resulting in pain or numbness of the affected limb;
- 5) Undesirable shortening or lengthening of the limb;
- 6) Traumatic arthrosis (joint disease caused by trauma) of the knee from intraoperative positioning of the extremity;
- 7) Cardiovascular disorders including venous thrombosis, pulmonary embolism, or myocardial infarction, or death;
- 8) Hematoma;
- 9) Delayed wound healing; and
- 10) Infection.

Late postoperative complications can include:

- 1) Trochanteric avulsion (tearing away of soft tissue from the bone) as a result of excess muscular tension, early weight bearing, or inadvertent intraoperative weakening;
- 2) Trochanteric nonunion due to inadequate reattachment and/or early weight bearing;
- 3) Aggravated problems of the affected limb or contralateral extremity caused by leg length discrepancy, excess femoral medialization, or muscle deficiency;
- 4) Femoral fracture by trauma or excessive loading, particularly in the presence of poor bone stock;
- 5) Periarticular calcification or ossification, with or without impediment to joint mobility; and
- 6) The following, singularly or in combination can lead to and/or cause decreased range of motion: improper selection or positioning of components, femoral impingement, and/or periarticular calcification.

How Supplied

- These products have been sterilized by either gamma radiation, hydrogen peroxide or ethylene oxide. Refer to the package label for the sterilization method.
- Do NOT resterilize.
- Take care to prevent contamination of any components.
- Inspect the packaging of sterile products for flaws before opening. In the presence of any flaws, assume that the product is nonsterile.
- Discard ALL nonsterile or contaminated products per applicable local regulations and hospital protocol.
- Device should not be used after the expiry date displayed on the label as packaging has not been validated beyond this date.

- Single use devices cannot be explanted and subsequently reimplanted as the physical forces exerted by these actions may compromise the physical integrity, dimensions and/or surface finishes of the devices. Also, sterility cannot be assured for reused devices as cleaning and re-sterilization procedures have not been verified.

Disposal

If a device is being returned for evaluation, please contact your local Stryker representative for shipping/handling information. If the device is not being returned to Stryker, implant components are to be disposed of in accordance with applicable laws, rules, and regulations for the disposal of biohazardous waste. Follow all guidelines for biohazardous waste in accordance with the Centers for Disease Control and Prevention guidelines as well as applicable federal/national, state and local regulations. As part of the disposal process, verify that the implant in its entirety has been explanted from the surgical site.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Caution: Federal Law in the USA restricts this device to sale by or on the order of a physician.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademark(s) or service mark(s): ADM, Howmedica, MDM, Osteonics, Restoration, Stryker, X3. All other trademarks are trademarks of their respective owners or holders.

Refer to product label for CE Mark Status, Notified Body Number and Legal Manufacturer. The CE mark is only valid if also found on the product label.

The following table contains a list of abbreviations that are used on Howmedica Osteonics Corp. product labeling:

| Term | Abbreviation | Term | Abbreviation |
|----------------|--------------|----------------|--------------|
| Alpha Code | ALPH CDE | Neck | NK |
| Angle | ANG | Offset | OFFST |
| Degree | DEG or ° | Outer Diameter | OD |
| Diameter | DIA | Right | RT ► |
| Extra Deep | XDP | Screw Holes | SCR HLS |
| Extra Large | XLGE | Side | SDE |
| Extra Small | XSM | Size | SZE |
| Head | HD | Small | SM |
| Height | HT | Standard | STD |
| Inner Diameter | ID | Taper | TPR |
| Insert | INSR | Thickness | THKNS |
| Large | LGE | Type | TYP |
| Left | ◀ LFT | With | W/ |
| Length | LNTH | Without | W/O |
| Medium | MED | | |

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