**stryker®**  
**Howmedica**  
**OSTEONICS**

## Constrained Acetabular Insert

Howmedica Osteonics Corp.  
325 Corporate Drive  
Mahwah, NJ 07430, USA  
A subsidiary of Stryker Corporation

Telephone #: +1 201-831-5000

**CE 0086**

©2016 Howmedica Osteonics Corp.  
QIN4357 Rev. AA

Refer to product label for CE mark status and Legal Manufacturer.  
The CE mark is only valid if also found on the product label.

## Labeling Symbols

Attention, See Instructions for Use

Do Not Reuse

Sterilized using Irradiation

Sterilized using Hydrogen Peroxide

Sterilized using Ethylene Oxide

Use by Date

Date of Manufacture

Legal Manufacturer

Authorized Representative in the European Community

Catalog Number

Batch Code

Serial Number

Constrained Acetabular Insert .....	1
Insert acétabulaire constraint .....	6
Constrained Acetabulumeinsatz .....	11
Inserto acetabolare vincolato .....	17
Inserto acetabular ligado .....	22
Inserção acetabular restrinvida .....	27
Trident begränsad acetabular insats .....	33
Rajoittava lonkkamajasisäke .....	38
Tvungen acetabulumindsats .....	43
Retentief acetabulair inlegdeel .....	48
Tvunget acetabuluminnlegg .....	54
Wkładka panewkowa Trident z mechanizmem blokującym .....	59
Περιορισμένη Κοτύλαια Πρόθεση .....	65
限制式髋臼内衬 .....	71
강제 관골구 삽입 .....	76
抑制寛骨臼ペアリングインサート .....	82

## English

### **CONSTRINED ACETABULAR INSERT**

#### **Description**

The HOWMEDICA OSTEONICS Constrained Acetabular Inserts are available in bipolar and unipolar designs. The TRIDENT Constrained Acetabular Insert is comprised of two pre-assembled components: an outer insert component and a captured UHR (Universal Head) component. The UHR component is comprised of an outer shell into which a bearing insert has been permanently assembled. The UHR bearing insert has a factory assembled UHMWPE retention ring. The outer acetabular insert has a Ti alloy retaining ring which retains the UHR head in the plastic portion of the insert. With the exception of the TRIDENT All-Poly Constrained Inserts, the TRIDENT Constrained Acetabular Inserts are designed to be assembled with TRIDENT metal acetabular shells. The TRIDENT All-Poly Constrained Inserts can be cemented directly into a GAP Cup, GAP Ring, TRIDENT metal acetabular shell, or directly into the acetabulum. The unipolar design consists of an UHMWPE insert, with a Ti alloy or CoCr locking ring/wire and a UHMWPE retention ring depending on the system. The assembled acetabular component is used in conjunction with any appropriately sized HOWMEDICA OSTEONICS stem of compatible head size, to achieve total reconstructive replacement of the hip joint.

#### **Materials:**

The devices are manufactured from materials that meet the following ASTM Standards:

- |  |  |
|--|--|
| • ASTM F-75 Cobalt Chromium Alloy                              | UHR Outer Shell  |
| • ASTM F-90 Cobalt Chromium Alloy                              | Locking Wire   |
| • ASTM F-648 Ultra-High Molecular Weight Polyethylene (UHMWPE) | UHR bearing insert, Acetabular bearing insert body, UHR retention ring |
| • ASTM F-136 Titanium 6Al-4V ELI Alloy                         | Retaining ring   |

#### **Compatibility**

- HOWMEDICA OSTEONICS Constrained Acetabular Inserts are compatible with all HOWMEDICA OSTEONICS Femoral Heads except V40 head catalog number series 6264-4-XXX and 6264-5-XXX and P.C.A. Head catalog number series 6284-0-XXX.

#### **Indications**

- A Constrained Acetabular Insert is indicated for use as a component of a total hip prosthesis in primary and revision patients at high risk of hip dislocation due to a history of prior dislocation, bone loss, joint soft tissue laxity, neuromuscular disease, or intraoperative instability.
- PCA, Series II, and System 12 Constrained Acetabular Inserts are indicated for use in revision patients only.

#### **Contraindications**

- Bone or musculature compromised by disease, infection or prior implantation, which cannot provide adequate support or fixation for the prosthesis.
- Infection in or about the hip joint.
- Skeletal immaturity.

#### **Warnings**

- Closed reduction of a dislocation of this device is not possible. Patients should be made aware that treatment of device dislocation will require additional surgery.
- Patients should be instructed on the impact of excessive loading that can result if the patient is involved in an occupation or activity which includes substantial walking, running, lifting, or excessive muscle loading due to patient weight causing extreme demands on the constrained insert can result in the failure of the device. Extreme demands on the device may also compromise the acetabular shell's fixation in the acetabulum.

- The UHMWPE UHR retention ring and Ti alloy retaining ring of the constrained insert should not be handled or removed as they are critical to the security of the assembly. Alteration of the factory preassembled device can result in improper function of the retaining mechanisms. Discard or return to manufacturer any constrained insert if the retaining mechanisms appear damaged or mishandled.
- Improper alignment of the acetabular insert within the acetabular shell prior to impaction may result in improper seating of the constrained acetabular insert.
- Removal of the constrained insert after its assembly into the metal shell results in the destruction of the insert. Discard any device removed after the locking mechanism has been engaged, do not reinsert the device.
- Care should be taken not to nick or notch the inner surface of the metal shell during insert removal, which could lead to premature wear of the UHMWPE.
- V40 femoral head catalog number series 6264-4-XXX and 6264-5-XXX and P.C.A. Femoral Head catalog number series 6284-0-XXX are not compatible with the UHR Bipolar and Constrained Acetabular Insert.
- Do not substitute another manufacturer's device for any component of the HOWMEDICA OSTEONICS Total Hip System. Any such use will negate the responsibility of Howmedica Osteonics Corp. for the performance of the resulting mixed component implant.
- Never reuse an implant, even though it may appear undamaged.
- Discard damaged, mishandled, or contaminated implants.
- See the Patient Counseling Information Section for more information.

#### **Precautions**

- Protect all components from contamination.
- Do not allow coated surfaces to contact cloth or other fiber-releasing materials.
- Protect polished bearing areas and machined taper surfaces from contact with hard or abrasive surfaces.
- Take care not to cut through surgical gloves when handling any sharp-edged orthopaedic device.
- In order to minimize the risks of dislocation and loosening of the shell-acetabular bone or shell-bone cement interface that may occur when using a metallic shell intended for biological fixation or cemented use only, surgeons should consider providing immediate resistance to tensile forces between the metallic shell and the acetabular bone or bone cement interface through the use of orthopedic bone fixation devices such as bone screws, spikes, screw threads, fins, or other bone fixation device.
- To correctly position the metallic locking ring, surgeons should consult the manufacturer's instructions for appropriate device assembly.
- Physicians should consider component malposition, component placement, and the effect on range of motion when using modular heads (with sleeves or skirts) and extended liners.

#### **Patient Selection**

- Proper implant selection is critical to the stability and longevity of the acetabular implant in hip arthroplasty. Proper implant selection must consider design, fixation, and environmental variables including: patient weight, age, bone quality and size, activity level and pre-operative level of health, as well as the surgeon's experience and familiarity with the device. Longevity and stability of the implant may be affected by these factors.
- Patients with high-activity level and/or higher weight patients are at greater risk for implant complications or failures. For patients with poor proximal bone quality, the use of supplemental adjunctive fixation/support is advised for implant stability.
- The surgeon must evaluate each situation carefully based upon the patient's clinical presentation before making any decisions regarding the selection of the implant.

#### **Adverse Effects**

- Dislocation of the hip prosthesis can occur due to inappropriate patient activity, trauma or other biomechanical considerations.

- Loosening of total hip components can occur. Early mechanical loosening may result from inadequate initial fixation, latent infection, premature loading of the prosthesis or trauma. Late loosening may result from trauma, infection, biological complications, including osteolysis, or mechanical problems, with the subsequent possibility of bone erosion and/or pain.
- Fatigue fracture of femoral stems has occurred in a small percentage of cases. Stem fracture is more likely to occur in heavy, physically active patients, or when contralateral joint disability results in a disproportionate distribution of weight on the reconstructed joint.
- Wear of polyethylene components may occur. Polyethylene wear has been associated with bone resorption, loosening, and infection.
- Peripheral neuropathies, nerve damage, circulatory compromise, and heterotopic bone formation may occur.
- Serious complications may be associated with any total joint replacement surgery. These complications include, but are not limited to: genitourinary disorders, gastrointestinal disorders, vascular disorders including thrombus, bronchopulmonary disorders including emboli, myocardial infarction, or death.
- Acetabular pain may occur due to loosening of the implant.
- Intraoperative fissure, fracture, or perforation of the femur, acetabulum, or trochanter can occur due to impaction of the component into the prepared femoral canal or acetabulum. Postoperative femoral or acetabular fracture can occur due to trauma, the presence of defects, or poor bone stock.
- Asymptomatic, localized progressive bone resorption (osteolysis) may occur around the prosthetic components as a consequence of foreign-body reaction to the particulate matter of cement, metal, ultra-high molecular weight polyethylene (UHMWPE) and/or ceramic. Particulate is generated by interaction between components, as well as between components and bone, primarily through wear mechanisms of adhesion, abrasion and fatigue. Secondarily, particulate can also be generated by third-body wear. Osteolysis can lead to future complications, including loosening, necessitating the removal and replacement of prosthetic components.
- Very small particles from metal and polyethylene components can be shed from the components during normal use and over time. Although most of this debris stays in the relevant joint (i.e. contained in the synovium) or is trapped by surrounding scar tissue, microscopic particles can migrate throughout the body and on occasions have been described as accumulating in lymph nodes and other parts of the body. Although no significant medical complications have been reported as a result of these particles, their migration and/or accumulation in the body have been described in the literature. Given the insufficient time period during which patients with these devices have been followed and the fact that these devices are currently being used in younger patients and remain in the body for increasingly longer periods of time, it should be said that the long-term effects, if any, from these particles, are unknown. The long-term effects have been theorized to include:
  - Cancer: There is presently no scientific evidence that links metallic or polyethylene debris with cancer. However, the possibility cannot be ruled out.
  - Lymphadenopathy and Accumulation in Other Tissues/Organs: There have been a few reports of the accumulation of wear debris in lymph nodes (proximate and distal). Although no medical complications or disease process has been reported as stemming from these accumulations, their existence should be recognized to facilitate diagnosis and avoid confusion with suspicious lesions, cancerous or otherwise.
  - Systemic Disease: There has been some speculation that there could be an association between migration of debris and as yet unidentified systemic effects. Long-term effects may be demonstrated at some point in the future, but because there is very little scientific data suggesting association between migration of debris and systemic disease, it is believed that the benefits of these devices clearly outweigh the potential risks for any such theoretical long-term effect.
- Metal sensitivity reactions have been reported following joint replacement.
- Undesirable shortening or lengthening of the limb.
- Infection can occur.
- Adverse effects may require medical intervention including reoperation, revision, arthrodesis of the involved joint, Girdlestone, or amputation of the limb.

## Use and Implantation

- Before clinical use, the surgeon should thoroughly understand all aspects of the surgical procedure and limitations of the device.
- The Surgical Protocol for the appropriate HOWMEDICA OSTEONICS hip system provides additional procedural information.
- Radiographic templates are available to assist in the preoperative prediction of component size and style.
- Specialized instruments are available and should be used to ensure accurate implantation of prosthetic components.
- To preserve the integrity of the actual implants and their sterile packaging, use the recommended trial components for size determination, trial reduction and range-of-motion evaluation.
- Proper selection, placement, and fixation of the implant components are critical factors affecting implant service life. The durability of prosthetic implants is affected by many biologic, biomechanical and other extrinsic factors that limit their service life. Accordingly, strict adherence to the indications, contraindications, precautions and warnings for this product is essential to potentially maximize service life.

## Information for the Patient

The surgeon should counsel the patient on the following information:

- Advise patients of the limitations of the reconstruction and the need to protect the implant from full weight bearing until adequate healing has occurred.
- Caution patients to protect the replaced joint from excessive loading, and to follow the physician's instructions regarding activity level, follow-up care, and treatment. Advise patients that the device cannot be expected to withstand indefinitely the same activity levels and loads as a normal healthy joint, that the implant can break or become damaged as a result of excessive loading, and that the device has a finite service life and may need to be replaced in the future.
- Warn the patient of surgical risks and possible adverse effects.
- Dental procedures, endoscopic examinations and other minor surgical procedures have been associated with transient bacteremia. Instruct the patient to inform their doctors that they have an artificial hip replacement, so that their doctors can decide whether to use antibiotic prophylaxis for such procedures.

## How Supplied

- These products have been sterilized by gamma radiation.
- Do NOT resterilize.
- Take care to prevent contamination of any components.
- Inspect the packaging of sterile products for flaws before opening. In the presence of any flaws, assume that the product is nonsterile.
- Discard ALL nonsterile or contaminated products.
- Single use devices cannot be explanted and subsequently reimplanted as the physical forces exerted by these actions may compromise the physical integrity, dimensions and/or surface finishes of the devices. Also, sterility cannot be assured for reused devices as cleaning and re-sterilization procedures have not been verified.

**Caution:      Federal Law in the USA restricts this device to sale by or on the order of a physician.**

**Stryker Corporation or its divisions or other corporate affiliated entities, own, use or have applied for the following trademarks: Howmedica, Osteonics, PCA, Stryker, System 12, Trident, UHR, V40. All other trademarks are trademarks of their respective owners or holders.**

Refer to product label for CE Mark Status and Legal Manufacturer. The CE mark is only valid if also found on the product label.

The following table contains a list of abbreviations that are used on Howmedica Osteonics Corp. product labeling:

Term	Abbreviation	Term	Abbreviation
Alpha Code	ALPH CDE	Neck	NK
Angle	ANG	Offset	OFFST
Degree	DEG or °	Outer Diameter	OD
Diameter	DIA	Right	RT ►
Extra Deep	XDP	Screw Holes	SCR HLS
Extra Large	XLGE	Side	SDE
Extra Small	XSM	Size	SZE
Head	HD	Small	SM
Height	HT	Standard	STD
Inner Diameter	ID	Taper	TPR
Insert	INSR	Thickness	THKNS
Large	LGE	Type	TYP
Left	◀ LFT	With	W/
Length	LNTH	Without	W/O
Medium	MED		

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