**stryker**  
**Howmedica**  
**OSTEONICS**

**HOWMEDICA OSTEONICS TRIDENT Acetabular Component System and  
HOWMEDICA OSTEONICS Alumina and BIOLOX *delta* Ceramic Heads**

Howmedica Osteonics Corp.  
325 Corporate Drive  
Mahwah, NJ USA 07430, USA  
A subsidiary of Stryker Corporation

Telephone #: +1 201-831-5000

**CE 2797**

©2020 Howmedica Osteonics Corp.

QIN 4350 Rev. AB

Refer to product label for CE mark status and Legal Manufacturer.  
The CE mark is only valid if also found on the product label.

## Labeling Symbols

Attention, See Instructions for Use

Do Not Reuse

Sterilized using Irradiation

Sterilized using Hydrogen Peroxide

Sterilized using Ethylene Oxide

Use by Date

Date of Manufacture

Legal Manufacturer

Authorized Representative in the European Community

Catalog Number

Batch Code

Serial Number

MR Safe

MR Conditional

MR Unsafe

HOWMEDICA OSTEONICS TRIDENT Acetabular Component System and HOWMEDICA OSTEONICS Alumina and BIOLOX <i>delta</i> Ceramic Heads .....	1
Système de composants acétabulaires TRIDENT de HOWMEDICA OSTEONICS et têtes en céramique d'alumine et BIOLOX <i>delta</i> de HOWMEDICA OSTEONICS .....	10
Trident Acetabulumkomponentensystem von Howmedica Osteonics und Aluminiumoxid- und BIOLOX <i>delta</i> -Keramikköpfe von Howmedica Osteonics .....	20
Sistema di componenti acetabolari TRIDENT di HOWMEDICA OSTEONICS e teste in allumina HOWMEDICA OSTEONICS e in ceramica BIOLOX <i>delta</i> .....	30
Sistema de componente acetabular TRIDENT HOWMEDICA OSTEONICS y cabezas de alumina HOWMEDICA OSTEONICS y de cerámica BIOLOX <i>delta</i> .....	40
Sistema de Componentes Acetabulares TRIDENT HOWMEDICA OSTEONICS e Cabeças de Alumínio HOWMEDICA OSTEONICS e Cabeças de Cerâmica BIOLOX <i>delta</i> .....	50
HOWMEDICA OSTEONICS TRIDENT acetabulart komponentsystem och HOWMEDICA OSTEONICS aluminiumoxid- och BIOLOX <i>delta</i> keramiska huvuden .....	60
HOWMEDICA OSTEONICS TRIDENT -lonkkamaljakomponenttijärjestelmä ja HOWMEDICA OSTEONICS -alumiinioksiidi- ja BIOLOX <i>delta</i> keraaminupit.....	69
HOWMEDICA OSTEONICS TRIDENT acetabulum komponentsystem og HOWMEDICA OSTEONICS alumina og BIOLOX <i>delta</i> keramikhoveder .....	78
HOWMEDICA OSTEONICS TRIDENT-acetabulumkomponentensysteem en HOWMEDICA OSTEONICS en BIOLOX <i>delta</i> -keramiekkopjes .....	87
HOWMEDICA OSTEONICS TRIDENT hoftleddskålskomponentsystem og HOWMEDICA OSTEONICS aluminium-oksidi og BIOLOX <i>delta</i> keramiske knokkelhoder .....	97
Komponenty panewkowe TRIDENT firmy HOWMEDICA OSTEONICS oraz glowy firmy HOWMEDICA OSTEONICS z ceramiki korundowej i ceramiki BIOLOX <i>delta</i> .....	106
Σύστημα Εξαρτήματος Κοτύλης HOWMEDICA OSTEONICS TRIDENT και HOWMEDICA OSTEONICS Alumina και Κεραμικές Κεφαλές BIOLOX <i>delta</i> .....	116
HOWMEDICA OSTEONICS TRIDENT 髌臼组件系统和 HOWMEDICA OSTEONICS 氧化铝和 BIOLOX <i>delta</i> 陶瓷股骨头 .....	126
HOWMEDICA OSTEONICS TRIDENT 관골구 콤포넌트 시스템과 Howmedica Osteonics 알루미늄 두부와 BIOLOX <i>delta</i> 세라믹 두부 .....	135
HOWMEDICA OSTEONICS TRIDENT 寛骨臼コンポーネントシステムおよび HOWMEDICA OSTEONICSアルミナ製およびBIOLOX <i>delta</i> セラミック製股骨头.....	146

## English

### HOWMEDICA OSTEONICS TRIDENT acetabular component system and HOWMEDICA OSTEONICS alumina and BIOLOX *delta* ceramic heads

#### Description

The HOWMEDICA OSTEONICS TRIDENT Acetabular Component System consists of a metal acetabular shell and TRIDENT acetabular bearing insert. The shells are available with a variety of surface enhancements including but not limited to Arc Deposition with or without Hydroxylapatite surface treatment. The shells are intended for cementless fixation within the prepared acetabulum. The HOWMEDICA OSTEONICS TRIDENT acetabular UHMWPE inserts may be used with any HOWMEDICA OSTEONICS Stem of compatible head size, with suitable stem size and style to achieve total reconstructive replacement of the hip joint. The dome hole plugs are optional devices which are available to seal the HOWMEDICA OSTEONICS Acetabular Shell. The plugs are to be threaded into the dome holes of the shell.

#### Compatibility

##### Shell-to-Insert

- TRIDENT Shells can be used with TRIDENT Polyethylene Inserts, TRIDENT Ceramic Inserts and MDM liners.

##### Insert-to-Head

- HOWMEDICA OSTEONICS Polyethylene Inserts can be used with all HOWMEDICA OSTEONICS Metal or Ceramic Heads.
- Outside the U.S., the HOWMEDICA OSTEONICS Ceramic Inserts can be used with either HOWMEDICA OSTEONICS Alumina or BIOLOX *delta* Ceramic Heads.
- Outside the U.S., the HOWMEDICA OSTEONICS BIOLOX *delta* Ceramic Heads can be used with Protheos XLfit BIOLOX *delta* Ceramic Cups.
- Within the U.S., the HOWMEDICA OSTEONICS Ceramic Inserts can be used **only** with the HOWMEDICA OSTEONICS Alumina Ceramic Heads.

##### Head-to-Stem

- HOWMEDICA OSTEONICS C-taper Alumina Ceramic Heads can be used with HOWMEDICA OSTEONICS C-taper titanium stems. When used with adaptor sleeve 17-0000E, the C-taper alumina heads and C-Taper BIOLOX *delta* Ceramic Heads can be used with HOWMEDICA OSTEONICS V40 Taper Titanium Stems and V40 Taper CoCr Stems. When used with adaptor sleeve 1034-0000J, the C-taper alumina heads can be used with HOWMEDICA OSTEONICS Morse Taper Titanium Stems and Morse Taper CoCr Stems.
- HOWMEDICA OSTEONICS V40 Alumina Ceramic Heads (series 6565-0-xxx) can be used with HOWMEDICA OSTEONICS V40 Titanium and V40 Stainless Steel Stems.
- HOWMEDICA OSTEONICS C-taper BIOLOX *delta* Ceramic Heads can be used with HOWMEDICA OSTEONICS C-taper Titanium or C-taper CoCr Stems.
- HOWMEDICA OSTEONICS V40 BIOLOX *delta* Ceramic Heads can be used with HOWMEDICA OSTEONICS V40 Titanium, CoCr or Stainless Steel Stems.
- HOWMEDICA OSTEONICS BIOLOX *delta* Ceramic Universal Heads must be used in conjunction with HOWMEDICA OSTEONICS Universal Adaptor Sleeves 6519-T-XXX or 19-0XXXX.
- HOWMEDICA OSTEONICS BIOLOX *delta* Ceramic Universal Heads when used with Universal Adaptor Sleeves 19-0XXXX can be used with HOWMEDICA OSTEONICS C-taper Titanium or CoCr Stems.
- HOWMEDICA OSTEONICS BIOLOX *delta* Ceramic Universal Heads when used with Universal Adaptor Sleeves 6519-T-XXX can be used with HOWMEDICA OSTEONICS V40 Titanium or CoCr Stems.

- HOWMEDICA OSTEONICS BIOLOX *delta* Ceramic Universal Heads when used with Universal Adaptor Sleeves 6519-T-XXX can be used with Stainless Steel Stems in the USA, European Union member states and Australia. In other international markets, please consult the local market's Stryker representative.

### Acetabular Bone Screws

- HOWMEDICA OSTEONICS 6.5mm or 5.5mm Bone Screws can be used with the dome screw holes of the acetabular shells.

### Materials:

- |  |   |
|--|---|
| • ASTM F-620, ASTM F-136<br>Ti 6Al -4V ELI Alloy               | Acetabular Shell, Acetabular Insert Sleeve,<br>Adaptor Sleeve |
| • ASTM F-67 CP Titanium  | Arc-Deposited Coating, Dome Hole Plugs                        |
| • ASTM F-1580  | CP Ti Powder for Tritanium Revision Coating                   |
| • ASTM F-1185 Hydroxylapatite                                  | Hydroxylapatite Powder  |
| • ASTM F-603 Aluminum Oxide (Al <sub>2</sub> O <sub>3</sub> )  | Alumina Ceramic Head, Acetabular Bearing Insert               |
| • ISO 6474-2 Zirconia Toughened Aluminum Oxide                 | <i>delta</i> Ceramic Head                                     |
| • ASTM F-648 Ultra-High Molecular Weight Polyethylene (UHMWPE) | Acetabular Bearing Insert                                     |

### Indications

- Painful, disabling joint disease of the hip resulting from: degenerative arthritis, rheumatoid arthritis, post-traumatic arthritis or late stage avascular necrosis.
- Revision of previous unsuccessful femoral head replacement, cup arthroplasty or other procedure.
- Clinical management problems where arthrodesis or alternative reconstructive techniques are less likely to achieve satisfactory results.
- Where bone stock is of poor quality or is inadequate for other reconstructive techniques as indicated by deficiencies of the acetabulum.

### Contraindications

- Any active or suspected latent infection in or about the hip joint.
- Any mental or neuromuscular disorder which would create an unacceptable risk of prosthesis instability, prosthesis fixation failure, or complications in postoperative care.
- Bone stock compromised by disease, infection or prior implantation which cannot provide adequate support and/or fixation to the prosthesis.
- Skeletal immaturity.

### Warnings

In using this system, the surgeon should be aware of the following:

- In selecting patients for total joint replacements, the following factor is of extreme importance to the eventual success of the procedure: The patient's weight. The heavier the patient, the greater the load on the prosthesis. As the loads on the prosthesis increase, the chance a patient will suffer adverse reactions increases, including but not limited to failure of fixation, loosening, fracture and dislocation of the device and can lead to a decreased service life. The effect of these loads will be accentuated when a small sized prosthesis is used in larger patients. Overweight or obese patients impose greater loads on the prosthesis.

As obesity is a clinical diagnosis, we leave it to the surgeon to make the diagnosis based on his/her own clinical judgment. However, the World Health Organization (WHO) defines "overweight" as a BMI equal to or more than 25, and "obesity" as a BMI equal to or more than 30.

- Do not reassemble a ceramic head and stem. Once a ceramic head has been assembled to a stem taper, it should never be reassembled to that stem or subsequently assembled to any other stem. In addition, a ceramic head should only be assembled to an unused stem taper. Once a stem taper has been assembled to any femoral head, it should never be subsequently assembled to any ceramic head component due to deformation of the stem's taper locking mechanism during initial stem/head assembly.
- Do not allow polished bearing areas and machined taper surfaces to come in contact with hard or abrasive surfaces, as scratching or in any way damaging these surfaces can significantly affect the structural integrity.
- Adaptor sleeve must be fully seated on the stem taper before the head is impacted. In no instance should any attempt be made to pre-assemble the adaptor sleeve to the ceramic bearing head.
- When using an adaptor sleeve, improper seating of the head and/or adaptor sleeve may result in a discrepancy in neck length, component disassociation and/or dislocation.
- Modular Juncions: Mate modular components firmly to prevent dissociation. Machined taper surfaces must be clean, dry and firmly mated to ensure proper seating and assembly. Repeated assembly/disassembly or failure to clean, dry and firmly mate the components could compromise the taper lock and contribute to wear/corrosion and significant clinical consequences to the patient. (See Adverse Effects below).
- With the exception of Protheos XLfit Cups, do not substitute another manufacturer's device for any of the Howmedica Osteonics TRIDENT System components because design, material, or tolerance differences may lead to premature device and/or functional failure. Components of these systems have been specifically designed to work together. Any such use will negate the responsibility of Howmedica Osteonics Corp. for the performance of the resulting mixed component implant. Outside the U.S., the HOWMEDICA OSTEONICS BIOLOX *delta* Ceramic Heads can be used with Protheos XLfit BIOLOX *delta* Ceramic Cups.
- Howmedica Osteonics Corp. strongly advises against the use of another manufacturer's bone screws with any HOWMEDICA OSTEONICS Acetabular System component, due to variations which exist between screw head and screw seat configurations.
- Do not handle the hydroxylapatite treated regions as it may compromise the sterility or cause failure under load.
- Do not use C-Taper alumina heads with CoCr stems.
- Do not use V40 alumina heads with CoCr stems without an adaptor sleeve (17-0000E).
- Do not use C-Taper alumina heads with Stainless Steel (ORTHINOX) stems.
- Avoid excessive verticalization of shell, which may accelerate bearing wear.
- Do not contour or bend an implant because it may reduce its fatigue strength and cause failure under load.
- Improper seating of the head may result in a discrepancy in neck length, component disassociation and/or dislocation.
- Ensure appropriate selection of bone screw length and location to avoid damage to underlying soft tissue structures. Perforation of the pelvic wall can result in internal bleeding and possible damage to vital organs.
- Discard all damaged or mishandled implants. Never reuse an implant, even though it may appear undamaged. It may have small defects and internal stress patterns which may lead to early failure of the device.
- **Do not resterilize.**
- Patient post-operative pain. Inherent to all joint replacement is the risk that a patient will develop post-operative pain; pain is a commonly reported symptom regardless of the device implanted. The clinical literature reveals numerous potential causes of pain not directly related to the implant performance including, but not limited to, prior history of trauma and natural disease progression.

For patients who present with pain following implantation of any orthopaedic implant system, the physician should consider all potential causes of the symptoms identified in the clinical literature, including infection, soft tissue impingement, and possible adverse local tissue reactions associated with wear debris, metal ions or corrosion. Accurate diagnosis of the source of pain and directed, timely intervention is essential to ensuring effective treatment of pain.

## Precautions

- Before clinical use, the surgeon should thoroughly understand all aspects of the surgical procedure and limitations of the device. Physicians must instruct patients in the limitations of the prosthesis, including, but not limited to, the impact of excessive loading through patient weight or activity, and be taught to govern their activities accordingly. If the patient is involved in an occupation or activity which includes substantial walking, running, lifting, or muscle strain, the resultant forces can cause failure of the fixation, the device, or both. The prosthesis will not restore function to the level expected with normal healthy bone, and the physician must advise the patient against having unrealistic functional expectations.
- Appropriate selection, placement and fixation of the total hip components are critical factors which affect implant service life. As in the case of all prosthetic implants, the durability of these components is affected by numerous biologic, biomechanic and other extrinsic factors, which limit their service life. Accordingly, strict adherence to the indications, contraindications, precautions and warnings for this product is essential to potentially maximize service life.
- If the ceramic component(s) fracture necessitating revision, take special care to remove all ceramic debris from the joint. Any remaining fragments could accelerate wear of the replacement components.
- Use caution when handling ceramic components during assembly because of the brittle nature of ceramic material.
- Intentional removal of an acetabular component can be accomplished by careful use of cutting burrs, thin and narrow osteotomes and cautious extraction forces. A threaded metal shell can be removed by carefully unscrewing the shell in a counterclockwise direction. If difficulty is encountered, the preceding techniques may be employed.
- Removal of an unloosened arc deposited or hydroxylapatite surface treated implant may require the use of special instruments to disrupt the interface at the implant surface.
- Care should be taken not to cut through surgical gloves when handling any sharp-edged orthopaedic device.
- Surgeons should warn patients with metallic implants of the potential risks of undergoing a Magnetic Resonance Imaging (MRI) scan. The electromagnetic field created by an MRI scanner can interact with the metallic implant, resulting in displacement of the implant, heating of the tissue near the implant, implant damage or malfunction, or other undesirable effects. In addition, the presence of a metallic implant can produce an image artifact that may appear as a void region or geometric distortion of the true image. If the image artifact is near the area of interest, it may make the MRI scan uninformative or may lead to inaccurate clinical diagnosis or treatment.

## MRI SAFETY INFORMATION

The following Stryker Hip Implants have been evaluated for safety and compatibility in the MR environment:

- Alumina C-taper heads (17-xxxx and 17-xx-xxE); Alumina V40 heads (6565-0-xxx) C-taper BIOLOX *delta* Ceramic Heads (18-xxxx, 18-xx-x); V40 BIOLOX *delta* Ceramic Heads (6570-0-xxx); Universal BIOLOX *delta* Ceramic Heads (6519-1-xxx); C-taper Universal Adaptor Sleeve (19-xxxxT); V40 to C-Taper Adapter Sleeve (17-0000E) and Universal V40 Adaptor Sleeve (6519-T-xxx)
- Trident Hemispherical Solid Back Shell (500-01-xxx/500-11-xxx); Trident Hemispherical Cluster Shell (502-01-xxx/502-11-xxx); Trident Hemispherical Multi-Hole (508-11-xxx); Tritanium Acetabular Shell (509-02-xxx); Trident AD Acetabular Shell- Solid Back AD Acetabular Shell (540-01-xxx); Trident PSL Solid Back Acetabular Shell (540-11-xxx); Trident AD Acetabular Shell – Cluster (542-01-xxx); Trident PSL Cluster Acetabular Shell (542-11-xxx); Trident PSL T Shell with HA (542T-11-xxx) and Trident AD Acetabular shell (543-01-xxx)

- Trident X3, N2Vac and Crossfire Polyethylene Liners (620-00-xxx; 620-10-xxx; 621-00-xxx; 621-10-xxx; 623-00-xxx; 623-10-xxx; 641-00-xxx; 643-00-xxx; 661-00-xxx; 661-10-xxx; 663-00-xxx; 663-10-xxx; 723-00-xxx; 723-10-xxx; 743-00-xxx; 763-00-xxx; 763-10-xxx)

Non-clinical testing has demonstrated that the devices listed above are MR Conditional. A patient with these listed devices can be safely scanned in an MR scanner meeting the following conditions:

- Static magnetic field of 1.5-Tesla and 3.0-Tesla, only.
- Maximum spatial gradient magnetic field of 2,310 gauss/cm (23 T/m).
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system.
- Evaluation was performed using a quadrature body coil only.

Under the scan conditions defined above, these devices are expected to produce a temperature rise of less than 4.1 °C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 84 mm from the device when imaged with a gradient echo pulse sequence using a 3.0 T/128 MHz MRI system.

This MRI information is also available at <https://www.stryker.com/us/en/joint-replacement/MRI.html>.

The Trident TC Threaded Cup with Dome Hole Plug (7849-1-xxxx – not available for sale in the United States) has not been evaluated for safety and compatibility in the MR environment. The Trident TC Threaded Cup with Dome Hole Plug has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Trident TC Threaded Cup with Dome Hole Plug in the MR environment is unknown. Scanning a patient who has this device may result in patient injury

### Utilization and Implantation

- The surgeon must be completely familiar with the implant system and surgical protocol, and complete preoperative planning should be carried out.
- The suggested surgical procedure should be strictly adhered to. Proper assembly of the ceramic inserts and the ceramic heads to their mating taper surfaces and proper assembly technique are critical to the success of ceramic hip systems.
- The recommended trial components should be used for size determination, trial reduction and range of motion evaluation, thus preserving the integrity of the actual implants and their sterile packaging.
- Radiographic templates are available to assist in the preoperative prediction of component size and style.
- The Surgical Protocol for the HOWMEDICA OSTEONICS TRIDENT Acetabular Component System provides additional procedural information.

### Information for patients

- The surgeon must warn patients of surgical risks, and inform them of possible adverse effects. The surgeon must warn patients that the implant does not replicate the flexibility, strength, reliability, or durability of a normal healthy joint, that the implant can break or become damaged for numerous reasons, including as a result of strenuous activity or trauma, and that the implant has a finite service life and may need to be replaced in the future.
- The surgeon must warn patients of the limitations of the reconstruction and the need to protect the implant from full weight bearing until adequate healing has occurred. The surgeon must advise the patient to limit activities and protect the implant from strenuous activity, trauma or impact loading, and to follow the surgeon's instructions regarding activity level, follow-up care, and treatment.
- The surgeon must advise patients that the implant cannot be expected to withstand the same activity levels and loads as a normal healthy joint, and that the implant will not restore function to the level expected with normal healthy bone. If the patient is involved in an occupation or activity which includes substantial walking, running, lifting, or muscle strain, the resultant forces can cause failure of the fixation, the implant, or both. The surgeon must advise the patient against having unrealistic functional expectations.

- The surgeon must warn patients that strenuous activity, trauma or impact loading affecting the implant have been implicated in failure of the implant by loosening, fracture and/or wear of the implants. Many factors, including loosening of the implant components can result in increased production of wear particles, as well as damage to the bone, making successful revision surgery more difficult.
- Transient bacteremia can occur in daily life. Dental manipulation, endoscopic examination and other surgical procedures have also been associated with transient bacteremia. To help minimize the risk of infection at the implant site, it may be advisable to use antibiotic prophylaxis before and after such procedures. Surgeons should advise the patient to inform their doctors/dentists if they have an artificial joint replacement so that a decision can be made regarding antibiotic prophylaxis for such procedures.

### Adverse Effects

- While the expected life of total hip replacement components is difficult to estimate, it is finite. These components are made of foreign materials which are placed within the body for the potential restoration of mobility or reduction of pain. However, due to the many biological, mechanical and physicochemical factors which affect these devices but cannot be evaluated *in vivo*, the components cannot be expected to indefinitely withstand the activity level and loads of normal healthy bone.
- Dislocation of the hip prosthesis can occur due to inappropriate patient activity, trauma or other biomechanical considerations. Muscle and fibrous tissue laxity can also contribute to these conditions.
- Loosening of total hip components can occur. Early mechanical loosening may result from inadequate initial fixation, latent infection, premature loading of the prosthesis or trauma. Late loosening may result from trauma, infection, biological complications, including osteolysis, or mechanical problems, with the subsequent possibility of bone erosion and/or pain.
- Fracture of ceramic components has been reported in a small percentage of cases.
- Intraoperative fissure, fracture, or perforation of the femur, acetabulum or trochanter can occur due to impaction of the component into the prepared femoral canal or acetabulum. Postoperative femoral or acetabular fracture can occur due to trauma, the presence of defects, or poor bone stock.
- If bone screws are used, appropriate selection of bone screw length and location is essential to avoid damage to underlying soft tissue structures. Perforation of the pelvic wall can result in internal bleeding and possible damage to vital organs.
- Peripheral neuropathies, nerve damage, circulatory compromise and heterotopic bone formation may occur.
- Serious complications may be associated with any total joint replacement surgery. These complications include, but are not limited to: genitourinary disorders; gastrointestinal disorders; vascular disorders, including thrombus; bronchopulmonary disorders, including emboli; myocardial infarction or death.
- Acetabular pain may occur due to loosening of the implant.
- Although rare, sensitivity/allergic reactions to the materials in the implant have occurred in patients following joint replacement. Implantation of foreign material in tissues can result in immune responses and in histological reactions involving macrophages and fibroblasts.
- Adverse effects may necessitate reoperation, revision, arthrodesis of the involved joint, Girdlestone and/or amputation of the limb. Surgeons should advise patients of these potential adverse effects.
- Polyethylene particles and metal particles from mechanisms other than wear. Very small particles from metal and polyethylene components can be shed from non-articulating surfaces during normal use and over time. Although most of these particles stay in the relevant joint (i.e. contained in the synovium) or are trapped by surrounding scar tissue, microscopic particles can migrate throughout the body and on occasions have been described as accumulating in lymph nodes and other parts of the body. Although no significant medical complications have been reported as a result of these particles, their migration and/or accumulation in the body have been described in the literature. The long-term effects, if any, from these particles, are unknown. The long-term effects have been theorized to include:
  - Cancer: There is presently no scientific evidence that links metallic or polyethylene particles with cancer. However, the possibility cannot be ruled out.

- Lymphadenopathy and Accumulation in Other Tissues/Organs: There have been a few reports of the accumulation of particles in lymph nodes (proximal and distal). Although no medical complications or disease process has been reported as stemming from these accumulations, their existence should be recognized to facilitate diagnosis and avoid confusion with suspicious lesions, cancerous or otherwise.
- Systemic Disease: It is possible that some long-term effects may be demonstrated at some point in the future, but because there is very little scientific data suggesting association between migration of particles and systemic disease, it is believed that the benefits of these devices clearly outweigh the potential risks for any such theoretical long-term effect.
- Non-metallic wear debris. Wear debris is generated by interaction between components, as well as between components and bone, primarily through wear mechanisms of adhesion, abrasion and fatigue. Secondly, particulate can also be generated by third-body wear. With all implant devices, asymptomatic, localized progressive bone resorption (osteolysis) may occur around the prosthetic components as a consequence of foreign-body reaction to the particulate matter of cement, ultra-high molecular weight polyethylene (UHMWPE) and/or ceramic. Osteolysis can lead to future complications, including loosening, necessitating the removal and replacement of prosthetic components.
- Metal wear debris. Metal wear debris, metal ions and corrosion of metal implants. Generation of metal wear debris, metal ions and/or corrosion occurs whenever two surfaces are in contact, and at least one surface is metal. There have been reports in the literature of cases of adverse local tissue reactions associated with wear and/or corrosion at modular junctions formed by modular heads (stem/head interface). Local joint chemistry and/or other patient-specific conditions including, but not limited to infection may affect the potential for *in vivo* corrosion and its possible clinical consequences.
- Metal debris and corrosion by-products. Modular junctions may release metal debris and/or metal ions due to fretting, galvanic corrosion, crevice corrosion, or other processes. There are several factors involved in these processes, including forces across a junction, which are not fully understood. These corrosion products or metal debris may affect the tissues surrounding the implant and adversely affect the duration of service life. There have been reports of patients developing adverse local tissue reactions (including, but not limited to, tissue necrosis, pseudotumors, cysts and fluid accumulations, metallosis and aseptic lymphocyte dominated vasculitis associated lesions), elevated metal ion levels in the blood and/or urine and hypersensitivity/allergic reactions associated with corrosion and/or wear-related debris in the implant vicinity. Affected patients may present with symptoms similar to those associated with infection, including pain (most likely during weight-bearing) and swelling at the local joint site. These reactions should be carefully monitored and may result in early revision surgery. Medical literature describes systemic reactions to by-products arising from contemporary metal on metal articulating bearing surfaces. While there are no metal on metal bearing interfaces in this system, it could be theorized that similar systemic reactions may arise from fretting and corrosion arising from any metal interface.
- Surgeons should warn patients of the above listed potential effects including the finite service life of the device and the need for post-operative protection of the implant.

## **Sterilization**

- These components have been sterilized by either gamma radiation, hydrogen peroxide, or ethylene oxide. Refer to the package label for the sterilization method.
- Do NOT re-sterilize.
- Inspect the packaging of ALL sterile products for flaws before opening. In the presence of any flaws, assume the product is not sterile.
- Take care to prevent contamination of ANY components.
- Discard ALL nonsterile or contaminated product.
- Device should not be used after the expiry date displayed on the label as packaging had not been validated beyond this date.
- Single use devices cannot be explanted and subsequently reimplanted as the physical forces exerted by these actions may compromise the physical integrity, dimensions and/or surface finishes of the devices.

Also, sterility cannot be assured for reused devices as cleaning and re-sterilization procedures have not been verified.


### Disposal

If a device is being returned for evaluation, please contact your local Stryker representative for shipping/handling information. If the device is not being returned to Stryker, implant components are to be disposed of in accordance with applicable laws, rules, and regulations for the disposal of biohazardous waste. Follow all guidelines for biohazardous waste in accordance with the Centers for Disease Control and Prevention guidelines as well as applicable federal/national, state and local regulations. As part of the disposal process, verify that the implant in its entirety has been explanted from the surgical site.

### DOME HOLE PLUG ASSEMBLY INSTRUCTIONS

#### INSERTION:

- Once the acetabular shell is seated in the acetabulum, the dome hole plug may be inserted. Place the dome hole plug onto the captive twist head of the driver (secure by tapping on a hard surface). Insert the dome hole plug into the threaded dome hole of the shell. Turn the driver clockwise until the plug is seated firmly. Extract driver from plug.

#### REMOVAL:

- Removal of the plug is the same as insertion, except the driver is turned counterclockwise.

**CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A LICENSED PHYSICIAN.**

The Trident TC Threaded Cup with Dome Hole Plug is not available for sale in the United States.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademark(s) or service mark(s): Crossfire, Howmedica, Orthinox, Osteonics, Stryker, Trident, V40 and X3. All other trademarks are trademarks of their respective owners or holders.

BILOX *delta* is a registered trademark of CeramTec AG.

Refer to product label for CE Mark Status and Legal Manufacturer. The CE mark is only valid if also found on the product label.

The following table contains a list of abbreviations that are used on Howmedica Osteonics Corp. product labeling:

Term	Abbreviation	Term	Abbreviation
Alpha Code	ALPH CDE	Neck	NK
Angle	ANG	Offset	OFFST
Degree	DEG or °	Outer Diameter	OD
Diameter	DIA	Right	RT ►
Extra Deep	XDP	Screw Holes	SCR HLS
Extra Large	XLGE	Side	SDE
Extra Small	XSM	Size	SZE
Head	HD	Small	SM
Height	HT	Standard	STD
Inner Diameter	ID	Taper	TPR
Insert	INSR	Thickness	THKNS
Large	LGE	Type	TYP
Left	◄ LFT	With	W/
Length	LNTH	Without	W/O
Medium	MED		

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- C-Konus-Köpfe aus Aluminiumoxid (17-xxxxx und 17-xx-xxE); V40-Köpfe aus Aluminiumoxid (6565-0-xxx) C-Konus BIOLOX *delta* Keramikköpfe (18-xxxx, 18-xx-x); V40 BIOLOX *delta* Keramikköpfe (6570-0-xxx); Universal BIOLOX *delta* Keramikköpfe (6519-1-xxx); C-Konus-Universal-Adapterhülse (19-xxxxT); Adapterhülse V40 zu C-Konus (17-0000E) und Universal-V40-Adapterhülse (6519-T-xxx)
- Trident halbkugelförmige, stabile Schale (500-01-xxx/500-11-xxx); Trident halbkugelförmige Cluster-Schale (502-01-xxx/502-11-xxx); TRIDENT halbkugelförmig mit Mehrfachbohrung (508-11-xxx); Tritanium Pfannenschale (509-02-xxx); Trident AD Pfannenschale- metallverstärkte AD Pfannenschale (540-01-xxx); TRIDENT PSL metallverstärkte Pfannenschale (540-11-xxx); Trident AD Pfannenschale – Cluster (542-01-xxx); Trident PSL Cluster-Pfannenschale (542-11-xxx); Trident PSL T-Schale mit HA (542T-11-xxx) und Trident AD Pfannenschale (543-01-xxx)
- Trident X3, N2Vac und Crossfire Polyethyleninlagen (620-00-xxx; 620-10-xxx; 621-00-xxx; 621-10-xxx; 623-00-xxx; 623-10-xxx; 641-00-xxx; 643-00-xxx; 661-00-xxx; 661-10-xxx; 663-00-xxx; 663-10-xxx; 723-00-xxx; 723-10-xxx; 743-00-xxx; 763-00-xxx; 763-10-xxx)

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- Teste in allumina a cono C (17-xxxxx e 17-xx-xxE); teste in allumina V40 (6565-0-xxx); teste in ceramica a cono C BIOLOX *delta* (18-xxxx, 18-xx-x); teste in ceramica V40 BIOLOX *delta* (6570-0-xxx); teste in ceramica universali BIOLOX *delta* (6519-1-xxx); manicotto adattatore universale a cono C (19-xxxxT); manicotto adattatore da V40 a cono C (17-0000E) e manicotto adattatore universale V40 (6519-T-xxx)
- Cotile emisferico solido Trident (500-01-xxx/500-11-xxx); cotile emisferico con fori a grappolo Trident (502-01-xxx/502-11-xxx); cotile emisferico multifori Trident (508-11-xxx); cotile acetabolare Tritanium (509-02-xxx); cotile acetabolare Trident AD - cotile acetabolare solido Trident AD (540-01-xxx); cotile acetabolare solido Trident PSL (540-11-xxx); cotile acetabolare Trident AD - con fori a grappolo (542-01-xxx); cotile acetabolare con fori a grappolo Trident PSL (542-11-xxx); cotile Trident PSL T con HA (542T-11-xxx) e cotile acetabolare Trident AD (543-01-xxx)
- Rivestimenti in polietilene Trident X3, N2Vac e Crossfire (620-00-xxx; 620-10-xxx; 621-00-xxx; 621-10-xxx; 623-00-xxx; 623-10-xxx; 641-00-xxx; 643-00-xxx; 661-00-xxx; 661-10-xxx; 663-00-xxx; 663-10-xxx; 723-00-xxx; 723-10-xxx; 743-00-xxx; 763-00-xyy; 763-10-xxx)

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- Aluminiumoxide C-conuskoppen (17-xxxx en 17-xx-xxE); aluminiumoxide V40 koppen (6565-0-xxx) C-conus BIOLOX *delta* keramische koppen (18-xxxx, 18-xx-x); V40 BIOLOX *delta* keramische koppen (6570-0-xxx); universele BIOLOX *delta* keramische koppen (6519-1-xxx); de C-conus universele adapterhuls (19-xxxxT); de V40 naar C-conus adapterhuls (17-0000E) en de universele V40 adapterhuls (6519-T-xxx)
- De Trident halfronde buitencup met dichte achterzijde (500-01-xxx/500-11-xxx); de Trident halfronde cluster-buitencup (502-01-xxx/502-11-xxx); de Trident halfronde buitencup met meerdere gaten (508-11-xxx); de Tritanium acetabulum-buitencup (509-02-xxx); de Trident AD acetabulum-buitencup-met dichte achterzijde AD acetabulum-buitencup (540-01-xxx); de Trident PSL acetabulum-buitencup met dichte achterzijde (540-11-xxx); de Trident AD acetabulum-buitencup-cluster (542-01-xxx); de Trident PSL cluster acetabulum-buitencup (542-11-xxx); de Trident PSL T buitencup met HA (542T-11-xxx) en de Trident AD acetabulum-buitencup (543-01-xxx)
- Trident X3, N2Vac en Crossfire polyethyleen-liners (620-00-xxx; 620-10-xxx; 621-00-xxx; 621-10-xxx; 623-00-xxx; 623-10-xxx; 641-00-xxx; 643-00-xxx; 661-00-xxx; 661-10-xxx; 663-00-xxx; 663-10-xxx; 723-00-xxx; 723-10-xxx; 743-00-xxx; 763-00-xxx; 763-10-xxx)

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