# MCIS Vol.5

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
Fig.2: A com AN rec 3x5mm、Neck 5mm

## Enterprise VRD (no tip) 22mm

Coils: Cashmere 14 3x6cm,  
Deltaplush 2x4, 2x2cm  
total 12cm VER 55.1%

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Fig.3: BA top unruptured AN rec 7.5x5mm, Neck 4.5mm

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- 1) Durst CR, Khan P, Gaughen J, et al. Direct comparison of Neuroform and Enterprise stents in the treatment of wide-necked intracranial aneurysms Clin Radiol.69(12):e471-6, 2014
- 2) Heller RS, Dandamudi V, Calnan D, Malek AM. Neuroform intracranial stenting for aneurysms using simple and multi-stent technique is associated with low risk of magnetic resonance diffusion-weighted imaging lesions. Neurosurgery 73(4): 582-90, 2013
- 3) Kadkhodayan Y, Rhodes N, Blackburn S, et al. Comparison of Enterprise with Neuroform stent-assisted coiling of intracranial aneurysms. Am J Roentgenol. 200(4): 872-8, 2013
- 4) King B, Vaziri S, Singla A, Fargen KM, Mocco J. Clinical and angiographic outcomes after stent-assisted coiling of cerebral aneurysms with Enterprise and Neuroform stents: a comparative analysis of the literature. J Neurointerv Surg. Dec; 7(12):905-9, 2015
- 5) Feng Z, Zhang L, Li Q, et al. Endovascular treatment of wide-neck anterior communicating artery aneurysms using the LVIS junior stent. J Clin Neurosci. 22(8): 1288-91, 2015
- 6) Feng Z, Li Q, Zhao R, et al. Endovascular treatment of middle cerebral artery aneurysms with the LVIS junior stent. J Stroke Cerebrovasc Dis. 24(6): 1357-62, 2015

All photographs taken by National Hospital Organization Mito Medical Center.  
Results from case studies are not predictive of results in other cases. Results in other cases may vary.

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Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Excelsior XT-27, Neuroform, Neuroform EZ, Target. All other trademarks are trademarks of their respective owners or holders.

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Copyright © 2016 Stryker