TACT for Peripheral Vol.3

stryker®
Neurovascular

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1st Coil : Interlock 3mm x 12cm
 2nd Coil : Interlock 3mm x 12cm
 3rd Coil : Target Helical Ultra 3mm x 4cm
 4th Coil : Target Helical Ultra 3mm x 4cm
 5th Coil : Target Helical Ultra 2.5mm x 6cm
 6th Coil : Target Helical Ultra 2.5mm x 6cm
 7th Coil : Interlock 3mm x 12cm
 8th Coil : Interlock 3mm x 12cm
 9th Coil : Target Helical Ultra 2mm x 8cm
 10th Coil : Target Helical Ultra 2mm x 8cm
 11th Coil : Target Helical Ultra 2mm x 8cm
 12th Coil : Target Helical Ultra 2mm x 6cm

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- 1) Guglielmi G, et al: Electrothrombosis of saccular aneurysms via endovascular approach. Part 2: Preliminary clinical experience. J Neuro-surg 75:8-14, 1991.
- 2) Sluzewski M, et al: Relation between Aneurysm Volume, Packing, and Compaction in 145 Cerebral Aneurysms Treated with Coils. Radiology 231:653-658, 2004.
- 3) Kawanabe Y, et al: Endovascular occlusion of intracranial aneurysms with Guglielmi detachable coils: correlation between coil packing density and coil compaction. Acta Neurochir 143:451-455, 2001.
- 4) Pollak JS, et al: Clinical and Anatomic Outcomes after Embolotherapy of Pulmonary Arteriovenous Malformations. JVIR 17:35-45, 2006.

All photographs taken by Fujita Health University Hospital.
 Results from case studies are not predictive of results in other cases. Results in other cases may vary.

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Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Excelsior 1018, Target, Transend. All other trademarks are trademarks of their respective owners or holders.

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