TACT Vol.11

Target Advanced Clinical Tactics

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
Fig.1: Preoperative MR Imaging showing enlarged left superior ophthalmic vein (arrow)

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

Fig.2: Preoperative digital subtraction angiography showing left carotid cavernous fistula fed by multiple meningeal feeders. Note shunting point is located in the posterior-internal portion of the cavernous sinus (arrow) and drains into left superior ophthalmic vein(double arrows).

A: Right external carotid angiogram (A-P view), B,C: left external carotid angiogram (A-P view and lateral view)


Fig.3: Intraoperative skull X ray showing the micro catheter located into the shunting point from left inferior petrosal sinus through origin of the superior ophthalmic vein.

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- Hayashi K, Horie N, Morofuji Y, Fukuda S, Yamaguchi S, Morikawa M. Quadruple coaxial catheter system on transvenous embolization for dural arteriovenous fistula. Neurol Res. 2015 Apr;37(4):328-31
- Fig 1: Preoperative MR Imaging showing enlarged left superior ophthalmic vein (arrow)
- Fig 2: Preoperative digital subtraction angiography showing left carotid cavernous fistula fed by multiple meningeal feeders. Note shunting point is located in the posterior-internal portion of the cavernous sinus (arrow) and drains into left superior onthalmic vein(double arrows)
 - ophthalmic vein(double arrows). **A**: Right external carotid angiogram (A-P view), **B,C**: left external carotid angiogram (A-P view and lateral view)
- Fig 3: Intraoperative skull X ray showing the micro catheter located into the shunting point from left inferior petrosal sinus through origin of the superior ophthalmic vein.

 A: A-P view, B: lateral view
- Fig 4: Intraoperative digital subtraction angiography showing target embolization with Target XL Mini $(A \to B \to C)$
- Fig 5: Postoperative digital subtraction angiography showing complete obliteration of the CCF $\bf A$: A-P view, $\bf B$: lateral view


Fig.4: Intraoperative digital subtraction angiography showing target embolization with Target XL Mini ($\mathbf{A} \rightarrow \mathbf{B} \rightarrow \mathbf{C}$)


Fig.5: Postoperative digital subtraction angiography showing complete obliteration of the CCF A: A-P view, B: lateral view

All photographs taken by Nagasaki University. Results from case studies are not predictive of results in other cases. Results in other cases may vary.

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Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Target, Target XL. All other trademarks are trademarks of their respective owners or holders.

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