
# Triathlon<sup>®</sup> TS

knee system

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## References

1. Stryker Test Report RD-07-025.
2. D'Lima, D., Chen, P., Kester, M., Impact of Patellofemoral design on Patellofemoral Forces and Polyethylene Stresses JBJS 2003: 85:85-93.
3. Hitt K, Shurman JR 2nd, Greene K, McCarthy J, Moskal J, Hoeman T, Mont MA. Anthropometric measurements of the human knee: correlation to the sizing of current knee arthroplasty systems. J Bone Joint Surg Am. 2003;85-A Suppl 4:115-22.
4. Hollisater AM, Jatana S, Singh AK, Sullivan WW, Lupichuk AG. The axes of rotation of the knee. Orthopaedic Laboratory, Harbor-UCLA Medical Center, Torrance. Clin Orthop. 1993 May;(290):259-68.
5. Freeman MAR, Pinskerova V, Hill PF, Vedi V, Williams A, Iwaki H, Kadoya Y, Kobayashi A, Nakagawa S, Sakamoto H, Todo S, Yamano Y, Brandsson S, Karrholm J. Tibiofemoral movement - Studies in the cadaveric and living knee using MRI, dissection and RSA; British Editorial Society of Bone and Joint Surgery, 2001.
6. Churchill DL, Incavo SJ, Johnson CC, Beynon BD. The transepicondylar axis approximates the optimal flexion axis of the knee. Clin Orthop. 1998 Nov;(356):111-8.
7. Robertson O, Knutson K, Lewold S, Lidgren L. The Swedish Knee Arthroplasty Register, Outcome with special emphasis on 1988-1997. Dept. of Orthopedics, University Hospital.
8. Engh GA, Lounici S, Rao AR, Collier MB. In vivo deterioration of tibial baseplate locking mechanisms in contemporary modular total knee components. JBJS 83A: 1660-1665, 2001.
9. Parks NL, Engh GA, Topoleski LD, Emperado J. Modular tibial insert micromotion: A concern with contemporary knee implants. Clin Orthop 356:10-15, 1998.
10. Series 7000 locking mechanism.
11. Rao AR, Engh GA, Collier MB, Lounici S. Tibial interface wear in retrieved total knee components and correlations with modular insert motion. J Bone Joint Surg Am. 2002 Oct;84-A(10):1849-55.
12. Tsao, A., Mintz, L., et al. Failure of the Porous-Coated Anatomic Prosthesis due to severe polyethylene wear JBJS vol. 75-A, No. 1 January 1993.
13. Haas, S., Insall, J.N. et al. Revision Total Knee Arthroplasty with the use of modular components with stems without cement.
14. Mahoney, O., Kinsey, T., Modular femoral Offset Stems Facilitate Joint Line Restoration in Revision Knee Arthroplasty Clin. Orthop. #446, pp. 93-98.
15. Stryker Test Report RD-06-013.
16. Stryker Test Data RD-03-041 and RD-04-027.
17. 8th Annual EFORT Congress Florence Italy 2007 - Dr. Stuckenburg Presentation.
18. Schoderbek, R.J., Brown, T., et al. Extensor Mechanism Disruption after Total Knee Arthroplasty Clin. Orth. Num. 446, pp. 176-185.
19. Rand, J., Extensor Mechanism Complications following Total Knee Arthroplasty JBJS Vol. 86-A, Num. 9, Sept. 2004.
20. Stryker® Orthopaedics Triathlon® CR Tibial Inserts made from X3® UHMWPE, 5530-G-409 show a 68% reduction in volumetric wear rate versus the same insert fabricated from N2\Vac™ gamma sterilized UHMWPE, 5530-P-409. The insert tested was Size 4, 9 mm thick. Testing was conducted under multiaxial knee simulator (multi-station MTS knee joint simulator) for five million cycles using appropriate size CoCr counterfaces, a specific type of diluted calf serum lubricant and the motion and loading conditions, representing normal walking, outlined in ISO/DIS 14243-3. Volumetric wear rates were  $17.7 \pm 2.2$  mm<sup>3</sup>/106 cycles for standard polyethylene inserts and  $5.7 \pm 1.5$  mm<sup>3</sup>/106 cycles for test samples. Test inserts were exposed to a gas plasma sterilization process. In vitro knee wear simulator tests have not been shown to quantitatively predict clinical wear performance.
21. Stryker® Orthopaedics Triathlon® PS Tibial Inserts made of X3® UHMWPE, 5532-G-409 show a 64% reduction in volumetric wear rate versus the same insert fabricated from N2\Vac™ gamma sterilized UHMWPE, 5532-P-409. The insert tested was Size 4, 9 mm thick. Testing was conducted under multiaxial knee simulator (multi-station MTS knee joint simulator) for five million cycles using a size 7 CoCr counterfaces, a specific type of diluted calf serum lubricant and literature or fluoroscopy based motion and loading conditions representing stair climbing. b, c Volumetric wear rates were  $3.6 \pm 0.61$  mm<sup>3</sup>/106 cycles for standard polyethylene inserts and were  $1.3 \pm 0.44$  mm<sup>3</sup>/106 cycles for test samples. Test inserts were exposed to a gas plasma sterilization process. In vitro knee wear simulator tests have not been shown to quantitatively predict clinical wear performance.
22. Stryker® Orthopaedics Scorpio® CR Tibial Inserts made from X3® UHMWPE, 72-22-0708, show a 79% reduction in volumetric wear rate versus the same insert fabricated from N2\Vac™ gamma sterilized UHMWPE, 72-2-0708. The insert tested was Size 7, 8 mm thick. Testing was conducted under multiaxial knee simulator (multi-station MTS knee joint simulator) for five million cycles using appropriate size CoCr counterfaces, a specific type of diluted calf serum lubricant and the motion and loading conditions, representing normal walking, outlined in ISO/DIS 14243-3. Volumetric wear rates were  $34.6 \pm 1.5$  mm<sup>3</sup>/106 cycles for standard polyethylene inserts and  $7.3 \pm 0.7$  mm<sup>3</sup>/106 cycles for test samples. Test inserts were exposed to a gas plasma sterilization process. In vitro knee wear simulator tests have not been shown to quantitatively predict clinical wear performance.
23. Stryker® Orthopaedics Scorpio® PS Tibial Inserts made from X3® UHMWPE, 72-23-0708, show a 77% reduction in volumetric wear rate versus the same tibial insert fabricated from N2\Vac™ gamma sterilized UHMWPE, 72-3-0708. The insert tested was Size 7, 8 mm thick. Testing was conducted under multiaxial knee simulator (multi-station MTS knee joint simulator) for five million cycles using appropriate size CoCr counterfaces, a specific type of diluted calf serum lubricant and literature or fluoroscopy based motion and loading conditions representing stair climbing. b, c Volumetric wear rates were  $35.8 \pm 1.7$  mm<sup>3</sup>/106 cycles for standard polyethylene inserts and were  $8.2 \pm 0.7$  mm<sup>3</sup>/106 cycles for test samples. Test inserts were exposed to a gas plasma sterilization process. In vitro knee wear simulator tests have not been shown to quantitatively predict clinical wear performance.
24. Stryker® Orthopaedics Trident® Acetabular Inserts made of X3® UHMWPE (unsterilized), 721-00-32E, show a 97% reduction in volumetric wear rate versus the same insert fabricated from N2\Vac™ gamma sterilized UHMWPE, 620-00-25E. The insert tested was 7.5 mm thick with an inner diameter of 32 mm. Testing was conducted under multiaxial hip joint simulation for 5 million cycles using a 32 mm CoCr articulating counterface and calf serum lubricant. X3® UHMWPE Trident® acetabular inserts showed a net weight gain due to fluid absorption phenomena but yielded a positive slope and wear rate in linear regression analysis. Volumetric wear rates were  $46.39 \pm 11.42$  mm<sup>3</sup>/106 cycles for N2\Vac™ gamma sterilized UHMWPE inserts and  $1.35 \pm 0.68$  mm<sup>3</sup>/106 cycles for X3® UHMWPE (unsterilized) Trident® Acetabular Inserts. Although in-vitro hip wear simulation methods have not been shown to quantitatively predict clinical wear performance, the current model has been able to reproduce correct wear resistance rankings for some materials with documented clinical results. d, e, f.
25. Yau SS, Wang A, Essner A, Manley M, Dumbleton J. "Sequential Irradiation and Annealing of Highly Crosslinked Polyethylenes: Resist Oxidation without Sacrificing Physical/Mechanical Properties," Transactions of the 51st Annual Meeting of the Orthopaedic Research Society; Washington, DC, 2005: 1670.
26. D'Antonio J, Manley M, Capello W, Bierbaum B, et al. "Five-year Experience with Crossfire® Highly Cross-linked Polyethylene". CORR, No 441, 2005, pp. 143-150.
27. Essner A, Yau SS, Schmidig G, Wang A, Dumbleton J, Manley M, Serekian P. "Reducing Hip Wear Without Compromising Mechanical Strength: A Next Generation Crosslinked and Annealed Polyethylene," Transactions of the 5th Combined Meeting of the Orthopaedic Research Societies of Canada, USA, Japan, and Europe; Calgary, Canada, 2004:80.
28. Greene, K.A. Range of Motion: Early Results from the Triathlon® Knee System, Stryker Literature Ref #LSA56., 2005.
29. Hemida, J et al. The Effect of Oxidative Aging on the Wear Performance of Highly Crosslinked Polyethylene Knee Inserts under Conditions of Severe Malalignment. Journal of Orthopaedic Research, 26.12, 2008: 1585-1590

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Copyright © 2021 Stryker