



Accolade® II

Femoral Hip Stem

Design rationale



**Designed to fit more patients,
designed to fit your approach.**



Addressing modern demands with novel technology.

The global THA population is evolving to include a younger¹, more active², and more demanding³ patient. Many femoral stem designs on the market today predate the emergence of this novel demand. Subsequently, an opportunity to enhance the conventional femoral stem design exists.

Conventional tapered wedge femoral stems have achieved popularity due to their simplicity and excellent clinical results^{4,5}. Despite these results, literature indicates that there are still unmet clinical needs^{6,7,8}. Incidence of subsidence⁶, distal-only implant engagement⁷, and peri-prosthetic fracture⁸ suggest a clinical need for an improved implant fit for this novel patient population.


Stryker, along with key industry leaders, embarked to develop a novel femoral stem. This design would build from the sound principles of tapered wedge philosophy to meet the unique needs of the current patient population. At the head of this development was a unique technology called Stryker Orthopaedic Modeling and Analytics, or SOMA. As a system that enables population-based design, SOMA has powerful functionality with which to design, model, and analyze novel orthopaedic devices.

Stryker utilized SOMA technology to design a novel stem building upon the conventional tapered wedge femoral design, incorporating unique features to allow for an enhanced implant fit in today's patient population⁹. By establishing an increased canal fit and fill⁹, Accolade II has been shown to allow for improved stability¹⁰, decreased intraoperative femoral fractures¹¹, as well as excellent survivorship and functional outcomes^{12,13}, ultimately leading to satisfied patients^{13,14}.



SOMA technology

Utilizing the proprietary SOMA technology, Stryker was able to complete one of the largest proximal femoral bone morphology studies ever undertaken¹³. An illustrated look at the process by which SOMA technology is employed in implant design is described below.




CT acquisition

The SOMA database continues to acquire new CT scans, and currently contains over 16,500 bones¹⁴




Segmentation

Once acquired, all bones are segmented into inner and outer cortices



Analysis

Using SOMA tools, bone morphology can be studied in a highly accurate and reproducible manner.



Design input

The results of these studies, such as the population Canal Flare Index, can be utilized in implant design



Validation

The resulting implant design can then be validated using SOMA fitting tools.

Three key SOMA-designed features

Bone morphology data allowed Stryker an unprecedented look at femoral anatomy, and assisted in the design of a novel femoral stem. The SOMA input* was instrumental in the establishment of three key design features of Accolade II:



Unique size-specific medial curvature

increasing proximal conformity to improve primary stability^{14,15}



Enhanced proximal-distal proportions¹⁶

shown to mimic canal anatomy to avoid distal only engagement and achieve cortical fit^{8,16}



Optimized stem length

enables muscle-sparing approaches without sacrificing stability^{15,17}

*SOMA-design of Accolade II based on 556 CT scans.



Unique size-specific medial curvature

Initial stability is critical to long-term implant performance⁶. Early subsidence and micromotion has been established as a strong indicator for implant failure⁶. Initial stability may be increased by creating a higher conformity between the implant and the femoral cortices, leading to a larger area of contact¹⁹.

Analyzing the SOMA morphology study data*, it was observed that a constant medial curvature may not allow for a conforming canal fit throughout varying femoral sizes.

This population-based input influenced Stryker to incorporate the market's first size-specific medial curvature into Accolade II. This feature was designed to enable a more conforming proximal cortical fit¹⁵, which has been shown to allow for improved implant stability¹⁹.

Using femurs from the SOMA database*, fit patterns of Accolade II can be compared to conventional tapered wedge designs. The three examples below illustrate how Accolade II achieves a more conforming canal fit throughout varying bone sizes.

Small femur



Accolade II

Conventional tapered wedge

Medium femur



Accolade II

Conventional tapered wedge

Large femur




Accolade II

Conventional tapered wedge



Enhanced proximal-distal proportions

Distal-only engaged femoral stems can experience stress shielding¹⁹ and consequently may lead to elevated failure rates due to loosening and migration⁹. In order to better mimic the femoral anatomy and avoid distal-only engagement, a more anatomic implant growth rate is needed.

Utilizing the SOMA femoral morphology study*, a more anthropomorphic proximal-distal stem growth rate was identified. This rate led to enhanced implant proportions¹³, as the distal geometry of Accolade II increases in size less than the proximal geometry. These proportions enable Accolade II to achieve a significantly better canal fit and fill,⁹ and Accolade II has shown a decreased incidence of distal-only engagement⁹.



Comparing implant fit




In the graph above, the proximal (P) and distal (D) measurements of a population of 556 femurs were plotted (black dots) against the corresponding stem diameters of Accolade II (gold) and a conventional tapered wedge design (grey). Accolade II achieved more fully-conforming and proximal-only fit types compared to the conventional design, while subsequently **reducing distal-only fit by 14%**.¹⁶

Optimized stem length

Popularity of muscle-sparing approaches and bone conserving fundamentals have led to a trend in shortening of stem length²⁰. However, there exists a complex relationship between stem length and implant stability²⁰. Shortening stem length without geometry optimization has been shown to increase the potential for micromotion²⁰, which is a strong indicator for implant failure⁶.

Accolade II utilized the SOMA database* and stability analyses to establish an optimized length for each stem size which not only accommodates muscle-sparing approaches²⁰, but demonstrated improved initial stability¹⁰.

“simply shortening a standard tapered wedge design may reduce the primary stability”¹⁷



Designed to fit your approach


Muscle-sparing surgical approaches continue to gain popularity, due to the potential patient benefits of faster recovery^{21,22}, less pain^{21,22} and greater satisfaction²³.

Stryker's portfolio of muscle-sparing techniques features modern instrumentation and dynamic Medical Education programs to support the **Direct Anterior Approach** and the **Direct Superior Approach**.

The new Direct Superior Approach was designed for surgeons who prefer the fundamentals and familiarity of the posterior approach, but seek to provide the next evolution in muscle-sparing THA surgical techniques for their patients.

Stryker has developed a training platform to help surgeons during the learning curve of a new surgical approach, Stryker's Training Academy.

For access to Stryker's Training Academy, contact your local Stryker sales representative.



DIRECT SUPERIOR APPROACH stryker

Toolkit Return to menu

Scenario: Femoral Neck Resection

You are doing your first DSA case. The procedure has gone well so far. You have dissected and retracted the confluence of the piriformis and obturator internus tendons and performed the capsulectomy, but you run into difficulty during the osteotomy. While the hip may be dislocated, you still cannot adequately visualize if the resection has been made at the templated level.

Patient profile:
• 58-year-old female • BMI of 32

What should you do? Make the best selection.

A Decision A
Internally rotate the femur to 90 degrees.

B Decision B
Perform a subcapital osteotomy. Then replace the femoral retractors to visualize the remaining neck and make a second neck cut. Assess the combined neck resections to the pre-op templated neck resection height.

C Decision C
Reduce the head back into the acetabulum and make an in situ neck resection.



Clinical performance

“Significantly better overall canal **fit**”⁹
than conventional tapered wedge design


over
400,000
implanted worldwide²⁴

over
16,500
bones*
in SOMA database¹⁴

99.2%
survivorship
for Accolade II pursuant
to a 3.5-year mean study¹²

less than
0.1 mm
subsidence pursuant to
2-year RSA study¹⁰

5x
less intraoperative
fractures
observed compared to
conventional tapered wedge¹¹



Accolade II Implant catalog numbers

Part number	Size	Neck angle
6720-0027	0	132°
6720-0127	1	
6720-0230	2	
6720-0330	3	
6720-0435	4	
6720-0535	5	
6720-0635	6	
6720-0737	7	
6720-0837	8	
6720-0937	9	
6720-1040	10	
6720-1140	11	

Part number	Size	Neck angle
6721-0027	0	127°
6721-0127	1	
6721-0230	2	
6721-0330	3	
6721-0435	4	
6721-0535	5	
6721-0635	6	
6721-0737	7	
6721-0837	8	
6721-0937	9	
6721-1040	10	
6721-1140	11	

References:

- National Center for Health Statistics. Health, United States, 2010: With Special Feature on Death and Dying, pgs. 17 & 60.
- Kurtz S, Lau E, Ong K, Zhao K, Kelly MP, Bozic KJ. Future Young Patient Demand for Primary and Revision Joint Replacement. CORR 2009; 467:2606-2612.
- Haddad F. A Positive End to 2016. Bone Joint J. 2016;98-B:1569-70.
- Casper D, Kim G, Restrepo C, Parvizi J, Rothman R. Primary Total Hip Arthroplasty With an Uncemented Femoral Component. The Journal of Arthroplasty. Vol 26, Issue 6, 2011; 838-841.
- McLaughlin JR, Lee KR. Total hip arthroplasty with an uncemented femoral component. J Bone Joint Surg 2008; 90:1290.
- Paul van der Voort, Bart G Pijls, Marc J Nieuwenhuijse, Jorrit Jasper, Marta Fiocco, Josepha W M Plevier, Saskia Middeldorp, Edward R Valstar & Rob G H H Nelissen (2015) Early subsidence of shape-closed hip arthroplasty stems is associated with late revision, Acta Orthopaedica, 86:5, 575-585, DOI: 10.3109/17453674.2015.1043832.
- Cooper J, et al. Distal Fixation of Proximally Coated Tapered Stems May Predispose to a Failure of Osteointegration. Journal of Arthroplasty. 2011. 10.1016/j.arth.2011.04.003.
- Abdel MP, et al. Epidemiology of periprosthetic fracture of the femur in 32 644 primary total hip arthroplasties. Bone Joint J 2016;98-B:461-7.
- Issa, K. et al. Radiographic Fit and Fill Analysis of a New Second-Generation Proximally Coated Cementless Stem Compared to its Predicate Design. Journal of Arthroplasty (2013). <http://dx.doi.org/10.1016/j.arth.2013.04.029>.
- Collopy, D. et al. A Prospective Roentgen Stereophotogrammetric Analysis (RSA) Study of the Stryker Accolade II Cementless Hipo Stem to 2 Years Post-Surgery. Abstracts from 2016 Annual AOA Meeting—Cairns, QLD, October 9-13, 2016.
- Bono, J. et al. Total Hip Intraoperative Femur Fracture: Does a Tapered Wedge Design with a Size Specific Medial Curvature Reduce the Incidence of Intraoperative Femur Fracture? Presentation at 47th Annual EOA Meeting -New Orleans, LA, October 19-22, 2016.
- Pierce T. et al. Second-Generation Versus First-Generation Cementless Tapered Wedge Femoral Stems. Orthopedics. 2015 38(9): 550-554.
- Kolisek, F. et al. Outcomes of Second-Generation Tapered Wedge Design. SurgTechnolInt. 2016 Apr 4;XXVIII. pii: sti28/702.
- Stryker Internal presentation 'Global Bone Morphology Study -Proximal Femur'. November 10, 2016.
- Narzikul, BS. et al. Alteration in Geometry of Femoral Stem Results in Better Fit and Fill: Comparison of Accolade I vs. Accolade II. Presentation at 47th Annual EOA Meeting -New Orleans, LA, October 19-22, 2016.
- Wuestemann T, Bastian A, Parvizi J, Nessler J, Kolisek F. A novel tapered hip stem design optimized for femoral fit in a wide array of bone types. 2011 EFFORT.
- Boucher, F. et al. Preclinical Finite Element Analysis of a Novel Tapered Wedge Stem Optimized for Stability at a Reduced Length. 2011 EFFORT Poster Presentation.
- Lovell T, Hozack W, Kreuzer S, Merritt P, Nogler M, Puri L, Wuestemann T, Bastian A. Influence of Stem Length on the Insertion Path in THR. 2009 ORS Poster #2040.
- Khanuja H, et al. Current Concepts Review: Cementless Femoral Fixation in Total Hip Arthroplasty. J Bone Joint Surg Am. 2011;93:500-9.
- Faizan A, et al. Development and Verification of a Cementless Novel Tapered Wedge Stem for Total Hip Arthroplasty, J Arthroplasty (2014), <http://dx.doi.org/10.1016/j.arth.2014.09.023>.
- Restrepo C, Parvizi J, Eslam Pour A, et al. Prospective Randomized Study of Two Surgical Approaches for Total Hip Arthroplasty. The Journal of Arthroplasty 2010; 25: 671-670.
- Vail T, Mariani E, Bourne M, Berger R, Meneghini M. Approaches in Primary Total Hip Arthroplasty. J Bone Joint Surgery 2009; 91:10.
- M.W. Zawadsky et al./ Early Outcome Comparison Between the Direct Anterior Approach and the Mini-Incision Posterior Approach for Primary Total Hip Arthroplasty: 150 Consecutive Cases, The Journal of Arthroplasty 29 (2014) 1256-126.
- Internal Sales Data (2017).

*SOMA-design of Accolade II based on 556 CT scans.

stryker

325 Corporate Drive
Mahwah, NJ 07430
t: 201.831.5000

www.stryker.com

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker's product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your sales representative if you have questions about the availability of any of Stryker's products in your area. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Accolade, SOMA, Stryker. All other trademarks are trademarks of their respective owners or holders.