

# Mako Total Hip arthroplasty: clinical summary



Mako clinical evidence

#### 1. Introduction

Total hip arthroplasty (THA) has been one of the most successful procedures within the field of orthopaedics since the late 1960s.1 The short- and long-term outcomes of THA may be influenced by several factors, including patient demographics, surgical technique and implant features.<sup>2</sup> One of the most important surgeon-controlled factors is component positioning.2 Component malposition has been linked to poor biomechanics, accelerated wear, leg length discrepancy (LLD), and revision surgeries.<sup>2</sup> In addition, component malposition is directly associated with dislocations and mechanical loosening, which account for approximately 40% of THA revisions.3 The Mako System was introduced with a goal of providing more accurate implant positioning and alignment to plan, to help restore anatomy and enhance patient outcomes. This document summarizes the evidence to date that supports the use of Robotic-Arm Assisted Surgery for total hip arthroplasty.

# 2. What is the evidence that Mako Total Hip works?


Successful clinical outcomes following total joint replacement are dependent on component placement and on restoring the natural joint anatomy of the hip.<sup>2</sup> Instability, early mechanical failures and dislocation in hip arthroplasty continue to be primary reasons for revision.<sup>2</sup> The Mako System is designed to minimize the margin of error associated with component placement, and to enhance the accuracy and reproducibility of THA.

#### 2.1 Accuracy and reproducibility in THA

In a multicenter clinical trial including 110 patients, acetabular cup position was compared between preoperative plan, assessment, and achieved radiographic measure.<sup>3</sup> Results confirmed that intraoperative robotic-arm assistance achieved greater accuracy in preparation and position of the acetabular cup during THA (Table 1).<sup>4</sup>

	Preoperative plan	Intraoperative robotic-arm measurements	Martell radiographic measurement
Inclination	$40.0^{\circ} \pm 1.2^{\circ}$	39.9°±2.0°	40.0°±4.1°
Version	18.7°±3.1°	18.6°±3.9°	21.5°±6.1°
Count (n)	119	119	110

Table 1. The average inclination and anteversion values of the acetabular components in the study, showing the preoperative plan, measures recorded intraoperatively and those measured from plan radiographs using the Martell method.<sup>2</sup>




Figures 1a and 1b. Scatterplots of the (a) robotic-assisted and (b) conventional cups in the safe zones of Lewinnek et al. and Callanan et al. are shown.<sup>6</sup>

Domb et al. (2015) conducted a study involving six surgeons at a single institute, in which 1,980 THA surgeries were evaluated.<sup>5</sup> The aim of this study was to understand the influence of surgical approaches and modes of guidance on accuracy of acetabular component placement.<sup>5</sup> Robotic-arm assisted surgery resulted in a significantly greater percentage of components placed in Callanan's safe zones (30°-45° inclination and 5°-25° version) than all other modalities, including navigationand fluoroscopy-guided approaches (p<0.05).<sup>5</sup> This study highlighted the consistency of the robotic-arm assisted technology, based on a large patient series.<sup>5</sup>

In another clinical study that compared robotic-arm assisted THA to manual THA, 100% of robotic-arm assisted THAs were within the Lewinnek safe zone  $(30^{\circ}-45^{\circ}$  inclination and  $5^{\circ}-25^{\circ}$  version), compared with 80% of the conventional THAs (p=0.001).<sup>6</sup> A total of 92% of robotic-arm assisted THAs were in Callanan's modified safe zone, compared with 62% of conventional THAs (p=0.001).<sup>6</sup> Use of the Mako System allowed for more consistent placement of the cup in both safe zones (Figure 1a-b).<sup>6</sup>

Clinical evidence continues to build on the potential benefits of robotic-arm assisted THA. Investigations have demonstrated robotic-arm assisted surgery is accurate to  $1.0\pm0.7$ mm for leg length/offset. Compared to manual THA, robotic-arm assisted THA was five times more accurate to plan in cup inclination and 3.4 times more accurate to plan in cup anteversion. A recent publication highlighted the influence of head center of rotation (COR) on the risk of hip dislocation. A potential benefit of robotic-arm assisted THA is that it has been shown to be significantly more accurate in reproducing COR when compared to manual implantation, which may result in reduced incidence of hip dislocation.

#### **Results: Bone stock**



Robotic-arm assisted THA

of the femoral head

c-f = bone thickness lost over course of surgery (c-f)/f = bone thickness lost through surgery per width

 $(c^3-f^3)/f^3$  = volume of bone lost through surgery

Figure 2. Illustrates the Mako System's single reaming technique preserves bone as compared to conventional THA's sequential reaming technique.9

The previous studies highlighted all included results from experienced surgeons. In a publication by Smith et al., they investigated if the accuracy of robotic-arm assisted surgery was translatable to newer surgeons in their fellowship.<sup>9</sup> In this cadaver study, two adult reconstruction fellows halfway through their training year performed manual THA (n=6) and robotic-arm assisted THA (n=6). The robotically prepared hips demonstrated statistically significant greater accuracy and precision to plan compared to the manually prepared hips when considering shell version, shell inclination and LLD. Error in shell placement was reduced by up to 9mm and error in LLD was reduced by up to 8mm when using robotic-arm assisted THA. The authors concluded that these findings suggest that CT-based preoperative planning and intraoperative robotic technology, such as the Mako Total Hip application, can allow less experienced surgeons to place implant components more consistently in the desired orientation, with comparable accuracy to what has been reported by experienced surgeons.9

In addition to providing accuracy to plan for less experienced surgeons, Shapira et al. studied the use of the Mako Total Hip system as a learning tool for fellows training in hip arthroplasty. The study evaluated the accuracy of fellows' estimation of cup and broach positioning using the Mako system. They found that the mean difference between estimated and actual cup inclination and version was 7.24° (P = 0.060) and 4.81° (P = 0.031), respectively. The mean difference in broach version was 7.00° (P = 0.159). Shapria et al. concluded that the robotic system is a useful learning tool for fellows in training to help them understand their own inaccuracies in estimating implant position and hence may help refine their abilities.

Robotic-arm assisted THA has also been associated with more precise reaming, which can not only influence recreation of COR, but also impact preservation of bone stock.<sup>11</sup>

Suarez-Ahedo et al. (2017) studied bone preservation during primary THA and performed a matched pair control study which demonstrated that when compared to conventional THA (n=57), robotic-arm assisted THA (n=57) allowed for more precise reaming. This led to the use of smaller acetabular cups in relation to the patient's femoral head size. Using acetabular cup size relative to femoral head size as a surrogate measure of acetabular bone resection, these results suggested greater preservation of bone stock using robotic-arm assisted THA compared to conventional THA. This may reflect increased translational precision during the reaming process Figure 2).

The potential benefits of using CT-based robotic technology such as Mako to assess the influence of native femoral version on final stem version (SV) and combined anteversion when using a straight, uncemented stem, was




Figure 3. Scatter graph of SV in respect to FNV. The stem "safe zone" was highlighted in green. When FNV was  $<5^{\circ}$ , stem version was "increased" 3% of the time, "normal" 37% of the time, and "reduced" 60% of the time, meaning that the surgeon was not always able to correct femoral retroversion. Also with a "normal" FNV, the stem was positioned with a SV  $<5^{\circ}$  34% of the time.

researched by Marcovigi et al. <sup>12</sup> A total number of 362 patients who underwent Mako Total Hip were enrolled from three different orthopaedic centers. <sup>12</sup> All patients underwent CT planning with measurement of femoral neck version (FNV) and intraoperative measurement of SV, acetabular component version (AV), and combined version (CV), with robotic instrumentation. <sup>12</sup> Results showed that the mean FNV was  $5.0^{\circ} \pm 9.6^{\circ}$ , and SV was  $6.4^{\circ} \pm 9.7^{\circ}$ . <sup>12</sup> A strong correlation was found between SV and CV (R = 0.89, P < .001) and a significant difference in SV was found between the three centers (P < .001). CV was <25° in 109 patients (30.1%) with relative risk of CV < 25° being 8.6 times greater with SV < 5° (P < .001) (Figure 3). <sup>12</sup>

From this data, it is important to note that when using an uncemented, single-wedge, straight stem, SV is highly variable. Thus, the greater variability of FNV in patients with osteoarthritis is confirmed. The Despite being moderately correlated with native FNV, SV can be partially influenced by the surgeon. The authors concluded that knowledge of preoperative and intraoperative stem version is fundamental to avoid abnormal combined version and therefore reduce risk of impingement, dislocation or acetabular uncoverage. They also emphasized that CT-based planning and robotic technology may be useful tools to have in the operating room, combined with stem designs which facilitate the achievement of desired version angles.






Figure 4. Mako Total Hip has a virtual range-of-motion tool indicating risk of impingement (highlighted in red). The surgeon is able to change implant position and/or implant systems to address impingement.

## 2.2 Functional implant planning

Surgeons have been using the Lewinnek safe zone as a guide for cup placement for over 40 years. <sup>13</sup> However, studies have shown that greater than 50% of total hip arthroplasty dislocators have cups placed within this safe zone. <sup>14,15</sup> One shortcoming of the Lewinnek safe zone is that it generically applies to all patients, regardless of their individual bone morphology, kinematics, implant choices or placement. Mako Total Hip has integration of features that allow a surgeon to assess a patient's functional pelvic tilt and virtual range-ofmotion (vROM) to help achieve functional implant planning (Figure 4).

O'Conner et al. evaluated the dynamic ranges of sagittal pelvic motion in 100 patients undergoing total hip arthroplasty. They found that pelvic tilt (PT) can significantly change the functional orientation of the acetabular component and may differ markedly between patients (Figure 5). 16 They then used the vROM tool from the Mako Total Hip application to investigate whether there is an ideal functional combined anteversion to help reduce risk of impingement for the patient. The authors found that the vast majority of THAs planned with standing combined anteversion between 30 to 50° and sitting combined anteversion between 45 to 65° (Figure 5) had minimal impingement. Those planned outside that combined anteversion window were more likely to have impingement. The authors concluded functional combined anteversion, which considers both cup and stem position, should be used when identifying an ideal position for impingement-free ROM.<sup>16</sup>






Figure 5. O'Conner et al. found high variability in preoperative lying, standing, and sitting pelvic tilt (PT) in patients undergoing total hip arthroplasty (a). The Mako Total Hip virtual ROM tool was used on these patients with varying PT to find an optimal zone for standing and sitting combined anteversion (b). A standing combined anteversion between 30 to 50° and a sitting combined anteversion between 45 to 65° were defined as the optimal zone where less impingement was detected when using the virtual ROM tool. 65°

Scholl et al. have reported on early outcomes out to 6-months and complications for a single surgeon's first 100 patients who underwent a Mako Total Hip procedure with functional implant planning.<sup>17</sup> At sixweeks postoperative, patients achieved a statistically significant improvement (p<0.05) in their HOOS JR scores when compared to their preoperative scores. This indicated that Mako Total Hip 4.0 provided a substantial clinical benefit to the patients. At three- and six-months follow-up, patients continued to have statistically significant improvement in HOOS JR when compared to both preoperative and six-week scores (p<0.05).<sup>17</sup>

# 2.3 Surgical team learning curve

In a retrospective, single-surgeon review of 100 consecutive Mako Total Hips, Bukowski et al. studied the effects of learning curve on the outcome of three groups of patients: 1) the surgeon's first 100 manual THA cases (2000-2001); 2) the surgeon's last 100 manual THA cases (2010-2011); and 3) the surgeon's first 100 Mako Total Hip cases (2011-2012). <sup>18,19</sup> Dislocation was more frequent in group one (5/100, 5%) and group two (3/100, 3%) than in group three (0/100, 0%) (p<0.05) at the one year follow-up interval. <sup>19</sup>

Similarly, Redmond et al. researched the learning curve during the adoption of robotic-arm assisted THA as measured by component position, operative time, and complications.<sup>20</sup>

The first 105 robotic-arm assisted THAs performed by a single surgeon were divided into three groups based on the order of surgery: 1) Group A consisted of the first 35 patients who underwent Mako Total Hip by the senior surgeon, 2) Group B consisted of patients 36–70; and 3) Group C consisted of patients  $71-105.^{20}$  The authors reported a decreased risk of acetabular component malpositioning with Mako experience (P < 0.05). Operative time appeared to decrease with increasing surgical experience with the Mako System (P < 0.05). A learning curve of 35 cases was observed, as a decreased incidence of acetabular component outliers and decreased operative time were noted with increased surgical experience with Mako. Operative time were noted with increased surgical experience with Mako.

Heng et al. carried out a retrospective comparison of a single surgeon's last 45 conventional THAs performed prior to changing to the robotic-arm assisted system, and compared them with the first 45 robotic-arm assisted THAs.<sup>21</sup> When comparing surgical times between the two groups, they found that the average surgical time was 96.7 minutes for the robotic-arm assisted group and 84.9 minutes for conventional group.<sup>21</sup> Upon further analysis, the authors determined that each robotic-arm assisted operation was approximately one minute shorter than the previous robotic operation and the average time for the last ten cases was reduced to 82.9 minutes, which was quicker than the average time of the conventional group.<sup>21</sup> It was concluded that surgical time

is comparable with conventional techniques after the initial learning curve of approximately 35 cases.<sup>21</sup>

Kong et al. published a retrospective comparative cohort study of an experienced manual surgeon's first 100 robotic-assisted THAs compared to the surgeon's last 100 manual cases.<sup>22</sup> The average operating time of roboticassisted THA was 95.92  $\pm$  15.64 minutes, ranging from 68 to 145 minutes.<sup>22</sup> The learning curve was assessed using a cumulative summation test for learning curve analysis which demonstrated that after the 14th case, a downtrend in the surgeon's operative time began.<sup>22</sup> There was no statistical difference between the first 14 cases versus cases 15 to 100 when considering cup positioning, postoperative LLD, offset and Harris hip score (HHS). Results indicate that there was a 14-case learning curve when considering operative time; however, the authors observed this learning curve did not impact patient outcomes.22

#### 2.4 Impact on surgical team

The Mako system provides a stereotactic boundary that guides the alignment of the robotic arm during acetabular reaming and cup insertion, alleviating the need for the surgeon to ensure proper alignment. Additionally, the system provides a single ream option, eliminating the need for the surgeon to perform multiple reams to achieve final ream size. It has been reported that 66.1% of arthroplasty surgeons have had a workplace-related injury, with 31% requiring surgery. Assistance in performing reaming and cup insertion may enhance the ergonomic health and reduce the workload demand on the surgeon. <sup>23</sup>

A cadaveric study was performed to determine how the use of Mako Total Hip can influence a surgeon's energy expenditure as well as mental and physical demand compared to manual THA.<sup>24,25</sup> Twelve THAs were performed by two surgeons, with varying robotic experience, in their fellowship training. Each cadaveric specimen received a manual THA on one hip and a Mako Total Hip on the contralateral side. The surgeons wore a biometric shirts that collected data on caloric expenditure<sup>24</sup> and were administered a modified surgery task load index questionnaire after each procedure to evaluate perceived mental and physical demand.25 Surgeons who performed Mako Total Hip demonstrated reduced caloric expenditure during acetabular reaming and acetabular implantation.<sup>24</sup> With the Mako system assisting through use of a stereotactic boundary, the surgeons also reported less physical and mental demand during acetabular reaming and acetabular implantation, with statistically significantly less mental demand during acetabular reaming in the Mako Total Hip group.<sup>25</sup>

# 3. What are the clinical benefits of Mako Total Hip?

Clinical benefits resulting from increased accuracy and precision afforded by Mako Total Hip have been investigated, including functional outcomes and levels of patient satisfaction. Results of studies in this area are promising.




Figure 6. Statistically higher modified HHS were shown for Mako Total Hip patients. $^{11}$ 

#### 3.1 Clinical and functional outcomes in THA

In the research conducted by Bukowski et al., outcomes for three groups of 100 consecutive THAs (first 100 manual THAs; last 100 manual THAs; and first 100 Mako Total Hips), were reviewed. Mako Total Hip resulted in significantly higher modified HHSs (92.1  $\pm$  10.5 vs. 86.1  $\pm$  16.2, p = 0.002) and UCLA activity levels (6.3  $\pm$  1.8 vs. 5.8  $\pm$  1.7, p = 0.033) than manual THA at minimum one-year follow-up (Figure 6 and 7, Table 2).  $^{18}$ 




Figure 7. Statistically higher UCLA scores were shown for Mako Total Hip patients. $^{11}$ 

	Group (RATHA n=100, MTHA n=100)	Preoperative	Postoperative	PROMs (postoperative- preoperative)	p-value
mHHS (mean and standard deviation)	RATHA	49.6 (16.3)	92.1 (10.5)	43.0 (18.8)	< 0.001
	MTHA	49.2 (14.8)	86.1 (16.2)	37.4 (18.3)	< 0.001
	p-value	0.865	0.002	0.035	
SF12-MCS (mean and standard deviation)	RATHA	54.1 (10.4)	54.6 (9.1)	0.4 (9.7)	0.629
	MTHA	53.1 (9.6)	53.0 (10.2)	0.5 (11.5)	0.970
	p-value	0.459	0.245	0.962	
SF12-PCS (mean and standard deviation)	RATHA	33.5 (9.6)	46.0 10.5)	12.5 (11.8)	< 0.001
	MTHA	30.3 (8.0)	44.4 (11.0)	14.0 (11.9)	< 0.001
	p-value	0.010	0.282	0.404	
WOMAC (mean and standard deviation)	RATHA	45.6 (18.9)	16.0 (14.9)	-29.6 (21.4)	< 0.001
	MTHA	47.1 (14.7)	17.3 (15.5)	-28.5 (18.3)	< 0.001
sidiladia devidiloli,	p-value	0.536	0.538	0.618	
UCLA (mean and standard deviation)	RATHA	5.1 (1.9)	6.3 (1.8)	1.2 (1.7)	< 0.001
	MTHA	4.8 (1.8)	5.8 (1.7)	1.0 (1.9)	< 0.001
	p-value	0.227	0.033	0.429	
Categorical analysis of	modified Harris H	ip Score			
	rTHA		mTHA		
90-100	75.0% (75)		61.0% (61)		0.034
80-89	13.0% (13)		15.0% (15)		0.684

5.0% (5)

19.0% (19)

70-79

< 70

6.0% (6)

6.0% (6)

0.756

0.005

Perets et al. have compared minimum two-year outcomes and complications for patients who underwent a Mako Total Hip procedure to patients who underwent manual THA. <sup>26</sup> Eighty-five Mako Total Hip patients were pair-matched to manual controls based on patient demographics. The patients prospectively reported on Harris hip score (HHS), Forgotten Joint Score-12 (FJS-12), Visual analog scale for pain (VAS), and satisfaction (scale ranging 0-10). <sup>26</sup>

The FJS-12 questionnaire has evidence of low-ceiling effects and is suitable for assessing longer term outcomes in well-performing groups after THA.<sup>27</sup> The literature has reported an FJS-12 ranging from 50.9  $\pm$  25.3 to 80  $\pm$  24 for patients who received manual THA.<sup>26,27</sup> Perets et al. reported an FJS-12 for the Mako Total Hip patients of 80.2, which was significantly greater (p=0.003) than a FJS-12 of 68.6 reported by the manual THA control group.<sup>26</sup> The Mako group had a significantly higher HSS (p<0.001) and a trend towards having lower VAS scores (p=0.120) when compared to the control group. Additionally, at two years, the Mako group had less LLD (1.0mm, p=0.013), less global offset (0.9mm, p=0.31) and no dislocations reported in either group.<sup>26</sup> Postoperatively, one robotic-arm assisted patient required a revision at 8.7 months after primary THA due to femoral stem loosening and three manual THA patients required revision at a mean 25.1 months postoperatively, all for femoral stem loosening.<sup>26</sup>

This same group of patients continued to be followed, and Maldonado et al. published on minimum five-year outcomes of this patient cohort.<sup>28</sup> When compared to a manual THA control group, the Mako Total Hip cases reported significantly higher Harris hip scores (p<0.001), FJS-12 (p=0.002), Veterans RAND-12 physical component scores (p=0.002), and Short Form Health Questionnaire (SF)-12 physical component scores (p=0.001) (Table 3).<sup>28</sup> While revision rates between these cohorts were similar (p=0.479), the acetabular components for the Mako Total Hip cases were more consistently placed within the Lewinnek (p=0.002) and Callanan (p=0.001) safe zones.<sup>28</sup> In addition, Mako Total Hip recipients had lesser absolute values of leg length discrepancy and global offset (p = 0.091, p = 0.001). This study used multiple validated functional hip outcome scores to determine that patients who received Mako Total Hip reported favorable outcomes at a minimum five-year follow-up.28

A similar trend was observed in a prospective study of 40 patients undergoing Mako Total Hip that were propensity matched to 80 patients undergoing manual THA.<sup>29</sup> Groups were matched based on age, sex, and preoperative function. At 12-month follow-up, the Mako Total Hip group had improved Oxford Hip Score (OHS, p=0.038), FJS (p<0.001), and less dissatisfied patients (Mako 0 vs. Manual 6). The FJS in the Mako group was 21.2 points higher than the manual group which

represented a minimal clinically important difference. Based on radiographic analysis, the manual THA group had a decrease in the horizontal centre of rotation, which was also associated with a decrease in acetabular offset. This shift of centre of rotation was not found in the Mako group. The authors hypothesised that this difference in restoration of hip centre and leg length could have impacted the differences in clinical outcomes between the two groups.<sup>29</sup>

#### 3.2 Patient satisfaction

THA has been one of the most successful surgeries in medicine, having demonstrated favorable shortand long-term outcomes and resulting in more than 95% survivorship at ten years. In addition, patient satisfaction post-THA is high, as demonstrated by Perets et al., where patient satisfaction at a minimum of two-year follow-up was assessed. For the 162 Mako Total Hip cases considered in this study, mean patient satisfaction was a high 9.3 out of 10.26

## 3.3 Patient recovery

When exploring a patient's road to recovery, their length of stay in hospital after surgery is a key factor to consider. Heng et al. retrospectively compared the length of stay of 45 patients who underwent Mako Total Hip against those who received conventional THA (n=45).<sup>21</sup> They reported similar results in both groups, however once the patients who required inpatient rehabilitation were excluded, the robotic group had a shorter hospital stay (4.22 days vs. 5.93).<sup>21</sup>

This finding was further validated by another study conducted by Banchetti et al., who retrospectively

Patient- reported outcomes	Robotic- assisted THA	Manual THA	p-value
HHS	90.57±13.46	84.62±14.45	<0.001
FJS-12	82.69±21.53	70.61±26.74	0.002
VAS	1.27±2.20	1.07±1.87	0.45
Satisfaction	8.91±2.00	8.52±2.62	0.35
VR-12 mental	60.76±5.94	58.97±6.03	0.17
VR-12 physical	50.30±8.83	45.92±9.44	0.002
SF-12 mental	56.59±5.60	56.20±6.62	0.81
SF-12 physical	48.97±9.21	44.01±10.26	0.001

Table 3. Minimum five-year patient-reported outcomes for a Mako Total Hip and manual THA cohort?<sup>2</sup>

analyzed 107 patients at 24-months follow-up (Mako Total Hip, n=56; standard technique THA, n=51). They found a significant difference in the length of hospital stay, defined by number of days hospitalized, between the Mako group (M=5.14. SD=1.98) and the standard group (M=8.11, SD=1.64).

Overall, early data from these studies suggests that patients who undergo Mako Total Hip may be able, on average, to return home sooner after surgery than those who undergo conventional THA. This may pose a great advantage for the patients' well-being and offer financial benefits to healthcare institutions, since a reduction in length of hospitalization has the potential to reduce economic burden to hospitals.<sup>30</sup> Furthermore, these findings have the potential to offer financial benefits to healthcare institutions since a reduction in the length of stay post–Mako Total Hip surgery potentially reduces the economic burden to hospitals. This is a key area being investigated by various surgeons worldwide.

Rosinsky et al. performed an analysis focused specifically on comparing patients who underwent robotic-arm assisted THA either at an inpatient or outpatient facility.31 The first 100 consecutive patients who underwent outpatient THA were matched to 100 patient who underwent inpatient THA during the same time period. They compared perioperative variables including surgical time, blood loss, and length of stay as well as 90-day complication rates and 2-year patient reported outcomes. The outpatient group had an average length of stay of 6.8 hours compared to 43.2 hours for the inpatient group (P < 0.001). There were no significant differences between the groups regarding readmissions, emergency room visits, and unplanned clinic visits. Complications and revision rates were similar in both groups. The group also concluded that in appropriately selected, younger patients, outpatient robotic-arm assisted THA can achieve improved postoperative 2-year PROs compared to inpatient manual THA.31

# 3.4 Role of Mako in complex cases

Chai et al. carried out a case study that included three complex cases with hip dysplasia, ankylosing spondylolysis and post-traumatic arthritis, respectively.<sup>32</sup> In all three cases, the Mako System was utilized to help accurately implement the surgical plan.<sup>32</sup> Since there was an absence of conventional bony landmarks, the preoperative CT scan in these cases was instrumental in planning.<sup>32</sup> The hip dysplasia patient reported at three months postoperatively that they were able to walk without assistance, had no hip pain and were satisfied with their leg lengths.<sup>32</sup> The patient with ankylosing spondylolysis reported no hip pain and was able to walk with a walking frame at three months postoperatively.<sup>32</sup> The patient with post-traumatic arthritis reported no hip pain and was able to walk without assistance at three months postoperatively.<sup>32</sup> According to this study, the planning and accuracy of execution in Mako Total Hip allowed the surgeon to give the patients excellent reconstruction of their hip joints which substantially enhanced their quality of life. The authors went on to say that Mako Total Hip surgery may be considered for complex THA cases in order to achieve the desired accuracy of the reconstruction, especially in the absence of conventional bony landmarks.<sup>32</sup>

Mako Total Hip has been shown to be a useful tool for complex cases. Kuroda et al. analysed a consecutive series of 69 patients who underwent robotic-arm assisted THA, where 30 of those patients had development dysplasia of the hip (DDH) and were classified according to the Crowe type. 33 Using the patients' preoperative plan and comparing to postoperative CT data, accuracy of cup alignment and 3D placement were compared between DDH and non-DDH patients. No significant differences were found in cup placement between the two groups and excellent restoration of leg length and combined offset were achieved in both groups. This study concluded that robotic-arm assisted THA may achieve accuracy and reproducibility of cup placement in both non-DDH and DDH patients, even those with severe DDH.<sup>33</sup>

# 4. Is Mako Total Hip cost effective?

In assessing the potential effects of Mako Total Hip on costs to U.S.-based private payers and Medicare, Maldonado and colleagues evaluated the long-term cost effectiveness of robotic-arm assisted THA (RAATHA) vs. manual THA (mTHA) through a Markov model. 34 The potential outcomes of THA were categorized into the transition states: infection, dislocation, no major complications, or revision. Cumulative costs and utilities were assessed using a cycle length of one year over a time horizon of five years. They found that RAA THA cohort was cost effective relative to mTHA cohort for cumulative Medicare and cumulative private payer insurance costs over the 5-year period. RAA THA cost saving had an average differential of \$945 for Medicare and \$1,810 for private insurance relative to mTHA while generating slightly more utility (0.04 quality-adjusted life year). The preferred treatment was sensitive to the utilities generated by successful RAA THA and mTHA. Microsimulations indicated that RAA THA was cost effective in 99.4% of cases. In the U.S. Medicare and private paver scenarios. RAATHA was shown to be more cost effective than conventional mTHA when considering direct medical costs from a U.S. payer's perspective.<sup>34</sup>

A separate U.S.-based Medicare analysis of the 90-day episode of care (EOC) of 938 RATHAs propensity matched to 4,670 MTHAs found that RATHA patients were less likely to have post-index inpatient rehabilitation (IPR) or skilled nursing facility (SNF) admissions (0.64% vs. 2.68%; p<0.0001 and 20.79% vs. 24.99%; p=0.0041, respectively).35 RATHA patients used fewer days in postindex inpatient and SNF care (7.15 vs. 7.91; p=0.8029 and 17.98 vs. 19.64; p=0.5080, respectively) and used fewer home health aide (HHA) visits, (14.06 vs. 15.00; p=0.0006) compared to MTHA. RATHA had lower 90-day EOC costs for: IPR (\$11,490 vs. \$14,674; p=0.0470), SNF (\$9,184 vs. \$10,408, p=0.0598) and HHA (\$3,352 vs.\$3,496; p=0.0133) compared to MTHA.<sup>35</sup> Overall, RATHA patients had 12% (\$948) lower average post-index costs compared to MTHA patients (p=0.0004).26 Total 90-day episode-of-care costs for RATHA patients were found to be \$785 less than those of MTHA patients (\$19,734 vs. \$20,519, p=0.0095).35

#### 5. Conclusions

Mako Total Hip offers the potential for surgeons to achieve component placement and alignment accuracy, as well as to enhance clinical outcomes. 4-8,11-12,18-21,26,27,30 Patients have reported tangible benefits of Mako Robotic-Arm Assisted Surgery, including treatment satisfaction and return to activities of daily living. 18,28 Surgeons are empowered to achieve their target preoperative plans with precision. Ultimately, the benefits of Mako Robotic-Arm Assisted Total Hip arthroplasty may be experienced by all key players – patients, surgeons and health systems.

#### References

- Knight SR, Aujla R, Biswas SP. Total Hip Arthroplasty over 100 years of operative history. Orthop Rev (Pavia). 2011;3(2):e16. doi:10.4081/or.2011.e16
- Callanan MC, Jarrett B, Bragdon CR, et al. The John Charnley Award: risk factors for cup malpositioning: quality improvement through a joint registry at a tertiary hospital. Clin Orthop. 2011;469(2):319-329. doi:10.1007/s11999-010-1487-1
- Bozic KJ, Kurtz SM, Lau E, Ong K, Vail TP, Berry DJ. The epidemiology of revision total hip arthroplasty in the United States. J Bone Joint Surg Am. 2009;91(1):128-133. doi:10.2106/JBJS.H.00155
- Elson L, Dounchis J, Illgen R, et al. Precision of acetabular cup placement in robotic integrated total hip arthroplasty. Hip Int. 2015;25(6):531-536. doi:10.5301/hipint.5000289
- Domb BG, Redmond JM, Louis SS, et al. Accuracy of component positioning in 1980 total hip arthroplasties: a comparative analysis by surgical technique and mode of guidance. J Arthroplasty. 2015;30(12):2208-2218. doi:10.1016/j. arth.2015.06.059
- Domb BG, El Bitar YF, Sadik AY, Stake CE, Botser IB. Comparison of roboticassisted and conventional acetabular cup placement in THA: a matched-pair controlled study. Clin Orthop Relat Res. 2014;472(1):329-336. doi:10.1007/ s11999-013-3253-7
- Nawabi DH, Conditt MA, Ranawat AS, et al. Haptically guided robotic technology in total hip arthroplasty: a cadaveric investigation. Proc Inst Mech Eng H. 2013;227(3):302-309. doi:10.1177/0954411912468540
- Jauregui JJ, Banerjee S, Elmallah RK, et al. Radiographic evaluation of hip dislocations necessitating revision total hip arthroplasty. Orthopedics. 2016;39(5):e1011-e1018. doi:10.3928/01477447-20160616-02
- Smith R, Borukhov I, Hampp E et al. Comparison of precision for manual versus robotic-assisted total hip arthroplasty performed by fellows. J of Hip Surgery 2020;4(03):117-123.
- Shapira J, Diulus SC, Rosinsky PJ, Maldonado DR, Lall AC, Domb BG. Robotics and navigation as learning tools for fellows training in hip arthroplasty. J Am Acad Orthop Surg. 2021 Feb 15;29(4):176-181. doi: 10.5435/JAAOS-D-20-00357.
- 11. Suarez-Ahedo C, Gui C, Martin TJ, Chandrasekaran S, Lodhia P, Domb BG. Robotic-arm assisted total hip arthroplasty results in smaller acetabular cup size in relation to the femoral head size: a matched-pair controlled study. Hip Int. 2017;27(2):147-152. doi:10.5301/hipint.5000418
- Marcovigi A, Ciampalini L, Perazzini P, Caldora P, Grandi G, Catani F. Evaluation of Native Femoral Neck Version and Final Stem Version Variability in Patients With Osteoarthritis Undergoing Robotically Implanted Total Hip Arthroplasty. J Arthroplasty. 2019;34(1):108-115. doi:10.1016/j. arth.2018.06.027
- Lewinnek et al. Dislocations after total hip replacement arthroplasties. J Bone Joint Surg Am 1978;60:217–220.
- 14. Callanan, M.C., Jarrett, B., Bragdon, C.R. et al. The John Charnley Award: Risk Factors for Cup Malpositioning: Quality Improvement Through a Joint Registry at a Tertiary Hospital. Clin Orthop Relat Res 469, 319–329 (2011).
- 15. Abdel MP, von Roth P, Jennings MT, Hanssen AD, Pagnano MW. What Safe Zone? The Vast Majority of Dislocated THAs Are Within the Lewinnek Safe Zone for Acetabular Component Position. Clin Orthop Relat Res. 2016.
- O'Conner P, Thompson M, Esposito C, Poli N, McGree J, Donnelly T, Donnelly W. Change in pelvic position in total hip arthroplasty patients through functional range of motion. Bone Jt Open 2021;2-10:834-841.
- Scholl LY, Taylor KB, Erwin DJ, Collins S, Marchand KB, Marchand RC. Early
  outcomes in robotic assisted total hip arthroplasty patients with functional
  hip positioning. ORS Annual Meeting, Feb 12-16, 2021. Virtual conference.
  PS-20.
- Bukowski BR, Anderson P, Khlopas A, Chughtai M, Mont MA, Illgen RL 2nd. Improved Functional Outcomes with Robotic Compared with Manual Total Hip Arthroplasty. Surg Technol Int. 2016;29:303-308.
- 19. Illgen R. Robotic Assisted THA: Outcomes after primary THA manual compared with robotic assisted presentation.
- Redmond JM, Gupta A, Hammarstedt JE, Petrakos AE, Finch NA, Domb BG.
   The learning curve associated with robotic-assisted total hip arthroplasty. J Arthroplasty. 2015;30(1):50-54. doi:10.1016/j.arth.2014.08.003

- Heng YY, Gunaratne R, Ironside C, Taheri A. Conventional vs robotic-arm assisted total hip arthroplasty (THA) surgical time, transfusion rates, length of stay, complications and learning curve. J Arthritis. 2018;7(4):1-4. doi:10.4172/2167-7921.1000272
- 22. Kong X, Yang M, Jerabek S, Zhang G, Chen J, Chai W. A retrospective study comparing a single surgeon's experience on manual versus robot-assisted total hip arthroplasty after the learning curve of the latter procedure A cohort study. Int J Surg. 2020;77:174-180. doi:10.1016/j.ijsu.2020.03.067
- Alqahtani SM, Alzahrani MM, Tanzer M. Adult Reconstructive Surgery: a high-risk profession for work-related injuries. J Arthroplasty. 2016;31(6):1194-1198. doi:10.1016/j.arth.2015.12.025
- 24. Abbruzzese K, Byrd Z, Smith R, et al. Assessment of surgeon biometrics during manual and robotic THA. Presented at: International Society for Technology in Arthroplasty (ISTA) meeting, 32nd Annual Congress; October 2-5, 2019; Toronto, Canada.
- Valentino A, Scholl L, Hampp EL, Smith R, Byrd Z, Mont M. Physical and mental demand during total hip arthroplasty. Presented at: Orthopaedic Research Society (ORS) Annual Meeting; February 8-11, 2020; Phoenix, AZ.
- Perets I, Walsh JP, Mu BH, Mansor Y, Rosinsky PJ, Maldonado DR, Lall AG, Domb BG. Short-term clinical outcomes of robotic-arm assisted total hip arthroplasty: a pair-matched controlled study. Orthopedics. Mar-Apr 2021;44(2):e236-e242. doi: 10.3928/01477447-20201119-10.
- 27. Hamilton DF, Loth FL, Giesinger JM, et al. Validation of the English language Forgotten Joint Score-12 as an outcome measure for total hip and knee arthroplasty in a British population. Bone Joint J. 2017;99-B(2):218-224. doi:10.1302/0301-620X.99B2.BJJ-2016-0606.R1
- Domb BG, Chen JW, Lall AC, Perets I, Maldonado DR. Minimum 5-year outcomes of robotic-assisted primary total hip arthroplasty with a nested comparison against manual primary total hip arthroplasty: a propensity score-matched study. J Am Acad Orthop Surg. 2020 Oct 15;28(20):847-856. doi:10.5435/JAAOS-D-19-00328
- Clement ND, Gaston P, Bell A, Simpson P, Macpherson G, Hamilton DF, Patton JT. Robotic-arm assisted versus manual total hip arthroplasty: a propensity score matched cohort study. Bone Joint Res 2021 Jan;10(1):22-30.
- 30. Banchetti R, Dari S, Ricciarini ME, Lup D, Carpinteri F, Catani F, Caldora P. Comparison of conventional versus robotic-assisted total hip arthroplasty using the make system: an Italian retrospective study. J Health and Soc Sci. 2018; 3(1):37-48.
- 31. Rosinsky PJ, Chen SL, Yelton MJ, Lall AC, Maldonado DR, Shapira J, Meghpara MB, Domb BG. Outpatient vs. inpatient hip arthroplasty: a matched case-control study on a 90-day complication rate and 2-year patient reported outcomes. J Orthop Surg Res. 2020;15:367. doi: 10.1186/s13018-020-01871-8.
- 32. Chai W, Guo RW, Puah KL, Jerabek S, Chen JY, Tang PF. Use of robotic-arm assisted technique in complex primary total hip arthroplasty. Orthop Surg. 2020;12(2):686-691. doi:10.1111/os.12659
- 33. Hayashi S, Hashimoto S, Kuroda Y, Nakano N, Matsumoto T, Ishida K, Shibanuma N, Kuroda R. Robotic-arm assisted THA can achieve precise cup positioning in developmental dysplasia of the hip. Bon Joint Res 2021:10(10):629-638.
- 34. Maldonado DR, Go CC, Kyin C, Rosinsky PJ, Shapira J, Lall AC, Domb BG. Robotic arm-assisted total hip arthroplasty is more cost-effective than manual total hip arthroplasty: a Markov model analysis. J Am Acad Orthop Surg. 2021 Feb 15;29(4):e168-e177. doi 10.5435/JAAOS-D-20-00498.
- Pierce J, Needham K, Adams C, Coppolecchia A, Lavernia C. robotic assisted total hip arthroplasty a 90-day episode of care cost analysis. ISPOR Annual Meeting, virtual. May 18-20, 2020.



325 Corporate Drive Mahwah, NJ 07430 t: 201 831 5000

stryker.com

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker's product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any of Stryker's products. The products depicted are CE marked according to the Medical Device Regulation 2017/745 or the Medical Device Directive 93/42/EEC. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your sales representative if you have questions about the availability of products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Mako, Stryker. All other trademarks are trademarks of their respective owners or holders.

### $MKOTHA\text{-}BRO\text{-}4\_Rev\text{-}3\_32051$