

Simple. Versatile.
Brilliant.



Meet the new **Sonopet iQ**

Ultrasonic Aspirator System

Soft tissue resection / Fine bone dissection / All-in-one system



Sonopet

iQ

A class of **its own**

The Sonopet System is a market leader, and Sonopet iQ delivers even more of what our customers love. We've raised simplicity and reliability to an even higher level for consistently excellent O.R. surgical experiences. We've also added semi-autonomous functionality and unparalleled customization so you can do more with the Sonopet iQ with less time and effort.


Trusted resection and power

367%
faster fibrous
tissue resection¹

133%
faster bone
cutting¹

75%
more precise
suction control²

48%
faster setup time¹




**setup success rate
on first attempt by
untrained users¹**

**"Sonopet iQ is amazing. This changes
neurosurgery, and provides confidence when
resecting bone from the most critical areas"**

Neurosurgeon
Current Sonopet user

**"If the equipment is not easy to set up,
it won't be used. You guys really figured
this out, when can we use this?"**

O.R. Nurse



**of ACGME facilities
rely on Sonopet³**

**"Sonopet has changed
my whole business.
And the new system
looks even better."**

Neurosurgeon
Current Sonopet user



Rapid readiness for any O.R. team

iQ Control: On-demand customized performance

Gain a new level of control with our exclusive blend of RFID technology, algorithms and other elements. These smart features are always in play, automatically adjusting and customizing settings to optimize performance throughout your entire procedure.

iQ Connect: RFID-enabled intelligence transfer

Create and store your custom user profile once, and it's ready to engage when you are. With RFID technology, your preference data for peripheral devices (e.g., tips, handpiece and irrigation cassette) is instantly transmitted to the console, effortlessly activating the exact settings you want.

iQ control and iQ connect

RFID technology instantly changes settings based on connected tip or user profile

Faster resection and more power

Large touchscreen interface

Precise settings for power, suction and irrigation at lower levels

Unparalleled customization
One-step retrieval of surgeon profile
Can start the case using **surgeon's defaulted profile**

Independent or synchronized settings for power, suction and irrigation

Hand controller for remote control from sterile or non-sterile field



Designed for controlled suction pinch valve with ability for reduced suction around critical structures

Integrated irrigation cassette eliminates exterior wheel and tube threading

Soft tissue/bone cutting tips and tubing all compatible with the same handpiece

Front facing connections with color-changing feedback; port light signals proper setup

One-step connections for irrigation, suction, tubing and sleeve

Irrigation integrated into the handpiece to improve visibility and control




Three foot pedal options;
Basic Foot Pedal
Advanced Foot Pedal
Advanced Wireless Foot Pedal

Intuitive setup;
all components labeled or color coded, and steps are non-sequential




Gain Stryker's advantage




Proprietary Longitudinal Torsional (LT) Technology



Couples longitudinal with torsional vibration for a 25,000x per second non-rotational motion that cuts bone without affecting soft tissue.



Emulsifies hydrous soft tissue (tumor, fat, parenchyma) without impacting vital anhydrous tissue such as vessels, nerves and tendons.



Resects rigid bone or hard tissue while surrounding elastic soft tissue vibrates away from the cutting tip.



Helical reliefs redirect the ultrasonic energy from purely longitudinal to combined torsional component so tissue distal to the tip is unaffected by cavitation.

Unmatched tip offering

15 tips to handle soft tissue, fibrous tissue and bone cutting all in one handpiece

Various lengths, diameters, tip styles and frangible sections for adjustable surface exposure

RFID technology recognizes tip and adjusts to the ideal settings

Improved resonance tracking per tip type (better bone cutting and fibrous tissue resection)

Complete tip is visible inside the unopened packaging

New Serrated Large tip for even faster resection

Versatility and value

Versatility optimizes your opportunities and your return on investment. So we've designed Sonopet iQ as a multi-functional system to address clinical and administrative needs. In this manner, we strive to maintain your satisfaction far beyond the initial sale.

In the O.R.

- One console serves multiple specialties
- One handpiece compatible with soft tissue/bone tips and tubing
- Only one system to learn; reduces staff training and potential user errors
- Unmatched simplicity; readiness with any O.R. team
- New customization
- Proprietary Longitudinal Torsional (LT) Technology
- Exclusive, semi-autonomous RFID system optimization
- More features, speed and power than ever before
- Fast setup and resection for smooth case flow
- Consistently excellent surgical experiences

And beyond

- Flex Financial customizable acquisition options
- Onsite implementation and training
- Year-round education and training options
- Localized ProCare Services; 250+ technicians nationwide

We're in this together

We share a goal to make healthcare better. Pairing your clinical expertise with our innovative procedural solutions – plus supportive customer service, continuing education and financial acquisition options – can help us jointly achieve this goal.

Test your iQ

Experience Sonopet iQ system's features and benefits first hand. For a demo or more information, call your Neurosurgical sales representative or **800 253 3210**. You may also learn more at neurosurgical.stryker.com.

1. All comparative metrics and claims are Sonopet iQ compared to Sonopet 1. Stryker internal test data on file.
2. Sonopet iQ has improved suction control for delicate resection rates with 75% lower achievable pressure settings compared to the first generation Sonopet. Stryker internal data.
3. The Sonopet Ultrasonic Aspirator is currently being used in 104 out of 112 ACGME teaching institutions in the U.S. Stryker internal data on file.

Neurosurgical

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate a Stryker product. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: Barracuda, Flex Financial, ProCare, Sonopet iQ and Stryker. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

Ordering information

Capital

Part number	Product description
5500-050-000	Sonopet iQ Ultrasonic Aspirator Console (includes power cord, irrigation pole and hand controller)
5500-255-000	Sonopet iQ Universal Angled Handpiece
5500-800-278	Sonopet iQ Sterilization Tray
5500-402-000	Sonopet iQ Hand Controller
5500-006-000	Sonopet iQ Basic Foot Pedal - Wired
5500-007-000	Sonopet iQ Advanced Foot Pedal - Wired
5500-008-000	Sonopet iQ Advanced Foot Pedal - Wireless (includes wireless receiver)
5500-008-300	Sonopet iQ Wireless Receiver
5450-800-039	Sonopet Torque Wrench, 7mm
5450-850-410	Sonopet Cart

Disposables

Part number	Product description	Qty.
5450-851-012	Suction Canister	12/box
5500-573-000	Sonopet iQ Irrigation Suction Cassette	4/box
5500-472-000	Sterile Drape	4/box

Soft tissue tips

5500-25S-309	12cm iQ Micro	4/box
5500-25S-307	12cm iQ Standard	4/box
5500-25S-308	12cm iQ Large	4/box
5500-25S-301	20cm iQ Standard	4/box
5500-25S-310	12cm iQ Barracuda	4/box
5500-25S-316	12cm iQ Serrated Large	4/box

Bone cutting tips

5500-25B-315	12cm iQ Micro Claw	4/box
5500-25B-318	12cm iQ Claw	4/box
5500-25B-317	12cm iQ 120 Claw	4/box
5500-25B-312	12cm iQ Apex 360	4/box
5500-25B-302	20cm iQ Apex 360	4/box
5500-25B-303	18cm iQ Micro Claw	4/box
5500-25B-319	18cm iQ Apex Micro Claw	4/box
5500-25B-107	11cm iQ Knife	4/box
5500-25B-114	11cm iQ Apex Knife	4/box

Stryker
4100 East Milham Avenue
Kalamazoo, MI 49001 USA
t: 269 323 7700 f: 269 323 2742
toll free: 800 253 3210

neurosurgical.stryker.com