

Jump distance

Raising the bar with Modular Dual Mobility (MDM[®])



Key takeaways


- Dual mobility has been shown to increase hip stability by providing greater implant jump distances compared to conventional polyethylene liners.^{1,2}
- Stryker's MDM achieved greater jump distances when compared to Smith+Nephew OR30[™] and Zimmer Biomet G7[®] dual mobility with numerical methods (CAD analysis).^{3,4}
- Stryker's 32 mm X3 polyethylene liners also demonstrated greater jump distances when compared to many competitive sizes of dual mobility (see graph).^{3,4}

Why is jump distance important?

Implant instability resulting in hip dislocation is a known complication after THA. A potential advantage of dual mobility implants is that they can enhance hip stability by providing larger femoral head diameters that increase jump distance of the bearing, or the distance the femoral head must travel to dislocate from the implant (Fig.1).¹

The greater the implant jump distance, the lower the THA dislocation risk.²

Factors that may affect jump distance include implant head diameter, design geometry of the implant bearing surface, position of the acetabular shell and orientation of the pelvis.²




Jump distance comparison

Jump distances were compared for MDM and other competitive systems in addition to a 32 mm polyethylene liner. Implants were oriented at 45° inclination, and jump distances were calculated by numerical methods using CAD software. The results presented here may be different in clinical practice.

As a baseline for comparing the concept of jump distance, conventional polyethylene liners with a fixed femoral head increase jump distance by 18% when using a 32 mm femoral head versus 28 mm.^{3,4}

MDM vs. the competition



MDM achieved an average of:

- **37% greater** jump distance than OR30⁴
- **18% greater** jump distance than G7³

Raising the bar

Stryker's MDM versus the competition

- Through numerical methods (CAD analysis), MDM showed the greatest jump distance as compared to corresponding same-sized shells from two competitive dual mobility systems.^{3,4} On average, **MDM demonstrated a 37% and 18% increase in jump distance over the OR30 and G7 systems**, respectively.^{3,4}
- MDM achieved a larger jump distance by providing more than 180° of poly insert coverage via a 2.4 mm cylindrical buildup of the liner. This feature allows for greater jump distance for a given polyethylene head diameter (Fig. 2). In contrast, the OR30 and G7 liners provide 170° and 180° of poly insert coverage, respectively.^{3,4}

Stryker's X3 liners versus the competition

- **32 mm polyethylene liners demonstrated a larger jump distance** versus OR30 in shells sizes up to 62 mm and G7 sizes up to 52 mm.^{3,4}
- Similar to MDM, the X3 polyethylene liners provide more than 180° of femoral head coverage via a 2.7 mm cylindrical buildup (Fig. 3).³



References:

1. Dankert JF, Lygrisse K, Mont MA, Schwarzkopf R. Dual mobility total hip arthroplasty in the United States: a review of current and novel designs. *Surg Technol Int.* 2020;36:379-387.
2. Nevelos J, Johnson A, Heffernan C, Macintyre J, Markel DC, Mont MA. What factors affect posterior dislocation distance in THA? *Clin Orthop Relat Res.* 2013;471(2):519-526. doi:10.1007/s11999-012-2559-1
3. Stryker Internal Technical Memo A0012226. 2D Jump Distance comparison for MDM and G7. March 2016.
4. Stryker Internal Technical Memo D0000055054. 2D Jump Distance comparison for MDM and OR30. July 2020.

This document is intended solely for healthcare professionals.

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker's product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your sales representative if you have questions about the availability of products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Howmedica, MDM, Osteonics, PSL, Stryker, Trident, Tritanium, X3. All other trademarks are trademarks of their respective owners or holders.

Howmedica Osteonics Corp.

325 Corporate Drive
Mahwah, NJ 07430, USA
A subsidiary of Stryker Corporation
t: 201 831 5000

stryker.com