Reprocessed by

## Sustainability Solutions

### Instructions for Use

## Reprocessed Ethicon ENDOPATH® XCEL™ Universal Trocar Stability Sleeve

### Reprocessed Device for Single Use

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

- STERILE

### Explanation of Symbols

 Sterilized by Ethylene Oxide Gas

 Date of Reprocessing

 Use by Date

 Product Code

 Do Not Reuse

 See Instructions For Use

**Reprocessed Ethicon ENDOPATH® XCEL™ Universal Trocar Stability Sleeve**

The Reprocessed Ethicon ENDOPATH® XCEL™ Universal Trocar Stability Sleeve, (hereinafter Reprocessed Endoscopic Sleeves) is a sterile single patient use device with a radiolucent sleeve. The trocar sleeve contains two seals that accommodate instruments 5 mm in diameter. Together, these two seals minimize gas leakage when instruments are inserted or withdrawn through the trocar. A stopcock valve is compatible with standard luer lock fittings and provides attachment for gas insufflation and desufflation. The stopcock is in the closed position when it is parallel to the sleeve.

**Indications for Use**

The Reprocessed Endoscopic Trocars and Sleeves have applications in abdominal, thoracic, and gynecologic minimally invasive surgical procedures to establish a path of entry for endoscopic instruments. The Bladeless Trocar may be used with or without visualization for primary and secondary insertions.

**Contraindications for Use**

Reprocessed Endoscopic Sleeves are contraindicated for the following uses:

- Any uses generally contraindicated for minimally invasive techniques.

**Warnings and Precautions**

- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure. The presence of the shield on the Reprocessed Dilating Tip obturator is not a substitute for proper endoscopic techniques.
- Minimally invasive instruments may vary in diameter from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure.
- A thorough understanding of the principles and techniques involved in laser, ultrasonic and electrosurgical procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device or other medical instruments. Ensure that electrical insulation or grounding is not compromised. Do not immerse electrosurgical instruments in liquid unless the instruments are designed and labeled to be immersed.
- Using a minimally invasive instrument with a diameter smaller than specified for the Reprocessed Endoscopic Sleeves may result in desufflation of the abdominal cavity.
- During insertion of the obturator into the sleeve housing, the obturator should not be in the activated position.
- The incorporation of the shield feature in the Reprocessed Dilating Tip Trocar obturator design is intended to minimize the likelihood of penetrating injury to intra-abdominal or intra-thoracic structures. However, because the trocar tip will be briefly unprotected prior to the shield advancement, the standard precautionary measures employed in all trocar insertions must be observed.
- Adhesions, anatomical anomalies, or other obstructions, if present, may prevent or delay advancement of the shield, leaving the tip uncovered, exposing internal structures to injury.
- Although the Reprocessed Dilating Tip Trocar is design with a shield and the Reprocessed Bladeless Trocar has a blunt tip, care must still be taken, as with all trocars, to avoid damage to major vessels and other anatomic structures (such as bowel or mesentery).

To minimize the risk of such injury, be sure to:

- Establish adequate pneumoperitoneum;
- Properly position the patient to help displace organs out of the area of penetration;
- Note important anatomical landmarks;
- Direct the trocar tip away from major vessels and structures;
- Do not use excessive force.
- Do not attempt to insert the Reprocessed Dilating Tip Trocar if the red shield reset button does not stay in the activated position.
- Once partial entry has been accomplished, very little pressure may be required to complete entry of the Reprocessed Bladeless Trocar. Excessive pressure could cause injury to intra-abdominal or intra-thoracic structures.
- If reactivation is necessary after partial entry has been made, the instrument should be removed, reactivated and then reinserted. Since partial entry has been accomplished, very little pressure may be required to complete entry. Excessive pressure could cause injury to intra-abdominal or intra-thoracic structures.
- Once complete entry has been made into the abdominal or thoracic cavity, the Reprocessed Dilating Tip Trocar should not be reactivated or the Reprocessed Bladeless Trocar should not be advanced for additional penetration. Continued entry of the obturator device at this point could cause injury to intra-abdominal or intra-thoracic structures.

- Use caution when introducing or removing instruments through the trocar sleeve in order to prevent inadvertent damage to the seals which could result in loss of pneumoperitoneum. Special care should be used when inserting sharp or angled edged endoscopic instruments to prevent tearing the seal.
- When using the Reprocessed Endoscopic sleeve, additional stability devices should not be used.
- After removing the trocar from the abdominal or thoracic cavity, always inspect the site for hemostasis. If hemostasis is not present, the appropriate techniques should be used to achieve hemostasis.
- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.
- Dispose of all opened instruments whether used or unused.
- This device is packaged and sterilized for single use only.

### Adverse Reactions

None

### Directions for Use

1. Verify compatibility of all instrument and accessories prior to using the instrument (refer to **Warnings and Precautions**).
2. Prepare the patient in accordance with standard surgical techniques prior to insertion of the trocar.
3. The package label is detachable and may be affixed to the medical record of the patient.

### Reprocessed Universal Trocar Stability Sleeve for Use with Reprocessed Ethicon ENDOPATH® XCEL™ Dilating Tip Trocar with OPTIVIEW® Technology

NOTE: Refer to the Reprocessed Ethicon ENDOPATH® XCEL™ Dilating Tip Surgical Trocar with OPTIVIEW® Technology Instructions for Use for additional information.

1. Using sterile technique, remove the trocar sleeve from the package. To avoid damage, do not flip the sleeve into the sterile field.
2. The trocar obturator and sleeve are packaged unassembled. Assemble the trocar by inserting the obturator into the trocar sleeve until they lock securely together.  
NOTE: The trocar sleeve is packaged with the stopcock in the open position. Close the stopcock before use. The stopcock is closed when the stopcock lever is parallel to the trocar sleeve.  
**Important:** To allow retraction of the shield, push the reset button on the trocar obturator forward to the activated position until the audible click is heard. The trocar cannot be deactivated by forcing the red shield reset button to the original position.
3. When the trocar is in the abdominal or thoracic cavity, press the locking buttons to remove the obturator handle assembly, leaving the sleeve in place. The internal seal in the sleeve automatically closes as the obturator is withdrawn. The seal system maintains insufflation in the absence of an instrument in the sleeve. Instruments of appropriate size can be easily introduced through the sleeve of the trocar.  
**Important:** If entry into the abdominal or thoracic cavity is incomplete or the surgeon is uncertain whether entry is complete, the instrument must be reactivated. In order to reactivate the instrument, it should first be removed. After removing the instrument, push the shield reset button forward to the activated position. The shield will again be free to retract when pressure is applied. Reinsert the instrument to complete the entry.  
**WARNING:** Since partial entry has been accomplished, very little pressure may be required to complete the entry. Excessive pressure could cause injury to intra-abdominal or intra-thoracic structures.
4. To insufflate, attach a gas line to the stopcock on the trocar sleeve and open the stopcock. The seal system maintains insufflation in the absence of an instrument in the sleeve.
5. Upon completion of the procedure, remove the gas line. Open the stopcock to rapidly deflate the abdominal cavity.

### Reprocessed Universal Trocar Stability Sleeve for Use with Reprocessed Ethicon ENDOPATH® XCEL™ Bladeless Trocar with OPTIVIEW® Technology

NOTE: Refer to the Reprocessed Ethicon ENDOPATH® XCEL™ Bladeless Trocar with OPTIVIEW® Technology Instructions for Use for additional information.

1. Using sterile technique, remove the trocar sleeve from the package. To avoid damage, do not flip the sleeve into the sterile field.  
NOTE: The obturator can be used for initial puncture or after insufflation.
2. The trocar obturator and sleeve are packaged unassembled. Assemble the trocar by inserting the obturator into the trocar sleeve until they lock securely together.  
**NOTE:** The trocar is packaged with the stopcock in the open position. Close the stopcock before use. The stopcock is in the closed position when the stopcock lever is parallel to the sleeve.

3. Connect the appropriately sized, 0° endoscope to the light supply and monitor as directed in the manufacturer's instructions. Verify proper connection of the endoscope, and ensure the clarity of the picture on the monitor.
4. Insert the endoscope into the opening at the proximal end of the obturator until it reaches the distal tip of the obturator.  
NOTE: Follow these instructions for all Reprocessed Bladeless Trocar Obturators.
5. Rotate the endoscope as desired. Secure the endoscope in the obturator using the scope locking cam on the obturator.
6. Upon completion of the procedure, remove the gas line. Open the stopcock to rapidly deflate the abdominal cavity.

## Warranty

### Reprocessed Products

Stryker warrants all reprocessed products, subject to the exceptions provided herein, to be free from defects in reprocessing and to substantially conform to the product specifications contained in the documentation provided by Stryker with the products for one use in accordance with the instructions for use of such product.

**STRYKER SHALL NOT BE LIABLE FOR ANY DAMAGES TO THE EXTENT CAUSED BY ANY DEFECT IN MATERIAL, WORKMANSHIP OR DESIGN BY THE ORIGINAL MANUFACTURER OF THE PRODUCT OR ANY ACT OR OMISSION OF THE ORIGINAL MANUFACTURER OF THE PRODUCT.**

### Products for which Stryker is the Original Manufacturer

Stryker warrants all products for which it is the original manufacturer, subject to the exceptions provided herein, to be free from defects in design, materials and workmanship and to substantially conform to the product specifications contained in the documentation provided by Stryker with the products for a period of one year from the date of purchase.

### General Warranty Terms Applicable to All Products

**TO THE FULLEST EXTENT PERMITTED BY LAW, THE EXPRESS WARRANTY SET FORTH HEREIN IS THE ONLY WARRANTY APPLICABLE TO THE PRODUCTS AND IS EXPRESSLY IN LIEU OF ANY OTHER WARRANTY BY STRYKER, EXPRESSED OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OR MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL STRYKER'S LIABILITY ARISING IN CONNECTION WITH THE SALE OF THE PRODUCT (WHETHER UNDER THE THEORIES OF BREACH OF CONTRACT, TORT, MISREPRESENTATION, FRAUD, WARRANTY, NEGLIGENCE, STRICT LIABILITY OR ANY OTHER THEORY OF LAW) EXCEED THE PURCHASE PRICE, CURRENT MARKET VALUE OR RESIDUAL VALUE OF THE PRODUCTS, WHICHEVER IS LESS. STRYKER SHALL NOT BE LIABLE FOR INDIRECT, SPECIAL, INCIDENTAL, PUNITIVE OR CONSEQUENTIAL DAMAGES RESULTING FROM ANY BREACH OF WARRANTY OR UNDER ANY OTHER LEGAL THEORY.**

This warranty shall apply only to the original end-user purchaser of products directly from Stryker or a Stryker authorized distributor. This warranty may not be transferred or assigned without the express written consent of Stryker.

This warranty does not apply to: (1) products that have been misused, neglected, modified, altered, adjusted, tampered with, improperly installed or refurbished; (2) products that have been repaired by any person other than Stryker personnel without the prior written consent of Stryker; (3) products that have been subjected to unusual stress or have not been maintained in accordance with the instructions in the user manual or as demonstrated by a Stryker representative; (4) products on which any original serial numbers or other identification marks have been removed or destroyed; or (5) products that have been repaired with any unauthorized or non-Stryker components.

If a valid warranty claim is received within thirty (30) days of the expiration of the applicable warranty period, Stryker will, in its sole discretion: (1) replace the product at no charge with a product that is at least functionally equivalent to the original product or (2) refund the purchase price of the product. If a refund is provided by Stryker, the product for which the refund is provided must be returned to Stryker and will become Stryker's property. In any event, Stryker's liability for breach of warranty shall be limited to the replacement value of the defective or non-conforming part or component.

If Stryker determines in its reasonable discretion that the claimed defect or non-conformance in the product is excluded from warranty coverage as described hereunder, it will notify the customer of such determination and will provide an estimate of the cost of repair of the product. In such an event, any repair would be performed at Stryker's standard rates.

Products and product components repaired or replaced under this warranty continue to be warranted as described herein during the initial applicable warranty period or, if the initial warranty period has expired by the time the product is repaired or replaced, for thirty (30) days after delivery of the repaired or replaced product. When a product or component is replaced, the item provided in replacement will be the customer's property and the replaced item will be Stryker's property. If a refund is provided by Stryker, the product for which the refund is provided must be returned to Stryker and will become Stryker's property.

The OEM information listed on the label is provided as device ID prior to reprocessing and may contain the trademarks of unrelated third parties that do not sponsor this device.

**Sterilization:** This product and its packaging have been sterilized with ethylene oxide gas (EtO). Even though the product then is processed in compliance with all applicable laws and regulations relating to EtO exposure, Proposition 65, a State of California voter initiative, requires the following notice:

**Warning:** This product and its packaging have been sterilized with ethylene oxide. The packaging may expose you to ethylene oxide, a chemical known to the State of California to cause cancer or birth defects or other reproductive harm.

XCEL™ is a trademark of Ethicon Endo-Surgery, Inc.

ENDOPATH® and OPTIVIEW® are registered trademarks of Ethicon Endo-Surgery, Inc.

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XCEL™ est une marque de commerce d'Ethicon Endo-Surgery, Inc.
ENDOPATH® et OPTIVIEW® sont des marques déposées d'Ethicon Endo-Surgery, Inc.