


# T2 Alpha<sup>®</sup>

## Femur Retrograde Nailing System

### Design Rationale



# T2 Alpha® Femur Retrograde Nail

## Design Rationale

### Contents

Indications, contraindications and intended use .....	2
T2 Alpha Femur Retrograde Nail Design .....	4
Short Nail.....	8
T2 Alpha Femur Retrograde Instrumentation .....	9
T2 Alpha Femur Retrograde Proximal Targeting ....	10

Stryker is a market leader in intramedullary nailing. For decades, we've been creating solutions that improve treatment and set new standards. Now, with the T2 Alpha Femur Retrograde Nail, you have a new option to treat distal and diaphyseal femoral fractures.

This publication explains the design rationale and features of the T2 Alpha Femur Retrograde Nail. For detailed procedural steps please reference the T2 Alpha Femur Retrograde Nail operative technique.

### Indications for use (US and Canada)

The indications for use of these internal fixation devices include:

- Open and closed femoral fractures.
- Pseudoarthrosis and correction osteotomy.
- Pathologic fractures, impending pathologic fractures and tumor resections.
- Supracondylar fractures, including those with intra-articular extension.
- Fractures involving osteopenic and osteoporotic bone.
- Fractures distal to a total hip prosthesis.
- Periprosthetic fractures.
- Nonunions and malunions.

### Contraindications

The physician's education, training and professional judgement must be relied upon to choose the most appropriate device and treatment. Conditions presenting and increased risk of failure include:


- Any active or suspected latent infection or marked local inflammation in or about the affected area.
- Compromised vascularity that would inhibit adequate blood supply to the fracture or the operative site.
- Bone stock compromised by disease, infection or prior implantation that cannot provide adequate support and /or fixation of the devices.
- Material sensitivity documented or suspected.

# T2 Alpha Femur Retrograde Nail Design


The T2 Alpha Femur Retrograde Long Nail features 7 distal locking options, the 4 most distal screw options are placed from 6 – 32mm from the driving end of the nail in three planes of fixation. This is the most distal locking configuration on the market\*. This system may allow surgeons to use an IM nail in lieu of plating distal femoral fractures. As a load sharing device, surgeons may prefer intramedullary nails due to early weight bearing potential and their minimally invasive nature<sup>6,7</sup>.

The locking configuration was designed for distal and diaphyseal femur fracture fixation<sup>4</sup>.


Based on SOMA design we validated that the screw trajectory captures the densest bone.<sup>9,10</sup>



Examples of locking options for distal diaphyseal and supracondylar fractures



3 Planes of Fixation



# T2 Alpha Femur Retrograde Nail Design

Mismatch of IM nail geometry and the bow of the femur increases the risk of anterior cortex perforation. The T2 Alpha Femur Retrograde Nail was designed using the Stryker Orthopedics Modeling and Analytics (SOMA) software. The software showed improved anatomic compliance with a length dependent radius of curvature in the long nails<sup>9,10</sup>. The radius of curvature ranges from .70 to 1.35m\*.

The distal end of the T2 Alpha Femur Retrograde nail also provides a graduated, length dependent posterior bend which is designed to allow for either a standard retrograde entry point or for a more posterior entry point in periprosthetic fractures<sup>10,12,16</sup>.



**T2 Alpha Femur Retrograde Nail shows improved anatomic compliance when compared to nails with a radius of curvature of 1.0 and 2.0M.**

Radius of Curvature:  
2.0m



**Previous generation  
Femur Retrograde  
Nail**

Radius of Curvature:  
.70-1.35m




**T2 Alpha Femur  
Retrograde Long Nail**



\*Long nails only

# T2 Alpha Femur Retrograde Nail Design


The T2 Alpha Retrograde Femur Nail allows the user to apply up to 10mm of compression through the nail. The user may choose to apply compression using the internal or external compression method. This allows for a more controlled reduction of the fragments when compared to the traditional backslapping method.




**External compression\***



**External compression device\***



**Internal compression**



\* Please note when using external compression device it is required to add a secondary locking screw proximally



# T2 Alpha Femur Retrograde Nail Design

T2 Alpha Femur Retrograde includes three screw options: The IMN Locking Screw, the Advanced Locking Screw and the Condyle Bolts. The Advanced Locking Screw can be used when additional stability is needed.

## Advanced Locking Screws\*

When treating distal femur fractures, fixed angle constructs may be advantageous in providing stability in short condylar segments and osteoporotic bone<sup>11</sup>.




The Advanced Locking Screws are designed with oversized threads which engage with the internal threads of the T2 Alpha Femur Retrograde Nail, thereby limiting relative movement between the nail and the screw<sup>1,3</sup>.

The operative technique requires use of a 5.5mm over-drill, which widens the cortex to allow for insertion of the 5.3mm section of the screw. Closest to the head of the screw, the 5.7mm threads are designed to ensure purchase of the near cortex.

## Mechanical Performance<sup>1</sup>

The usage of Stryker's Advanced Locking Screws did not affect torsional and bending stiffness, whereas axial stiffness was increased compared to the fully threaded locking screws. [T2]\*\*.

The usage of Stryker's Advanced Locking Screws resulted in a negligibly low axial play and significantly decreased both torsional and bending play compared to the locking screws, fully threaded [T2]\*\*.



\*Advanced Locking Screws can only be used in 5mm round holes  
 \*\*Locking screws, fully threaded [T2, Stryker]


# T2 Alpha Femur Retrograde Nail Design

## IMN Locking Screws

The T2 Alpha IMN Locking Screws show an increased fatigue strength of 30% when compared to the legacy T2 Screws<sup>5,12</sup>.



Threaded head of both the IMN Locking Screw and Advanced Locking Screw engages with the self retaining screwdriver sleeve



## T2 Alpha Femur Retrograde Short Nail

To meet surgeon preference and diverse fracture patterns, the T2 Alpha Femur Retrograde Nailing System offers a 170mm and 200mm short nail option.

The short nail has 4 distal locking options 6 – 32mm from the driving end of the nail in three planes of fixation. Proximally, the nail has two lateral to medial locking options which are inserted using the short proximal targeting arm.



\*Plate and nail in the image shown is not a linked construct. The T2 Alpha Femur Retrograde nail is not indicated for linking with plates.


**Image of short nail used with a femoral stem and anatomic plate.\***



# T2 Alpha Femur Retrograde Instrumentation

All T2 Alpha Nailing Systems share a core set of general instruments to ensure consistency across all indications. The T2 Alpha basic instruments and the T2 Alpha Femur Retrograde indication instruments are designed to improve instrument continuity and functionality for a streamlined operative workflow<sup>19</sup>.

- Robust self-retaining screwdriver.
- Depth gauges with improved legibility.
- Universal modular interfaces for quick-lock assembly.
- Static distal targeting arm eliminates the need to rotate targeting arm for distal screw insertion.




## T2 Alpha Femur Retrograde Proximal Targeting

Based on the T2 Alpha GT/PF & Tibia Distal Targeting Technology, The T2 Alpha Proximal Targeting Device provides guided locking for the proximal AP holes of the long T2 Alpha Femur Retrograde nails. The guided proximal targeting device is designed to provide reproducible locking and reduce the total number of x-ray shots compared to freehand locking<sup>12,17,18</sup>.



### Proximal Targeter Positioning

To perform guided proximal locking, it is essential to place the X-ray beam of a C-arm approximately 30 degree oblique to the axis of the drill sleeve assembly. Please defer to the T2 Alpha Femur Retrograde operative technique for a procedural overview.



## Notes

## References

1. T2-WP-6\_Rev-1, 12-2019
2. T2-WP-8, 03-2019
3. Internal Test Report A0030043
4. T2-ST-37, 06-2021
5. D0000082802, KIE 914-4- 5550\_T2 Alpha Femur Retrograde\_DOF 30-502 SOMA Report\_NSTD\_2021\_04
6. G. T. Gänsslen A, Hildebrand F, Pape HC, Oestern HJ, "Femoral shaft fractures in adults: treatment options and controversies.," (in English), Acta Chir Orthop Traumatol Cech, vol. 81, no. 2, pp. 108-17, 2014,
7. AO Surgery Reference, Baumgaertner M. et al., Plate and retrograde nailing, Published , Accessed July 18, 2022, <https://surgeryreference.aofoundation.org/orthopedic-trauma/periprosthetic-fractures/knee/fracture-around-a-stable-femoral-component-with-good-bone-stock/plate-and-retrograde-nailing>
8. D0000082804, KIE 914-4- 5550\_T2 Alpha Femur Retrograde\_DOF 30-502 SOMA Report\_NSTD\_2021\_06
9. Stryker Internal Test Report D0000082804, Stryker, Trauma & Extremities, 2021
10. Stryker Internal Test Report D0000084065, Stryker, Trauma & Extremities, 2021
11. Ghandour A, Cosker TD, Kadambande SS, Pallister I, Russell ID. Experience of the T2 supracondylar nail in distal femoral fractures. Injury. 2006 Oct;37(10):1019-25. doi: 10.1016/injury.2006.07.001. Epub 2006 Aug 28. PMID: 16934811
12. Stryker White Paper, T2-WP-6, Stryker, Trauma & Extremities, 2019
13. Zimmer® Natural Nail® System Retrograde Femoral Nail Surgical Technique (zimmerbiomet.com)
14. 103661491 Rev C - RFNA Surgical Technique Guide.pdf (llnwd.net)
15. trigen-metanail-retro-fem-st-01362v3.pdf (smith-nephew.com)
16. Stryker Internal Test Report D0000082802, Stryker, Trauma & Extremities, 2021
17. Adam, P., Taglang, G., Brinkert, D., Bonnomet, F., Ehlinger, M. (2012) External distal targeting device for safe insertion of distal locking screws when performing long gamma nail: comparison with the freehand technique. Orthopaedic proceedings, 94-B (SUPP XXXVII), 435.
18. M.Ehlinger et al. Distal targeting device for long Gamma nail Monocentric observational study. Orthopaedics and Traumatology:Surgery and Research. Vol 99-Issue 7. November 2013, Pages 799-804
19. EPS Final Report T2 Alpha Femur Retrograde, 2021

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate a Stryker product. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker, T2 Alpha. All other trademarks are trademarks of their respective owners or holders.

Content ID: T2-SS-23, 07-2022

Copyright © 2018 Stryker



### Manufacturer:

Stryker Trauma GmbH  
Prof.-Küntschers-Str. 1-5  
24232 Schönkirchen  
Germany  
[stryker.com](http://stryker.com)