


T2

Recon Nailing System R1.5

Operative Technique



Introduction

Contributing Surgeons:

Kyle F. Dickson, MD, MBA

Professor and Chairman
University of Texas Medical School at Houston
Department of Orthopaedic Surgery
Houston, Texas
USA

We greatly acknowledge and appreciate the contributions to this operative technique made by:

Kevin W. Luke, M.D.

Parkview Orthopaedic Group
Assistant Clinical Professor
Department of Orthopaedic Surgery
University of Illinois at Chicago
Illinois, Chicago
USA

Anthony T. Sorkin, M.D.

Rockford Orthopaedic Associates, LLP
Clinical Instructor
Dep. of Surgery University of Illinois
College of Medicine
Director Orthopaedic Traumatology
Rockford Memorial Hospital
Rockford, Illinois
USA

Ariaan D.P. van Walsum, MD

Trauma surgeon
Medical Spectrum Twente
Enschede
Netherlands

Don Weber, MD, FRCSC

Associate Clinical Professor of
Orthopaedics
Chief of Orthopaedics
University of Alberta Hospital
Edmonton, Alberta
Canada

This publication sets forth detailed recommended procedures for using Stryker Trauma devices and instruments.

It offers guidance that you should heed, but, as with any such technical guide, each surgeon must consider the particular needs of each patient and make appropriate adjustments when and as required. A workshop training is required prior to first surgery.

See package insert (L22000007) for a complete list of potential adverse effects, contraindications, warnings and precautions. The surgeon must discuss all relevant risks, including the finite lifetime of the device, with the patient, when necessary.

Contents

| | Page |
|---|-------------|
| 1. Introduction | 4 |
| Implant Features | 4 |
| Technical Specifications | 5 |
| Instrument Features | 6 |
| 2. Indications and Contraindications | 7 |
| 3. Pre-operative Planning | 7 |
| 4. Locking Options | 8 |
| 5. Operative Technique | 9 |
| Patient Positioning and Fracture Reduction | 9 |
| Incision | 9 |
| Entry Point | 10 |
| Reaming | 12 |
| Nail Selection | 13 |
| Assembly of the Targeting Device and the Nail | 13 |
| Nail Insertion | 14 |
| Guided Locking for the Recon Mode | 15 |
| Guided Locking for Antegrade Femoral Mode | 23 |
| Freehand Distal Locking | 26 |
| Set Screw or End Cap Insertion | 27 |
| Nail Removal | 27 |
| <hr/> | |
| Ordering Information – Implants | 28 |
| Ordering Information – Instruments | 30 |

Introduction

1. Introduction

Over the past decades **antegrade femoral nailing** has become the treatment of choice for most femoral fractures.

As an addition to the **T2 Nailing System**, Stryker Trauma has created a new generation femoral implant: the **T2 Recon Nail** for the treatment of complex, as well as more common fractures.

The advantages of using intramedullary fixation for the treatment of proximal femur fractures include less soft tissue dissection and stable fracture fixation with a load sharing device.

Through the development of a common, streamlined instrument system and intuitive surgical approach, both in principle and in detail, the **T2 Recon Nail** offers **the opportunity for significantly increased speed and functionality** for the treatment of fractures and simplifies the training requirements for all personnel involved.

Furthermore, the **T2 Recon Nail** offers the following competitive advantages:

- **Versatility to switch from a standard antegrade femoral nailing to a recon option without changing the nail**
- **An optional Set Screw, Recon that can be tighten down onto the superior Lag Screw, thus minimizing the potential sliding of the proximal Lag Screw.**
- **Static or dynamic distal locking options for the antegrade femoral mode**
- **Trochanteric entry point**
- **Reduced proximal nail diameter allowing freehand placement of accessory K-Wires around the nail (anterior and posterior) for precise femoral neck fracture reduction**

Implant Features

The **T2 Recon Nail** is the realization of superior biomechanical intramedullary stabilization using strong, cannulated implants for the internal fixation of the Femur.

As with all other T2 Nails, the **T2 Recon Nail** is made of **Type II anodized Titanium Alloy (Ti6Al4V)** for **enhanced biomechanical and biomedical performance.**

The **T2 Recon Nail** features a **125° CCD angle and a 10° anteversion angle for the 2 proximal holes which utilize 6.5mm cannulated Lag Screws.** With this lower **CCD angle**, easy insertion of 2 screws into the femoral head can be achieved.

Alternatively a **proximal 70° Oblique hole with 7° retroversion** for a 5mm Fully Threaded Screw can target the Lesser Trochanter in the Femoral Antegrade mode.

The 6.5 mm cannulated Lag Screws have a unique thread design providing a better grip, improved front cutting flutes for a lower insertion torque and thinner flanks for less bone removal. Secure placement of the Lag Screws within even very small neck diameters can be achieved due to the design of a 10.5mm inner, therefore 17mm outer distance between the two 6.5mm Lag Screws.

Two Set Screws are available:

- a **Set Screw, Recon** to tighten down on the proximal Lag Screw (for the Recon Mode) and
- a **Set Screw, Antegrade** to tighten down on the oblique Fully Threaded Screw (for the Femoral Antegrade Mode).

Available as **left and right versions**, the **T2 Recon Nail** incorporates an antecurvature radius of **1.5M of the shaft**, as well as a **4° Medial Lateral bend** for trochanteric insertion.

The **distal locking configuration** features a round and an oblong hole to allow for **static or dynamic distal locking.**

Low profile 5mm cortical screws, common to the T2 Nailing System, are designed to simplify the surgical procedure and promote a minimally invasive approach.

- Fully Threaded Locking Screws are available for distal locking (Recon or Femoral Antegrade Mode) and for the proximal locking in Femoral Antegrade Mode.

End Caps are available in various length to provide an **improved fit** for every indication.

See the **detailed chart on the next page** for the design specifications and size offerings of the implants.

Features

Technical Specifications


Nail Diameter 9, 10, 11, 12 and 13mm
(Left and Right)

Sizes 260–420mm, in 20mm increments

Note:
The proximal diameter is 13mm for the 9, 10, 11 and 12mm Nails and 15mm for the 13mm Nail.


Set Screw, Recon

Set Screw, Antegrade



5.0mm Fully Threaded Locking Screws

L = 25–120mm



6.5mm cannulated Lag Screws

L = 65–130mm



Note:
Screw length is measured from top of head to tip.

End Caps



Features

Instrument Features

The major advantage of the instrument system is a breakthrough in the integration of a **core instrument platform** which can be used not only for the complete **T2 Nailing System**, but represents the platform for future Stryker Trauma nailing systems, reducing complexity and inventory.

The T2 instrument platform offers advanced precision and usability, and features ergonomically styled targeting devices.

Except for the addition of a small number of dedicated instruments, the T2 Femur instrument platform is used for the **T2 Recon Nail**.


Dedicated instruments for the **T2 Recon Nail** include the **Recon Targeting Device** which has one **Nail Adapter** and two Targeting Arms:

- **Targeting Arm, Recon** used for the placement of two 6.5mm cannulated Lag Screws into the femoral head in the Recon mode
- **Targeting Arm, Antegrade** used for insertion of the oblique screw in the Antegrade mode

In addition to the advanced precision and usability, the instrument tray is numbered and color coded to indicate their usage during the surgical procedure. The number coding indicates the step during the procedure in which the instrument is used.

With the exception of the carbon fiber targeting device, dedicated instruments for the recon mode are color coded with “bronze”.


This makes it easy to differentiate them from the core platform instruments.



Nail Adapter



Targeting Arm, Recon



Targeting Arm, Antegrade

Drills

Drills feature color coded rings:

4.2mm = Green

(Consistent with the Gamma3 and T2 Instrument Platform, this drill features a green color ring.)

The 4.2mm drills are used for 5.0mm Fully Threaded Locking Screws (either for distal locking or for proximal oblique locking).

6.5mm

The Solid Stepdrill for Lag Screw is color coded with “bronze”.

Indications & Contraindications

2. Indications and Contraindications

The **T2 Recon Nail** is indicated for:

- Open and closed femoral fractures
- Pseudarthrosis and correction osteotomy
- Pathologic fractures and impending pathologic fractures
- Intertrochanteric and Subtrochanteric fractures
- Ipsilateral neck/shaft fractures

Relative Contraindications:

The physician's education, training and professional judgement must be relied upon to choose the most appropriate device and treatment. Conditions presenting an increased risk of failure include:

- Any active or suspected latent infection or marked local inflammation in or about the affected area.
- Compromised vascularity that would inhibit adequate blood supply to the fracture or the operative site.
- Bone stock compromised by disease, infection or prior implantation that can not provide adequate support and/or fixation of the devices.
- Material sensitivity, documented or suspected.
- Obesity. An overweight or obese patient can produce loads on the implant that can lead to failure of the fixation of the device or to failure of the device itself.
- Patients having inadequate tissue coverage over the operative site.
- Implant utilization that would interfere with anatomical structures or physiological performance.
- Any mental or neuromuscular disorder which would create an unacceptable risk of fixation failure or complications in postoperative care.
- Other medical or surgical conditions which would preclude the potential benefit of surgery.

3. Pre-operative Planning

An X-Ray Template, Recon (1806-3082) is available for pre-operative planning. Thorough evaluation of pre-operative radiographs of the affected extremity is critical. Careful radiographic examination of the trochanteric region and neck regions can reduce the potential of intra-operative complications.

According to the fracture type, either Recon or Antegrade Femoral Mode can be chosen.


Evaluation of the femoral neck angle on the pre-operative X-Rays is mandatory as the **T2 Recon Nail** has a fixed 125° neck angle for the two Lag Screws and proper placement of both Lag Screws in the femoral head is essential.

If possible, X-Rays of the contralateral side should be used to determine the normal neck angle and length of the femur.

The proper nail length should extend from the Tip of the Greater Trochanter to the Epiphyseal Scar.

Note:

Check with local representative regarding availability of nail sizes.



Operative Technique

4. Locking Options

The T2 Recon Nail can be locked proximally either with two Lag Screws (Recon Mode, Fig. 1) or with one Fully Threaded Screw (Antegrade Femoral Mode, Fig. 2).

For both Recon and Antegrade Femur applications, depending on fracture pattern, either static or dynamic distal locking can be used.

Recon Mode




Fig. 1

Antegrade Femoral Mode




Fig. 2

Operative Technique

5. Operative Technique

Patient Positioning and Fracture Reduction

Patient positioning for **T2 Recon Nail** insertion is surgeon dependent. However, it is recommended to position the patient in supine or lateral position on a fracture table, to allow closed reduction of the fracture (Fig. 3).

Manipulate and reduce the fracture in the usual fashion, according to the fracture type. Reduction should be achieved as anatomically as possible. If this is not possible, reduction in one plane should be complete, leaving reduction in the other plane to be achieved prior to reaming and nail insertion.

The unaffected leg is abducted as far as possible to ease image intensifier positioning. This will also allow easier access to entry point.




Fig. 3

Incision

The design of the **T2 Recon Nail**, with a **4° Medial Lateral bend**, will **only** allow the insertion through the **Tip of the Greater Trochanter**.

With experience, the Tip of the Greater Trochanter can be identified by palpation (Fig. 4).

A longitudinal skin incision of approximately 3–5cm is made starting just above the Greater Trochanter to the Iliac Crest (Fig. 5). The incision is then deepened to expose the Tip of Greater Trochanter.




Fig. 4

Smaller or larger incisions may be used based on individual patients anatomy and at the surgeons discretion.

Note:

The targeting instruments of the T2 Recon Nail have been designed to allow for a more percutaneous approach.




Fig. 5

Operative Technique

Entry Point

• The Tip of the Greater Trochanter

The entry point is located at the junction of the anterior third and posterior two-thirds of the Greater Trochanter, on the medial edge of the tip itself (Fig. 6).

Note:
Before opening the Tip of Greater Trochanter, use image intensifier views (A/P and M/L) to confirm correct identification of the entry point.

The medullary canal can be opened with the

- Curved Awl/Curved Awl, 90° Handle or
- One Step Conical Reamer.

Note:

During opening of the entry portal with the Awl, dense cortex may block the tip of the Awl. An optional Awl Plug can first be inserted through the Awl to avoid penetration of bone debris into the cannulation of the Awl shaft, and then removed for Guide Wire insertion.

• Entry point with Curved Awl

Once the Tip of the Greater Trochanter has been opened, the Ø3 × 1000mm Ball Tip Guide Wire may be advanced through the cannulation of the Curved Awl with the Guide Wire Handle and Chuck (Fig. 7).

The proximal femur may then be prepared with the One Step Conical Reamer.




Fig. 6



Fig. 7

Operative Technique

• Entry point with One Step Conical Reamer

Alternatively, the 13mm diameter One-Step Conical Reamer for the 9, 10, 11, 12 and 13mm nails or the 15mm diameter Reamer for the 13mm nail may be used for opening the medullary canal and reaming of the trochanteric region.

Under image intensification control, the entry point is made with a $\text{Ø}3.2 \times 400\text{mm}$ K-Wire, Recon attached to the Guide Wire Handle and advanced into the medullary canal. Confirm its placement within the center of the medullary canal on A/P and lateral image intensifier views.

Note:

The K-Wire, Recon used for the entry point should not be used again for the Lag Screw insertion. It is recommended to utilize a new K-Wire.

The Protection Sleeve, Recon and Multi-hole Trocar are positioned with the central hole over the K-Wire.

Note:

The Multi-hole Trocar has a special design for more precise insertion of the $\text{Ø}3.2\text{mm}$ Recon K-Wire (Fig. 8). Beside the central hole, 4 other holes are located eccentrically at different distances from the center (Fig. 8a) to easily revise insertion of the guiding K-Wire in the proper position (entry point).

When correct placement of the guiding K-Wire is confirmed on image intensifier views (A/P and lateral), keep the Tissue Protection Sleeve in place and remove the Multi-hole Trocar.

The T-Handle is attached to the One-Step Conical Reamer and **hand reaming** is performed over the K-Wire through the Tissue Protection Sleeve (Fig. 9).

The K-Wire is then removed and replaced with the $\text{Ø}3 \times 1000\text{mm}$ Ball Tip Guide Wire.




Fig. 8




Fig. 9

Operative Technique

Reaming

The $\text{\O} 3 \times 1000\text{mm}$ Ball Tip Guide Wire is inserted with the Guide Wire Handle through the fracture site to the level of the Epiphyseal Scar.

The $\text{\O} 9\text{mm}$ Universal Rod with Reduction Spoon may be used as a fracture reduction tool to facilitate Guide Wire insertion through the fracture site (Fig. 10).

Note:

The Ball Tip at the end of the Guide Wire will stop the Bixcut reamer* head (Fig. 11).

Note:

Prior to reaming, it is important to check the centered intramedullary position of the Guide Wire with the image intensifier. Lateral displacement of the Guide Wire could lead to resection of more bone on the lateral side of the wire, which in turn will lead to an offset position of the nail and the risk of fracturing the shaft.

Note:

Make sure that the reduction is maintained throughout the reaming process.

Reaming is commenced in 0.5mm increments until cortical contact occurs (Fig. 12).

For easier nail insertion, the medullary canal should be reamed 2mm more than the diameter of selected nail.

Note:

The **T2 Recon Nail** may be inserted without reaming of the subtrochanteric and dyaphyseal region of the femur, particularly in elderly patients with wide medullary canals. If appropriate, after the trochanteric region is prepared with the One-Step Conical Reamer, the nail can be inserted without further reaming of the medullary canal.

Reaming of the trochanteric region is needed (Fig. 13) as the proximal nail diameter (driving end) is larger than the nail diameter (13mm for the 9, 10, 11, 12 and 13mm diameter nails and




Fig. 10




Fig. 11




Fig. 12




Fig. 13

15mm for the 13mm diameter nail). For both reamed or unreamed applications, the proximal 5cm of the trochanteric region must be opened to at least 13mm or 15mm, (depending on

the proximal diameter of the nail).

* see pages 32-33 for additional Bixcut Reamer system details

Operative Technique


Nail Selection

Diameter

The diameter of the selected nail should be 1.5–2.0mm smaller than that of the last reamer used.

Length

Nail length may be determined by measuring the remaining length of the Guide Wire. The Guide Wire Ruler may be used by placing it on the Guide Wire and reading the correct nail length at the end of the Guide Wire on the Guide Wire Ruler (Fig. 14 and 15).



End of Guide Wire Ruler equals Measurement Reference

Fig. 14




Fig. 15

Assembly of the Targeting Device and the Nail

Upon completion of reaming, the appropriate size nail is ready for insertion. A unique design feature of the T2 Recon Nail is that the $\text{Ø}3 \times 1000\text{mm}$ Ball Tip Guide Wire does not need to be exchanged.

The selected nail is assembled onto the Nail Adapter with the Nail Holding Screw (Fig. 16a). Be sure to securely tighten the Nail Holding Screw with the Screwdriver Shaft, Ball Tip and T-Handle so that it does not loosen during nail insertion.



Fig. 16a

Note:

Prior to the nail insertion, check the correct assembly by passing the Stepdrill for Lag Screw through the Tissue Protection Sleeve, Recon and/or Drill Sleeve, Recon placed in the corresponding hole of the Targeting Arm, Recon, into the holes of the nail (Fig. 16b).

For the Femoral Antegrade Mode, use the Targeting Arm, Antegrade with the Tissue Protection Sleeve and Drill Sleeve assembly to pass the $\text{Ø}4.2 \times 340\text{mm}$ Drill through the oblique hole of the nail.



Fig. 16b

Operative Technique

Nail Insertion

The nail is advanced through the entry point passing the fracture site to the appropriate level.

If dense bone is encountered, first re-evaluate that sufficient reaming has been achieved, then, if necessary, the Strike Plate, Recon can be attached to the Nail Adapter and the Slotted Hammer may be used to further insert the nail (Fig. 17).

Note:

The nail must progress smoothly, without excessive force. If too much resistance is encountered, removal of the nail and additional reaming is recommended.

Note:

Remove the Guide Wire prior to drilling or K-Wire insertion.



Fig. 17

Operative Technique

Guided Locking for the Recon Mode

Nail/Lag Screws Positioning

Drive the T2 Recon Nail to the depth that correctly aligns the proximal screw holes parallel with the femoral head and neck under fluoroscopic control (Fig. 18).

Two aspects regarding the Nail/Lag Screws position must be carefully checked with the image intensifier before drilling into the femoral head:

- Alignment of the anteversion (ML view)
- Depth of Nail insertion (A/P view).

The distal Lag Screw should run along the calcar region (on the A/P view) and centered into the femoral neck and head (on the ML view).

Note:

The use of the One Shot Device (1213-3010) is recommended to predetermine the optimal Lag Screw placement. Details are described on Page 18 to 19.

Attach the Targeting Arm, Recon to the Nail Adapter. The Fixation Screw can be firmly tightened by using T2 Recon/Femur Wrench (Fig. 19a).

Slide in the Targeting Arm, Recon as close to the limb as possible in order to obtain the highest mechanical stability (Fig. 19), then securely tighten the Fixation Screw.

Note:

Make sure that no pressure is applied on the Targeting Arm, Recon while tightening the Fixation Screw to avoid possible malalignment.




Fig. 18




Fig. 19



Fig. 19a

Operative Technique

Now attach the Paddle Trocar, Recon, to the T-Handle, AO Medium Coupling (Fig. 20). Then, advance them together with the Tissue Protect Sleeve, Recon to the skin through the distal (A) hole on the Targeting Arm, Recon by pressing the distal safety clip. The Targeting Arm has a locking mechanism and it opens and closes by pressing the distal safety clip. This mechanism will keep the sleeve in place and prevent it from sliding out. A small skin incision is made and the assembly is pushed through until it is in contact with the lateral cortex (Fig. 21).




Fig. 20

Then remove the Trocar and then insert the K-Wire Sleeve, Recon through the Targeting Arm. Place a K-Wire, Recon into the K-Wire Inserter and attach it to the T-Handle. The K-Wire is then manually advanced through the K-Wire Sleeve until it reaches the subchondral bone of the femoral head (Fig. 22). Alternatively, the K-Wire Inserter can be attached to a Power Tool and the K-Wire, Recon is inserted to the same depth.




Fig. 21




Fig. 22

Operative Technique

Note:

With the image intensifier, verify if the K-Wire is placed along the calcar region in the A/P view, and central on the lateral view (correct anteversion) (Fig. 23).

If the K-Wire is incorrectly positioned, the first step is to remove it and then to reposition the nail.

More commonly, the nail is positioned too proximal and repositioning should be carried out either by hand or by using the Strike Plate, Recon placed into the Nail Adapter, inserting the nail further. If a higher position is required, the Universal Rod and Slotted Hammer may then be attached to the Strike Plate, Recon to carefully and smoothly extract the assembly (Fig. 24).

The new position is checked again with the image intensifier as described above.




Fig. 23



Fig. 24

Operative Technique

Nail/Lag Screws Positioning with the One Shot Device

The use of the One Shot device (1213-3010) is recommended to predetermine the optimal Lag Screw placement* (Fig. 25). The One Shot Device is made of carbon fiber and works by providing a target to indicate the position of the K-wire on the fluoroscope screen. The target contains 3 radio-opaque wires embedded in the arm – a dashed inner wire and two solid outer wires. These wires work like a gun sight to indicate the position of the K-wire.

The One Shot Device is attached by slightly pressing the grip and releasing it when positioned onto the Tissue Protection Sleeve. To reposition or remove the device, the grip must be depressed.

Note:

The use of the One Shot Device should not interfere with or replace any steps in the **T2 Recon Nail Operative Technique**.

While depressing the attachment grip, the device is positioned between the anterior aspect of the patient's hip and the fluoroscope screen positioned for an A/P view of the hip (Fig. 25 and 26).

Note:


It is important to drape the patient such that the One Shot Device does not interfere with any drapes anterior to the patient's hip.

When positioned correctly, the target will appear in the fluoroscopic image (A/P view), with the dashed inner wire in the middle of the two solid outer wires (Fig. 26). If it does not, the One Shot Device should be moved towards or away from the patient by depressing the grip slightly until the target is seen as described above.

* Tokunaga et al, Correct lag screw positioning for the Gamma Nail: Development for the targeting device for insertion, *Osteo Trauma Care* 2005; 13:14-17



Fig. 25



A/P view

Fig. 26

Operative Technique

Note:

The **One Shot Device** cannot be rotated 90° with the Tissue Protection Sleeve because it will be stopped by the Nail Adapter. Therefore, the C-Arm must be turned more than 90°, in order to get a lateral view of the **One Shot Device**.

To identify the accurate positioning, the dashed wire of the target must appear between the two solid wires at the desired position. If the position is incorrect the T2 Recon Nail may be repositioned by either pulling backwards or pushing forwards. The K-wire can then be placed into the femur and the targeting arm is held in place until the K-wire's position in the Lateral view has been determined (Fig. 27). When positioned correctly, the target will appear in the fluoroscopic lateral views (Fig. 28).

If the dashed wire of the target appears between the two solid wires, then advance the Tissue Protection Sleeve and Trocar, Recon as shown in Fig. 20.

Note:


Prior to advancing the K-Wire, check the correct guidance through the K-Wire Sleeve. Do not use bent K-Wires.

Note:

The K-Wire inserted into the most distal Lag Screw hole of the nail helps in achieving the correct positioning of the nail (depth & rotation) with minimal resection of bone in case repositioning is needed.



Fig. 27



Lateral view

Fig. 28

Operative Technique

Solid Stepdrill Technique

For the insertion of proximal screws in Recon Mode, the Cannulated Drilling Technique was also mentioned in a previous version of the Operative Technique. Several years of T2 Recon Nailing experience has shown that the Solid Stepdrill Technique, which is mentioned in this chapter, is the recommended method to optimize the proximal targeting accuracy.

Attach the Paddle Trocar, Recon to the T-Handle, AO Medium Coupling. Then slide the Tissue Protection Sleeve, Recon together with the Paddle Trocar assembly to the skin through the proximal (B) hole on the Targeting Arm, Recon by pressing the safety clip. A small skin incision is made and the assembly is pushed through until it is in contact with the lateral cortex (Fig. 29).

Then remove the Trocar assembly and insert the Drill Sleeve for the Solid Stepdrill, Recon while the distal K-Wire, Recon and K-Wire Sleeve are still left in place. The Drill Sleeve for the Solid Stepdrill, Recon, is inserted through the proximal hole of the Targeting Arm, Recon. The Ø6.5mm Solid Stepdrill for Lag Screw, Recon, is forwarded through the Tissue Protection Sleeve and Drill Sleeve assembly and pushed onto the lateral cortex (Fig. 30). There is a dedicated Drill Sleeve for the Solid Stepdrill Technique. This Sleeve is marked "Use with Solid Step Drill" as shown (Fig. 30b).

Reaming is performed under fluoroscopic control just until the tip of the Solid Stepdrill for Lag Screw reaches the subchondral bone. The required length of the Lag Screw can be read directly off the Solid Stepdrill for Lag Screw, Recon at the end of the Drill Sleeve (Fig. 30a).




Fig. 29



Fig. 30




Fig. 30b

Operative Technique

Using the Screwdriver, Recon, the correct Lag Screw is inserted through the Tissue Protection Sleeve and threaded up to the sub- chondral part of the femoral head. The screw is near its proper seating position when the groove around the shaft of the screwdriver is approaching the end of the Tissue Protection Sleeve (Fig. 31 & 31a).

The required length of the second Lag Screw can be measured using the Lag Screw Gauge, Recon.

Remove the Distal K-Wire, Recon and K-Wire Sleeve and insert the Sleeve for the solid Stepdrill into the distal Tissue Protection Sleeve.

Repeat the same surgical steps for drilling and insertion of the distal Lag Screw without K-Wire guidance.




Fig. 31




Fig. 31a

Operative Technique

Alternatively, the K-Wire can be used prior to drilling with the Solid Drill.

Place a second Recon K-Wire into the K-Wire Inserter and attach it to the T-Handle or power tool. The K-Wire is then advanced through the K-Wire Sleeve until it penetrates the subchondral bone of the femoral head.

Note:

Correct placement of the K-Wire tip in subchondral bone must be checked with image intensifier in both A/P and Lateral views.




Fig. 30c

The required length of the Lag Screw is measured using the Lag Screw Gauge, Recon. Before starting to measure, ensure that the Tissue Protection Sleeve and K-Wire Sleeve assembly is firmly pressed against the lateral cortex of the femur (Fig. 30d).

Take the Lag Screw Gauge, Recon and place it directly under the distal K-Wire and against the K-Wire Sleeve (Fig. 30d). The correct Lag Screw length corresponds to the measurement indicated at the end of the K-Wire on the Lag Screw Gauge.

After the measurement, remove the K-Wire and drill the channel with the Solid Stepdrill according to the Solid Stepdrill technique described on page 20.



Fig. 30d

Operative Technique

Guided Locking for Antegrade Femoral Mode

For Antegrade Femoral Mode, attach the Targeting Arm, Antegrade onto the Nail Adapter. The Targeting Arm will sit in the groove positioned on the upper surface of the Nail Adapter because of an integrated spring locking mechanism which prevents sliding. A “click” will confirm correct positioning of the Targeting Arm (Fig. 32). Tightening of the Locking Knob is mandatory and it can be firmly retightened by using the T2 Recon/Femur Wrench in order to achieve a precise proximal locking (Fig. 32a).

Now attach the Paddle Trocar, Antegrade and the T-Handle, AO Medium Coupling (Fig. 33). Then, advance them together with the Tissue Protection Sleeve, Long, to the corresponding hole of the Targeting Arm, Antegrade (for left or right) by pressing the safety clip (Fig. 34). The mechanism will keep the sleeve in place and prevent it from falling out. It will also prevent the sleeve from sliding during screw measurement. To release the Tissue Protection Sleeve, the safety clip must be pressed again.

A small skin incision is made and the assembly is pushed through by manipulating the T-Handle until the Tissue Protection Sleeve is in contact with the lateral cortex (Fig. 34).



Fig. 32



Fig. 32a




Fig. 33



Fig. 34

Operative Technique

Pre-drilling the lateral cortex

Pre-drilling offers a possibility to open the lateral cortex for the drill entry. Pre-drilling helps to prevent a possible slipping of the drill on the cortex and may avoid deflection within the cancellous bone. This helps to perform the following drilling procedure without nail contact.

The Paddle Trocar Assembly is then removed and the Drill Sleeve is inserted through the Tissue Protection Sleeve, Long, (Fig. 35). With the Tissue Protection Sleeve, Long firmly engaged in the cortex, the lateral cortex should be opened using the centered tip green coded 4.2mm Drill.

The Drill can be connected with the Teardrop Handle, AO Coupling allowing pre-drilling by hand (Fig. 36). It also can be done by power.

Note:

For optimal stability, the tip of the oblique screw should be positioned at the level of the Lesser Trochanter (Fig. 37).

Then use the center-tipped, calibrated $\text{Ø}4.2 \times 340\text{mm}$ Drill and drill through both cortices (Fig. 38).

The screw length may be read directly from the Calibrated Drill, at the end of the Drill Sleeve (Fig. 38a).

Note:

Start the drill before touching the bone and then keep gentle pressure on the pre-drilled cortex to ensure accurate drilling.

Note:

The position of the end of the Drill, as it relates to the far cortex, is the same as where the end of the screw will be.




Fig. 35




Fig. 36




Fig. 37




Fig. 38



Fig. 38a

Operative Technique

Therefore, if the end of the Drill is 3mm beyond the far cortex, the end of the screw will also be 3mm beyond (Fig. 39). **Check the position of the end of the Drill with image intensification before measuring the screw length.** If Screw measurement with the Screw Gauge, Long, is preferred, first remove the Drill Sleeve, Long and read the screw length directly at the end of the Tissue Protection Sleeve, Long.

Note:

The Screw Gauge, Long is calibrated so that with the bend at the end pulled back flush with the far cortex, the screw tip will end 3mm beyond the far cortex (Fig. 40).

When the Drill Sleeve is removed, the correct Locking Screw is inserted through the Tissue Protection Sleeve using the Long Screwdriver Shaft with Teardrop Handle (Fig. 41). The screw is advanced through both cortices. The screw is near its proper seating position when the groove around the shaft of the screwdriver is approaching the end of the Tissue Protection Sleeve (Fig. 41a).



Fig. 39



Fig. 40




Fig. 41a




Fig. 41

Operative Technique

Freehand Distal Locking

The freehand technique is used to insert Fully Threaded Locking Screws into both distal transverse holes in the nail.

Rotational alignment must be checked prior to locking the nail. This is performed by checking a lateral view at the hip and a lateral view at the knee. The anteversion should be the same as on the controlateral side.

Multiple locking techniques and radiolucent drill devices are available for freehand locking. The critical step with any freehand locking technique, proximal or distal, is to visualize a perfectly round locking hole with the C-Arm.

The center-tipped $\text{ø}4.2 \times 180$ Drill is held at an oblique angle to the center of the locking hole (Fig. 43). Upon X-Ray verification, the Drill is placed perpendicular to the nail and drilled through the lateral and medial cortex (Fig. 44). Confirm in both the A/P and Lateral views by X-Ray that the Drill passes through the hole in the nail.

After drilling both cortices, the screw length may be read directly off of the Long Screw Scale at the green ring on the center-tipped 4.2×180 Drill (Fig. 45).

Alternatively, the Screw Gauge for Freehand technique can be used instead of the Screw Scale, Long, to determine the screw length.

Routine Locking Screw insertion is employed with the assembled Long Screwdriver Shaft and Teardrop Handle.

Note:

The Screwdriver Shaft can be used in conjunction with the Long Screw Capture Sleeve.

Repeat the locking procedure for the insertion of the second 5mm Fully Threaded Locking Screw into the oblong hole, in a static position (Fig. 47).

The **T2 Recon Nail** may be used in the **dynamic locking mode**. Only when the fracture pattern permits, dynamic




Fig. 43




Fig. 44




Fig. 45




Fig. 47

locking may be utilized for transverse, rotationally stable fractures. While dynamic locking can only be performed at the end of the nail, this will require a freehand distal targeting of the oblong hole in a dynamic position. This allows the nail to move and the

fracture to settle while torsional stability is maintained.

Operative Technique

Set Screw or End Cap Insertion

After removal of the Target Device, a Set Screw or End Cap can be used.

Two different Set Screws are available (Fig.48a):

- a **Set Screw, Recon** to tighten down on the Proximal Lag Screw for the Recon Mode
- a **Set Screw, Antegrade** to tighten down on the oblique Fully Threaded Screw for the Femoral Antegrade Mode

Note:

If a Set Screw is used, an End Cap can no longer be inserted.

Four different sizes of End Caps are available to adjust nail length and to reduce the potential for bony ingrowth into the proximal thread of the nail (Fig. 48b).

The Set Screw or End Cap is inserted with the Long Screwdriver Shaft and Teardrop Handle after intra-operative radiographs confirm satisfactory reduction and hardware implantation (Fig. 49). Be sure to fully seat the End Cap or Set Screw to minimize the potential risk for loosening.

Nail Removal

Nail removal is an elective procedure. The Set Screw or End Cap is removed with the Long Screwdriver Shaft and Teardrop Handle (Fig. 50).

The Universal Rod is inserted into the driving end of the nail. Alternatively, the Extraction Rod, conical, can be attached to the Universal Rod to facilitate extraction of the nail. All Locking Screws are removed with the Long Screwdriver Shaft and Teardrop Handle. The optional Long Screw Capture Sleeve may be used on the Screwdriver Shaft. For removal of the Lag Screws, the Recon Screwdriver or the Recon Screwdriver Shaft and T-Handle are to be used.

The Slotted Hammer is used to extract the nail in a controlled manner (Fig.51).



Fig. 48a



Fig. 48b




Fig. 49




Fig. 50




Fig. 51

A captured Sliding Hammer is available as an optional addition to the dedicated instrument set.


Ordering Information - Implants

T2 Recon Nail R1.5, Left



| Titanium REF | Diameter mm | Length mm |
|--------------|-------------|-----------|
| 1843-0926S | 9.0 | 260 |
| 1843-0928S | 9.0 | 280 |
| 1843-0930S | 9.0 | 300 |
| 1843-0932S | 9.0 | 320 |
| 1843-0934S | 9.0 | 340 |
| 1843-0936S | 9.0 | 360 |
| 1843-0938S | 9.0 | 380 |
| 1843-0940S | 9.0 | 400 |
| 1843-0942S | 9.0 | 420 |
| 1843-1026S | 10.0 | 260 |
| 1843-1028S | 10.0 | 280 |
| 1843-1030S | 10.0 | 300 |
| 1843-1032S | 10.0 | 320 |
| 1843-1034S | 10.0 | 340 |
| 1843-1036S | 10.0 | 360 |
| 1843-1038S | 10.0 | 380 |
| 1843-1040S | 10.0 | 400 |
| 1843-1042S | 10.0 | 420 |
| 1843-1126S | 11.0 | 260 |
| 1843-1128S | 11.0 | 280 |
| 1843-1130S | 11.0 | 300 |
| 1843-1132S | 11.0 | 320 |
| 1843-1134S | 11.0 | 340 |
| 1843-1136S | 11.0 | 360 |
| 1843-1138S | 11.0 | 380 |
| 1843-1140S | 11.0 | 400 |
| 1843-1142S | 11.0 | 420 |
| 1843-1226S | 12.0 | 260 |
| 1843-1228S | 12.0 | 280 |
| 1843-1230S | 12.0 | 300 |
| 1843-1232S | 12.0 | 320 |
| 1843-1234S | 12.0 | 340 |
| 1843-1236S | 12.0 | 360 |
| 1843-1238S | 12.0 | 380 |
| 1843-1240S | 12.0 | 400 |
| 1843-1242S | 12.0 | 420 |
| 1843-1326S | 13.0 | 260 |
| 1843-1328S | 13.0 | 280 |
| 1843-1330S | 13.0 | 300 |
| 1843-1332S | 13.0 | 320 |
| 1843-1334S | 13.0 | 340 |
| 1843-1336S | 13.0 | 360 |
| 1843-1338S | 13.0 | 380 |
| 1843-1340S | 13.0 | 400 |
| 1843-1342S | 13.0 | 420 |

T2 Recon Nail R1.5, Right



| Titanium REF | Diameter mm | Length mm |
|--------------|-------------|-----------|
| 1845-0926S | 9.0 | 260 |
| 1845-0928S | 9.0 | 280 |
| 1845-0930S | 9.0 | 300 |
| 1845-0932S | 9.0 | 320 |
| 1845-0934S | 9.0 | 340 |
| 1845-0936S | 9.0 | 360 |
| 1845-0938S | 9.0 | 380 |
| 1845-0940S | 9.0 | 400 |
| 1845-0942S | 9.0 | 420 |
| 1845-1026S | 10.0 | 260 |
| 1845-1028S | 10.0 | 280 |
| 1845-1030S | 10.0 | 300 |
| 1845-1032S | 10.0 | 320 |
| 1845-1034S | 10.0 | 340 |
| 1845-1036S | 10.0 | 360 |
| 1845-1038S | 10.0 | 380 |
| 1845-1040S | 10.0 | 400 |
| 1845-1042S | 10.0 | 420 |
| 1845-1126S | 11.0 | 260 |
| 1845-1128S | 11.0 | 280 |
| 1845-1130S | 11.0 | 300 |
| 1845-1132S | 11.0 | 320 |
| 1845-1134S | 11.0 | 340 |
| 1845-1136S | 11.0 | 360 |
| 1845-1138S | 11.0 | 380 |
| 1845-1140S | 11.0 | 400 |
| 1845-1142S | 11.0 | 420 |
| 1845-1226S | 12.0 | 260 |
| 1845-1228S | 12.0 | 280 |
| 1845-1230S | 12.0 | 300 |
| 1845-1232S | 12.0 | 320 |
| 1845-1234S | 12.0 | 340 |
| 1845-1236S | 12.0 | 360 |
| 1845-1238S | 12.0 | 380 |
| 1845-1240S | 12.0 | 400 |
| 1845-1242S | 12.0 | 420 |
| 1845-1326S | 13.0 | 260 |
| 1845-1328S | 13.0 | 280 |
| 1845-1330S | 13.0 | 300 |
| 1845-1332S | 13.0 | 320 |
| 1845-1334S | 13.0 | 340 |
| 1845-1336S | 13.0 | 360 |
| 1845-1338S | 13.0 | 380 |
| 1845-1340S | 13.0 | 400 |
| 1845-1342S | 13.0 | 420 |

Implants are packed sterile.

Ordering Information - Implants

5mm Fully Threaded Locking Screws



| Titanium REF | Diameter mm | Length mm |
|--------------|-------------|-----------|
| 1896-5025S | 5.0 | 25.0 |
| 1896-5030S | 5.0 | 30.0 |
| 1896-5035S | 5.0 | 35.0 |
| 1896-5040S | 5.0 | 40.0 |
| 1896-5045S | 5.0 | 45.0 |
| 1896-5050S | 5.0 | 50.0 |
| 1896-5055S | 5.0 | 55.0 |
| 1896-5060S | 5.0 | 60.0 |
| 1896-5065S | 5.0 | 65.0 |
| 1896-5070S | 5.0 | 70.0 |
| 1896-5075S | 5.0 | 75.0 |
| 1896-5080S | 5.0 | 80.0 |
| 1896-5085S | 5.0 | 85.0 |
| 1896-5090S | 5.0 | 90.0 |
| 1896-5095S | 5.0 | 95.0 |
| 1896-5100S | 5.0 | 100.0 |
| 1896-5105S | 5.0 | 105.0 |
| 1896-5110S | 5.0 | 110.0 |
| 1896-5115S | 5.0 | 115.0 |
| 1896-5120S | 5.0 | 120.0 |

6.5mm Lag Screws



| Titanium REF | Diameter mm | Length mm |
|--------------|-------------|-----------|
| 1897-6065S | 6.5 | 65 |
| 1897-6070S | 6.5 | 70 |
| 1897-6075S | 6.5 | 75 |
| 1897-6080S | 6.5 | 80 |
| 1897-6085S | 6.5 | 85 |
| 1897-6090S | 6.5 | 90 |
| 1897-6095S | 6.5 | 95 |
| 1897-6100S | 6.5 | 100 |
| 1897-6105S | 6.5 | 105 |
| 1897-6110S | 6.5 | 110 |
| 1897-6115S | 6.5 | 115 |
| 1897-6120S | 6.5 | 120 |
| 1897-6125S | 6.5 | 125 |
| 1897-6130S | 6.5 | 130 |

End Caps



| Titanium REF | Diameter mm | Length mm |
|--------------|-------------|-----------|
| 1822-0003S | 8.0 | Standard |
| 1847-0005S | 13.0 | + 5mm |
| 1847-0010S | 13.0 | +10mm |
| 1847-0015S | 13.0 | +15mm |

Set Screws



| Titanium REF | Diameter mm | |
|--------------|-------------|----------------------|
| 1847-0003S | 8.0 | Set Screw, Antegrade |
| 1847-0001S | 8.0 | Set Screw, Recon |


Implants are packed sterile.

Ordering Information - Instruments

| | REF | Description |
|---|------------|---|
|  | 1806-0020 | Guide Wire Ruler |
|  | 1806-0040 | Awl Curved |
|  | 1806-0073 | Universal Rod |
|  | 1806-0085S | Guide Wire, Ball Tip, Ø3 x 1000mm, sterile* |
|  | 1806-1095 | Guide Wire Handle, blue coded |
|  | 1806-1096 | Guide Wire Handle Chuck, blue coded |
|  | 1806-0125 | Reduction Spoon |
|  | 1806-0130 | Wrench |
|  | 1806-0170 | Slotted Hammer |
|  | 1806-0185 | Tissue Protection Sleeve, Long |
|  | 1806-0215 | Drill Sleeve, Long |
|  | 1806-0227 | Screwdriver Shaft AO, Long |
|  | 1806-0232 | Screwdriver, Long |
|  | 1806-0240 | Screw Capture Sleeve, Long |
|  | 1806-0292 | Screw Driver Shaft, 3.5 x 85mm |
|  | 1806-0325 | Screw Gauge, Long |
|  | 1806-0350 | Extraction Rod, conical |
|  | 1806-0365 | Screw Scale, Long |
|  | 1806-0480 | Screw Gauge, Femur |
|  | 1806-3047 | T2 Paddle Trocar Recon Mode |
|  | 1806-3048 | T2 Paddle Trocar Antegrade Mode |
|  | 1806-4260S | Drill Ø4.2 x 340mm, AO* |
|  | 1806-4270S | Drill Ø4.2 x 180mm, AO* |
|  | 702429 | Teardrop Handle, AO Coupling |
|  | 702628 | Teardrop Handle, AO Medium Coupling |

* For non-sterile, leave "S" off the REF number when ordering.


Ordering Information - Instruments



| REF | Description |
|------------|------------------------------------|
| 1806-1007 | Fixation Screw |
| 1806-3001 | Nail Adapter, Recon |
| 1806-3002 | Targeting Arm, Recon |
| 1806-3003 | Targeting Arm, antegrade |
| 1806-3005 | Nail Holding Screw, Recon |
| 1806-3010 | One Step Conical Reamer Ø13, Recon |
| 1806-3015 | One Step Conical Reamer Ø15, Recon |
| 1806-3026S | Solid Stepdrill for Lag Screw* |
| 1806-3030S | ReconK-Wire, Recon* |
| 1806-3031S | K-Wire, Recon, CoCr |
| 1806-3035 | Lag Screw Gauge, Recon |
| 1806-3040 | K-Wire Sleeve, Recon |
| 1806-3041 | Drill Sleeve for Solid Stepdrill |
| 1806-3045 | Tissue Protection Sleeve, Recon |
| 1806-3050 | Screw Driver Shaft, Recon |
| 1806-3055 | Multi-hole Trocar |
| 1806-3057 | Protection Sleeve, Antegrade |
| 1806-3060 | Screw Driver, Recon |
| 1806-3070 | K-Wire Inserter |
| 1806-3090 | Screw Driver Shaft, AO, Ball Tip |
| 1806-3096 | Strike Plate, Recon |
| 1806-8018 | Drill Ø4.2×250mm, oblique, AO |
| 1806-3082 | X-Ray Template, Recon R1.5 |
| 1806-9400 | Instrument Tray, empty |
| 1806-9410 | Add-on Instrument Tray, empty |

* For non-sterile, leave "S" off the REF number when ordering.


Ordering Information - Instruments



| REF | Description |
|-----------------------------|--|
| Optional Instruments | |
| 1213-3010 | One Shot Device |
| 1806-0032 | Awl Plug |
| 1806-0041 | Awl, Curved, Ø10mm, 90° Handle (optional,not stored on tray) |
| 1806-0315 | Trocar, Long |
| 1806-3020 | T2 Fixation Wrench Recon |
| 1806-3025 | Stepdrill for Lag Screw, Recon |
| 1806-3065 | Extraction Screw Driver |
| 1806-0233 | Selfholding Screwdriver, long |
| 1806-0203 | Selfholding Screwdriver, extra short |

Ordering Information - Instruments

Bixcut



Complete range of modular and fixed-head reamers to match surgeon preference and optimize O.R. efficiency, presented in fully sterilizable cases.

Large clearance rate resulting from reduced number of reamer blades coupled with reduced length of reamer head to give effective relief of pressure and efficient removal of material.


Cutting flute geometry optimized to lower pressure generation.

Forward- and side-cutting face combination produces efficient material removal and rapid clearance.

Double-wound shaft transmits torque effectively and with high reliability. Low-friction surface finish aids rapid debris clearance.


Smaller, 6 and 8mm shaft diameters significantly reduce IM pressure.

Typical Standard
Reamer Ø14mm



Clearance area:
32% of cross section


Bixcut
Reamer Ø14mm



Clearance area:
59% of cross section

Recent studies¹ have demonstrated that the pressures developed within the medullary cavity through the introduction of unreamed IMnails can be far greater than those developed during reaming – but this depends very much upon the design of the reamer.

After a three year development study² involving several universities, the factors that determine the pressures and temperatures developed during reaming were clearly established. These factors were applied to the development of advanced reamers that demonstrate significantly better performance than the best of previous designs.



¹ Jan Paul M. Frolke, et al.; Intramedullary Pressure in Reamed Femoral Nailing with Two Different Reamer Designs, Eur. J. of Trauma, 2001 #5

² Mehdi Mousavi, et al.; Pressure Changes During Reaming with Different Parameters and Reamer Designs, Clinical Orthopaedics and Related Research Number 373, pp. 295-303, 2000

Ordering Information – Instruments

Bixcut Modular Head

| REF | Description | Diameter mm |
|-----------|-------------|-------------|
| 0226-3090 | Bixcut Head | 9.0 |
| 0226-3095 | Bixcut Head | 9.5 |
| 0226-3100 | Bixcut Head | 10.0 |
| 0226-3105 | Bixcut Head | 10.5 |
| 0226-3110 | Bixcut Head | 11.0 |
| 0226-3115 | Bixcut Head | 11.5 |
| 0226-3120 | Bixcut Head | 12.0 |
| 0226-3125 | Bixcut Head | 12.5 |
| 0226-3130 | Bixcut Head | 13.0 |
| 0226-3135 | Bixcut Head | 13.5 |
| 0226-3140 | Bixcut Head | 14.0 |
| 0226-3145 | Bixcut Head | 14.5 |
| 0226-3150 | Bixcut Head | 15.0 |
| 0226-3155 | Bixcut Head | 15.5 |
| 0226-3160 | Bixcut Head | 16.0 |
| 0226-3165 | Bixcut Head | 16.5 |
| 0226-3170 | Bixcut Head | 17.0 |
| 0226-3175 | Bixcut Head | 17.5 |
| 0226-3180 | Bixcut Head | 18.0 |
| 0226-4185 | Bixcut Head | 18.5 |
| 0226-4190 | Bixcut Head | 19.0 |
| 0226-4195 | Bixcut Head | 19.5 |
| 0226-4200 | Bixcut Head | 20.0 |
| 0226-4205 | Bixcut Head | 20.5 |
| 0226-4210 | Bixcut Head | 21.0 |
| 0226-4215 | Bixcut Head | 21.5 |
| 0226-4220 | Bixcut Head | 22.0 |
| 0226-4225 | Bixcut Head | 22.5 |
| 0226-4230 | Bixcut Head | 23.0 |
| 0226-4235 | Bixcut Head | 23.5 |
| 0226-4240 | Bixcut Head | 24.0 |
| 0226-4245 | Bixcut Head | 24.5 |
| 0226-4250 | Bixcut Head | 25.0 |
| 0226-4255 | Bixcut Head | 25.5 |
| 0226-4260 | Bixcut Head | 26.0 |
| 0226-4265 | Bixcut Head | 26.5 |
| 0226-4270 | Bixcut Head | 27.0 |
| 0226-4275 | Bixcut Head | 27.5 |
| 0226-4280 | Bixcut Head | 28.0 |

Bixcut Fixed Head – AO Fitting**

| REF | Diameter mm | Length mm |
|-----------|-------------|-----------|
| 0225-5060 | 6.0* | 400 |
| 0225-5065 | 6.5* | 400 |
| 0225-5070 | 7.0* | 400 |
| 0225-6075 | 7.5 | 480 |
| 0225-6080 | 8.0 | 480 |
| 0225-6085 | 8.5 | 480 |
| 0225-6090 | 9.0 | 480 |
| 0225-6095 | 9.5 | 480 |
| 0225-6100 | 10.0 | 480 |
| 0225-6105 | 10.5 | 480 |
| 0225-6110 | 11.0 | 480 |
| 0225-8115 | 11.5 | 480 |
| 0225-8120 | 12.0 | 480 |
| 0225-8125 | 12.5 | 480 |
| 0225-8130 | 13.0 | 480 |
| 0225-8135 | 13.5 | 480 |
| 0225-8140 | 14.0 | 480 |
| 0225-8145 | 14.5 | 480 |
| 0225-8150 | 15.0 | 480 |
| 0225-8155 | 15.5 | 480 |
| 0225-8160 | 16.0 | 480 |
| 0225-8165 | 16.5 | 480 |
| 0225-8170 | 17.0 | 480 |
| 0225-8175 | 17.5 | 480 |
| 0225-8180 | 18.0 | 480 |

Bixcut Fixed Head – Modified Trinkle fitting+

| REF | Diameter mm | Length mm |
|-----------|-------------|-----------|
| 0227-5060 | 6.0* | 400 |
| 0227-5065 | 6.5* | 400 |
| 0227-5070 | 7.0* | 400 |
| 0227-6075 | 7.5 | 480 |
| 0227-6080 | 8.0 | 480 |
| 0227-6085 | 8.5 | 480 |
| 0227-6090 | 9.0 | 480 |
| 0227-6095 | 9.5 | 480 |
| 0227-6100 | 10.0 | 480 |
| 0227-6105 | 10.5 | 480 |
| 0227-6110 | 11.0 | 480 |
| 0227-8115 | 11.5 | 480 |
| 0227-8120 | 12.0 | 480 |
| 0227-8125 | 12.5 | 480 |
| 0227-8130 | 13.0 | 480 |
| 0227-8135 | 13.5 | 480 |
| 0227-8140 | 14.0 | 480 |
| 0227-8145 | 14.5 | 480 |
| 0227-8150 | 15.0 | 480 |
| 0227-8155 | 15.5 | 480 |
| 0227-8160 | 16.0 | 480 |
| 0227-8165 | 16.5 | 480 |
| 0227-8170 | 17.0 | 480 |
| 0227-8175 | 17.5 | 480 |
| 0227-8180 | 18.0 | 480 |

Bixcut Shafts (Sterile)^{1,2,3,4}, Shaft Accessories and Tray, empty

| REF | Description | Length mm |
|------------|---|-----------|
| 0227-8240S | Mod. Trinkle | 284 |
| 0227-3000S | Mod. Trinkle | 448 |
| 0227-8510S | Mod. Trinkle | 510 |
| 0227-8885S | Mod. Trinkle | 885 |
| 0226-8240S | AO | 284 |
| 0226-3000S | AO | 448 |
| 3212-0-210 | Grommet (pack of 25) | |
| 3212-0-220 | Grommet inserter/extractor | |
| 0225-6010 | Grommet Case | |
| 0225-6000 | Tray, Modular Head (up to size 22.0mm) | |
| 0225-6001 | Tray, Modular Head (up to size 28.0mm) | |
| 0225-8000 | Tray, Fixed Head (up to size 18.0mm) | |
| 0225-6040 | Mini Trauma Tray (for modular heads 9-18) | |
| 0225-6050 | Mini Revision Tray (for modular heads 9-28) | |

Optional Instruments

| REF | Description |
|------------|--|
| 5235-6-606 | Hand Reamer 6 mm w/T-Handle |
| 5235-6-607 | Hand Reamer 7 mm w/T-Handle |
| 5235-6-608 | Hand Reamer 8 mm w/T-Handle |
| 5235-6-609 | Hand Reamer 9 mm w/T-Handle |
| 0227-0060 | Hand Reamer 6 mm w/Mod Trinkle connection |
| 0227-0070 | Hand Reamer 7 mm w/Mod Trinkle connection |
| 0227-0080 | Hand Reamer 8 mm w/Mod Trinkle connection |
| 0227-0090 | Hand Reamer 9 mm w/Mod Trinkle connection |
| 1806-6520 | Curved Reduction Rod 8.5 mm w/Mod Trinkle connection |
| 1806-6500 | T-Handle w/Mod Trinkle connection |

* Use with 2.2mm × 800mm Smooth Tip and 2.5mm × 800mm Ball Tip Guide Wires only.

** Use with Stryker Power Equipment

1.Non-Sterile shafts supplied without grommet. Use new grommet for each surgery. See Shaft Accessories.

2.Sterile shafts supplied with grommet pre-assembled.

3.For Non-Sterile leave “S” off the REF Number when ordering (510 and 885mm available only sterile Modified Trinkle Fitting).

4.Non-Sterile, AO Fitting Shafts in 510 and 885mm are available as build to order items:

- CM810921 AO Fitting Shaft, length 510mm
- CM810923 AO Fitting Shaft, length 885mm

Notes

Joint Replacements

Trauma, Extremities & Deformities

Craniomaxillofacial

Spine

Biologics

Surgical Products

Neuro & ENT

Interventional Spine

Navigation

Endoscopy

Communications

Imaging

Patient Care & Handling Equipment

EMS Equipment

Stryker Trauma GmbH
Prof.-Küntschers-Strasse 1-5
D-24232 Schönkirchen
Germany

www.osteosynthesis.stryker.com

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery. The information presented in this brochure is intended to demonstrate a Stryker product. Always refer to the package insert, product label and/or user instructions including the instructions for Cleaning and Sterilization (if applicable) before using any Stryker products. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker, Bixcut, T2, Gamma3, One Shot Device.

All other trademarks are trademarks of their respective owners or holders.
The products listed above are CE marked.

Literature Number: B1000029
LOT C1509

Copyright © 2009 Stryker

