Reprocessing Guide

stryker[®]

Shaver Handpiece Tray

REF 0272700000

CE R_X ONLY

Contents

User Guide	EN-1
Guide de l'utilisateur	FR-13
Benutzerhandbuch	DE-25
Manuale d'uso	IT-39
Manual do Utilizador	PT-51
Guía del usuario	ES-63
Gebruikershandleiding	NL-75
Brugervejledning	DA-87
Käyttöohje	FI-99
Brukerhåndbok	NO-111
Användarhandbok	SV-123
Podręcznik użytkownika	PL-135
Οδηγός χρήσης	ES-147
用户指南	CHS-159
ユーザーガイド	JA-169
사용 설명서	KO-181

Contents

Introduction EN-2

Intended Use of Sterilization Trays EN-2

Warnings EN-3

Cautions EN-3

Instructions EN-4

Sterilization Tray Setup EN-10

References EN-11

Introduction

This reprocessing guide provides instructions for the proper cleaning and sterilization of the following sterilization tray:

272-700-000 Shaver Handpiece Tray

This sterilization tray is intended for use **only** with the devices listed below. For tray/device configuration, see the “Sterilization Tray Setup” section of this guide.

Part Number	Description
375-708-500	Formula 180 Shaver Handpiece (with buttons)
375-704-500	Formula Shaver Handpiece (with buttons)
375-701-500	Formula Shaver Handpiece (no buttons)
275-601-500	TPS Small Joint Shaver

While these instructions have been validated by the manufacturer as being capable of preparing the tray for re-use, it remains the responsibility of the processor to ensure that the reprocessing (as actually performed, using equipment, materials, and personnel in the reprocessing facility) achieves the desired result. This normally requires validation and routine monitoring of the process.

Intended Use of Sterilization Trays

Sterilization trays are plastic and/or metal containers used to hold and protect surgical devices during the sterilization process. They consist of an interlocking tray and lid, which are both perforated to allow the passage of sterilizing agent from outside the tray to the devices placed inside.

Sterilization trays typically feature a silicon finger mat or group of device holders that secure devices during the sterilization process. Some models feature stacking internal trays to allow the segregation of devices.

Warnings

- These instructions are validated only for sterilization of the tray and device(s) identified herein. Using combinations or parameters not described in this manual may result in incomplete sterilization.
- These instructions do not replace the cleaning instructions provided with individual devices. Prior to sterilization, clean all devices as specified in their respective user manuals.
- Wear appropriate protective equipment (gloves, eye protection, etc.) when reprocessing any medical device.
- Both the device and the tray must be cleaned prior to sterilization, or incomplete sterilization will result.
- The sterilization tray, its lid, and any internal components have been designed and validated for use as a single system. Do not separate components from the system, for use individually or in combination, or incomplete sterilization may result.
- Inspect the tray and its components for visible damage, such as cracking or chipping, prior to use. Do not use the tray if it is damaged.

Cautions

- Before lifting the tray assembly, verify that the latches connecting the lid to the tray are secure.
- The tray is not designed for use as a shipping container. To avoid damage, remove all devices and pack them separately.

Limitations on Reprocessing

- Proper processing has a minimal effect on the tray. End of life is normally determined by wear and damage due to use.
- Do not cross-sterilize the tray. Using multiple sterilization methods may significantly reduce the performance of the tray.
- Do not leave the tray in solutions longer than necessary. This may accelerate normal product aging.
- Damage incurred by improper processing will not be covered by the warranty.

Instructions

Instructions apply to the **sterilization tray only**. For instructions on how to reprocess the devices, consult their respective instructions for use.

Point of Use	• Wipe excess soil from the tray using disposable paper towels. • If an automated cleaning method will be used, rinse the tray in sterile distilled water immediately after use.
Containment and Transportation	Reprocess device(s) / tray as soon as reasonably practical following use.
Preparation for Manual Cleaning ¹ Cleaning was validated using Enzo [®] at 1 oz/gal. at 35–40°C.	1 Disassemble the tray into its individual components: base, lid, and (if applicable) silicon mat. 2 Prepare an enzymatic detergent¹ according to the manufacturer's recommendations. 3 Wipe the entire tray with a clean cloth dipped in detergent. 4 Immerse the tray in the detergent. 5 Soak the tray in the detergent for 15 minutes.
Manual Cleaning	1 Brush • Thoroughly brush the exterior of the tray with a soft-bristled brush, focusing on any mated or rough surfaces. • Using a pipe cleaner and brush, brush and flush the hard-to-reach areas of the tray. • Actuate the tray, brushing any movable parts in their extreme open and closed positions.

2 Rinse

- Rinse the tray with reverse osmosis/deionized (RO/DI) water at ambient temperature to remove all detergent residues. After all detergent residues have been removed, continue to rinse for a minimum of 30 seconds.
- Drain excess water from the tray and dry it using a clean cloth and filtered pressurized air (40 psi).
- Visually inspect the tray for cleanliness, paying close attention to hard-to-reach areas. If visible soil remains, repeat steps 1 and 2.

3 Soak

- Prepare a non-enzymatic detergent² according to the manufacturer's recommendations.
- Fully immerse the tray and inject any mated surfaces with at least 50mL of the detergent.
- Soak the tray for a minimum of 15 minutes.

4 Brush

- Thoroughly brush the exterior of the tray using a soft-bristled brush.
- Inject the prepared detergent into any mated surfaces a minimum of five times.
- Actuate the tray, brushing any movable parts in their extreme open and closed positions.
- Using a pipe cleaner and syringe, brush and flush the hard-to-reach areas of the tray.

5 Rinse

- Thoroughly rinse the tray with ambient RO/DI water until all detergent residue is removed. Flush any crevices five times with a syringe. After the detergent residue is removed, continue rinsing for a minimum of 30 seconds.
- Drain the excess water from the tray and dry it with a clean cloth and filtered pressurized air (40 psi).
- Visually inspect the tray for cleanliness, paying close attention to hard-to-reach areas.

² Cleaning was validated using Renu-Klenz™ at ¼ oz/gal. at 35-40°C.

Automated Cleaning

1 Rinse

- Disassemble the tray into its individual components: the base, lid, and (if applicable) silicon mat.
- Rinse the tray with RO/DI water at ambient temperature until there is no visible detergent residue. Use a syringe to assist in rinsing. After all detergent residues have been removed, continue to rinse for a minimum of 30 seconds.
- Place the tray in the washer on an incline to facilitate drainage.

2 Automated wash

- Program the washer using the settings below.

	Recirculation Time	Water Temperature	Detergent (type and concentration)
Pre Wash	2 min.	Cold tap water	—
Enzyme Wash	2 min.	Hot tap water	Enzymatic Detergent ³
Wash 1	2 min.	Set point 60°C (140°F)	Non-enzymatic Detergent ⁴
Rinse 1	2 min.	Heated 60°C (140°F)	—
Dry Phase	7 min.	Drying Temperature 115°C (239°F)	—

3 Drying

- Remove the tray from the washer after the completion of the dry phase.
- Dry the tray using a clean soft cloth. Use filtered pressurized air (40 psi) to assist in drying.
- Visually inspect the tray for cleanliness, paying close attention to hard-to-reach areas.

³ Cleaning was validated using Enzol® at 1 oz/gal.

⁴ Cleaning was validated using Renu-Klenz™ at ¼ oz/gal.

Maintenance	No routine maintenance is required for the sterilization tray.		
Inspection and Testing	- Inspect the tray for dents and cracks. If a problem is observed or suspected, the tray should be returned for repair. - Inspect all components for cleanliness. If fluid or tissue buildup is present, repeat the above cleaning and disinfection procedures.		
Packaging for Sterilization	Double-wrap the tray unless otherwise noted below.		
Sterilization	- Prepare the device(s) and the sterilization tray as indicated in the “Sterilization Tray Setup” section below. - Sterilize the device(s) inside the tray according to the parameters indicated below.		
	Ethylene oxide (EtO)		
	Preconditioning	Sterilization	Aeration
	Temperature	55°C (131°F)	55°C (131°F)
	Relative Humidity	70%	—
	Vacuum Set Points	1.30 psia	—
	EtO Concentration	—	600 mg/L (100% EtO)
	Time	60 minutes	240 minutes
			12 hours

STERRAD®

Follow the manufacturer's instructions to sterilize the device(s) and tray. Select the appropriate cycle.

STERRAD® 100S and STERRAD® NX

375-708-500 Formula 180
Shaver Handpiece

375-704-500 Formula
Shaver Handpiece

375-701-500 Formula
Shaver Handpiece

STERRAD® Advanced NX

275-601-500 TPS Small Joint
Shaver

Autoclave (Steam)

	Gravity	Prevacuum	"Flash" Prevacuum
Wrapping	double	double	—
Temperature	132°C (270°F)	132°C (270°F)	132°C (270°F)
Time	20 minutes	4 minutes	4 minutes
Dry Time	60 minutes	50 minutes	—

⚠ Warning: Drying time depends on several variables, including: altitude, humidity, type of wrap, preconditioning, size of chamber, mass of load, material of load, and placement in chamber. Users must verify that drying time set in their autoclave yields dry surgical equipment.

⚠ Warning: These instructions were developed using the guidance from AAMI TIR 12, ISO 17665 and AAMI ST79 and Stryker recommends users observe these standards.

Storage	Trays with non-absorbent liners (such as plastic or silicone-fingered organizing mats) can cause condensate to pool. To preserve sterility, remove devices from standing solution.
Additional Information	N/A

Sterilization Tray Setup

Maximum Weight Load (devices and tray combined)	1.52 kg (3.35 lbs)
Internal Stacking	No internal stacking is permitted with this tray.
External Stacking	Do not stack other trays or devices on or below this tray.
Accessories	There are no accessories available for use with this tray.
Device Distribution	Devices should be placed in the tray as illustrated below.
Chemical Indicator* 	Place the chemical indicator () in the bottom-left corner of the tray as illustrated in the following picture.
Biological Indicator* 	Place the biological indicator () on device tip as illustrated in the following picture.

*Note: Placement of the biological or chemical indicator is not required during sterilization of this tray. If, per hospital procedure, placement of an indicator is desired, the recommended placement sites are as follows.

References

- **ISO 17664**—Sterilization of medical devices—Information to be provided by the manufacturer for the processing of resterilizable medical devices
- **ISO 17665-1**—Sterilization of health care products—Moist heat—Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices
- **ANSI/AAMI ST77**—Containment devices for reusable medical device sterilization
- **ANSI/AAMI ST79**—A comprehensive guide to steam sterilization and sterility assurances in health care facilities
- **ANSI/AAMI ST81**—Sterilization of medical devices—Information to be provided by the manufacturer for the processing of resterilizable medical devices
- **AAMI TIR12**—Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: a guide for device manufacturers

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European Representative
Regulatory Manager, Stryker France
ZAC Satolas Green Pusignan
Ave. De Satolas Green
69881 MEYZIEU Cedex, France

P10299 B
2014/02

Stryker® is a registered trademark of Stryker Corporation.

ENZOL® Enzymatic Detergent is a registered trademark of Advanced Sterilization Products, Division of Ethicon Inc., a Johnson and Johnson company.

Renu-Klenz™ is a trademark of STERIS corporation.

STERRAD® is a registered trademark of Advanced Sterilization Products.