

Integrations

IHE Integration Statement

Vendor	Product Name	Version	Date
Stryker Imaging	OfficePACS® Power™	4.0	31 Jan 2007
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Radiology			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Image Manager		
	Image Archive		
	Performed Procedure Step Manager		
Patient Information Reconciliation	Image Manager		
	Image Archive		
	Performed Procedure Step Manager		
Access to Radiology Information	Image Manager		
	Image Archive		
IT Infrastructure			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Consistent Time	Time Server		
	Time Client		
Personnel White Pages	DNS Server		
Internet address for vendor's IHE information:			
Copies of Standards Conformance Statements for the Implementation			
HL7	Please contact Stryker Imaging for copies of HL7 and DICOM Conformance Statements		
DICOM			
Links to general information on IHE			
In North America www.rsna.org/IHE	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	