

InTouch™

Critical Care Bed

Medical

Multi-function foot end touch screen provides ability to collect, display and process patient and bed status information

iBed™ Awareness monitors smart bed positions, providing local or remote alarms if pre-set bed configurations change

Only-in-class, three-zone Chaperone® Bed Exit system monitors patient position, alerting changes via local or remote alarm, and includes spoken patient commands

Protocol Reminders™ help organize and promote compliance with proven best practices

Spoken Language Translations help improve communication with Limited English Proficiency (LEP) patients

On-board Sound Therapy helps promote a soothing environment for patients

Support surfaces can be changed as therapies dictate, and Stryker's powered XPRT™ and PositionPRO® surfaces can be directly controlled by the touch screen

Industry's first USB upgradable bed allows InTouch to be upgraded as best practices evolve over time



A Higher Standard of Care Through Protocol Compliance

InTouch represents a higher standard of care in the ICU, helping to make protocol compliance easier for caregivers. All bed functions are controlled with an easy, intuitive foot end touch screen interface using a powerful on-bed computer. The *InTouch* bed allows you to set Protocol Reminders™ to help make your busy schedule simple, and document each therapy you perform to help reduce paperwork. With a touch, you can see a patient weight history graph of the last 50 weights, and *InTouch* even leads you through each step to assign a Braden Scale assessment to your patient. Pre-recorded clinical phrases in 12 languages help you communicate with Limited English Proficiency patients while Sound Therapy allows you to play music or nature sounds as part of a soothing environment. Real-time information is available from all *InTouch* beds in the ICU, as part of the Connected Hospital® from Stryker.



iBed Awareness alerts you to changes in preset smart bed positions while Protocol Reminders help ensure that best practices, such as patient turns, are done promptly and consistently.



Zoom on *InTouch* features front wheel drive for better maneuverability and a 3.6 mph maximum speed for efficient transport.



Through the intuitive touch-screen interface, caregivers can monitor bed status over time, as well as operate Stryker powered support surfaces.

InTouch™ Critical Care Bed

stryker®



Medical

3800 E. Centre Ave.
Portage, MI 49002 U.S.A.
t: 269 329 2100 f: 269 329 2311
toll free: 800 787 9537

www.stryker.com

Standard Features

- Touch screen interface
- Bed history screen
- Smart maintenance screen
- Event Manager™ screen
- Smart scale
- BackSmart Pivot in-bed patient positioning
- Chaperone three-zone bed exit alarm
- Electric brake with mechanical back-up
- CPR foot pedal and CPR-ready headboard
- Full functionality battery back-up
- Enhanced cardiac chair
- One-handed siderail release
- Siderail-storable patient control pendant
- Radiolucent litter top
- Patient lock-out controls
- Integrated pump holder
- Single 110-volt auxiliary outlet
- LED nightlight
- DC actuators
- Foley bag hooks

Specifications

• Model Number	2131 (2141 Zoom)
• Overall Length	90" (229 cm)
• Overall Width	
Siderails Down	40" (102 cm)
Siderails Up	42" (107 cm)
• Safe Working Load including all accessories	550 lbs (249 kg)
• Static Load Capacity (per IEC 60601-2-38)	1,100 lbs (499 kg)
• Height Range	
High	33" (84 cm)
Low	16" (41 cm)
• Litter Positioning	
Backrest	0° – 73°
Knee Gatch	0° – 17°
Trend./Reverse Trend.	±12°
• Patient Surface	35" x 84" (89 x 213 cm)
• Caster Diameter	5.5" (14 cm)

Electronics

• Agency Approvals	CSA
• Volts	120 VAC
• Ampere Ratings	10.0 A
• Current Leakage	< 300 microamperes
• Frequency	60 Hz

Optional Features

- iBed Awareness
- Protocol Reminders
- Calculator
- Converter
- Translator
- USB port
- Audio
- Documentation
- Zoom Motorized Drive System
- Oxygen bottle holders
 - Headboard-mounted
 - Upright
- Dual 110-volt auxiliary outlets
- Four-sided controls
- Bed extender and pad
- Pendant clip
- Transport tray/patient helper
- IV poles
 - 1 or 2 permanent
- Smart TV
- Room light controls
- Nurse call
- Pillow speaker
- X-ray cassette holder
- WallSaver™ Breakaway Wall Connector
- Foot end roller bumpers
- Line management accessory
- Service integration port

Copyright © 2008 Stryker
PC/SCG-7.5M-0508, 08-SMD-2262
Mkt Lit-180-042707 Rev D



Mixed Sources
Product group from well-managed
forests, controlled sources and
recycled wood or fibre
Cert no. SW-COC-1530
www.fsc.org
© 1996 Forest Stewardship Council