


stryker®

Medical



BackSmart™
Bariatric Solutions

One of a Kind

Stryker's BackSmart Bariatric Bed

Hospitals told us they need a bariatric bed that offers advanced ergonomics, is exceptionally easy to operate and is aesthetically pleasing. Stryker has the solution.



The Stryker Bari 10A offers the perfect blend of aesthetics, versatility, and simple operation.

Our new best-in-class Bari™ 10A bariatric bed meets and exceeds those requirements, including your number one need in a bariatric bed — the capacity to easily expand and retract without caregiver back strain.

With the simple touch of a lever, the Bari 10A bed fits through spaces as narrow as 39-inches yet expands to 48-inches to comfortably accommodate larger patients. And when combined with our other BackSmart mobility solutions, make bariatric patient handling and caregiving faster, safer, and easier.

Aesthetics Practically designed and quiet, the Bari 10A compares to a standard bed in appearance, meeting the aesthetic needs of modern health-care facilities and respecting the dignity of the patient.

Versatility The Bari 10A bed's easy expansion/retraction capabilities, in-bed scale, battery back-up and BackSmart ergonomics make it the easiest bariatric bed to use on the market. Coupled with Stryker's open architecture platform that accepts multiple surfaces, you'll find the Bari 10A enhances your ability to care for a challenging patient population.

Mobility The Bari 10A delivers the kind of superior mobility you've come to expect from Stryker. Its fifth wheel offers excellent traction and cornering, resulting in reduced back strain, improved efficiencies and increased throughput. Paired with Stryker's other bariatric solutions like our M-Series Stretchers, Stryker Glide™ transfer device and patient room seating, bariatric patient handling just got faster, safer, and easier.



SM304

With motorized Zoom® drive, Stryker's SM304 Stretcher helps reduce the high risk of caregiver back strain that comes with bariatric patient transport.



Bariatric Chair

Transitional styling and clean lines offer a classic look in bariatric seating support.



Stryker Glide

BackSmart design offers a smart solution to bariatric patient transfer with Stryker's lateral air transfer system.

Meet all of your bariatric patient mobility needs with Stryker BackSmart solutions.



Narrow Doorways Are No Problem

The Stryker Bariatric bed has a generous width for patient comfort, yet easily contracts to 39-inches to fit through standard, 48-inch doorways. Easy-to-use color-coded levers release the bed in sections so that even a single caregiver can change the width.



Bed Contracted



Left Side Extended



Both Sides Extended

Stryker Bariatric innovations help increase patient comfort and safety while making everyday tasks easier for caregivers.

Fifth Wheel Mobility

A center-mounted fifth wheel provides an additional pivot point to turn the bed in a tighter arc. The fifth wheel retracts to clear the under-bed space.

In-bed Scale

A built-in scale provides accurate patient weight on an easy-to-read LCD screen on the foot-board.

Battery Back-up

When the bed is unplugged for transport or in the event of power loss, motors and display continue to function.

Stryker Arise® 1000 EX

Made especially for Stryker's Bari 10A, the new Arise 1000 EX features custom sizing, low-air loss pulsation and simple operation. The Arise 1000 EX expands and contracts in length and width with inflatable air bladders.



Bariatric Patient Helper

Made with heavy-gauge, tubular steel, the Bariatric Patient Helper adds to patient dignity by allowing the patient to lift the upper body with arm strength. The trolley can travel the length of the cantilevered arms as well as side to side.



Fold-down Siderails

Large, robust hinges support the siderails that provide safety when up and fold down for easy access to the patient.

Lockout Control

For bariatric patients, easy breathing can depend on proper bed positioning. The pendant can be locked out to prevent changes in bed posture.

Flexible Width Adjustment

Color-coded grip-levers on each side of the bed allow sections of the bed to be easily expanded or contracted with a pull or push.

Joint Replacements

Trauma, Extremities & Deformities

Craniomaxillofacial

Spine

Biologics

Surgical Products

Neuro & ENT

Interventional Pain

Navigation

Endoscopy

Communications

Imaging

Patient Care & Handling Equipment

EMS Equipment

3800 E. Centre Ave.
Portage, MI 49002 USA
t: 269 329 2100 f: 269 329 2311
toll free: 800 787 9537

www.stryker.com

The information presented in this brochure is intended to demonstrate a Stryker product. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker Account Manager if you have questions about the availability of Stryker products in your area.

Products referenced with [™] designation are trademarks of Stryker.
Products referenced with [®] designation are registered trademarks of Stryker.

PC/SCG-7.5M-0109, 07-SMD-2096
Mkt Lit-128-102105 Rev C

Copyright © 2008 Stryker
Printed in USA

