

# Electric & Manual Beds



Electric

**Choice of full- or half-length siderails in painted finish**

**Centrally-located lock and steer system**

**Optional fixed patient siderail controls**

**Fifth wheel helps improve maneuverability and ease transport**

**Optional removable crank handles for manual bed positioning**

## Designed to Simplify Patient Care


Stryker's Electric and Manual Med/Surg Beds address the demands put on today's nursing professional. Keeping our product design focused on key needs ensures a practical and user-friendly patient platform. Safety features such as low bed height, and convenience features, which include a centrally-located lock and steer system, increase comfort for patients and help relieve strain on caregivers.



Manual



Optional Complete Siderail Coverage



Low Bed Height



Retractable Fifth Wheel

# Electric & Manual Beds

**stryker**<sup>®</sup>



Electric



Manual

## Medical

3800 E. Centre Ave.  
Portage, MI 49002  
t: 269 329 2100 f: 269 329 2311  
toll free: 800 787 9537

[www.stryker.com](http://www.stryker.com)

### Standard Features

- 6" (15 cm) casters
- Centrally-located lock and steer system
- Two drainage bag hooks
- Eight IV pole and traction equipment sockets
- Patient restraint locations
- Lateral and foot-end mattress retainers
- Trend./reverse Trend. controls
- Roller bumpers at head-end of bed
- Half-length siderails
- Painted siderail finish

### Specifications

	Electric	Manual
• <b>Model Number</b>	FL23SE	FL23SM
• <b>Overall Length</b>	88.5" (225 cm)	88.5" (225 cm)
• <b>Overall Width</b>	41.5" (105 cm)	41.5" (105 cm)
• <b>Weight Capacity</b>	500 lbs (226 kg)	500 lbs (226 kg)
• <b>Height Range (to litter top)</b>		
High	30" (76 cm)	29" (74 cm)
Low	13" (33 cm)	15" (38 cm)
• <b>Litter Positioning</b>		
Backrest	0° - 60°	0° - 62°
Knee Gatch	0° - 30°	0° - 32°
Trend./reverse Trend.	±17°	±14°
• <b>Patient Surface</b>	35" x 78" (89 x 198 cm) 35" x 80" (89 x 203 cm)	35" x 78" (89 x 198 cm) 35" x 80" (89 x 203 cm)
• <b>Caster Diameter</b>		
Standard	6" (15 cm)	6" (15 cm)

### Electronics

• <b>Agency Approvals</b>	C-US CSA
• <b>Volts</b>	120 VAC (Domestic)/240 VAC (International)
• <b>Ampere Rating</b>	4A (Domestic)/2.7A (International) 5A (rating for 120 volt outlet option)
• <b>Current leakage</b>	<300 microamperes
• <b>Frequency</b>	50/60 Hz
• <b>Hospital Grade Plug</b>	Yes

*Stryker reserves the right to change specifications without notice.*

### Optional Electric

- PVC siderail protector
- Removable crank handle
- Oxygen bottle holder
- Patient helper
- Monitor tray
- Retractable fifth-wheel steering
- Auto contour
- 120 volt outlet
- CPR release
- Night light
- Nurse call

### Optional Manual

- PVC siderail protector
- Removable crank handle
- Oxygen bottle holder
- Patient helper
- Monitor tray
- Retractable fifth-wheel steering
- Stabilizer Legs
- Permanent head and footboards

Copyright © 2005 Stryker  
PC/SCG-7.5M-0109, 05-SMD-1009  
Mkt Lit-134-11032005 Rev C



**Mixed Sources**  
Product group from well-managed  
forests, controlled sources and  
recycled wood or fibre  
Cert no. SW-COC-1530  
[www.fsc.org](http://www.fsc.org)  
© 1996 Forest Stewardship Council