

Symmetry® Plus

Treatment Recliner



500 lb weight capacity

Unique recline mechanism offers infinite number of recline positions

Posture-correct lumbar and lateral cushioning

Sinuuous wire-spring construction provides premium comfort over time

BackSmart™ Trendelenburg positioning

Wide seat comfortably accommodates larger patients

Comfortable full flat-bed position

Extended, wide arm caps facilitate easy entry and exit

Limited lifetime warranty

Maximum Comfort, Clinical Support

Perfect for patient ambulation, recliners promote recovery by helping patients get out of bed faster. Stryker's Treatment Recliner features an infinite number of recline positions and full flat-bed position.



Full Flat-bed Position



Step-lock Brake



Trendelenburg Position

Symmetry Plus

Treatment Recliner

stryker®



Medical

3800 E. Centre Ave.
Portage, MI 49002 USA
t: 269 329 2100 f: 269 329 2311
toll free: 800 787 9537

www.stryker.com

Standard Features

- Independent ottoman allows feet to be raised without moving the back
- Mid-ottoman cushion provides extra support when feet are elevated
- Chair construction includes metal seat frame and sinuous wire springs
- Highly-resilient polyurethane foam cushions
- Seamless upholstery and removable back cushion enables thorough cleaning
- Pull-out foot tray gives patients added comfort while in transit
- Rigid push bar for easy maneuverability
- Fixed sides provide stability for electromechanical mounts
- Egress assist polyurethane arm caps
- Meets flammability standards of BIFMA F1-1978 - revised 1980 and CAL 117
- Meets ANSI//BIFMA X5.4 2005 per tests conducted in a certified test facility
- Oversized tracking, lockable casters provide agile maneuverability (front rigid, rear swivel)

Specifications

• Style	Mobile Base Part #3500-000-710
• Overall Width	36" (91.4 cm)
• Overall Height	47" (119.4 cm)
• Depth	
Overall	38" (96.5 cm)
Seat	20.5" (52.1 cm)
Recline/Legrest Up	74" (188 cm)
• Weight	165 lbs (74.8 kg)
• Seat Width	26" (66 cm)
• Back Height From Seat	31" (78.7 cm)
• Seat Height From Floor	20.5" (52.1 cm)
• Arm Height From Floor	28" (71.1 cm)

*Stryker reserves the right to change specifications without notice.
Furniture is manufactured exclusively for Stryker Medical.*

Warranty

- Limited lifetime warranty
- Warrantied 500 lbs (226.8 kg) weight capacity

Optional Features

- Fold-down side tray provides prep area for caregivers and serves patients' ancillary needs (permanently mounted to right side, left side or both)
- Rear-mounted IV pole and brackets save valuable floor space and allow for easier transport
- Egress assist wood arm caps
- Hinged push bar for maneuverability
- Four-wheel swivel with locking casters
- Choose from a variety of stains and fabrics
- CAL 133

Copyright © 2007 Stryker
LSS-457-01E
PC/SCG-7.5M-1207, 07-SMD-2148
Mkt Lit-84-122007 Rev E