

Symmetry® Plus

Patient Room Recliners



Closed-Arm Recliner

Open-Arm Recliner

Rocker or Wall-Saver models available

Exclusive recline mechanism offers infinite number of recline positions

Posture-correct lumbar and lateral cushioning

Sinuous wire-spring construction provides premium comfort over time

Limited lifetime warranty

Versatile Recliners for a Healing Environment

Stryker's stylish Patient Room Recliners maximize patient and guest comfort, while providing uncompromised durability. Choose from rocker or Wall-Saver models, with two arm options — open arm or upholstered arm.



Mobile Base



Stationary Base

Symmetry Plus

Patient Room Recliners

stryker[®]



Closed-Arm
Recliner

Open-Arm Recliner

Medical

3800 E. Centre Ave.
Portage, MI 49002 USA
t: 269 329 2100 f: 269 329 2311
toll free: 800 787 9537

www.stryker.com

Standard Features

- Stationary or mobile
- Mobile with 3 inch front rigid/rear swivel locking casters
- Mid-ottoman cushion gives extra support when feet are elevated
- Chair construction includes metal seat frame, sinuous wire springs and 7-ply engineered plywood frame
- Seamless upholstery and removable back cushion enables thorough cleaning
- Highly-resilient polyurethane foam cushions
- Stationary arms and fixed sides provide stability for drainage catheters or electromechanical mounts (Open-Arm only)
- Egress assist polyurethane arm caps
- Meets flammability standards of BIFMA F1-1978 – revised 1980 and CAL 117
- Meets ANSI/BIFMA X5.4 2005 per tests conducted in a certified test facility
- Oversized tracking, lockable casters provide agile maneuverability
- Rear push bar allows for easy patient transport
- Pull-out rolling footrest gives patient added comfort while in transit

Optional Features

- Open arm or upholstered arm
- Egress assist wood arm caps
- Choose from a variety of stains and fabrics
- CAL 133 in selected fabrics
- Fold-away push bar
- 3 inch front/rear swivel locking casters

Specifications

• Style	Rocker Recliner Part #3500-000-750	Wall-Saver Part #3500-000-730
• Overall Width	31" (78.7 cm)	31" (78.7 cm)
• Overall Height	45.5" (115.6 cm)	47.5" (120.7 cm)
• Depth		
Overall	41.5" (105.4 cm)	40.5" (102.9 cm)
Seat	19.5" (49.5 cm)	19.5" (49.5 cm)
Recline/Legrest Up	70" (177.8 cm)	70" (177.8 cm)
• Weight	138.2 lbs (62.7 kg)	130.5 lbs (331.5 kg)
• Seat Width	23" (58.4 cm)	23" (58.4 cm)
• Back Height From Seat	32.5" (82.6 cm)	33.5" (85.1 cm)
• Seat Height From Floor	20.5" (52.1 cm)	20.5" (52.1 cm)
• Arm Height From Floor	27" (68.6 cm)	27" (68.6 cm)

*Stryker reserves the right to change specifications without notice.
Furniture is manufactured exclusively for Stryker Medical.*

Warranty

- Limited lifetime warranty
- Warrantied 350 lbs (158.8 kg) weight capacity

Rocker Recliner

- Rear lock lever allows chair to remain stationary

Wall-Saver

- Reclines away from the wall to save space