

Evacuation Chair

civilian use stair chair

stryker[®]



EMS Equipment

3800 E. Centre Ave.
Portage, MI 49002 USA
t: 269 329 2100 f: 269 329 2311
toll free: 800 784 9537

www.ems.stryker.com

Standard Features

- Innovative Stair-TREAD™ system
- Easy to understand user instructions incorporated on seatback
- Light-weight, rugged aluminum construction
- Extendable head and foot end lift handles
- Molded hand grips
- Positive action locking mechanism
- Compact storage size
- Oversized rear wheels with sealed bearings
- Color-coded controls
- Extendable upper control handle
- Dual wheel locks
- Power washable
- Grease-free maintenance
- Passenger restraint set¹
- 4" (10 cm) front caster wheels
- Head support
- Two-piece molded ABS seat
- In-service training video

Optional Features

- Wall storage bracket
- Wall storage cover

Specifications²

• Model	6253	
• Height	45"	(114 cm)
• Width	20.5"	(52 cm)
• Depth³	28"	(71 cm)
• Folded Depth	11"	(28 cm)
• Weight	34 lbs	(15 kg)
• Maximum Load⁴	500 lbs	(227 kg)

Warranty

- 15 YEAR WARRANTY

¹Passenger restraint set includes lap, chest and ankle restraints.


²Dimensions are measured from the outermost edges of the main frame. Specifications are rounded to the nearest whole number. Conversions are calculated before rounding.

³Depth dimensions are measured with extendable handles retracted.

⁴Maximum load capacity is total weight distributed in accordance to basic human anatomy. Operators must consider the weight of the passenger, equipment and accessories when determining the total weight on the product.

In-service video included with every order.

Stryker reserves the right to change specifications without notice.



Instructions on back of seat