

# M-1®

## roll-in ambulance cot



**stryker®**

### EMS Equipment

3800 E. Centre Ave  
 Portage, MI 49002 USA  
 t: 269 329 2100 f: 269 329 2311  
 toll free: 800 784 4336

www.ems.stryker.com

### Standard Features

- One-person operation
- Roll-in loading with two-stage load release safety mechanism
- Oversized wheels with sealed caster and wheel bearings
- Color-coded controls
- High visibility powder-coated frame
- Light-weight, rugged aluminum construction
- Shock, knee gatch or flat leg positioning
- Head-end steer lock system
- Ergonomic lift-and-grip design
- Seven height positions
- Foot-end wheel locks
- Reflective labeling
- Release lever guards
- Removable litter
- Extendable one-hand release lift handles
- Integrated bumper system
- One-hand release, infinite positioning, lift assist, pneumatic backrest
- One four-point shoulder harness restraint
- Two lap-belt restraints
- Lubrication-free maintenance
- Pressure-washable

### Optional Features

- Premium sealed bolster mattress that is antimicrobial
- Sealed bolster mattress
- Flat mattress
- Head extension with pillow
- Pillow (premium mattress only)
- One-hand release, fold-down siderails
- Heavy-duty, three-stage IV pole (patient right or patient left)
- DIN lock-down pins
- Head-end retractable lift handles (base)
- Foot-end cot retaining loop
- Foot-end cot retaining post
- Removable oxygen bottle holder

### Specifications

<b>• Model Numbers</b>	<b>6100</b>	
<b>• Height<sup>1</sup></b>		
Position 1	14"	(37 cm)
Position 2	24"	(62 cm)
Position 3	30"	(76 cm)
Position 4	35"	(88 cm)
Position 5	38"	(96 cm)
Position 6	39"	(100 cm)
Position 7	13"	(33 cm)
<b>• Length</b>	77"	(197 cm)
<b>• Width</b>	22"	(56 cm)
<b>• Weight<sup>2</sup></b>	74 lbs	(33 kg)
<b>• Wheels</b>		
Diameter	6"	(15 cm)
Width	2"	(5 cm)
<b>• Maximum weight capacity</b>	500 lbs	(228 kg)
<b>• Recommended loading height</b>	27" - 31"	(69 - 79 cm)
<b>• Minimum operators required</b>		
Occupied cot	1	
Unoccupied cot	1	
<b>Recommended fastener system</b>		
Model 6373	center mount, DIN guide	
Model 6376	center mount, straight guide w/o belts	

### Litter Specifications

	Flat Foot		Knee Gatch	
<b>• Length</b>	75"	(190 cm)	75"	(190 cm)
<b>• Width</b>	22"	(56 cm)	22"	(56 cm)
<b>• Height</b>	7"	(18 cm)	7"	(18 cm)
<b>• Weight<sup>2</sup></b>	43.5 lbs	(20 kg)	49 lbs	(22 kg)
<b>• Maximum weight capacity</b>	500 lbs	(228 kg)	500 lbs	(228 kg)
<b>• Handle extension</b>	8"	(20 cm)	8"	(20 cm)
<b>• Articulation</b>				
Backrest	0° - 75°		0° - 75°	
Shock position	N/A		+17°	
Knee gatch	N/A		30°	

### Warranty

- One year parts, labor and travel or two years parts

<sup>1</sup>Height measured from bottom of mattress at seat section to ground level.

<sup>2</sup>Cot is weighed without mattress and restraints.

Stryker reserves the right to change specifications without notice.